
Data Ally™

Chromatography Manager

User's Manual



Contents

Figures	7
License Agreement	11
1.0 Introduction: Data Ally Chromatography Manager	13
2.0 Data Ally Specifications	14
2.1.0 Number of instruments.....	14
2.2.0 Signal Inputs/Acquisition.....	14
2.3.0 Graphics.....	14
2.4.0 Methods	14
2.5.0 Integration/Peak Identification.....	15
2.6.0 Calibration	15
2.7.0 Batch Operations.....	15
2.8.0 Reports.....	16
2.9.0 Files Handling.....	16
3.0 Installing Data Ally Chromatography Manager	17
3.1.0. Contents Listing	17
3.2.0 Installation	17
3.2.1. Install the Data Ally Card	17
3.2.2 Install the Data Ally Software.....	18
3.3.0. Starting Data Ally	20
3.4.0 Installing Data Ally Software Upgrades.....	20
4.0 The Data Ally Desktop: Basics.....	21
4.1.0 Menus and Commands.....	21
4.1.1 Dialog Boxes	22
4.2.0 Data Ally Windows.....	25
4.2.1 Switching Between Windows	26
4.2.2 Making a Window "Active".....	26
4.2.3 Moving Windows and Dialog Boxes	26
4.2.4 Expanding Windows to Maximum Size.....	27
4.2.5 Resizing Windows	27
4.2.6 Closing a Window.....	28
4.2.7 Minimizing a Window	28
4.2.8 Splitting a Window	28
4.3.0 Mouse Pointers and Techniques	29
4.3.1 Pointer Shapes	29
4.3.2 Basic Mouse Techniques	31
4.4.0 Data Ally Tables	31
4.4.1 Moving Around in Tables.....	31
4.4.2 Adding, Moving, Copying, and Changing Rows	33
5.0 Configuring Data Ally.....	35
5.1.0 Configuration Files	35
5.2.0 Performing a Startup Configuration.....	35
5.2.1 Selecting Configured Components.....	37
5.2.2 Configuring Signal Inputs.....	38
5.2.3 Configuring Relay Flags	41
5.2.4 Deleting Components.....	42
5.2.5 Setting Controlled System Instrument Components.....	43
5.2.6 Program Screen Options, Menus, and Defaults	45
5.3.0 Configuration Files	89
5.3.1 Saving Configuration Files	89
5.3.2 Loading Configuration Files	90
5.3.3 Deleting Configuration Files.....	90
6.0 Getting Started	91
6.1.0. Operating Structure and Programs	91

6.2.0 Files and File Manipulations	92
6.3.0 Entering the Data Ally Software.....	93
6.4.0 Using the Default System.....	93
6.4.1 Starting a Run.....	93
6.4.2 Analyzing Run Data Manually	94
6.4.3 Completing the Integration Peak Table	94
6.4.4 Setting Up a Calibration File.....	95
6.4.5 Programming a Simple Sequence.....	97
6.5.0 Exiting the Data Ally System	98
7.0 Chromatogram Display, Review, and Editing: Graph Screen	100
7.1.0 GRAPH Screen Commands.....	101
7.1.1 Application and Program Control Menu Commands	101
7.1.2 Screen-Based Commands	101
7.1.3 Graph Command Menu Bar	104
7.2.0 Working With Chromatogram Files	105
7.2.1 Files Commands	105
7.3.0 Editing Graph Files	108
7.3.1 Editing Commands	108
7.3.2 Integrating Chromatograms Manually	109
7.3.3 Manual Chromatogram Baseline Editing	111
7.3.4 Manual Zone Code Editing	114
7.3.5 Manual Chromatogram Analysis	117
7.3.6 Manual Calibration.....	119
7.3.7 Manual System Suitability Calculations.....	121
7.3.8 Cutting the Current Display Image.....	121
7.4.0 Math Functions.....	122
7.4.1 Adding/Subtracting/Ratioing Chromatograms	122
7.5.0 View Commands	123
7.5.1 View Commands	123
7.6.0 Changing Graph Display Characteristics.....	125
7.6.1 Display Attributes.....	125
7.6.2 Screen Display Characteristics Settings	127
7.6.3 Status Box	131
7.6.4 Monitor Time Bases Function.....	132
7.7.0 Manual Printing of Graph Documents.....	133
7.7.1 Printing Chromatograms.....	133
7.7.2 Selecting Information Included On Chromatogram Printouts	135
7.7.3 Specifying Report Types	136
7.8.0 Running From the Graph Screen	137
8.0 Programming and Executing Runs: The Method Screen.....	139
8.1.0 Method Functions.....	139
8.2.0 The Top Level Method Dialog.....	140
8.2.1 Method File Identification.....	140
8.2.2 Setting Inputs Criteria	140
8.2.3 Setting Run Criteria.....	141
8.2.4 Setting Results Criteria.....	142
8.3.0 Events Table.....	144
8.3.1 Events Table Functions	144
8.3.2 Events Table Column Parameters.....	145
8.3.3 Entering Events Table Parameters.....	148
8.4.0 Manual Operation.....	153
8.4.1 Hold Function.....	153
8.5.0 Working with Method Files.....	154
8.5.1 Loading Methods.....	154
8.5.2 Clearing Methods	155

8.5.3 Saving Methods	156
8.5.4 Deleting Methods.....	156
8.5.5 Printing Methods.....	157
9.0 Integration Screen: Finding and Identifying Peaks	158
9.1.0 Top-Level Integration Screen	158
9.1.1 File Identification.....	159
9.1.2 Integration Algorithm	159
9.1.3 System Suitability	163
9.1.4 Peak Identification	165
9.2.0 Integration Peak Table	165
9.2.1 Peak Table Functions.....	166
9.2.2 Peak Table Components	167
9.2.3 Entering Peak Table Values.....	169
9.3.0 Working with Integration Files	173
9.3.1 Loading Integration Files.....	173
9.3.2 Clearing Integration Peak Table	174
9.3.3 Saving Integration Files	175
9.3.4 Deleting Integration Files.....	175
9.3.5 Printing The Integration File.....	175
10.0 Reviewing and Editing Calibrations: Calibration Screen.....	176
10.1.0 Calibration File Definition	176
10.2.0 The Calibration Display	177
10.2.1 The Reference Peak Table	177
10.2.2 The Calibration Point Table.....	187
10.2.3 The Calibration Plot.....	188
10.3.0 Theory of Calibration.....	190
10.3.1 External Standard Calibration.....	190
10.3.2 Internal Standard Calibration	193
10.4.0 Building, Reviewing and Modifying Calibrations.....	196
10.3.1 Building a New Calibration Table	196
10.3.2 Reviewing the Calibration Point by Point	197
10.3.3 Editing Individual Standard Points	198
10.3.4 Changing the Calibration Curve Fit	199
10.3.5 Changing the Calibration Basis.....	200
10.5.0 Working With Calibration Files.....	200
10.5.1 Loading Calibration Files.....	200
10.5.2 Clearing Calibrations	201
10.5.3 Saving Calibrations.....	202
10.5.4 Deleting Calibrations	203
10.5.5. Printing Calibrations	203
11.0 Running and Reprocessing Batches: The Sequence Screen.....	204
11.1.0 Sequence Screen Functions.....	204
11.2.0 Top-Level Sequence	204
11.2.1 Operating Modes.....	205
11.3.0 Sequence Table.....	206
11.3.1 Sequence Table Functions	206
11.3.2 Sequence Table Columns.....	208
11.3.3 Fill Functions	210
11.4.0 Running a Sequence.....	213
11.4.1 Modifying a Sequence While Executing.....	214
11.5.0 Working with Sequence Files	214
11.5.1 Loading Sequence Files	215
11.5.2 Clearing Sequence Files.....	215
11.5.3 Saving Sequence Files	216
11.5.4 Deleting Sequence Files.....	217

11.5.5	Printing Sequence Files	217
11.5.6	Automatically Starting Sequence Files	217
12.0	Utilities	218
12.1.0	Utilities Files Commands	218
12.1.1	Resetting Configuration	218
12.1.2	Translate Functions	219
12.2.0	Configure	222
12.2.1	Exiting the Data Ally System	222
13.0	Data Ally Help Functions	223
13.1.0	Using the HELP System	223
13.2.0	Marking Help Topics	225
13.3.0	Annotating Help Topics	226
14.0	Reporting with Microsoft Excel	228
14.1.0	Using Excel to Review and Edit Report Templates	228
14.1.1	Loading Excel	228
14.1.2	Opening An Existing Report Template File	228
14.1.3	Elements of a Report Template File	229
14.2.0	Excel Basics	230
14.2.1	Changing Column Widths and Row Heights	231
14.2.2	Changing Contents of a Single Cell	231
14.2.3	Moving Data From One Cell To Another	231
14.2.4	Using The Name Box	231
14.2.5	Using the Define Name Dialog Box	232
14.2.6	Inspecting Defined Names Used in Report Templates	232
14.3.0	Modifying Existing Report Template Styles	233
14.3.1	Modifying A Template Paper Size	233
14.3.2	Placing And Formatting Report Headers	234
14.3.3	Placing And Formatting Page Headers Using Print_Titles	234
14.3.4	Placing and Formatting Data/Result Tables	234
14.4.0	Data Ally DDE Names For Use On Reports	235
14.4.1	Current Context Sensitive DDE Names	236
14.4.2	Calibration DDE Names	236
14.4.3	Integration DDE Names	238
14.4.4	Method DDE Names	238
14.4.5	PDF DDE Names	240
14.4.6	Sequence DDE Names	241
15	Hardware Specific Interface Information	243
15.1	Autosamplers	243
15.1.1	Marathon Autosampler by Spark-Holland	243
15.1.2	Parameters which can be controlled on all Models	244
15.1.3	Parameters which apply to XT and Basic-Plus Models Only	244
15.1.4	Connections Between Data Ally and Marathon	244
15.1.5	Placing the Marathon in Remote Control Mode	244
15.1.6	Autosampler Sequence Control Issues	244
15.1.7	Injection Modes	245
15.1.8	Dispenser Controlled Injection Mode	245
15.1.9	Pressure Controlled Injection	245
15.1.9	Marathon Error Codes	245
15.2	Kontron 460 Autosampler	245
15.2.1	Setting up the Kontron for Remote Control	246
15.2.2	Connections Between Data Ally and Kontron 460	246
15.2.3	Working with Kontron 460 Series Autosamplers	246
15.2.4	Kontron Autosampler Sequence Control Issues	247
15.3	Detectors	247
15.3.1	General Instructions for connection of detectors to Data Ally	247

15.3.2 UV and Fluorescence Detector Control	248
15.3.3 Controlling the Antec Decade Electrochemical Detector.....	248
15.3.4 Controlling the ESA Electrochemical Detector	249
15.4 Photo Diode Array Detectors.....	253
15.4.1 General PDA Formulas used for Spectra, Similarity, and Purity Plots	254
15.5 LabAlliance Model 506 or Linear 206 UV/Vis Control Installation Instructions.....	255
16 Data Ally System Security & Log On Feature	257
16.1.0 System Log On.....	257
16.2.0 Adding Users	257
16.2.2 Setting Users Role.....	259
Access Roles.....	260
17 Installation and Routine Troubleshooting.....	262
17.1.0 Serial Mismatch Errors	262
17.1.2 Other Startup Error Messages	263
17.1.3 Errors During Data Ally Operation.....	263

Figures

Figure 3.2 Data Ally Application Group and Icon in Windows Program Manager	19
Figure 3.8 Exiting Data Ally In Utilities Screen.....	19
Figure 4.1 Command Menu Bar for Graph Screen.....	21
Figure 4.2 Pull-Down Display Menu for Graph Screen	21
Figure 4.3 Files/Load Dialog Box for Graph Screen.....	22
Figure 4.4 Control Application Menu for Graph Screen Window	24
Figure 4.5 Maximize and Minimize Buttons for Time Base (top) and Program Windows	27
Figure 4.6 Graph Screen Program Control Menu.....	29
Figure 5.1. Selecting Configuration Mode From Utilities Screen.	36
Figure 5.2 Configuration Window.....	36
Figure 5.3 Component Set Dialog Box	37
Figure 5.4 Time Base Window System Setup Dialog	37
Figure 5.5 Entering Time Base Description	38
Figure 5.6 Analog Input Setup Dialog Box.....	38
Figure 5.7 Entering a Description For Detector Input.....	39
Figure 5.8 Add Component Dialog Box.....	40
Figure 5.9 Component Setup Dialog After Adding Second Analog Detector Input.....	40
Figure 5.10 Adding Relay Flags.....	41
Figure 5.11 Input Flags Setup Dialog.....	42
Figure 5.12 Removing a Defined Component	43
Figure 5.13 Selecting Desired Detector Type	43
Figure 5.14 ESA 5200 Detector Configuration Setup Dialog	44
Figure 5.15 Configuration Graph Window.....	45
Figure 5.16 Loading a Chromatogram In Configuration Graph Screen	46
Figure 5.17 Graph Edit Default Submenu	47
Figure 5.18 Analyze Default Dialog.....	47
Figure 5.19 Selecting Default Report From Existing Template Files.....	48
Figure 5.20 Calibrate Default Dialog	48
Figure 5.21 Integrate Default Dialog.....	49
Figure 5.22 Configuration Graph Display Submenu	50
Figure 5.23 Setting Default Display Attributes	50
Figure 5.24 Setting Default Display Colors	51
The Component Set dialog box (Figure 5.3) will appear.	51
Figure 5.26 Offset Default Dialog.....	52
Figure 5.27 Pumps Offset Display Dialog.....	52
Figure 5.28 Display Scaling Default Dialog.....	53
Figure 5.29 Pump Scaling Default Dialog.....	54
Figure 5.30 Status Box Default Dialog	54
Figure 5.31 Graph View Default Submenu	55
Figure 5.32 Zoom Center Default Submenu.....	56
Figure 5.33 Zoom Scale Factors Default Dialog	56
Figure 5.34 Print Command Default Submenu.....	57
Figure 5.35 Printer Setup Default Dialog.....	57
Figure 5.36 Printout Info Default Dialog	58
Figure 5.37 Chromatogram Report Default Dialog.....	58
Figure 5.38 Graph Screen Button Configuration Pop-up Submenu	59
Figure 5.39 Set Button Color Dialog.....	60
Figure 5.40 Set Button Title Dialog	60
Figure 5.41 Method Configuration Dialog	61
Figure 5.41 Method Configuration Dialog	62
Figure 5.43 Method Reports Default Dialog.....	65
Figure 5.44 Method Files Configuration Submenu	66

Figure 5.45	Default Method Description Dialog	66
Figure 5.46	Assign Default Method Dialog	67
Figure 5.47	Method Events Table Configuration Window	68
Figure 5.48	Table Button Edit Pop-up	68
Figure 5.49	Configure Events Table Columns Dialog	69
Figure 5.50	Events Table With Inputs Column	69
Figure 5.51	Set Button Title Dialog	70
Figure 5.52	Button Title Color Select Dialog	71
Figure 5.53	Integration Configuration Default Window	72
Figure 5.54	Integration Default File Submenu	73
Figure 5.55	Integration Default Description Dialog	73
Figure 5.56	Integration Default Assignment Dialog	74
Figure 5.57	Configuration Peak Table Window	74
Figure 5.58	Configure Peak Table Columns Dialog	75
Figure 5.59	Configuration Calibration Window	76
Figure 5.60	Calibration Default File Submenu	77
Figure 5.61	Assign Default Calibration File Dialog	77
Figure 5.62	Calibration Tables Configuration Window	78
Figure 5.63	Calibration Default File Submenu	78
Figure 5.64	Calibration Default Setup Dialog	79
Figure 5.65	Configure Calibration Reference Table Columns Dialog	80
Figure 5.66	Column Selector Dialog for Calibration Data Table	81
Figure 5.67	Configuration Sequence Window	82
Figure 5.68	Sequence File Default Submenu	83
Figure 5.69	Sequence Default File Description Dialog	84
Figure 5.70	Assigning Sequence Default File	84
Figure 5.71	Configure Sequence Table Columns Dialog	85
Figure 5.72	Utilities/Files Configuration Submenu	87
Figure 5.73	Save As Configuration Dialog	87
Figure 5.74	Load Configuration Dialog	88
Figure 5.75	Exiting Configuration Dialog	88
Figure 7.1	Graph Screen	100
Figure 7.2	Graph Application Control Menu	101
Figure 7.7	Graph Files Submenu	105
Figure 7.8	Files/Clear Dialog	106
Figure 7.9	Files/Load Dialog	107
Figure 7.10	Files/Save As Dialog	108
Figure 7.11	File Overwrite Trap Dialog	108
Figure 7.12	Graph Edit Command Submenu	109
Figure 7.13	Integration Dialog	110
Figure 7.14	Baseline Editing Submenu Commands	111
Figure 7.15	Baseline Edit Pop-up Selector	112
Figure 7.16	Placing a Drop Line With Pointer	113
Figure 7.17	Erase Baseline Code Commands	114
Figure 7.18	Define Zone Codes Submenu	115
Figure 7.19	Zone Codes Editor Pop-up	116
Figure 7.20	Zone Codes Placed Along Graph Time Axis	116
Figure 7.21	Manual Analyze Dialog	118
Figure 7.22	Example of Typical Area PerCent Report via Excel	119
Figure 7.23	Manual Calibration Dialog	119
Figure 7.24	Example of Calibration Report	120
Figure 7.25	Cut Image Command	121
Figure 7.26	Math Functions Menu	122
Figure 7.27	View Submenu Commands	123
Figure 7.28	Zoom Box Enabled in Graph Window	124

Figure 7.29	Zoom Center Commands	124
Figure 7.30	Zoom Center Scale Factors Dialog	125
Figure 7.31	Graph/Display Submenu Commands	125
Figure 7.32	Graph Display/Attributes Dialog	126
Figure 7.33	Display/Colors Dialog.....	127
Figure 7.34	Display/Offset Dialog	128
Figure 7.35	Display/Scaling Dialog	129
Figure 7.36	Pumps/Scaling Dialog.....	129
Figure 7.37	Display Submenu with HPLC Pumps	130
Figure 7.38	Setting Pump Trace Plotting	130
Figure 7.39	Display/Status Dialog.....	131
Figure 7.40	Status Box Enabled in Graph Window	132
Figure 7.41	Monitor Time Bases Screen Display.....	132
Figure 7.42	Print Commands Submenu.....	133
Figure 7.43	Chromatograms Print Dialog	134
Figure 7.44	Printer Setup Dialog.....	134
Figure 7.45	Print Options Dialog	135
Figure 7.46	Printout Info Dialog	136
Figure 7.47	Print Report Dialog.....	136
Figure 8.1	Top-Level Method Screen	140
Figure 8.2	Detectors Setup Fields	141
Figure 8.3	Sample Information Dialog.....	142
Figure 8.4	Result Files Definition Fields.....	142
Figure 8.5	Report Set Dialog.....	143
Figure 8.6	Method Events Table	144
Figure 8.7	AutoFill Dialog for Events Table Column.....	149
Figure 8.8	Method Files Commands Submenu.....	154
Figure 8.9	Method Files/Load Dialog	155
Figure 8.10	Clearing Method Events Table.....	156
Figure 8.11	Method File/SaveAs Dialog.....	156
Figure 8.12	Reporting Status Window	157
Figure 9.1	Top Level Integration Screen.....	158
Figure 9.2	Integration Algorithm Parameter Fields.....	159
Figure 9.3	Use of Slope Sensitivity Parameter.....	160
Figure 9.4	Use of Baseline Drift Parameter	160
Figure 9.5	Integration Peak Table Window.....	166
Figure 9.6	Automatically Filling Retention Times	170
Figure 9.7	Integration Files Command Submenu	173
Figure 9.8	Integration Files/Load Dialog	174
Figure 9.9	Clearing Peak Table.....	174
Figure 9.10	Integration Files/Save As Dialog	175
Figure 10.1	Calibration Screen.....	177
Figure 10.2	Reference Peak Table Window.....	178
Figure 10.3	Calibration Files Submenu	180
Figure 10.4	Calibration Setup Dialog.....	180
Figure 10.5	AutoFilling Curve Fit Type.....	182
Figure 10.6	Changing Basis for Existing Peak in Reference Table.....	183
Figure 10.6	Calibration Table	187
Figure 10.8	Calibration Plot Window	189
Figure 10.9	Using Manual Calibrate Function to Process a Standard Chromatogram	197
Figure 10.10	Automatic Indexing of Calibration Table After Plot Point Selection.....	198
Figure 10.11	Failing a Calibration Data Point.....	199
Figure 10.12	Calibration Files/Load Dialog.....	201
Figure 10.13	Clear Submenu for Calibration Files.....	201
Figure 10.14	Calibration Files/SaveAs Dialog.....	202

Figure 10.15 Set Calibration Description Dialog.....	202
Figure 11.1 Sequence Screen	204
Figure 11.2 Auto Filling the METH Column.....	211
Figure 11.3 Split-Screen Sequence Table with Graph	214
Figure 11.4 Sequence Files Menu.....	215
Figure 11.5 Sequence Files/Load Dialog.....	215
Figure 11.6 Sequence Clear Function	216
Figure 11.7 Sequence Files/Save As Dialog	216
Figure 12.1 UTILITIES Screen	218
Figure 12.2 Utilities Files Submenu.....	218
Figure 12.3 Load Configuration Dialog.....	219
Figure 12.4 Translate Submenu Commands	219
Figure 12.5 Axxiom 700 Series File Conversion Dialog	220
Figure 12.6 File Export Dialog	221
Figure 12.7 Export/Files Dialog	221
Figure 13.1 Data Ally Help Window	223
Figure 13.2 Help Contents Index Window.....	224
Figure 3.3 Help Search Dialog.....	224
Figure 13.4 Help History Window.....	225
Figure 13.5 Bookmark Define Dialog.....	226
Figure 13.6 Help Annotation Dialog.....	226
Figure 14.1 Excel Files/Open Dialog Box	229
Figure 14.2 Examining an Excel Worksheet Report Template	229
figure 15.1 Marathon Autosampler Setup Dialog	243
Figure 15.2 Typical Sequence Table For Autosampler Operation.....	243
figure 15.3 Kontron Autosampler Setup Dialog	246
figure 15.4 Typical Sequence Table For Autosampler Operation.....	247
Figure 15.5 Typical Data Ally PDA Graph Screen.....	253
Figure 15.6 Integration Peak Table with Peak Spectra Column.....	255
Figure 15.7 Top Level Method Screen for PDA with four individual channels available	255
Figure 1 System Security Startup Greeting Dialog	257
Figure 2 System Security menu items Edit Users & Logon	257
Figure 3 Set Users Dialog	258
Figure 4 Edit Users Dialog	258
Figure 5 User Access Rights Dialog	259
Figure 6 Set Access Role Dialog	259
Figure 17.9 Communications Error Information Message	264
Figure 17.10 Removing Driver Lines in SYSTEM.INI File.....	265

License Agreement

You have purchased Data Ally™ Chromatography Manager and agree to the following terms.

The Data Ally software is owned by SSI Distribution, Inc. and is protected by United States copyright laws and international treaty provisions. You are permitted to make backup copies of the software for your own use. You may not copy the written materials accompanying the software.

This License Agreement is your proof of license to exercise the rights granted herein and must be retained by you. You may not rent or lease the Data Ally software, but you may transfer your rights under this Agreement on a permanent basis provided you transfer this License Agreement, the software, its accompanying hardware, and all accompanying written materials and retain no copies, and the recipient agrees to the terms of this Agreement. You may not reverse engineer, decompile, or disassemble the software, or modify the program code in any way, without prior written consent from SSI Distribution. Any transfer of the software must include the most recent update and all prior versions.

This License Agreement shall be construed, interpreted, and governed by the laws of the Commonwealth of Pennsylvania. If any provision of this Agreement is found void or unenforceable, it will not affect the validity of the balance of the Agreement, which will remain valid and enforceable according to its terms. If any remedy provided is determined to have failed its essential purpose, all limitations of liability and exclusions of damages set forth in the Agreement shall remain in full force and effect. This Agreement may only be modified in writing with consent of an authorized officer of SSI Distribution. Use, duplication, or disclosure by the United States Government of the software and documentation in this product shall be subject to the restricted rights applicable to commercial computer software (under DARS 52.227-7013).

All rights not specifically granted in this Agreement are reserved by SSI Distribution.

1.0 Introduction: Data Ally Chromatography Manager

Data Ally™ Chromatography Manager is a system of hardware and software for IBM-compatible personal computers which acquires, manages, and analyzes signal information from any type of chromatography instrumentation. It can also control the operation of those instruments.

Data Ally can accommodate any types of chromatography instruments (HPLC, GC, IC, SFC, CZE, or other techniques). The instrument is linked to the computer via the Data Ally card installed in one of the computer's bus slots.

The Data Ally software programs Methods for acquiring chromatography signal data and for controlling the desired instruments. The card at the start of each run sends all required commands to the instrumentation to execute the Method. At the same time, the card captures incoming analog signals from detectors through its high performance analog to digital converters, or acquires digitized signals directly. Each captured and digitized signal is then saved in real time to the Data Ally computer's hard disk drive, and also displayed, if desired, for monitoring purposes.

At the end of each run, the chromatogram signal data from each detector is saved as a unique DOS file on the computer's hard disk. These files can be recalled, displayed, and edited at any time, individually or in batches. Data Ally provides a complete set of analytical and editing tools for integration, baseline modification, comparisons, and graphical manipulations of raw data files. The Data Ally Calibration Module allows any types of calibrations to be performed for quantitative methods. Reports can be created using Microsoft Excel, with transparent communication via Direct Data Exchange (DDE). Data Ally reports can include any information known to the system, including results, comments, text, graphics, and other file-specific data.

The Data Ally software requires and operates within the Microsoft Windows 3.1 or Microsoft Windows 95 environment. It provides many graphical tools and capabilities for programming, editing files, and creating a fully customized user interface, consistent with other Windows applications. It also supports Windows-based multitasking of DOS or Windows applications in 386-Enhanced Mode simultaneously with all Data Ally operations.

Setting Up Your Computer For Data Ally


It is suggested that you perform the following procedures prior to installing the Data Ally hardware PCB in the computer. This will avoid any hardware conflicts that may prevent the computer from booting up. These conflicts should be eliminated by the procedures outlined in this document.

Reserving Interrupt Resources in Windows 95 and the computer's CMOS Setup Utility

In order to for your Data Ally Hardware to function properly, you must inform your computer's BIOS and Windows 95 operating system that certain resources are required by Data Ally. This is a two step process. The first step is to reserve the IRQs 10 & 11 in Windows 95. The second step requires you to reserve IRQs (Interrupts) 10 & 11 in the CMOS Setup Utility. The instructions that follow are intended to be a general guide as to how to reserve these resources. CMOS Setup Utilities vary somewhat in how the user must proceed in order to reserve resources so the guide lines provided here should may not be exact but they should be close.


Reserving IRQs in Windows 95

Go into control panel (click on Start/Settings/Control Panel) and click on the SYSTEM icon and then settings



When you see the above dialog box, click on Properties Button

Make sure the IRQ radio button is selected and you will see the listing for those that exist.



Next Click on the Reserve Resources Tab at the top

Click on the Add button and type in “10” with out the quotes
Repeat and Add “11”

Then click OK and close up all the dialogs.

Re-boot the computer. When you re-boot, Windows 95 will reassign any devices that are using the reserved IRQs thus freeing them up for use with Our ISA Card.

Reserving IRQs for ISA Legacy Cards in the CMOS BIOS Setup Utility

1. You must enter the Setup Utility during the boot up process when your computer is first powered up or reset. Usually this is done by pressing the Delete Key quickly when the message appears on the screen. Some computers may require a different Key press in order to enter the Setup program. Most all will post a message on the screen during boot up informing you as to what key press will invoke the Setup Utility.
2. Once you have entered the Setup Utility, you usually must proceed to the Advanced Setup Menu and then locate the Reserve Resources or PnP PCI Setup menu.
3. When you have found the PnP Setup menu, you may have to select Manual Setup in order to see the choices for reserving you ISA IRQs.
4. Once you have found the ISA Legacy IRQ menu, proceed to select and reserve IRQs 10 and 11.
5. Finally you will want to exit and Save your changes. This usually is done by pressing F10 key but there will be instructions displayed on the screen as to how to do it. Save and Exit the BIOS Setup and then proceed to the Windows 95 instructions shown below.


Checking on I/O Port Address Conflicts

Data Ally uses only 2 locations in the PC I/O address space. Those locations are **300hex & 301hex** . These two locations are not normally used by other devices and are designated for the “prototype card” by IBM.

Some devices that may be using these address will show up in the Windows 95 System Resources list but any devices installed via MS-DOS configuration or autoexec.bat files are not likely to show up. You may examine the autoexec.bat and config.sys files and look for device driver lines that reference I/O port addresses 300hex and 301hex. If you find any you will have to consult the manuals for those devices and relocate them to a non-conflicting I/O address.

To examine the I/O port usage in Windows 95..

1. Go to the Control Panel (click on Start/Settings/Control Panel) and double click the System Icon.
2. Next click on the Device Manager Tab
3. Click on the Properties Button
4. Click on the input/output radio button. You should see a dialog box like the one shown below...
5. You may scroll the list up and down to see if any devices are shown as using 300hex or 301hex.



Important Note Regarding Windows 95^a Installation Error #70

When installing The Data System software onto a Windows 95[®] you must run the Setup program on Disk #1 after starting Windows. It is important that you not run any previous version of the Data System software after re-booting the computer and before installing the new version of software.. If this procedure is not followed an **“Error 70”** message will be encountered when the installation program attempts to install the file named PYRXFR32.EXE.

The reason for the error is whenever you run the Data System software, Windows 95 determines that Excel owns the current copy of PYRXFR32.EXE and will not allow it to be replaced until the system is re-booted and prior to the Data System software being executed.

2.0 Data Ally Specifications

2.1.0 Number of instruments

Available: One Active/Edit Pair

Analog Inputs per : 2

Digital Inputs per : 2

Contact Closure Relays Out per : 4

Contact Closure Relays In per : 4

Serial I/O Ports per : 4

Analog Pump I/O Ports per : 2

2.2.0 Signal Inputs/Acquisition

Analog: Effective Sampling Rate per : 0.0001 Hz - 400 Hz

Maximum Input Range: 1 V, 10 V

Signal-to-Noise at 1 Hz: 25 bits

Signal-to-Noise at 200 Hz: 20 bits

Digital: Effective Sampling Rate per Input:

Maximum Input Range: 1 V, 10 V

Maximum Run Time: Based on available disk memory

Raw data file size: 4 - 10 bytes/point (depending upon file configuration)

2.3.0 Graphics

Real-time Displays: Up to four signal inputs per plus up to 12 saved reference files
Hardware control profiles
Configurable status information box
Split or multi-pane displays
Monitor function for display of any combination
Autosizing or fixed-scale scrolling
Unlimited zoom, rescale, and offset in real time

Editing/Comparison: Up to 16 files at once
Reposition, stack swap functions
Switchable colors, attributes, filling on display
Baseline and zone code markers
Peak names and numbers

2.4.0 Methods

Number of Method programs: Unlimited, any path combination, with one default

Maximum number of timed event steps per Method: Unlimited

Maximum Number of zone codes per Method: Unlimited

Automatable Method Functions: File saves for all inputs
Integration for all inputs using unique programs
Calibration for all inputs using unique programs
Report Generation and saving

2.5.0 Integration/Peak Identification

Number of Integration Programs: Unlimited, any path combination, with one default

Maximum Number Peaks in Peak Table: Unlimited

Automatable Integration Functions: Slope Sensitivity
Slope Calculation Interval
Baseline Drift Sensitivity
Area Threshold
Height Threshold
Peak Width Threshold

Automatable Peak Finding Functions: Absolute Retentions
Absolute Retention Windows

2.6.0 Calibration

Number of Calibration Files: Unlimited, any path combination, with one default

Maximum Number of Entries in Calibration File: Unlimited

Maximum Number of Levels per Calibration: 16

Calibration Type Options: External Standard
Internal Standard (up to 5 different IS)
Area or Height Basis for any peak

Curve Fitting Types Available: Linear
Linear Force Through Origin
Quadratic
Quadratic Force Through Origin
Interpolated
Interpolate Extrapolate Through Origin

Other Calibration Options: Manual Pass/Fail of Individual Calibrations
Omit/Replace Entire Level

2.7.0 Batch Operations

Number of Sequence Programs: Unlimited, any path combination, with one default

Maximum Number of Injections per Sequence: Unlimited

Processing Modes: Acquire, Reprocess and Summary

2.8.0 Reports

Number of Reports per Method, each Input: 1 to 4

Maximum Number of Reports in Template Library: Unlimited

Report Generation Modes: via interface with Microsoft Excel
 Manual Chromatogram Graphics Print functions

2.9.0 Files Handling

System File Types: Chromatogram - *.PDF
 Method - *.MET
 Integration - *.INT
 Calibration - *.CAL
 Sequence - *.SEQ
 Report - *.RPT/*.XLS
 Excel Report Template - *.XLS
 Excel Report Macro - *.XLM

3.1.0. CONTENTS LISTING

Data Ally consists of the following:

- Data Ally Software CD-ROM Disk
- Data Ally Interface P1 Module
- Detector Input Cables (mini-phono plugs) - 2 cables, 2 meter length
- Relay Flag Cables (Twisted Pair) - 2 cables, 2 meter length
- User's Manual

Before beginning any installation procedure, verify that all the components above are present and undamaged.

The computer requirement for installing and operating Data Ally is as follows:


- Pentium processor or equivalent, 400 MHz or faster
- At least 64 MB of system RAM memory
- At least 50 MB free hard disk space recommended
- One Free Serial Port (16550 Uart)
- CD-ROM Drive
- Microsoft Windows 95, 98 or NT 4.0
- Microsoft Excel Version 7.0 or latter (required for report generation)

Data Ally is compatible with any printers or see your Windows documentation for a list plotters for which Windows of supported devices. drivers are available. Please

3.2.0 INSTALLATION

3.2.1. INSTALL THE DATA ALLY HARDWARE

.



3.2.2 INSTALL THE DATA ALLY SOFTWARE.

Before beginning the software installation, be certain that the Data Ally that all PCPs are connected to the Data Ally and powered ON.

Reboot the PC and start Microsoft Windows. Insert the Data Ally CD-ROM in the CD-ROM drive. The Setup program will load and appear onscreen. Follow the instructions given to complete the installation process. If, after you begin installation, the Install program determines that there is not sufficient free hard disk space to install Data Ally, a message box will appear alerting you to this condition and requesting that disk space be freed before re-attempting the installation. Follow all instructions given on screen while using Setup. The Setup program will add a new "Data Ally" application icon to the Windows Program Manager, which can be used to start the Data Ally software directly from Windows.




Figure 3.2 Data Ally Application Group and Icon in Windows Program Manager

NOTE: MICROSOFT EXCEL, V. 7.0 OR OFFICE 2000, MUST BE INSTALLED

in order for reports to be generated.

➡ To exit Data Ally, click on the UTIL left-screen SoftButton to see the Utilities screen, then click the left mouse button on Quit on the top right command line. The Quit submenu will appear. If you wish to exit Data Ally and remain in Windows, move the highlight to Exit Application and click. To abort the quit, click on the CANCEL button. To quit, click on QUIT.

Figure 3.8 Exiting Data Ally In Utilities Screen

When the desired entries appear in the fields, and you have checked them carefully, click on OK. A prompt will appear to confirm that you want to change the current settings, and to remind you that you must remember these new settings in order to change anything in the Configure submenu later on. Click on OK if you are certain you wish to finalize the changes and that the changes shown are correct. If not, click on CANCEL to return to the previous dialog box for further changes.

After completing your installation, secure the Installation disks in a safe place where they can be retrieved if the system must be reloaded in the future. Please complete the user registration card supplied with the disks and return it to Your Distributor. Completing this card will ensure that you receive all future technical bulletins and all software updates.

3.3.0. STARTING DATA ALLY

Depending upon how you installed Data Ally, you can start Data Ally in any of the following ways:

Turn the PC power off and then on again. start Data Ally automatically, the Data Ally the "Log On" screen.

If you have installed AUTOEXEC.BAT and WIN.INI files to software will be initialized. After a short delay you will see

If the PC is powered on and you see the DOS prompt (C>), either:

- Start Windows in the normal manner.
- Click on the Data Ally Icon. After a short delay you will see the "Log On" screen.
- Log on to the system.

When the Data Ally Graph screen appears, the system is ready to use.

Next configure the system to your exact requirements (Section 6), create programs immediately using the factory-supplied configurations, or run Methods. A set of "default" files for each program type are provided in the software with the names DEFAULT.MET, DEFAULT.INT, DEFAULT.CAL, and DEFAULT.SEQ, and are ready to be modified or run as desired.

Note: To using the take full advantage of Data instructions in Section 5. Ally and your system you can custom configure the software

3.4.0 INSTALLING DATA ALLY SOFTWARE UPGRADES

Periodically, upgraded versions of Data Ally software will be improved performance and features.

When re-installing an upgrade for Data Ally, you should first ensure that all existing Method, Integration, Calibration, Sequence, and Configuration files are backed up on floppy disks, tape, optical media, or on the hard drive of another PC. Although re-installation of the same or a newer version will not overwrite any raw data, program, or report files, it is always good practice to maintain separate copies of such files.

➡ *To re- install a new version of Data Ally, Follow the instructions for a first time installation.*

The new software version will be installed into the C:\Data Ally or other designated directory on the hard disk in the same manner as the original installation. When the installation is complete, you must exit and restart Windows as advised by the completion message screen in order to load the new Data Ally version.

In some cases, you may be able to use existing configuration and/or program files - documentation provided with the upgrade software package or in the Setup program message screens will advise you as to the limitations for each type of upgrade.

4.0 The Data Ally Desktop: Basics

As a Microsoft Windows application, Data Ally Chromatography Manager employs many of the basic Windows functions, concepts, and commands found in other Windows software programs. Moreover, Data Ally conforms completely to all Windows rules, so that it can operate properly in the Windows environment alongside other running applications.

This section reviews basics: how to use Data Ally menus and commands, how to work with various program "windows" and dialog boxes, how to move and resize objects on screen, and how to use the mouse.

4.1.0 Menus and Commands

To tell Data Ally what to do, you can choose various types of commands from menus or by using "buttons."

You will use three types of command "menus:" the "Command Menu Bar" commands listed horizontally along the top line of each screen; "pull-down" menus which are usually hidden from view but are made accessible by selecting the "root command" item from a Menu Bar; and the vertical row of "SoftButtons" at the left edge of all channel display screens..




Figure 4.1 Command Menu Bar for Graph Screen

Command items found in top-of-screen Command Menu bars (Figure 4.1) generally perform "global" functions related to the type of screen being displayed. Pull-down menus or submenus are associated with each Command Menu Bar item, are displayed by selecting a Command Menu item (click on it with the left mouse button) and usually contain a set of related functions (Figure 4.2).




Figure 4.2 Pull-Down Display Menu for Graph Screen

The left-screen SoftButtons, which are always visible and accessible from any part of the software, are used to switch between "time bases" or between "program windows" within one time base.

In some pull-down menus, some commands may be "grayed out" at certain times and will not be useable. If a command is "grayed out", it is either inappropriate at the present time, or its function has been defeated or disabled by some part of the software. It is also possible that certain Command Menu items, or even SoftButtons, will not be shown at all in channels in which those commands have been "hidden" by the active Configuration.

Some command items in pull-down menus are followed by a "..." or ellipsis, which indicates that a dialog box will be displayed if that command is selected. The options you choose in the dialog box will determine how that command is executed. Command items may be preceded or followed by a "check mark" which indicates that the particular command is active or enabled.

All Data Ally commands can be chosen directly from their menus by using the mouse pointer or by using a "shortcut" key assignment. Where a set of keystrokes can be used to run a command, the command will show a notation listing the appropriate keys to its right. Not all Data Ally commands can be executed via keystrokes.

➤ *To choose a command from a menu*, point to the command item name with the pointer cursor and click the left mouse button. From the keyboard, you can also press the ALT key to activate the menu, and then press the underlined letter for the menu item you wish to select.

➤ *To cancel a menu without using any command*, click the left mouse button anywhere on screen outside the menu, or press the <ESC> key.

4.1.1 Dialog Boxes

Data Ally utilizes a large number of different dialog boxes to help you choose how commands will be carried out. Each dialog box is a special type of "window" which appears in a predetermined position on screen, and remains displayed as long as you are working with it.

An example of a dialog box is the Files/Load dialog from the Data Ally Graph Screen (Figure 4.3).




Figure 4.3 Files/Load Dialog Box for Graph Screen

Dialog boxes can contain several types of "objects:"

List Boxes: List Boxes are areas which show a "directory list" of alternative choices for some entry - a list box for files appears in the upper left corner of Figure 4.3. List Boxes usually include a vertical "scroll bar" which allows you to move quickly upward or downward through the listing. Often, a list box will be used to "highlight" or select one or more items from the listing, using the mouse. If the listed items fill up all the empty lines inside a List Box, there may be more items which are "hidden" and which can be viewed by scrolling the box down with the lower scroll bar arrow. If one or more items in the box must be highlighted to make a selection, you can click the left mouse button when pointing at a desired item to "toggle" the highlight on or off for that item.

Entry Fields: Entry fields are intended for entering (typing) values for various parameters. An entry field for filename to be loaded appears above the list box in top left of Figure 4.3. To use an entry field box, click the left mouse button with the pointer cursor anywhere inside the box. The flashing "vertical-bar" symbol will appear inside the box. Now type in the desired entry. You can use the Delete key to delete the entire entry in the box, or use the Backspace key to backspace out incorrect characters entered. If you double-click inside the entry box, the entire contents of the box will become highlighted - typing any characters now replaces the old contents. Note that you can also position the "I-bar" entry point at any character within the field for editing. Any entered characters, or deletions, will appear immediately after the entry point. It is also possible to move the entry point from the keyboard within a field, one character at a time, by pressing the left or right arrow keys. As many characters as are necessary can be typed into any entry field, regardless of the physical size of the entry box.

All Data Ally entry fields, as well as all table cells, can be made to show either a listing of all possible entry values for each particular field or an instructional message regarding the type of entry that can be made. To see the "menu listing" for any entry field, click the cursor anywhere inside that field with the left mouse button, and then click the RIGHT mouse button - a "pop-up" menu will appear listing all selections possible for that field. To select one of the entry options, highlight that option with the pointer and click on it with the left button; this will cause the selected value to appear in the field, replacing any previous value. For fields in which filenames must be entered, clicking the right button in this manner will display the normal Files/Load dialog for the type of file required for that field, which can then be used to select a file in the normal manner.

In Data Ally dialog box entry fields, it is usually necessary to move the entry point out of a field into which an entry or modification has just been made, in order for that entry or modification to be "recognized" by the system. As a rule, it is advisable to always click the mouse pointer inside at least one other entry field in each dialog after making or editing entries, to assure that all the new entries will be accepted, before closing or executing the dialog.

Data Ally will verify that entries made in any entry fields fall within the ranges of values or types of values defined in Configuration, to ensure that incorrect or invalid data will not corrupt system operation. In some instances, making an inappropriate entry will immediately display an error message advising the operator to retry; in other cases, an error message may not appear until some activity, such as clicking the RUN button to begin a Method, occurs.

In those fields which specify disk path names (directory locations) for various types of files, any filename can be entered by typing the complete path name (e.g. C:\Data Ally\SUB1\data). All normal DOS file naming rules and conventions must be followed. Filenames cannot include spaces, colons, \, or / characters. Any desired number of subdirectory levels, and any real or virtual disk or memory device, can be specified in any filename field. Note that if any program file type (e.g. Method) is selected from a path other than the "default" path, the specified path will be retained for that program file type until it is changed again or until the default path is re-selected.

Check Boxes: Check boxes are used to enable or disable certain items, independently from all other items. If a check box for a particular item or option in a dialog box is "checked", that item is selected or enabled when you exit the dialog box. If the item is "unchecked", it is inactive or disabled. To

check or uncheck any check box, click the left mouse button with the pointer on the desired checkbox or on the item label for that checkbox. Each click will "toggle" the box to its alternate entry.

Radio Buttons: "Radio Buttons" are like check boxes, but differ in that they select only one from a group of mutually exclusive options or items in a dialog box. For example, if you select the radio button for the second item in a group of three items, in which the first item was previously selected, the second item's button will be "dark" or active, and the first item's button will become "inactive". To set a radio button, point at the desired button with the pointer and click the left mouse button. The setting for the button will toggle, along with the setting of any other active button. A group of radio buttons for selecting the display register for file loading appears at the middle right in Figure 4.3.

Example Areas: These are areas in a dialog box which are used to show an example of some function or attribute, such as a picture of a font you have just selected, or a color. Examples are provided to make it easier to choose the most desirable setting for certain parameters.

Command Buttons: Command Buttons are single-function buttons such as "OK" or "CANCEL" in Figure 4.3. Usually, command buttons cause another dialog box to appear (if the command on the button is followed by an ellipsis) or some other action to occur. If certain command buttons are not functional or available at certain times, they will be "grayed out" and cannot be used. To activate a command button, click on it with the left mouse button.

Control-Menu Box: The Control Menu box (the small "horizontal dash" button) at the upper left corner of certain dialog boxes has the same function as for any Window. Clicking the left mouse button with the pointer on the Control Menu box displays the Control "pulldown" submenu common to all Windows, with commands for Restore, Move, Resize, Minimize, Maximize, and Close (Figure 4.4). Not all dialog boxes in Data Ally include the Control Menu Box; those that do include it are effectively "windows" which can be moved, resized, minimized or "iconized", maximized, or closed. A "Switch-to" command is often provided which allows direct switching to a different application in Windows. In some cases only certain commands are available in the Control Menu for certain dialogs. For example, the Control Menu pulldown for the Files/Load dialog in Figure 4.3 only includes commands for closing and moving the dialog box.




Figure 4.4 Control Application Menu for Graph Screen Window

There are two ways to "close" a dialog box. If a Control-Menu Box is present at its upper left corner, click on the Box and select Close. If no Control Menu Box is present, generally clicking on either the OK or CANCEL button will close the dialog box.

The OK button in any dialog box is used to close that box and accept all the current settings in the box at the time OK is clicked. Pressing <ENTER> on the keyboard when viewing a dialog box has the same effect as clicking on the OK button.

The CANCEL button discards any options you have selected in the dialog box, and closes the box without changing anything in the system controlled by that dialog box. CANCEL differs from Close in that Close will close a dialog box without canceling any of the changes made in that box before it is closed. Instead of pressing the CANCEL button, you can use the <ESC> key. If the Control Menu Box exists, you can click on the Control Menu box and press either <ESC> or <ENTER> to close the dialog box.

4.2.0 Data Ally Windows

Like other Windows applications, Data Ally displays information inside various "windows" on the screen. Each window is surrounded by a double-line "box".

Data Ally organizes its windows into two classes.

- "Time base" windows are present for each "system", with only one time base window normally drawn on the full screen at any time. Separate and totally independent windows are provided for each defined "system" for both active data collection and file editing, as well as for individual "systems". You can switch between different time base windows by using the appropriate SoftButtons on the left screen. Each time base window contains one or more "program windows".
- "Program" windows correspond to different sections of the Data Ally program which is running in the time base being observed in the current time base window. Once inside a time base window, individual program windows are selected using the blue SoftButtons along the left edge of the screen. In certain cases it is possible to view more than one program window at a time on a single screen.

All windows of both types have some common features:

- Title Bar: The top, colored bar on the window displays the name of the Data Ally time base, if a time base window, or the name of the program window type, if a program window. Other information, such as loaded chromatogram or file name in program windows or the running mode indicator in a time base window, may also be found in the title bar.
- Application Control Menu: The "horizontal bar" symbol inside the button at the top left corner of all time base windows displays a menu for sizing and moving the time base window, or switching to other applications.
- Program Control Menu: The "horizontal bar" symbol inside the button at the top left corner of all program windows displays a menu for sizing, moving, or splitting the program window.
- Maximize Button: The Maximize button (upward arrow) at the top right of each time base and program window causes either the time base or program window to fill the screen. When you click on the Maximize button, it is replaced by the Restore button (double up/down arrow) which returns the window to its previous size.
- Minimize Button: The Minimize button (downward arrow) at the top right of each time base window causes the time base to be "iconized" without affecting its current function. Minimize does not operate in program windows.
- Restore Button: The Restore button (double up/down arrow) which replaces the Maximize button in the top right corner of time base or program windows returns the window to its previous size.
- Mouse Pointer: In program windows the "arrow" pointer or cursor is used to select commands and manipulate graphic images on screen. The appearance of the pointer may change

depending on where it is placed on screen (see section below on Mouse Pointers and Techniques).

- **Insertion Point:** In certain program windows only (Method, Integration, Calibration, Sequence) this blinking vertical bar indicates where the next character typed will be entered. The insertion point can be moved at will with the mouse pointer before making any entries.
- **Scroll Bars:** In program windows and "subwindows" such as Tables, a set of horizontal and vertical scroll bars can be used to change the current view of the contents of the screen area bordered by the scroll bars. To use scroll bars, you can either click the left mouse button while pointing at the arrow at either end of a bar, to move the display in small increments, or you can click on top of the "slider" button inside each scroll bar and "drag" it across the bar while holding down the left mouse button. Also, clicking the left button with the pointer anywhere in the bar will cause the scrolling to "jump" to that new position.
- **Help Message:** Most Data Ally screens have, in addition to the built-in HELP functions, a short "help message" area to the right of the Command Menu bar. This area will usually contain a brief phrase describing the current feature selected or the state of the system.

4.2.1 Switching Between Windows

➤ To switch from the current Data Ally "time base window" to see the current status of other time bases in the software, click on the SoftButton corresponding to the desired time base. Each time you change to a view of a new time base window, you will see the last program window which was viewed in that time base, in the exact state in which it was left when the time base was switched. It is also possible to switch directly to a different time base window (or to a different program running in Windows) by selecting the Switch-To command in any channel window's Application Control Menu, and then choosing the name of the desired Data Ally time base or program from the resulting listing. A third way to switch between various running windows is to press Alt + Tab simultaneously, which will cause an immediate jump to the last open program or time base window.

➤ To see a graphic display of all Data Ally time bases on a single screen, click on the MONITOR button at bottom right of any GRAPH screen, which has the special function of putting views of all the "active" time base windows present on screen at once. When the Monitor screen is displayed, clicking on any of the individual time base SoftButtons in any of the separate time base windows will return the full screen view to the selected time base only.

➤ To switch the program window being viewed in the current time base, click on the SoftButton corresponding to the screen type (e.g. Method, Sequence, etc.) you wish to see.

4.2.2 Making a Window "Active"

When more than one window is displayed on a screen, such as the Method Events Table or Integration Peak Table displays, only one window is "active" at any time, and only one window can accept entries. The active window is denoted by the "active" color scheme you designate when setting up Microsoft Windows.

➤ To select a window and make it "active", click anywhere inside that window. Alternatively, you can click on the Program Control button in the upper left corner of the main window on the screen, and choose "Next Window" or press CTRL+F6 until the window you want is active.

➤ To make the previous window active again, click inside it or press CTRL+SHIFT+F6.

4.2.3 Moving Windows and Dialog Boxes

If you need to move an open window on the display in order to see something behind it, you can "drag" it to any new position.

➡ To move a window, click the left mouse button with the pointer anywhere in the title bar of the window to be moved, and hold the mouse button down while you drag the window to the new position. At least some of the title bar must be visible after a window is moved.

You can also use the Move command from the Program Control Menu by clicking on its button.

➡ To move a dialog box displayed on screen, click the left mouse button with the pointer on its title bar and dragging.

A window cannot be moved if it is currently "maximized" and fills the screen.

4.2.4 Expanding Windows to Maximum Size

The maximum size of any window is the entire screen display. Channel windows in Data Ally are normally maximized.

➡ To maximize a program or other window, click the Maximize (up arrow) button in the upper right corner of the window's top bar (Figure 4.5). The window will expand to fill the screen. You can also use Maximize in the Program Control menu in any window by clicking on its button, by pressing CTRL+F10, or by double-clicking with the left mouse button on the top bar of the window.



Figure 4.5 Maximize and Minimize Buttons for Time Base (top) and Program Windows

4.2.5 Resizing Windows

Any window can be resized to fill any desired space on the screen. If you are using other software programs at the same time as Data Ally, you may wish to resize the Data Ally time base windows so that they appear beside windows in which other programs are running.

You can also enlarge windows which show, for example, a Peak Table or a Calibration Plot, for added convenience or to see more information in a single view.

➡ To resize a window which has not been maximized, move the pointer to any border or corner of that window until the pointer arrow changes to a "two-headed" arrow. With the "two-headed" arrow displayed, depress the left mouse button and begin to drag the border or corner until the new outline for the window reaches the desired size. Release the mouse button, and the window will be redrawn to the new size.

A maximized window cannot be resized - you must use the Restore function to "un-maximize" it first, by clicking on the "up/down arrow" Restore button at its top right corner.

You can also resize any window from its Program Control menu by clicking on the menu button at the top left corner of the window (Figure 4.4), and selecting the Size function (ALT, HYPHEN, R or CTRL+F8). The pointer will change to a "four-headed" arrow - press one of the arrow keys to indicate which border you want to move, and the pointer will change to a "two-headed" arrow.

Use the appropriate arrow keys to move the border as desired, then press <ENTER>. Pressing <ESC> will return the window to its previous size.

4.2.6 Closing a Window

➤ *To close a window*, double-click on the Application or Program Control Menu button at the upper left corner of the current window (Figure 4.4). The window will close and be removed from the display.

Alternatively, click on the Application or Program Control Menu button, and select the Close option by clicking on it or by pressing CTRL+F4.

4.2.7 Minimizing a Window

In Data Ally, you can "minimize" any of the available time base windows so that they become "iconized", using the Minimize function. Iconizing a time base will remove that time base's display from the screen and place an icon representing that time base at the bottom edge of the screen. Iconized time bases remain fully functional and will continue to execute any ongoing Methods, Sequences, or other activities while iconized, even though the executing functions cannot be observed while a window is iconized. Note that iconizing a time base actually iconizes both the "active" and "edit" time bases associated with the selected system.

➤ *To minimize a time base window*, click the left mouse button with the pointer on the Minimize button (down arrow) at top right of the screen. Alternatively, click on the Application or Program Control menu button and select the Minimize command.

➤ *To restore an iconized time base*, double-click the left mouse button with the pointer on its icon, which will place the time base window back on the screen in its last size and orientation before it was iconized.

If time base or other icons are not visible on screen after minimization, you can display them at any time by pressing the ALT+ESC keys. This will cause all icons for all current minimized application programs in Windows to be shown at the bottom of the display.

4.2.8 Splitting a Window

Data Ally Graph windows can be split into as many desired "panes" as you like, by using either the Split command in the Program Control Menu or the Display/Split command from the Command Menu bar. Splitting a Graph window lets you see different parts of the same chromatogram in different "panes" on screen at the same time, or puts different chromatogram traces in different windows in real time for differentiation. When Graph windows are split into more panes, all panes use the identical commands from the Command Menu bar on top of the screen, but only one pane is active at any one time - you select the active pane by clicking anywhere inside that pane. Each time a Graph window is split, it retains all the display characteristics of its "parent" window, but any such characteristics can be changed as desired in the "child" window.

➤ *To use the Program Control Menu to split Graph windows*, click on its button at the top left corner of the current Graph window, and select Split (ALT, HYPHEN, T) (Figure 4.6). The space filled by the current window will be split into equally sized panes. Each pane can be resized and repositioned by dragging, or panes can be subsequently split into additional panes.




Figure 4.6 Graph Screen Program Control Menu

To use the Display/Split command, see the section below on Graph screen.

4.2.9 Switching Between Programs in Windows

Since Data Ally runs in Windows "386 Enhanced" mode, other Windows or DOS applications can be running on the PC simultaneously. These applications can either be maximized to full screen, placed in windows of any desired size, or minimized as icons. Windows corresponding to non-Data Ally application programs can be on screen at any time along with Data Ally channel windows, in any relative sizes or orientations.

If you are running Data Ally as a full screen application, you may wish to switch to other programs.

➤ *To switch from Data Ally to another application window*, press the ALT + TAB keys simultaneously. Pressing ALT + TAB repeatedly will toggle through all non-iconized software applications currently in Windows.

➤ *To switch from Data Ally to another iconized application*, click on the Program Control button and select the Switch To command from the menu. A listing of all loaded software applications will appear. Click on the name of the application you wish to see and then click on the Switch To button. The display will change to show the newly selected application.

4.3.0 Mouse Pointers and Techniques

Your mouse will work with Data Ally in mostly the same ways as with other Windows applications.

4.3.1 Pointer Shapes

When you move the mouse pointer to different areas on the Data Ally screens, which correspond to different functions, the shape and appearance of the pointer may change, to help you recognize which functions are available.

If the pointer assumes a shape representing a function you do not wish to use, you can either press <ESC> to restore the pointer to its normal (arrow) shape, or move the pointer to a different part of the screen so that its shape changes automatically.

The table below shows the common pointer shapes and their meanings:

"Arrow:" The pointer is on any of the menus, inactive windows, scroll bars, Graph screen display, SoftButtons, or in areas of any dialog box which include buttons, check boxes, or radio buttons.

| "Vertical Bar:" The pointer is inside any text or numerical entry field into which values can be typed. The position of the vertical bar before or after any already-existing entry in a field box shows exactly where new entered characters will begin to be inserted.

←|→ "Left/Right Arrows Around Vertical Bar:" The pointer is on the edge of a table column header button, and can be used to resize that column. To resize, hold down the left mouse button and "drag" the edge of the column button left or right to change its size - release the left button when the desired size is achieved.

➡ "Large Left, Right, Up, or Down Arrow:" The pointer is in the furthestmost top, bottom, left, or right edge of a Graph window, and can be used to "pan" the graph display in the direction of the arrow by clicking the left mouse button while the large arrow is displayed.

⇅ "Double-Headed Vertical Arrow:" The pointer is in a window border area, and you can change the upper or lower horizontal border position of that window by clicking the left mouse button and dragging the existing border up or down.

⇆ "Double-Headed Horizontal Arrow:" The pointer is in a window border area, and you can change the left or right vertical border position of that window by clicking the left mouse button and dragging the existing border left or right.

"Double-Headed Diagonal Arrow:" The pointer is in a window corner border area, and you can change both its vertical and horizontal dimensions by clicking the left mouse button and dragging both borders to the new position.

"Four-Headed Arrow (Heavy):" You have selected the Move or Size command from the Application or Program Control Menu button. You can drag the entire window to a new position or move any borders of that window to resize it.

⇄ "Two-Headed Horizontal Arrow (Light):" The graphical editing function for timed events or peak windows is enabled. You can move the position of the marked point or code by depressing the left mouse button and dragging it to the desired location. This cursor is also used for graphic moving or offsetting of chromatograms in the Graph window.

| B "Baseline Code/Vert. Cursor:" The graphical baseline editing function is active. You can move the enabled code marker anywhere in the graph window and click the left mouse button to set it.

☞ "Hand Symbol:" The baseline or zone codes move function is enabled in the graph window. Click on the left mouse button to "grab" the code beneath the cursor symbol, and hold the button down while you drag it to the desired position. The "hand" symbol cursor is also used to drag chromatogram filename icons in the register stack, and to move other objects such as annotations or peak names in the Graph window.

✱ "Bull's Eye Symbol:" The baseline or zone code clear function is enabled in the graph window. Position the circular "bull's eye" over the code to be cleared and click the left mouse button to clear it.

☐ ↓ "Zone Code Box/Arrow:" You are in zone code placement mode in the graph window. Move the box containing the zone code symbol so that the arrow indicates its desired placement, and then click the left mouse button to place it.

- ↑ "Double Vertical Arrow Offset Cursor:" Graphical chromatogram offset has been activated. Moving the cursor up or down in the graph window will pull the selected chromatogram trace to a new position.

4.3.2 Basic Mouse Techniques

There are several basic mouse operations needed to make use of all Data Ally graphical functions.

<u>To</u>	<u>Action Required</u>
Point:	Move the mouse pointer until it is on or next to the desired object.
Click:	With the pointer in position, press and release the left (or right) mouse button.
Double-Click:	With the pointer in position, quickly press and release the left mouse button twice (there is no double-click function for the right button).
Drag:	With the pointer in position, press and hold down the left mouse button while you move the mouse. The object you are moving on screen will follow the motion. Release the button when the object is in the desired position.

4.4.0 Data Ally Tables

The Data Ally system utilizes "tables" in order to program method events for execution (Events Table), peak identifications, referencing, and handling (Peak Table), calibration standards and results (Calibration Reference Peak and Calibration Tables), and series of batched methods (Sequence Table). All the tables incorporate the same commands and techniques for making entries, moving through the table to review, adding or deleting lines, and filling or graphically programming the contents of columns.

All Data Ally tables are essentially "spreadsheets" consisting of a series of vertical columns, defined in Configuration, and at least one horizontal row. Rows are defined and added as tables are programmed. Columns can be moved around by direct "drag and drop", and their individual widths can easily be modified.

Any entry in any table column can be modified by pointing to that cell/entry, clicking the left mouse button to move the Edit cell highlight and vertical bar cursor inside that cell/entry, and typing in the new entry or change. Double-clicking with the pointer inside any cell will cause the existing entry in that cell to be completely cleared as soon as any new character is typed.

All table cells with predefined sets of allowed entries can be edited by first clicking inside the desired cell with the left mouse button, and then clicking the RIGHT mouse button to see a "pop-up" listing of all possible entries. Highlight the desired entry and click it with the left mouse button to place it in the cell - this will overwrite any previous value found in that cell. If a table cell does not have completely predefined entries (such as PeakName column in the Peak Table), a message will appear if the right button is clicked on that cell advising you of the type of entry that is required.

4.4.1 Moving Around in Tables

Whenever you make any Data Ally table window active by clicking inside it, a blue "highlight" appears around all four edges of the upper left cell of the table, called the "Edit cell". This highlight can be moved around the table in any direction to select any cell for making an entry, or to select a column or row for other operations.

You can simply point at any cell in the table with the cursor pointer, and click with the left button, to move the Edit cell immediately to the new cell location. The following keys and commands can be used to move around in any table:

Right and Bottom Scroll Bars	Use the vertical right side scroll bars to move up or down in the table, if it has too many rows to show inside the existing table window. Use the bottom scroll bar to view columns not observed in the existing table window.
ENTER Key	Press the ENTER key to move the Edit cell highlight one column position to the right. At the end of each row (last column in that row), pressing ENTER again will move the cursor to the first column in the next row.
TAB Key	Press TAB to move the Edit cell highlight one column position to the right, in the same manner as ENTER.
CTRL + ⇨	Press the CTRL and right arrow keys simultaneously to move the Edit cell highlight one cell to the right.
CTRL + ⇩	Press the CTRL and left arrow keys simultaneously to move the Edit cell highlight one cell to the left.
CTRL + ⬆	Press the CTRL and up arrow keys simultaneously to move the Edit cell highlight one row upward in the table.
CTRL + ⬇	Press the CTRL and down arrow keys simultaneously to move the Edit cell highlight one row downward in the table.
CTRL + HOME	Press the CTRL and HOME keys simultaneously to move the Edit cell highlight to the first (leftmost) cell in the current row.
CTRL + END	Press the CTRL and END keys simultaneously to move the Edit cell highlight to the last (rightmost) cell in the current row.
PgUp	Press the PgUp key to move the Edit cell highlight upward in the current column by the number of rows shown in the current table window.
PgDn	Press the PgDn key to move the Edit cell highlight downward in the current column by the number of rows shown in the current table window.
CTRL + PgUp	Press the CTRL and PgUp keys simultaneously to move the Edit cell highlight upward to the first (top) row in the current column.
CTRL + PgDn	Press the CTRL and PgDn keys simultaneously to move the Edit cell highlight downward to the last (bottom) row in the current column.

In addition to moving the Edit highlight from cell to cell, there are also a number of ways to move the "vertical bar" editing point to various positions inside the current entry in the Edit cell, after clicking inside the Edit cell to position the vertical bar cursor. Besides simply pointing to a certain position in

the entry inside the current Edit cell and clicking the left mouse button to fix the vertical bar at that position, you can use the following keys:

⇒	Press the right arrow key to move the vertical editing point one character to the right inside the Edit cell.
⇐	Press the left arrow key to move the vertical editing point one character to the left inside the Edit cell.
⇧	Press the up arrow key to move the vertical editing point to the first (leftmost) character inside the Edit cell.
⇩	Press the down arrow key to move the vertical editing point to the last (rightmost) character inside the Edit cell.
HOME	Press the HOME key to move the vertical editing point to the first character inside the Edit cell (same as ⇧).
END	Press the END key to move the vertical editing point to the last character inside the Edit cell (same as ⇩).

4.4.2 Adding, Moving, Copying, and Changing Rows

All Data Ally tables (with the exception of the Calibration Table) can be modified using the four yellow command buttons found beneath the left portion of each table: INSERT, CUT, PASTE, and COPY. These commands are analogous to the similarly-named editing commands found in many Windows applications. They are used as follows:

INSERT	Click the left mouse button on INSERT to insert a new blank row immediately above the currently selected row in which the Edit cell highlight appears. This is used to place new rows into a table at any desired position.
DELETE	Click the left mouse button on DELETE to delete the current row of the table in which the Edit cell highlight appears. Unlike CUT (below) DELETE commands are irreversible and cannot be used to PASTE rows back into the table.
CUT	Click the left mouse button on CUT to "remove" the current row in the table in which the Edit cell highlight appears. After cutting, the contents of this row will be placed in a special "buffer" in the Data Ally software, similar to the Windows Clipboard (but not actually using Clipboard), and can thereafter be PASTEd into any other row in the table. The CUT buffer can only hold one row's contents at a time, so CUTting a second row will cause the information from the first CUT row to be lost. CUT is used to either delete a row entirely or to remove it (CUT) and PASTE it into the same table in a different position.
COPY	Click the left mouse button on COPY to copy the contents of the current row in the table in which the Edit cell highlight appears into the table buffer. The current row is unaffected by a COPY command. COPY is used to insert a copy of

any existing row's cell contents into another blank row via a PASTE command.

PASTE

Click the left mouse button on PASTE after using CUT or COPY to transfer a row's contents into the table buffer. Before using PASTE, you must select a new blank row by either moving the Edit cell highlight to the target blank row, or by using INSERT to insert a new blank row and then moving the Edit cell into that blank row. PASTE causes the contents of the table buffer to be placed into the row currently containing the Edit cell highlight.

Note that the "spreadsheet editing" command buttons operate on only one row at a time.

5.0 Configuring Data Ally

Data Ally's Configuration capabilities constitute one of its most powerful and unique features. With Configuration, the system software can be personalized for any type(s) of applications or instrumentation. Configurability makes Data Ally the first truly "personal" data system since it can be customized in any manner to the preferences of its operator(s).

Use of configuration is normally recommended at the time of system installation

5.1.0 Configuration Files

Every time the Data Ally software is loaded, a "Configuration file" is automatically recalled from memory and used to structure the visual display screens of the user interface. Each Configuration file is essentially a "database" filled with user-specified settings and values for all parameters in the system. These settings include names of instruments, detectors, program files, and report templates, as well as default values for numeric entry fields such as run time or slope sensitivity. Configuration also includes any special software drivers needed to directly control specific instrumentation connected to each individual time base's Communications Processor.

At the time of installation, "generic" Configuration files for the intended system types (e.g. HPLC or GC) will be automatically loaded, so that you will not need to perform a configuration right away. This allows you to use the system immediately, and make detailed configuration changes later to adjust your system to fit your personal preference. The "generic" configurations represent the "default" settings for each time base.

The Configuration file loaded at bootup can be modified at any time, or replaced by a new Configuration. You can always use any existing Configuration file as a "template" which is modified and resaved to construct a new configuration. Although only one Configuration file can be loaded and used by Data Ally at any time, you can create as many such files as you desire for a "library" of instantly available potential software configurations.

Configuration files are designated by the file extension ".CFG". They are normally saved by Data Ally in a special directory called \Data Ally\CONFIG. In most cases, you will only need to configure your Data Ally system once, at the time you first install the system and connect your instrumentation. If you modify or add to your initial analytical system, or your requirements change, you may wish to reconfigure Data Ally at that time.

5.2.0 Performing a Startup Configuration

The following describes the initial configuration of a Data Ally system.

After the Data Ally software has been loaded, click on the UTIL SoftButton (F7) in any time base window to display the Utilities access screen. Click on Configure on the Command Menu Bar (Figure 5.1) to enter Configuration Mode. Note that the active window is now labeled 'Configuration' (Figure 5.2) and the top line Command Menu Bar has changed.




Figure 5.1. Selecting Configuration Mode From Utilities Screen.

When you are in Configuration mode and the Configuration window is visible, all Data Ally time bases are disabled **and any runs in progress will be terminated**. Therefore, you should be certain that all runs (acquire or reprocess) have been stopped before beginning reconfiguration.

In the Configuration Window, the name of the current time base view appears in the center top bar followed by <<Configure>>; the name of the currently-loaded Configuration file appears in the subordinate window bar below.

The acquisition time base has an associated "Edit" time base with its own dedicated SoftButton. The Edit time base, which cannot support live data acquisition from detectors but is used for off-line results review and processing, can be configured independently of the "live" time. You can use this capability, as an example, to set up special features in an "edit" time base which you do not wish to be present in the associated "live" time base used to capture data.




Figure 5.2 Configuration Window

5.2.1 Selecting Configured Components

The current Configuration screen is for the "live" time base now named TB1. To set up its parameters, click on Component and then on Component Setup to begin configuration. The Component Set dialog box (Figure 5.3) will appear.




Figure 5.3 Component Set Dialog Box

A list will appear in the box showing all the currently enabled objects in the current time base. For a new Data Ally system, only two "objects" will be shown in the listing: the current time base, called "TB1" or "TB2", and one analog detector input, called "Det1". If instrument control options have been included, the drivers for each instrument will also be shown. Begin by selecting the top item, the current time base object, by clicking on it with the left mouse button to "highlight" its name. Click the SETUP button to view the Setup Dialog for the System or time base window (Figure 5.4).




Figure 5.4 Time Base Window System Setup Dialog

This dialog provides information describing the current time base, including the application type, the desired name for labeling the window, the identity of the Communications Processor connected to this time base, serial and instrument numbers assigned, and any desired text description of the time base itself.

Entering System Descriptors

Click inside the "System Name" entry field and enter the desired name for labeling the time base, up to twelve characters long - this name will be shown in the top bar of the time base window. This name will be used to label both this "live" time base and its companion edit time base. Click inside the "System S/N" field and enter a specific serial or instrument number, if desired - this information will be permanently stored in Data Ally as part of this Configuration and can be recalled any time for reporting or labeling purposes.

To enter an additional text description of the current time base for documentation purposes, click the "DESC" button and type the description into the entry field provided (Figure 5.5). Click OK when the entry is completed (the text can be as long as desired) or on CANCEL to escape without changing the current entry. Note that you can also import text or graphics from another file into the Description box to save it on Data Ally with the current time base information simply by copying the text into Windows Clipboard and then Pasting it into the Description entry box.




Figure 5.5 Entering Time Base Description

5.2.2 Configuring Signal Inputs

If you will be connecting analog detector signal inputs to the Communications Processor for this time base, click on "Det1 Input" in the Component List dialog (Figure 5.3) to highlight it. Click again on the SETUP button to display the Analog Input Setup dialog (Figure 5.6).




Figure 5.6 Analog Input Setup Dialog Box

The Analog Inputs Setup box is used to set up each detector input. You can define a maximum of two different analog inputs for a single time base. Again, this setup can be done while configuring either the live time base or its associated editing time base - all entries will appear the same in both time bases.

Click inside the "Name" entry field and type the desired name for this detector, up to eight characters long (only six characters maximum will be shown in some locations in the software). You must enter at least a single character name for any defined input. This name will be used for labeling all Data Ally functions, displays, and reports which reference this detector. Normally, the name selected should clearly describe the type or model of the detector to help the operator immediately recognize the nature and source of each signal input.

At the bottom of the dialog box in the area marked "Hardware" are entry fields designating the number of the currently assigned Processor (Interface Module) and the identity of the input on that Module (A or B) assigned to this analog signal. You must return to the System Setup dialog to change the Interface Module number. Click inside the "Input Signal" field and enter the correct input identity (A or B) for this detector connection if this is not already correct.

Click on the DESC button to enter a text description for this analog detector input - the text can be of any desired length or format (Figure 5.7).




Figure 5.7 Entering a Description For Detector Input

When you are finished, click OK or CANCEL to return to the Analog Input Setup dialog.

If all entries are correct, click OK. Click CANCEL to avoid any changes to the previous entries for this detector.

Adding a Second Detector

To specify two or more detectors for the current time base, click on the ADD button in the Component Setup dialog (Figure 5.3) to see the Add Component dialog box (Figure 5.8).




Figure 5.8 Add Component Dialog Box

This dialog presents a listing of all possible system components which are possible to install via Configuration given the set of drivers supplied with this version of the Data Ally software. This listing appears in the menu box on the left side. Each type of component (e.g. detector, pump, autosampler, GC, CE, etc.) appears as a separate listing.

To install a component, first click the left mouse button on the desired item category on the left listing box to highlight it - the right box marked "Drivers" will immediately show a listing of all items in that category.

Begin installing a second analog input by clicking on "Input Signal" in the left list box, and then on the "Analog" item which appears in the right box to highlight it. When both items are highlighted, click OK to add the new item to the existing Components list, or click OK to CANCEL. The dialog will close and the Components Setup dialog will be displayed again (Figure 5.3). Note that a new item - a second "Analog Input" - now appears in the Components list. You can proceed to use the SETUP button to name and describe the second analog input in the same manner used to configure the first input above.




Figure 5.9 Component Setup Dialog After Adding Second Analog Detector Input

This procedure is used repeatedly to install any number of additional items to the existing components configuration for this time base. You can quickly review all the available items,

component types, and drivers available by clicking the BROWSE button, which will show a listing of all such items as you page through the list.

You can proceed to assign a third or even fourth analog detector via the same method, using the "Additional Inputs" function in the Detector Setup dialog to borrow analog inputs from different Processor Modules.

If digital detector inputs are to be used, they will be programmed in conjunction with a serial or other dedicated interface which is specified by selecting a detector driver using the Add Components function. See Section 5.2.4 below for a complete explanation.

5.2.3 Configuring Relay Flags

When all the desired analog inputs have been defined, the next configuration step is setting and enabling any relay or contact closure flags which will be used in the Communications Processor assigned to the current time base. The electronic relay flag contacts, although always provided in each Processor, must be configured in order to be functional and programmable while running the system.

Click on the ADD button in the Components Setup dialog menu (Figure 5.9) to see the Add Components dialog (Figure 5.8) again. Now click on "Flags" to highlight it in the left list box - two new items, Input and Output, will immediately appear in the right list box (Figure 5.10).




Figure 5.10 Adding Relay Flags

To set an Input Relay flag for starting Data Ally methods in response to an Injection contact closure from an autosampler, reed switch, or other device, click on Input and then click OK to add Input Flags to the Components listing (or click CANCEL to cancel the addition). The Components dialog will redisplay with the new Input Flags item listed, exactly as you added the second analog detector input above.

Now click on Input Flags in the Components Setup dialog to highlight it and then on the SETUP button - the Input Flags Setup dialog will appear (Figure 5.11).




Figure 5.11 Input Flags Setup Dialog

Each Communications Processor has four input relay flags or contact closures, identified on its connections panel by code letters. One of the input flags is marked "S" and is normally reserved for use with injection markers for triggering the start of a Method. However, any of the four inputs available can be used for this or any other purpose. If you have already wired one or more of the input relays to your installed instruments, you will now simply specify which inputs are used for what purposes, and give them appropriate names to help you when programming methods in this time base.

To assign a remote start function, click on the "Enable" field for the "S" input to check-mark it. Now click inside the Name field (on the same row) and type a name such as "Inject/Start". The name can be up to sixteen characters. Note this name will appear when you select the pop-up listing showing choices for input flag entries while programming method events.

Repeat the same procedure to enable, disable, or name any or all of the other input flags. Once you have configured flags in this dialog, the Input Flags column will be available for construction of the Method Events Table (below).

When you are finished specifying all input flags, click on OK, or CANCEL to avoid changing the previous entries. The Components Setup dialog will be redisplayed.

You can repeat the same procedure to enable and define any of all of the eight available output relay flags on each Processor Module. Use the ADD button to select Output flags, and then use SETUP in the Components Setup dialog to specify and name each desired flag.

5.2.4 Deleting Components

Once a Components listing for the current time base is generated, you can also delete or remove items previously added for the current Configuration. From the Components Setup menu (Figure 5.12), click on the desired item to be removed to highlight it. For example, to delete a second analog input from the current Configuration, click on Analog Input as shown.




Figure 5.12 Removing a Defined Component

Now click the REMOVE button - the highlighted item will immediately disappear from the listing. You can use REMOVE repeatedly to delete several items. You can replace each item by repeating the ADD command series as detailed above.

5.2.5 Setting Controlled System Instrument Components

Data Ally allows you to identify any and all instrumentation and accessories which form your chromatographic system. This is particularly important when instrument control options are present.

To select an instrument control driver from the Components Setup dialog, click the ADD button to see the Add Components dialog (Figure 5.8). Now click to highlight the item in the left list box corresponding to the type of instrument/component you wish to add (You must have Purchased the Control Package for that Component). For this example, click Detector to see the list in the right box of all available detector interfaces (Figure 5.13).




Figure 5.13 Selecting Desired Detector Type

begin doing so, click the left mouse button on Components in the Command Menu bar, to display the Components submenu.

Depending upon the type of system (e.g. HPLC, GC, etc.) you specified in the Set System command above, certain choices for component type may be "grayed out" in the selection listing which are not appropriate for that type.

As an example, assume your system is configured as HPLC. To set up this Data Ally time base to use the ESA 5200 electrochemical detector, click on ESA5200 in the right "Drivers" list box to highlight it and then click OK. The Detector Setup dialog box specific for this detector driver will appear (Figure 5.14).




Figure 5.14 ESA 5200 Detector Configuration Setup Dialog

You can now select the mode in which the ESA 5200 will be used for programmed control by clicking the appropriate selector button (Pulsed or DC) and indicating the number of analytical cells to be used if DC mode is chosen. In the lower area marked "Hardware" the currently assigned Processor/Interface Module number is listed and the serial port number on that Processor assigned to the detector is given - you can change this number if a different port is to be used. Ports on the Processor are automatically assigned during configuration on each time base in the order in which devices are defined and configured.

To enable the actual programmed control via Data Ally, click the Enable Control field at lower right to check mark it. If this field is not checked, the ESA detector in this example is still configured but the actual control mechanism within Method programs will be inactive.

Click the DESC button to enter a text description for the detector of any length, and click OK when done to return to the Setup dialog.

When the setup appears correct, click OK to finish the detector assignment. This detector is now set until the current configuration is changed or a new Configuration file is loaded to replace it.

You can proceed to set up other components in the same manner, if their drivers are supplied with your Data Ally software.

All instrument control interfaces ordered as options with each Data Ally unit should be listed as Drivers by brand and model under each component type. If control interfaces you have ordered do not appear in the listing, contact your distributor for further information. In some cases, a particular instrument can use a driver protocol originally intended for a different instrument.

Clicking SETUP from the Components Setup screen for specific instrument brands or models added to Configuration will often show unique Setup dialog boxes reflecting the independent parameters required to define and operate each instrument interface. All setup dialogs permit the recording of basic items such as brand, model, S/N, instrument number, and description.

You are not required to select nor complete the setup dialogs for any specific instruments you may be using with Data Ally. It is best to configure only those for which you wish to perform direct programmed control and those whose function and presence in the analytical system must be thoroughly documented. Even if your Data Ally software includes control drivers for one or more specific instrument types, you do not need to configure and use them in order to actually run Data Ally.

When you are completely finished selecting and setting all configured elements, objects, and instruments, click OK in the Components Setup dialog to return to the main Configuration Window (Figure 5.2).

5.2.6 Program Screen Options, Menus, and Defaults

The next step in Configuration is to specify the contents of each of the primary programming and display screens, select which options are to be used, and set all relevant defaults for variables in each screen. By permitting you to customize every aspect of the user interface in this way, Data Ally provides for the most perfect fit possible between the features you can see and use and the nature of your application(s).

Each of the program screens can be configured differently for any "live" time base and its associated editing time base, for maximum convenience and to fit your applications. The only parameters which will remain constant are the time base name, detector input names, other instrument component names, and flags.

You can also set names and colors for all screen elements, such as buttons, fields, and labels, which will appear in the operating software.

You can set up each of the screens which are normally displayed by pressing their respective SoftButtons, by pressing the usual button for the desired screen. The HELP SoftButton provides help during this process, and the "channel selector" SoftButtons are inactive, except for the Edit channel SoftButton which contains configure functions specific to editing.

5.2.6.1 Graph Default Screen Setup

Begin by clicking on the Graph SoftButton with the left mouse button to set the normal Graph display configuration. The Configuration Graph window will appear (Figure 5.15).




Figure 5.15 Configuration Graph Window

This Graph window has nearly all the functionality of the normal Data Ally Graph screen (Section 6 below) so that you can try each feature and option available in order to determine which Graph configuration is best as a default for your applications. Note that the RUN SoftButton is non-

functional since no live data acquisition or reprocessing is possible from Configuration Mode - clicking the RUN button has no effect.

Configuring Chromatogram File Handling Functions

The Graph Files commands are available during Configuration both to load chromatograms for experimentation and to set defaults for file loading and saving.

If you wish to load a chromatogram file, assuming at least one such file is available, in order to experiment with configuring display options, click the Files command on the top menu bar and select Load from the submenu (all other commands in the submenu are inactive) - the Load dialog will appear (Figure 5.16).




Figure 5.16 Loading a Chromatogram In Configuration Graph Screen

Select a *.pdf (Data Ally data file) file from any directory or disk using the dialog, and click OK to load this chromatogram file into the Graph window. The chromatogram will appear labeled with any integration codes, peak names, and other annotations.

The normal default directory specified for the Data Ally software is a directory named C:\Data Ally. This directory path will be selected first when you utilize the File Load dialog. By using this dialog to load a file from a different directory or disk path, you can reset the default path for loading chromatograms. WHEN YOU EXIT AND RESAVE THE CONFIGURATION FILE FOR A TIME BASE, THE FILE LOAD AND SAVE DIALOGS WILL BE DEFAULTED TO THE PATH/DISK/DIRECTORY FROM OR TO WHICH YOU LAST LOADED OR SAVED A CHROMATOGRAM FILE. This is a convenient means of resetting the default file handling commands if you wish to primarily use a directory other than C:\Data Ally. Note that you can always select any path for file transfers using the Graph Files commands in the Data Ally software, regardless of the default path settings.

Configuring Default Edit Functions

Click on the Edit command on the Graph top menu bar to see the Edit submenu in configuration. Although all the commands are available, including Cut, Define Baseline, and Define Zone Codes, only the Analyze, Calibrate, and Integrate functions can be configured for specific requirements in the current time base. See Section 7.3 below for a complete explanation of the use of these commands for analyzing and manipulating chromatogram data; the following section deals only with configuration related changes to the Analyze, Calibrate, and Integrate dialogs.




Figure 5.17 Graph Edit Default Submenu

Click Analyze in the Edit submenu to see the Analyze Default dialog (Figure 5.18). This dialog is used to produce result reports manually from loaded and integrated chromatograms.

Figure 5.18 Analyze Default Dialog

The Sample related information fields can be completed to set the defaults for sample name, number, and description. These should be set only if your applications in this time base will require frequent use of the same sample descriptors - otherwise, it will be easier to make these entries as needed while processing chromatograms.

Any entries made here in the Report Template Name or Report Filename fields or the Print/Show checkboxes will be used as the defaults for this dialog in the current time base when the default Method is loaded. This allows you to specify the type of reporting to be done for any chromatogram or sample as the default. To select a default report template from any selection of *.XLS report template files in the \Data Ally\REPORTS subdirectory, click on the arrow button selector at the right side of the Template field to see a listing of all such template files (Figure 5.19) - click on the filename desired as the default, if any.




Figure 5.19 Selecting Default Report From Existing Template Files

During actual operation, Data Ally will capture all other sample-specific information in this dialog from the Method program itself or from individual chromatogram file records, if and only if such information exists. You can of course change any of the default reporting settings when running Data Ally, but it is usually good practice to specify the most commonly used reporting type as the default. Do not specify a default report file name if you plan to frequently use different report file specifications.

Note that the contents of the boxed areas marked "Correction Factors" and "Internal Standard Amounts" are actually determined by how the master Calibration screen setup is performed (below, Section 5.2.5.4). If you wish to specify defaults for different types of correction factors from those actually shown here, or for a different number of internal standard peaks, you should first click the Calib SoftButton in the left column and reset these items in Calibration screen configuration.

When all desired defaults are entered correctly, click OK to save them - if you have specified a report in the default settings, the system will attempt to produce that report with any currently loaded chromatogram.

Now, click Edit in the Graph top bar menu and then on Calibrate in the submenu to see the Calibrate Default dialog (Figure 5.20).




Figure 5.20 Calibrate Default Dialog

Again, as with the Analyze dialog above, if you wish to add or remove correction factors or internal standards from this dialog, you must first modify the Calibration screen configuration by clicking the Calib SoftButton. If defaults values for any factor(s) are entered here, those defaults will appear in the actual Edit/Calibrate dialog in this time base only if information specifying those factors is not acquired from the loaded chromatogram raw data file.

Just as with the Analyze dialog above, the items in the "Report" area at the bottom can be defaulted during configuration in this dialog. All other items will be automatically entered from the current Method file or will be loaded from chromatogram records when processing files. Enter any defaults as desired, and click OK when done, or Cancel to retain the previous settings.

To review and reset defaults for manual integration, click the Graph Edit command and then Integrate in the submenu to see the Integrate Default dialog (Figure 5.21).




Figure 5.21 Integrate Default Dialog

This dialog can be used for manually integrating loaded chromatograms, once the required algorithm parameters have been entered in the appropriate fields. You can click inside any of the parameter fields in the dialog and type in a suitable value - these entries will also appear in the Integration screen configuration dialog (below). If you have already completed the Integration Configure dialog, the corresponding algorithm values entered there will appear in the same field locations here. See Section 8. below for a complete explanation of the meaning and interpretation of these parameters.

Checking the "View Only" checkbox will default to manual integration of the current graph view of the foreground chromatogram only. If you wish to perform primarily global integrations (entire chromatograms), leave this box unchecked in this dialog.

When all settings are as desired, click OK to close. You can experiment with trial parameter settings to determine the best values if you have loaded a typical chromatogram in the Graph window. Each time you click OK, the system will attempt to integrate the current foreground chromatogram as specified.

Configuring Default Graph Display Functions

Now, set the Graph Display default parameters by clicking Display on the top menu bar to view the Display submenu (Figure 5.22).




Figure 5.22 Configuration Graph Display Submenu

You can select any of the submenu items to set the defaults associated with the assigned default Method file.

Click on Attributes to set default display attributes (Figure 5.23).




Figure 5.23 Setting Default Display Attributes

In the Attributes dialog, click to check or uncheck the boxes to enable or disable all the desired default options. These are divided into sections - Fill entries determine color filling, Show entries determine display of codes or annotations, Baseline entries determine what baseline parameters will be shown, and Peaks determines what peak descriptors or peaks will be displayed. You will see the effects of changing any attributes on the loaded chromatogram in the Graph Window. When all the desired attributes are set, click on OK, or on Cancel to use the previous settings.

See Section 7.6.1 below for a more detailed explanation of the various attribute functions. Note that regardless of your settings, you will have complete access to changing any attributes at any time in the Data Ally software after configuration is completed.

Now click the Colors command in the Graph Display submenu to see the default colors setup dialog (Figure 5.24). This dialog selects the display colors associated with the default method.





Figure 5.24 Setting Default Display Colors

To change any of the colors, click the left mouse button with the pointer inside the "color box" corresponding to the parameter or detector trace you wish to reset. The series of color boxes at the top left marked "Trace Colors" will correspond to the detector inputs you have just specified for this time base as you configured system components and inputs. Grid colors correspond to background colors for the Graph window, Fill colors to chromatogram filling, and Attribute colors to baseline and peak descriptors.



The Component Set dialog box (Figure 5.3) will appear.

To modify any default color, click inside the color box corresponding to the desired parameter to see the Color Selector Dialog (Figure 5.25). The color palette can be used to change to the exact color you prefer - click on the desired color in the palette to highlight it and then click OK, or click Cancel to retain the current color setting. The color box on the Colors Dialog will change accordingly. When all colors are adjusted, close the Default Colors Dialog (Figure 5.24) to return to the Graph Window.




Figure 5.26 Offset Default Dialog

Click the Graph Display top menu bar command and then click Offset to see the Offset Default Set Dialog (Figure 5.26). This allows you to specify the signal offset level (either positive or negative) from the true baseline level for display of chromatograms in any of the defined signal input and display registers. If one or more chromatogram files are loaded, their current offsets in μV will appear adjacent to their names in fields labeled with the detector input names you have configured. If you wish to specify a permanent offset on the signal axis associated with the default Method file, click inside the offset value field corresponding to that input register and type in the desired offset. When all offsets have been entered as needed, click OK to close the dialog or Cancel to keep the previous settings. Note that an entry of "0" (zero) will select no offset from the true signal value. You have complete access to setting offsets either graphically or numerically at any time in the Graph screen - entries of values other than zero are suggested primarily in cases in which you are aware that the overlap of multiple detector signals will make them difficult to distinguish unless one or more signal traces are offset.




Figure 5.27 Pumps Offset Display Dialog

If the current time base is configured for HPLC and pump control is defined in the Components Configuration, a button will appear in the Set Offset dialog to set the default pump trace display offsets. Click this button to specify these offsets with the Pumps Offset Display Dialog (Figure 5.27). Fields are provided in this dialog for up to four individual pump gradient percentages (A through D) and for a Flow value in mL/min depending upon the actual configuration. Click inside any of the fields and enter an offset value from zero in the appropriate units. When finished, click OK to close the

dialog or Cancel to retain its previous entries. Again, this dialog sets the defaults only; the actual offsets can always be changed during actual Data Ally operation.




Figure 5.28 Display Scaling Default Dialog

Click on the Graph Display top menu bar command and then on Scaling in the submenu to set the graph scaling associated with the default Method file. The Scaling Default dialog (Figure 5.28) will appear.

This dialog specifies the type of default scaling for chromatograms in the Graph window in this time base. Scaling is set separately for the Signal and Time axes, either by checking the "Auto" box adjacent to either field to specify "autoscaling" or by clicking inside the minimum and maximum value fields and typing in actual scaling display limits in the appropriate units. If "Auto" is checked for either Signal or Time, that axis will be adjusted automatically according to the magnitude of the chromatogram(s) loaded or captured in the time base. See Section 7.6.3 below for a more complete explanation of Scaling functions. As with all other parameters, you can change scaling at any time during actual operation, including during real time data processing. It is recommended that you set the defaults for what is normally the most convenient display presentation for your chromatograms and applications. When scaling defaults are set, click OK to save them or Cancel to retain the current settings.

If the current time base is configured with HPLC pump control, a Pumps button will appear in the dialog for setting the pump trace scaling defaults. To adjust these settings, click on Pumps to view the Pump Scaling Default dialog (Figure 5.29).




Figure 5.29 Pump Scaling Default Dialog

You can specify either autoscaling independently for the Flow and Gradient Per Cent traces, by checking the "Auto" check boxes corresponding to either or both, or you can click inside the Minimum and Maximum value entry fields to set the default scaling limits. Of course, it is possible to rescale pump parameter traces at any time while running Data Ally regardless of this configuration. When setup is as desired, click OK, or Cancel to retain the previous settings.

If you wish to set the default Graph display with more than one window "pane" when the default Method is loaded, click Display in the top menu bar and then click Split in the submenu - a "child" window which is an exact duplicate of the existing window will appear. You can repeat using the Split command to create up to four separate panes as the default scenario. Each can have its own Graph defaults set independently; simply click anywhere inside any window to make it "active", and then select the Graph Display Attributes, Colors, Offsets, or Scaling dialog to specify the exact defaults for that particular window. You can also close any window by using the Window control button at upper left and selecting Close in the submenu.

To set a Status Window as the default, click on Graph Display in the top menu bar and then on Status to see the Status Box Default Dialog (Figure 5.30).




Figure 5.30 Status Box Default Dialog

This dialog allows you to select a status window which will appear in the current time base when the default Method file is loaded, containing specified parameters regarding chromatogram or trace signals, run timing, and sample or sequence related information. The left portion of the dialog contains a list of parameters which can be included in the status display window. To select the desired parameters, click on each one to highlight its name. When all the desired items are highlighted, be sure the Enable check box is checked and click OK to enable the defaulted Status Box, which will appear on the current Graph window as specified. If you do not wish to have the Status Box enabled, click the Enable checkbox to "uncheck" it.

The Status Box is the primary means of reporting real time information in Data Ally concerning the signals being captured and displayed, and the progress of the running method. You should enable the Status Box as default if you will normally use this function to track the operation of your methods. Like all other display functions, you can switch the Status Box on or off at any time during actual operation, or move or resize it at will.

The Swap function in the Graph Display configuration submenu has the same purpose as in live operation, to move the chromatogram file in the "bottom" display register to the "top" register. This command is active during configuration for use in manipulating files but has no other setup or defaulting requirements. See Section 7.6.4 below for a more complete explanation of the Swap command.

Configuring Default Graph View Functions

Next, configure the defaults for the View functions, the Zoom Box and Zoom Center mode.

Click the Graph View command in the top menu bar to see the submenu (Figure 5.31).




Figure 5.31 Graph View Default Submenu

To specify the Zoom Box to appear when the default Method is loaded, click on Zoom Box - the Zoom Box window will appear. The Zoom Box is a "range-finding" tool which is most useful when editing long chromatograms, to help you locate your current position when zooming various views. Normally, the Zoom Box will not be enabled as a default. If the Zoom Box is already enabled as the default, simply close it to defeat this setting.

Click on Zoom Center to set the Zoom/Center magnification parameters for both the time and signal axes - the Zoom Center submenu (Figure 5.32) will appear.




Figure 5.32 Zoom Center Default Submenu

To activate the Zoom Center function as default for this time base, click on Enable Zoom Center Scaling to "check mark" it. When Zoom/Center is activated, double-clicking anywhere on the current Graph window will "zoom" the display centered around the cursor position by preset expansion factors.

To reset the expansion factors as default, click on Set Zoom Scale Factors to see the Zoom Center Scale Factors Default dialog (Figure 5.33).




Figure 5.33 Zoom Scale Factors Default Dialog

You can now click the left mouse button inside each of the set field boxes and type in the desired magnification factor. Note that a factor greater than one will cause the display to "expand", while a factor less than one will cause it to "contract". Recommended settings for are factors of between two and five for both axes. When the factors are set as desired, click OK to save them, or Cancel to retain the existing factors. You can of course reset the default factors using this dialog without actually enabling Zoom Center as a default function.

Configuring Print Functions

You can set the default configuration of the manual printing functions in the Graph Print submenu by clicking the Print top bar menu command to view the Print Default submenu (Figure 5.34).




Figure 5.34 Print Command Default Submenu

Click on Chromatograms in the submenu to see the printer setup dialog (Figure 5.35). This dialog permits specification of the manual chromatogram print output associated with the default Method file.




Figure 5.35 Printer Setup Default Dialog

This dialog selects the current printer type, the layout of the default chromatogram print, the number of copies printed, whether the supplementary information is included, and allows access to additional setup functions. A more detailed and complete explanation of the use of this dialog is provided in Section 7.7.1 below. It is advisable to configure the default for the actual printer to be used with this time base; remember that only a single printer can be specified for Data Ally.

When all printer and layout information is completed as desired for the default, click OK, or Cancel to keep the existing settings.

Click the Graph Print top bar menu command and then Printout Info to default the desired chromatogram information to be included on printouts in this time base. The Printout Info Default dialog will appear (Figure 5.36).




Figure 5.36 Printout Info Default Dialog

This dialog specifies information specific to the foreground chromatogram which will be included in chromatogram graphic printouts in the current time base when the default Method is loaded. Note that unlike the printer setup dialog, the contents of this configuration dialog can be completely independent for each time base. Check all the checkboxes corresponding to the printout information to be included as default. Click OK when finished. Note that the printout information can be modified prior to any printing activity while running Data Ally no matter how the defaults are set; it is recommended that you specify only those information elements as default which you will always want to include.

Now click Graph Print on the top bar menu and click Report in the submenu to see the Chromatogram Report Default dialog (Figure 5.37).




Figure 5.37 Chromatogram Report Default Dialog

This dialog sets the default Excel template file for printing chromatogram reports and the file name used for saving such reports, which can be specific to each independent time base. Use these settings only if you intend to print out special report formats routinely when manually producing hard copy chromatograms.

Click the "down arrow" button at the right of the Report Template field to see a complete listing of all the report template *.XLS files available in the current \Data Ally\REPORTS subdirectory. If you wish to select one of these reports as the default, highlight it and click the left mouse button. If you want to omit any default for the report template type, simply leave this field blank.

When you have made the desired entries, click OK or click Cancel to escape without changes.

You have now completed initial configuration of the complete Graph menu which will be associated with the default Method file in this time base. There is one other default setup task relative to the Graph screen - the labeling and color of the Hold, Clear, and Monitor buttons at bottom right.

All Data Ally Softbuttons and display buttons can be reset individually and renamed as desired. To configure the "Hold" button, place the cursor on the button and click the RIGHT mouse button - the button configure pop-up submenu will appear (Figure 5.38).




Figure 5.38 Graph Screen Button Configuration Pop-up Submenu

There are two items in the poppa - Set Color and Set Title. To change the color of this button, click on Set Color to see the Set Button Color Dialog (Figure 5.39).




Figure 5.39 Set Button Color Dialog

The Set Button Color Dialog contains a color palette from which you can select the desired button face color by clicking on any color field. When the selected field is highlighted, click OK to recolor the button or Cancel to escape. The button color will automatically be reset.

To change the title of the button, click the right mouse button again with the pointer on the button - the Set Button Title dialog will appear (Figure 5.40).




Figure 5.40 Set Button Title Dialog

Click inside the entry field and type in the new name for this button, up to eight characters per line. You can also enter a second line - press the CTRL and ENTER keys together to move to the second line and type characters. Note that the "nominal" title of each button is always listed above the entry field (in this example, Hold) to help you identify the function of each button you are renaming no matter what actual name you assign. If you wish to change the text color for a better match to or contrast against the selected background color, click inside the Title Color field. When the button name is indicated correctly, with the desired color, click OK to rename the button or Cancel to escape.

This same mechanism is used to rename and recolor all the Data Ally left column SoftButtons at any time, and with all other buttons found at the bottom of the Method, Events, Integration, Peak Table, Calibration Tables, and Sequence screens. You can select any combination of colors and names desired. You should assign the names which are the most familiar to the system operators and the most appropriate for your specific applications.

You can also rearrange the order of the three buttons at bottom right of the Graph configuration screen, within the "solid border" that surrounds them. To move a button to a different position, simply move the mouse pointer onto that button, hold down the left mouse button, and "drag" the entire button to the new location - then "drop" it by releasing the mouse button. The other buttons will be displaced and the "dragged" button inserted in the "dropped" location. You can also drag-and-drop in this manner to rearrange the order of the SoftButtons in the left screen margin or any other sets of screen-specific buttons.

Now proceed to configuring defaults for the Method program screen.

5.2.6.2 Method Default Screen Setup

Click on the METHOD SoftButton to set the Method screen defaults. A screen similar to the Graph default screen will appear (Figure 5.41).




Figure 5.41 Method Configuration Dialog

There are two sections to the Method screen which are configured separately: the top level or "global" method dialog, and the Method Events table window. The contents of both are completely controlled by the manner in which you have specified Components Configuration for this time base above.

Top Level Method Setup

The Method Configuration Dialog (Figure 5.41) contains information regarding the default Method file name, the detector signal inputs, the handling of chromatogram data files, and the processing of runs, in separate "blocked" sections of the screen. Since two independent detector inputs have been configured for this time base in this example above, the Method dialog shown has parameter entries for two inputs named "UV-Vis" and "ECD" - these could be renamed or one or the other removed by clicking the Utilities SoftButton and selecting the Components configuration commands again as in Section 5.2.1 above.

Default File Name: When Data Ally is first initialized all default program files will appear in their Configuration windows with the name "Default.xxx" where xxx is the correct name extension, in this case, .met. There are three ways to assign the default method file which will automatically be loaded each time the current Data Ally time base is started, and whenever another method is cleared from the time base; you can accept the name "default.met" and proceed, you can click inside the Name field and type in any desired name (up to eight characters), or you can click on Files in the top command bar to select a Method file which already exists on disk and "assign" it as the default with all its existing parameters (see below). It is recommended that all default files be given the name best suited to describe the nature of the method itself - if you perform mostly one type of method on this time base, you may not need to program any Data Ally methods other than the one defined as the default to do all your applications.

Once the default filename is shown correctly, click inside the Run Duration default field to set the default run time in any desired time units. Type the value and units designator into the field and either press ENTER or click inside another field on the screen. The new units will appear in the gray box at right of the DURATION field. If the value already present for DURATION is acceptable, you do not need to make any entry in the field.

Click the SAMPLE button to define the sample-specific entry fields associated with methods in the current time base. The Sample Information Configuration dialog will appear (Figure 5.42).




Figure 5.41 Method Configuration Dialog

This dialog defines what information fields will be available for entering sample-specific data for processing while executing individual methods, and sets the default values for those parameters, if any.

In the top area of the dialog marked "Sample", click inside the check boxes in front of each sample parameter listed to check those parameters you wish to be visible and accessible for entering sample information. Do not checkmark those items you do not wish to use in your method files. You can click inside the entry field corresponding to any parameter(s) and make entries which will

automatically appear in the default method. Such entries can of course be changed manually when running the default method or when programming other methods; however, this is recommended only if the default method will be used for nearly all applications in this time base and common sample-related information is desired for output to typical reports.

NOTE THAT ITEMS WHICH ARE NOT CHECKMARKED WILL NEVER APPEAR IN ANY METHOD PROGRAMS CREATED IN THIS TIME BASE.

In the center portion of the dialog marked "correction factors" are fields with check boxes which correspond to the currently-enabled and available correction factors specified in the Calibration screen configuration for this time base (Section 5.2.5.4 below). If the factors shown are sufficient or appropriate for your applications, you can enter any desired default values for any such factors. If you wish to add or delete factors for this time base associated with the method files, including Dilution, Weight, and Injection Volume, click the Calib SoftButton and enable the desired factors, and then return to this screen by clicking the Method SoftButton followed by the Sample button to complete your factor assignments.

You can click inside the respective entry fields for each type of factor and type in a default entry value, which will always appear when the default method is loaded. These entries can be modified before executing the method.

You can also "lock" any or all factor values as desired simply by typing in a suitable value for any factor's field and then un-checking its check box. In this manner, you can still use individual factors for custom computations in reports, etc., but you will not be able to access them at all using the SAMPLE button from methods in this time base.

In the bottom portion of the screen are several fields and check boxes corresponding to Internal Standards (1 to 5). The number of internal standard peaks available for methods in this time base is also specified by the Calibration configuration - to change the number of such internal standards, click the Calibrate SoftButton and enable or disable internal standards as needed before returning to this screen to set appropriate values.

When the correct number of available internal standard peaks is defined for the default method, you can click inside any or all of the enabled value entry fields to insert IS concentration values as defaults, if desired.

Click on the Sample Information Enable check box at right center to check it if the SAMPLE button is to be visible and accessible with the defined information in your method files in the current time base. If the Enable box is not checked, the SAMPLE button will not be displayed and it will not be possible to enter sample-specific information for method processing in any method file.

When all the desired entries have been made, click OK to return to the Method Configuration Dialog or CANCEL.

In the boxed area at top right marked "Detectors", separate rows and fields will appear for all the detector inputs you have chosen in the current Components list for this time base. Each can be enabled or disabled by default independently of all others by checking/unchecking the box before each detector's name. Be sure all the detectors you wish to default as "active" are checked.

Each detector can also have its resolution or sampling rate set independently. Click inside the RESOLUTION field for any detector and type in the desired value followed by a designator for any acceptable units. If only a number is entered in the field, the existing units shown in the gray box at right of the field will be used - if you have also typed a new units designator, the new units will be shown in the gray box after you confirm your entry by pressing ENTER or clicking in any other field on the screen. Again, if you wish to remove or rename a detector signal input, you must return to the Utilities screen in Configuration and modify the current Components configuration settings.

The boxed area marked "Results" lists separate rows for each defined detector input, showing the default chromatogram raw data file name to be used, a button accessing Reports for each chromatogram, and fields to select the Integration and Calibration files to be used, if any. If you wish to default to a specific filename for any or all detectors, click inside their respective FILENAME fields and type the desired filename(s) up to eight characters. If no filename is given for any detector input, the FILENAME field will come up blank by default and you will need to enter or select a filename for each input before starting any runs using the default method.

If you wish to have any or all detector input signals automatically integrated at the end of your methods, including the default method, click the check box in front of the INTEGRATION file entry field for each such detector input to check it. **IF YOU DO NOT CHECK THIS BOX FOR ANY DETECTOR, YOU WILL NOT BE ABLE TO DEFINE AUTOMATIC INTEGRATION AND PEAK FINDING FOR THAT INPUT AS PART OF ANY METHOD.**

Enter an Integration program file name, if desired, as the default for any enabled/checked INTEGRATION fields for detectors. This field supplies the program filename to be used for handling the data from each detector independently in the default method. If you do not enter a filename here for any detector, you will need to enter such a name before actually running any Data Ally method in this time base if you want automatic analysis of that detector's data file during that method. Clicking on the "down arrow" button at the right of this entry field will show a current listing of all available *.INT files which you can associate with this default method - you can highlight the desired filename, if any, and click to enter it.

If you wish to perform automatic calibrations using standard materials with the default method in this time base for any detector inputs, click on the check box in front of each such detector's CALIBRATION filename entry field to check it. **IF YOU DO NOT CHECK THIS BOX FOR ANY DETECTOR, YOU WILL NOT BE ABLE TO DEFINE AUTOMATIC CALIBRATION FOR THAT INPUT AS PART OF ANY METHOD.** Use the down arrow button at right of this field to choose from the current selection of *.CAL files available as defaults.

For each check-marked detector, you can now enter a default calibration filename in the CALIBRATION field, if desired, and/or designate a calibration TYPE or LEVEL number in the appropriate fields. If you leave any of these fields blank for any detector, you will need to make entries to use any calibration functions before actually running the default method in this time base. It is best to leave these fields blank unless you will be using applications in which very frequent calibrations are performed.

To set the default reports for chromatograms from any detector input, click on the REPORTS button corresponding to that detector to view the Method Reports Default dialog (Figure 5.43).




Figure 5.43 Method Reports Default Dialog

To set reporting capability for any individual detector input, first click the ENABLE check box to check it. The "number of reports" selector buttons will immediately become accessible as shown in Figure 5.43.

Now click on the button corresponding to the maximum number of reports you wish to be able to produce during method execution for the current detector input (1 to 4). Immediately after clicking on the desired number, that number of individual report fields will become enabled in the bottom portion of the dialog.

Now click inside each of the check boxes corresponding to the individual report designator fields to check the number of reports/rows you wish to be active for this detector in the default method. You can also click inside the Report Template and/or Filename fields to enter Excel template filenames and/or saved report filenames, respectively, as defaults. In the Template field, you can click the right mouse button to see a "pop-up listing" of all report templates provided with your Data Ally system, from which you can select a template for any row - of course, you can add new or modified report templates to the default method at any later date once they have been created. If you leave these fields blank for any or all lines, you will need to make entries before you can use the automatic reporting functions when actually running the default method.

When all reports are shown as desired, click OK to return to the Configuration Method Window or click CANCEL to avoid making any changes.

To enter a description for the default method file, click on Files in the top command bar menu to see the Method Files Configuration submenu (Figure 5.44).




Figure 5.44 Method Files Configuration Submenu

Click on Method Description in the submenu to see the Default Method Description dialog (Figure 5.45).




Figure 5.45 Default Method Description Dialog

Type in a text description of any length for the default Method file - you can begin new lines by pressing CTRL and ENTER simultaneously. This description will be available in the Data Ally Method files directory to assist in selecting Methods to be loaded. Click OK when done or CANCEL.

It is recommended that, as with all other default settings, the default Method screen be configured for the most common routine situation for this channel, which is most convenient since no additional programming will need to be entered whenever that situation is present. The default Method can always be modified in any manner to create custom Methods, which can be loaded for use at any time.

Assigning an Existing Method File As Default: Instead of creating a new name and completely entering a default method file, you can "re-assign" an existing Data Ally method file and its file configuration as the default/configuration for the current time base. An assigned file will replace any existing file shown in the Configuration Method Window and will substitute the new file's parameters and configuration in place of the old.

To do this, click on File in the top command menu bar and select Assign Method from the pulldown submenu to view the Assign Default Method dialog (Figure 5.46).




Figure 5.46 Assign Default Method Dialog

Use the path selector and file listing functions to find and highlight the desired method file (in any path accessible) to use as the new default file. When the desired method name appears, click OK to select it, or CANCEL to abort the assignment.

If a valid method file is chosen, it will be loaded into the Configuration Method Window and all its settings and parameter values will be used to modify the current default configuration. You can then proceed to edit this new default configuration item by item in the same manner as described above, if desired, or you can accept this complete method as the default and proceed.

At this point all the global Method defaults have been defined - the next step is to set the default Method Events table. Before proceeding, you may wish to reset the Events button color or label at bottom right of the top level Method dialog by clicking the right mouse button for the pop-up button edit menu as already described. When finished, click Events to see the default Events Table screen (Figure 5.47).

Bottom Level Method (Events Table) Setup

The Events Table determines what chromatography system parameters you will be able to control from Data Ally in real time while executing your Methods. You must include all columns in the table which correspond to features or elements you will wish to control directly.

The normal Events Table is a split-screen window with a table in its bottom portion and the Graph window visible in the top portion (Figure 5.47). This arrangement permits you to actually load a chromatogram file (one file only) into the Graph window to assist you in properly programming the default method, using all the same table editing tools normally available in Data Ally. Click anywhere inside the Events Table window to make it active.




Figure 5.47 Method Events Table Configuration Window

Only a few columns are shown as "preliminary defaults" in the table, including the "Time" column and "Comment". At this stage, you can add or remove columns from the table, you can change the sizes/widths of individual columns as they are added, and you can change the order of the individual columns in the table.

To begin, you can review the available selection of column types based on the current Components configuration by placing the cursor pointer on any of the Events table column header buttons and clicking the RIGHT mouse button - the table button edit pop-up submenu will appear (Figure 5.48).




Figure 5.48 Table Button Edit Pop-up

You can use this button to reset the color or button label/title for each specific button in the same manner as previously described. To view the complete selection of available buttons, click on Enable/Disable in the poppa. The Configure Events Table Columns dialog will appear (Figure 5.49).




Figure 5.49 Configure Events Table Columns Dialog

The dialog shows a listing of all the elements which have been set in the current Components configuration which are appropriate for inclusion in the Method Events Table. Those currently in the table are highlighted. To add new columns, click on the names of the desired column titles to highlight them in the listing; to remove now-highlighted column titles, click them again to un-highlight. When only the desired column items are highlighted, click OK, or Cancel to retain the previous columns. Since at least one Input Flag for remote injection start was configured for this time base, you would click to highlight "Input Flags" in the list to add this parameter as a new column for timed event control. Upon clicking OK, the table will be restructured and the new columns will appear (Figure 5.50).




Figure 5.50 Events Table With Inputs Column

You can repeat using the Enable/Disable command to add or delete other columns, or return to Components configuration to make new options available for events, such as pump control, depending upon the overall Data Ally configuration in this time base.

Next, you should adjust the size (width) and order of any column(s) in the table in a manner which is optimized for your applications. Unlike the addition or removal of columns, which can be performed only in configuration mode, column sizes and order can be modified at any time while running Data Ally. If you will generally perform the same types of applications in this time base, it may be convenient to set the best column size and order to avoid the need for further adjustments later on - column order and width settings associated with the default method will always be retained when Data Ally is started.

To change column order, place the mouse pointer on any column header button to be re-positioned in the table and hold down the left mouse button; now "drag" the column to its new position and "drop" it by releasing the left mouse button. The table will immediately be redrawn with the column in its new position. Note that the pointer cursor will change while a column is being dragged to help you determine in what position you will be "dropping" the column.

To resize any column, move the pointer cursor to the vertical boundary between any two column header buttons until the pointer changes to a "left/right arrow" cursor. Now hold down the left mouse button and move the special resize cursor to the left or right to expand or reduce the width of the column selected - when the desired size is reached, release the left mouse button and the table will be redrawn showing the column with its new width setting. You can use these tools repeatedly to set the order and width of any columns in the default method. Note that the reorder and resize tools are also available outside of Configuration mode for modifying the table columns at any time during operation.

If the existing names shown on all the column buttons are not satisfactory, or you would like to assign alternative names or descriptions on each button, click the RIGHT mouse button with the pointer on any button to see the button edit pop-up (Figure 5.48), and then click on Set Title in the pop-up menu. The Set Button Title dialog will appear (Figure 5.51).




Figure 5.51 Set Button Title Dialog

Enter the modified title on one or two lines as done previously in the Graph screen. Press CTRL and ENTER at any time to start a second line. The title text can be of any length up to 256 characters, but only as many characters will be shown in the table as can fit inside the width allotted for each button. If you wish to decrease the size of a given column, for example, so that more columns will fit on the usual single-page Events Table window, you may want to use retitling to enter a more abbreviated title appropriate for the narrower column needed.

If you want to change the title color, click the TITLE COLOR button at the bottom of the dialog to see the Button Title Color Select dialog (Figure 5.52).




Figure 5.52 Button Title Color Select Dialog

Again, a new color can be applied to any individual button as above.

Each column header button can be titled and colored independently of all other buttons using the button edit pop-up.

If you wish, you can also change the colors and titles of the COPY, CUT, PASTE, INSERT, DELETE, FILL, HOLD, CLEAR, and METHOD buttons below the table, by using the RIGHT mouse button to display the same Button Edit Pop-up dialog. You can change the relative order of these buttons within each group of three on the screen by "drag and drop" in the same manner as reordering the column header buttons.

You can now proceed, if you wish, to enter an actual Events Table program for the default Method file. This could be as simple as a one-line program with hold-for-inject input flag S programmed at 1 sec, or as elaborate as a complete "standard" pump and detector control protocol with timed integration codes. You can type directly into table cells or use the spreadsheet or FILL functions for cells or ranges. The Clear Method command in the Default Method File submenu can also be used to clear all existing Events Table entries. See Section 8.3 below for a complete explanation of Events Table programming.

Any program entered in the default method will always appear as part of this file automatically when the current time base is initialized or when another loaded method file is cleared. Therefore, it is recommended that only Events table parameters which are very commonly used in applications be entered as defaults, since they will have to be cleared for programming or running alternative methods. For most installations, it will be most appropriate to leave the default method events table blank (empty) after its configuration has been set.

You have now completed basic method default configuration for this time base, and can proceed to configuring the Integration programming screen.

5.2.6.3 Integration Default Screen Setup

Click on the INTEGRATION SoftButton to set the Integration Configuration Default Window (Figure 5.53).




Figure 5.53 Integration Configuration Default Window

This program screen is used to set the default or basic peak integration algorithm and to specify peak identification options and the format of the peak identification table. Like the Method Default Configuration screen, it contains two levels - a "top" or "global" window, and a Peak Table window, which are set up independently.

Top Level Screen Setup

Default Filename: At top left, click inside the NAME field box and enter the filename (eight characters maximum) for the default Integration program file, if you do not wish to use the already assigned name "default.int". The name extension ".int" will be added automatically, which designates all Data Ally Integration files. Alternatively, you can use the File/Assign Integration command in the top command menu bar to use an existing Integration file already on disk as the new default file in this time base (see below).

"Algorithm Parameters:" This boxed area at left contains a series of fields comprising the peak-finding algorithms used for automatic chromatogram integration. When shipped, the initial "default.int" file includes values in each of these fields which may or may not be suitable for the types of chromatograms or applications to be run in this time base. You can click inside any of these fields and type in a new value with or without a new units designator - if the units designator is different from that shown in the gray box to the right of each field, the gray box will change to reflect the new units when you press ENTER or click inside another field after completing your entry. Note that if you have already entered or modified the integration algorithm fields in the Graph Edit/Integrate dialog in Section 5.3 above, those default values will automatically appear in this configuration screen as well; the two locations in the software are interconnected.

It is normally suggested that you use the default values supplied for these fields when beginning to use Data Ally, until you acquire some representative chromatograms and gain a better understanding of the operation of the integration algorithm. However, if you have already used Data Ally or have other Data Ally units on hand, you may wish to transfer a known integration algorithm parameter set from another existing program file by using the Assign command (see below).

To enter a description for the default Integration program file, click the File command on the top command menu bar to see the File submenu (Figure 5.54).




Figure 5.54 Integration Default File Submenu

Click on Integration Description in the submenu. A Description Entry dialog box will appear (Figure 5.55); type in the desired description of any length on one or more lines and click OK when finished.




Figure 5.55 Integration Default Description Dialog


Assigning an Existing Integration File As Default: Instead of creating a new name and completely entering a default integration file, you can "re-assign" an existing Data Ally integration file and its file configuration as the default/configuration for the current time base. An assigned file will replace any existing file shown in the Configuration Integration Window and will substitute the new file's parameters and configuration in place of the old.

To do this, click on File in the top command menu bar for the File submenu (Figure 5.54) and select Assign Integration from the pulldown submenu in the same manner as described for Method files above - the Integration Default Assignment dialog will appear (Figure 5.56).



The color and/or title of the PEAK TABLE button at bottom right can be modified by using the RIGHT mouse button to display the Button Edit Pop-up.

When the "top level" integration configuration is completed, proceed to the Peak Table by clicking the button marked "Peak Table" at lower right. The Configuration Peak Table Window will be displayed (Figure 5.57).



In the same manner as for the Events Table, you should begin by determining if all the column items needed for your applications in this time base are present. With the pointer on any table header button, click the RIGHT mouse button to see the Table Button Edit Pop-up, exactly the same as that for the Events Table. Click Enable/Disable to view the Configure Peak Table Columns dialog (Figure 5.58).




Figure 5.58 Configure Peak Table Columns Dialog

The list box in the dialog shows the complete set of column elements available for this table - click to highlight all the desired items and to un-highlight/defeat the unneeded items. When all the required items are highlighted, click OK, or click CANCEL to restore the previous selection.

The table will immediately be redrawn showing only the columns selected. You can now use the resize and drag-and-drop reordering tools to reformat the column buttons as you wish, in the same manner as for the Method Events Table.

Now, use the Button Edit Pop-up to reset the text labels, text colors, or button colors for any or all column buttons, in the same manner as for the Events Table described above.

If you wish to color or rename the spreadsheet command buttons (COPY, CUT, PASTE, INSERT, DELETE, FILL), HOLD/CLEAR, or INTEGRATION buttons beneath the Peak Table, use the Button Edit Pop-up to do so by placing the pointer on the button and clicking the RIGHT mouse button to access the editing poppa menu.

If you wish, you can now enter peak data into the peak table either by manual entry or by automatic filling from a chromatogram loaded into the Graph window above the Peak Table (see Section 8 for instructions on creating a Peak Table). The spreadsheet or FILL functions can be used to edit or insert new table entries. Any such entries will be saved as the default and will always appear in the peak table when the time base is initialized. Therefore, if you do not perform applications that nearly always use the same Peak Table information, it is suggested that you leave the Peak Table blank/empty in the default Integration file. The Clear Integration command in the Default Integration File submenu can also be used to clear any existing entries from the Peak Table.

When all Integration screen and Peak Table elements are set as desired, proceed with configuring the next program screen.

5.2.6.4 Calibration Default Setup

Click on the CALIBRATE SoftButton to display the Configuration Calibration Window (Figure 5.59).




Figure 5.59 Configuration Calibration Window

The Calibration Default screen contains two general levels, a "top" or "global" window which sets overall options and defaults, and a lower "table" level which accesses and sets formats for both types of tables incorporated into the normal Calibration window.

"Top Level" Calibration Functions

Default Calibration Filename: Click inside the NAME field and type in a new filename to be assigned to the default Calibration file which is loaded when the current time base is initialized, or use the standard "default.cal" name. The name extension ".cal" will be added automatically. Alternatively, you can use the File/Assign function to denote an existing calibration file as the default file and format.

On the left side of the screen, in the area marked "Standard Type", the types of calibrations available for the current time base are specified by clicking on the check boxes to check the options which are to be enabled. Click also on the "Def" button corresponding to the option to be defined as the default. If no internal standard calibrations are to be performed for any applications in this time base, simply leave the "Internal" type box unchecked and this option will never appear in any calibration screens.

In the area marked "Basis" are the two options - area or height - for quantitating peaks. Click to check which options are to be available, and on the "Def" button to assign the default condition if more than one selection is indicated.

In the area marked "Curve Type" click on all the curve fitting options to be accessible in the current time base for multiple-point standard computations. Also click on the single option to be set as the default (most often Linear) which will be used in the absence of any manual modifications.

In the area marked "Correction Factors", check the checkboxes of each of the factor types you wish to be active in this time base. These checkboxes will enable the use of each type of factor in setup of your Method and Sequence programs, as well as the manual Analyze and Calibrate dialogs in the Graph screen, for this time base. IF A FACTOR IS NOT ENABLED HERE, IT WILL NOT BE AVAILABLE FOR USE IN ANY MANNER IN THE CURRENT TIME BASE.

NOTE, AS FOR THE OTHER DEFAULT SETUP SCREENS, THAT THOSE OPTIONS NOT CHECKMARKED FOR CALIBRATION WILL NOT BE DISPLAYED NOR AVAILABLE FOR USE WITH ANY CALIBRATION PERFORMED USING THIS TIME BASE CONFIGURATION.

To enter a description for the default calibration file, click File in the top command menu bar to see the File submenu (figure 5.60).




Figure 5.60 Calibration Default File Submenu

Click on Calibration Description and enter the desired text description of any length into the Description dialog entry field - click OK when complete.

To assign an existing Calibration file on disk as the default file for the current time base, click File in the top command menu bar, and select Assign File to see the Assign Default Calibration File dialog (Figure 5.61).



Figure 5.61 Assign Default Calibration File Dialog

Choose a Calibration file from any disk or directory using the selector dialog, and click OK - the selected file will be loaded into the default register and its configuration will update the previous configuration for calibration.

If you like, you can use the Button Edit pop-up accessed with the RIGHT mouse button to change the label or color of the Calib. Tables button at bottom right - however, this button will not appear in the operating software and exists in Configuration mode only.

Calibration Tables

Now proceed to set up the appropriate calibration tables by clicking the CALIBRATION TABLES button at bottom right to view the Calibration Tables Configuration Window (Figure 5.62)




Figure 5.62 Calibration Tables Configuration Window

This screen is very similar in appearance to the actual "Calibration Screen" in the Data Ally operating software - see Section 9 below - and is used to specify the contents of the Calibration Reference table (top right), which sets standard concentrations, basis, curve fitting, and other parameters for each peak, and of the Calibration Data Table (bottom), which displays desired information about all standards which have been run and processed. Before structuring both tables, it is necessary to perform the Default Calibration Setup using the File/Setup command.

Click File in the top menu bar to see the Calibration Default File Submenu (Figure 5.63).




Figure 5.63 Calibration Default File Submenu

Now click Setup to view the Calibration Default Setup Dialog (Figure 5.64).




Figure 5.64 Calibration Default Setup Dialog

To create the Default calibration Reference Table, a valid Integration program file saved on disk - usually the file designated as the default integration above - is designated as the "model" for defining the required Peak Table. As for all calibrations, defining the relationship between the default calibration file and an integration file establishes a permanent link between the two files, so that whenever the "model" Integration file's peak table is modified, the same modification will automatically occur in the default Calibration file. Click on the arrow selector to see a listing of all existing Integration files, which could be used as models, and select one from the list, or simply type the name of the selected model file.

Now set the number of calibration concentration or standard levels for the default by either clicking inside the "Standard Levels" field and typing the desired number (between 1 and 16) or clicking the arrow selector and choosing the desired levels from the listing. When both the Integration file reference and number of standards are shown as desired, click OK to close the dialog and proceed. the Data Ally Configuration System now has sufficient information to allow you to fully define the Reference Table.

If the Integration file referenced in Setup actually contains a defined Peak Table, the data from this table will now appear in the Calibration Reference Table. It is perfectly legitimate to have actual calibration data saved in the default Calibration file for general use, although it will be more common to have only the format of a calibration saved as the default file, and to create new Calibration files with specified names by using the Default Calibration as a "template" only.

Calibration Reference Table: To review and configure this table, first click anywhere inside this table window to make it active. Then place the pointer on any of the column header buttons in the table and click the RIGHT mouse button to display the Header Button Edit Pop-up, which is employed in exactly the same manner as the pop-ups in the Events and Peak Tables to Enable/Disable table column items, and change button text labels and colors. Click on Enable/Disable in the pop-up to see the list selector dialog for the Calibration Reference Table (Figure 5.65).




Figure 5.65 Configure Calibration Reference Table Columns Dialog

The RT, Peak No., and Peak Name columns are "imported" from a corresponding Peak Table screen at the time a calibration file is established. You can select which of these items to show in the Reference Table by enabling (highlighting) or disabling them. A series of "Level" columns are also available, according to the maximum number of levels you have permitted in Setup.

You should select only those column functions in the Reference Table that you will wish to use to change the way calibration data are interpreted for your applications. Columns representing those parameters which will always remain fixed should normally be omitted by disabling them; they can always be restored later on by re-configuring. As an example, if your quantitative methods will always employ peak area rather than height measurements, and if linear regression curve fitting will always be used for single level calibration, you should disable the Basis and Fit columns and all the Level columns except for Level 1; in the top level calibration configuration you should check Area Basis and Linear Regression fit as the only options. This will "lock" all calibrations performed with this configuration for these settings.

When selecting the number of Calibration Level columns, set the maximum number of different standard concentration levels you ever anticipate for any application using this time base. Up to 16 different level columns can be chosen. The only drawback to configuring many level columns is that this will tend to slow down redraws of the Reference table window. Again, if you configure too few levels now to accommodate a future application, you can always re-configure the time base later when necessary.

NOTE: IN ORDER TO COMPLETE CONFIGURATION OF THE REFERENCE TABLE YOU MUST SELECT AT LEAST ONE LEVEL COLUMN. Otherwise an error message will appear when you attempt to proceed with further configuration or to exit the configuration mode.

Highlight all the column items to be enabled and click OK to redraw the Reference Table. You can now edit the column header button titles, colors, and column order using the Button Edit Pop-up dialog and the button drag-and-drop tools in exactly the same manner as for the other tables.

You can type information directly into the Reference Table cells, modify peak table information already present there, or use the spreadsheet or FILL functions to complete cell entries. See Section 9.0 below for a complete explanation of how calibration reference tables are created and utilized.

Calibration Data Point Table:

The final step in configuring Calibration files is to set the contents of the Data table. Click anywhere inside the Calibration Data Table window in the bottom portion of the screen to make this window

active. To review the current contents and available fields for the table, put the pointer on any of the column header buttons, click the RIGHT mouse button to see the Button Edit Pop-up, and click Enable/Disable to view the Configuration Calibration Point Table dialog (Figure 5.66).




Figure 5.66 Column Selector Dialog for Calibration Data Table

This table contains selected information about each standard sample used to calibrate the system, providing a mechanism for reviewing pertinent data during routine calibrations. You should select any or all of the parameters in the listing which will be of value in monitoring the quality of calibrations for your applications, and which will help you establish an adequate historical record of series of calibrations as they are performed and logged. There are a large number of descriptive parameters available and the table can be as wide or as large as needed or desired.

Parameter items which can be included in this table display include:

- Area % - Area Percent
- Height % - Height Percent
- Level - Calibration Level
- Dil - Dilution
- Int. File - Integration Filename
- IS Area - Internal Standard Area
- Method - Method
- IS Conc. - Internal Standard Concentration
- Cg Filename - Chromatogram Data Filename
- IS Height - Internal Standard Height
- Date - Date of Calibration
- Time - Time of Calibration
- Area - Area
- KC - Polynomial Coefficients
- KB - Polynomial Coefficients
- KA - Polynomial Coefficients
- Conc. - Concentration
- Height - Height
- RF - Response Factor
- RT - Retention Time
- Run Date - Run Date
- Run Time - Run Time
- Smp. Desc. - Sample Description
- Smp. Name - Sample Name
- Smp. No. - Sample Number

Click on as many items as you wish to include in the table to highlight them, disable all unwanted items by un-highlighting them, and click OK when finished to redraw the table. You can now proceed to set new titles and colors for any or all of the column header buttons using the Button Edit Pop-up tools, and to reorder the columns if you wish using the drag-and-drop tool.

You can also change the colors and titles of the spreadsheet function buttons (CUT, COPY, PASTE, FILL, INS, DEL, CALIB) beneath the Reference and Data Tables using the Button Edit Pop-up, accessed with the RIGHT mouse button. Within each group of these buttons, you can also use the "drag-and-drop" tools to reorder the individual buttons.

If you wish, you can now load one or more chromatograms, integrate them in the Graph window, and use them to create a default calibration data table for one or more peaks. See Section 9 below for a detailed description of how calibrations are performed.

You can use the Clear Calibration command in the Default Calibration File submenu to remove all existing entries from the Reference or Data tables at any time.

The Calibration defaults should now be completely configured - click on OK to leave the Calibration defaults screen.

5.2.6.5 Sequence Default Screen Setup

Click on the SEQUENCE SoftButton to view the Configuration Sequence Window (Figure 5.67). This screen configures all parameters related to batch analysis and reprocessing, and the use of an autosampler. The contents of this screen may vary depending upon which, if any, autosampler has been designated in the Components list.




Figure 5.67 Configuration Sequence Window

Default Filename: Click inside the NAME field and type the desired filename for the default Sequence file, if "default.seq" is not satisfactory. The name extension ".seq" will be added automatically. Alternatively, use the File/Assign function in the top command bar menu to assign an existing Sequence file on disk as the new default file and format, in the same manner as described above for the other default program file types.

Inside the area marked "Process Mode", click on all the options to be available in Sequences in this time base, and click the "Def" box to mark which mode will be the default. The Acquire mode is for capturing live run data, and is necessary for connection to chromatographic instruments. If this

channel is to be used as an edit-only channel, Acquire mode should not be checked, since no data acquisition can be performed - only Reprocess mode can be used.

In the area marked "Position", you can set the defaults for the number of rows to be processed in the batch table by clicking inside the Start or End entry fields and typing a row position number. Regardless of how many rows are actually entered into the batch table below, these entries will force execution of the batch within the defined range of rows unless changed manually before starting the Sequence.

In the area marked "Options", click on the check boxes for "Suppress Integration" or "Suppress Calibration" to set these functions on or off for the default sequence file. These check boxes override any specific integration or calibration instructions entered in the methods being run in a sequence or in the sequence table itself.

Click the "Validate" checkbox to specify sequence table validation as the default. With Validation enabled, all the parameters specified in the sequence table (existence of methods, integration files, calibration files, etc.) will be verified by Data Ally when a Sequence is started before actually attempting to perform the first method, to be sure all the required information is available to perform the entire sequence of operations as stated.

To enter a description for the default Sequence file, click File on the top command menu bar to view the File Default submenu (Figure 5.68).




Figure 5.68 Sequence File Default Submenu

Click Sequence Description to view the Description dialog box (Figure 5.69). Type the text description into the entry box in the Description dialog and click OK when finished.




Figure 5.69 Sequence Default File Description Dialog

To assign an existing Sequence file as the default file, click File in the top command menu bar and click Assign File. The selector dialog will appear (Figure 5.70); choose a Sequence file from disk or other device and click OK. The indicated Sequence file will be loaded into the default register and will be used to reset the default configuration.




Figure 5.70 Assigning Sequence Default File

Sequence Table:

The Sequence or Batch table at the bottom of the screen can include any series of column functions for complete automated processing of single sample injections or batches of injections.

Begin by reviewing the current contents and available columns for this table; click inside the table window to make it active, place the pointer on any column header button, and click the RIGHT mouse button to see the Button Edit Pop-up. Click Enable/Disable to view the Configure Sequence Table Columns Dialog (Figure 5.71).




Figure 5.71 Configure Sequence Table Columns Dialog

You can now choose all desired column elements for this table by highlighting those entries and un-highlighting unwanted items. Note that any columns placed in the Sequence table can contain run instructions which will override the same parameters specified in Method files called by the sequence. You should include all parameters as sequence table columns you will need if you are routinely setting up batches of runs and entering sample-specific data, random combinations of methods, integration files, and calibration instructions, or sample correction factors such as IS amounts, dilutions, or injection volumes. As with all Data Ally tables, it is advised that only those columns which will be useful for your anticipated applications be included in the table configuration.

The following column parameters are available for the Sequence Table:

- # Injs. - Number of Injections
- No. - Vial Number
- Method - Method Filename
- [Det. A] Cg File - Chromatogram Filename from Detector A input
- [Det. A] Int. File - Integration Filename for Detector A input
- [Det. A] Cal. File - Calibration Filename for Detector A input
- [Det. A] Rpt. Templ. - Sample Report Template File for Detector A input
- [Det. A] Cal. Hist. Templ. - Calibration History Report Template for Detector A input
- Cal. Level - Calibration Level
- Cal. Code - Calibration Code
- IS1 - Internal Standard Amount
- [Det. A] Rpt. File - Sample Report Filename for Detector A input
- PrtCg - Print Chromatogram Command
- Dil. - Dilution
- Smp. Name - Sample Name
- Smp. No. - Sample Number
- Smp. Desc. - Sample Description
- Comment - Comment

When all desired items are highlighted, click OK to proceed and redraw the table.

You can now use the Button Edit Pop-up tools to retile or recolor any of the column buttons, or the drag-and-drop tool to change the column order and relative widths.

You can also use the Button Edit Pop-up to change the spreadsheet function buttons (COPY, PASTE, CUT, FILL, INSERT, DELETE), HOLD/CLEAR, and screen selector button (GRAPH) below the table. The order of these buttons can be changed within each of the three functional groupings by dragging-and-dropping.

If appropriate, you can now enter information directly into the Sequence table to create a default table. Entries can be made by typing directly into cells or by using the spreadsheet or FILL functions. This is recommended only if a "generic" batch listing will be used to run most applications on this time base. See Section 11 below for a complete explanation of the use of Sequence functions. You can use the Clear Sequence command in the Sequence Default File submenu to clear all existing entries from the default sequence table.

The Sequence defaults should now be completely defined. Click on OK to finish.

5.2.6.6 Retitling, Reordering, and Coloring Program SoftButtons

If you have not already done so, like all other Data Ally screen elements, the SoftButtons displayed on the left screen margin button bar can now be individually renamed and recolored and their order changed as desired. With the mouse pointer on the desired SoftButton, click the RIGHT mouse button to see the Button Edit Pop-up; make changes as for any other system buttons above. You can also drag and drop the SoftButtons to change their relative positions along the height of the button bar.

In addition to changing the SoftButtons, it is also possible to change the information in and format of the "Last Changed" field at the bottom of each screen type. To do this, place the pointer on the "Last Changed" button and click the RIGHT mouse button - the Button Edit Pop-up will appear.

5.2.6.7 Completing Configuration

If all Components specifications are correct and all the default program screens are set up as desired, you are ready to finish Configuration of the current time base. You may wish to use the Print Configuration commands for any or all of the program default screens (found in their respective File submenus) to print a hard copy record of each part of the overall default configuration.

If you wish to set up the configuration for the edit time base associated with the live time base you have just completed, click on the "Edit TB" SoftButton for the edit time base and review all of its program screens, making appropriate changes. **Note that the Method screen parameters and default values must be identical for each pair of live and edit time bases - however, all other screens can be varied markedly by means of configuration.**

Once you have reviewed the edit time base configuration, proceed to any other live or edit time bases available by clicking their respective left-column SoftButtons and then clicking the Util SoftButton to begin setup of each new time base as described above.

ALL TIME BASES PRESENT MUST BE COMPLETELY CONFIGURED BEFORE EXITING AND SAVING THE CONFIGURATION FILE. Only one Configuration File sets the characteristics of the entire Data Ally software regardless of the number of time bases present.

To finish configuration, click the UTIL SoftButton in any time base window to view the Utilities Configuration screen. Click Files in its Command Menu bar to see the submenu (Figure 5.72).




Figure 5.72 Utilities/Files Configuration Submenu

To save the new configuration as part of the current configuration file (the name of which is shown on the top bar) click Save Configuration. To save the new configuration as a new and unique configuration file, click Save Configuration As to view the Save Configuration dialog (Figure 5.73).




Figure 5.73 Save As Configuration Dialog

Enter a description of the current configuration by clicking inside the Description box and typing the desired entry.

5.2.6.8 Completing Time Base Configuration

When all Configuration settings and default files are shown as desired for the current time base (in other words, the desired configuration file is currently loaded), you can exit Configuration to begin using Data Ally with any new configuration changes, or you can restore the previous or any other saved configuration.

If you wish to load a previously saved Configuration file, either to edit it or to set up the Data Ally software as it specifies, click on Files/Load in the UTIL screen to view the Load Configuration dialog (Figure 5.74)




Figure 5.74 Load Configuration Dialog

Select the desired configuration file and click OK to load it, or click CANCEL. You can now edit this file in the same manner as you have defined a configuration above.

To exit Configuration Mode:

In the UTIL screen, click the File command on the top menu bar to see the Files submenu (Figure 5.72). Click on Exit Configuration - the Exiting Configuration dialog will appear (Figure 5.75).




Figure 5.75 Exiting Configuration Dialog

This dialog specifies whether or not the CURRENTLY LOADED Configuration file will become the new "default" configuration for Data Ally, or whether the last-used configuration will remain as the default.

Click "Yes" to make the current configuration file - which you have loaded and/or edited - the new "default" configuration and to initialize Data Ally using that configuration.

Click "No" to revert to the previously-used configuration file, located in any disk directory path, and initialize Data Ally using that configuration.

If you need to make additional configuration changes before initializing Data Ally, click CANCEL to return to configuration mode.

After clicking Yes or No, the Data Ally software will be restarted according to the instructions in the selected configuration file. When the hourglass cursor disappears, all time bases will be ready for use. You can proceed to write programs, or run methods directly using the default program files you have created.

Any Configuration file can be used at any time to reset and reinitialize the Data Ally software simply by entering configuration mode, loading the desired configuration file or modifying any configuration file as needed, and using the Exit Configuration command followed by a "Yes" response in the Exiting Configuration dialog.

5.3.0 Configuration Files

When you enter Configuration Mode, the Data Ally software changes so that only the contents of Configuration files are displayed and available for editing. Each time the Configure command is used in the Utilities screen, the currently-applied Configuration file is loaded for viewing in Configuration Mode. If you wish to review and edit a Configuration file which is not currently in use, you can easily load any *.CFG file on any existing disk or device. You can also save the complete contents of Configuration Mode to overwrite an existing configuration file with new information, or as a completely new file, in any location on disk.

It is recommended that you create and archive as many different Configuration files as you will need to properly set up all applications you may perform with your Data Ally system. These files should be given unique names and descriptions so that they can be located and identified easily, and should be copied onto floppy disks or some other storage media for secure backup archiving purposes.

REMEMBER THAT THE ONLY WAY THE SOFTWARE CONFIGURATION CAN BE CHANGED IS BY ENTERING CONFIGURATION MODE, SELECTING A CONFIGURATION SO THAT IT APPEARS IN CONFIGURATION MODE, AND EXITING CONFIGURATION MODE. The act of exiting Configuration Mode always automatically applies the Configuration currently loaded and viewed in Configuration Mode. To change configurations, you must first load the desired new Configuration file in Configuration mode and then exit Configuration to apply that file to the system. The Configuration file used to set up the system can be located in any physical path on disk.

5.3.1 Saving Configuration Files

To save the current Configuration on which you have been working under its previous name, click on Files on the top Command Menu bar, to see the Files submenu (Figure 5.72). Click on Save Config on the submenu to resave the current configuration, overwriting the old file.

To save the current Configuration under a new name, click on Files on the top command menu bar, and click on Save As - the Save As dialog box will appear (Figure 5.73).

Enter the desired name, or modify an existing filename, in the dialog box. When the desired name is shown, click OK to save a new file with that name.

5.3.2 Loading Configuration Files

The last Configuration file loaded in Data Ally is assumed to be the file which should be loaded the next time the software is initialized.

In Configuration Mode, you can load any existing *.cfg file for review or modification, or to implement that file to change the system configuration.

To load a configuration file, click File on the top command menu bar, and then click Load Config in the submenu. The Files/Load dialog will appear (Figure 5.74).

Select a drive, path, and filename for a *.cfg file to load, and click OK. The selected file will be loaded and displayed in all the screens in Configuration Mode.

5.3.3 Deleting Configuration Files

You can delete Configuration files which are no longer useful from Configuration Mode by using the Windows File Manager - there is no direct Delete command within the configuration files dialog.

6.0 Getting Started

Data Ally Chromatography Manager is a multitasking data acquisition, management, and control system for use with computers running the Microsoft Windows 3.1 or Windows 95 graphical user interface. One time base is available with two inputs for detectors. Each detector can have its own unique "Configuration" based on the analytical instrumentation assigned to that channel and its specific applications.

6.1.0. Operating Structure and Programs

Data Ally runs one time base with two detectors simultaneously. Each uses the same standard operating structure, customizing changes made via the unique Configurations of the separate detectors. The row of SoftButtons at the left of all Data Ally display includes a number of selector buttons used to set which time base is currently being observed and programmed.

Data Ally user display interface utilizes a set of six basic, distinct screens for visual presentation and editing of all raw data, calibration, report, and program files. Each of these screens shows a "snapshot" of information currently in memory which is relevant to the particular part of the system that screen represents.

The best way to understand the Data Ally user interface is to think of the display as a "viewer" which focuses on only one part of the overall program. Using the mouse to click on any of the buttons in the left screen global column selects an alternate view of a current program.

In addition to the program screen and time base selector buttons, a HELP button is available to show on-screen help messages, and a RUN/STOP button is used to start or terminate runs in the currently-observed time base.

Data Ally organizes information into the following screen display types:

GRAPH Screen:	Shows plots of chromatogram files while system is running or for editing purposes. Provides all raw data file handling, display, and comparison functions.
METHOD Screen:	Shows system instructions for running a Method to capture chromatograms, save and analyze them, and make reports. Also contains all events instructions for instrument control during the course of the Method.
INTEGRATION Screen:	Shows system instructions for integrating a chromatogram and identifying its peaks. Also contains calibration information regarding each defined peak.
CALIBRATION Screen:	Shows all current calibration information for a loaded chromatogram, including a summary of selected data from all standard injections and a graphic plot of the calibration response curve.
SEQUENCE Screen:	Shows system instructions for running a batch analysis using an autosampler, or for reprocessing batches of chromatograms. Selects order of injections and sets all control parameters for an autosampler.
UTILITIES Screen:	Accesses configuration settings.

Each screen type has its own associated blue "SoftButton" drawn in the left edge area of every Data Ally display. Clicking on a SoftButton at any point in the software will display the selected screen type in the current time base.

The contents of each of the principal display screens can be modified by using Data Ally's Configuration functions, for complete customization of each time base in the system.

Data Ally is organized into "full screen" or "maximized" time base windows, each of which contains a family of "program windows" at all times. Changing time base windows by means of the time base SoftButtons makes the contents of all of the program windows in the new time base immediately accessible to the operator by clicking their respective SoftButtons.

Each Data Ally time base window can be "minimized" or "iconized" while it is idle or actively running. Iconizing a channel window captures all its currently loaded program windows, and effectively removes the time base from the current display, without interfering in any activities proceeding in that time base. It is possible to iconize all or several time base windows at once and to "remaximize" only the time base window you wish to review. Whether any time base windows are iconized or not, as a Microsoft Windows application, Data Ally permits any other software applications to be running concurrently.

You cannot "iconize" individual program windows, since each is associated with a time base window at all times. However, you can "maximize" program windows in the same manner as time base windows, if you wish to see a full-screen display of the information in any such window.

In program windows where "subwindows" appear on the display, such as the Method Events Table/Graph screen, you can change the relative sizes and positions of any of the subwindows in any way you like, within the space defined by the time base window.

6.2.0 Files and File Manipulations

The Data Ally system utilizes and creates DOS-compatible files. The following file types can be produced, and edited, within the data system:

- **Chromatogram Raw Data Files (*.PDF):** These files contain the set of detector input signal versus time information corresponding to a chromatogram run, one chromatogram per file. Each raw data file also contains header information regarding the date, time, and operator for the file, the time base and programs with which the file was created, all the current integration data, and, in some cases, calibration information.

Raw data files can be loaded into any Data Ally time base and processed at any time. Raw data files can also be as long as necessary. The only limitation to file length and complexity is the available disk space when a file is originally saved.

- **Program Files (Method - *.MET, Integration - *.INT, and Sequence - *.SEQ):** These files contain all the information describing each program type. Program files can be of variable size depending on the length of the Table portion of each file type.
- **Calibration Files (*.CAL):** These files contain all the calibration data associated with a specific set of standards, referencing a particular Integration program file and its Peak Table. Calibration files can be of variable length depending upon the number of standards processed.
- **Report Template Files (*.XLS):** These files are produced with Microsoft Excel (in Data Ally v. 1.00), and contain a set of instructions for formatting and creating a single report type. Template files do not actually include any report or results data.

- **Report Files:** These files are produced using Report Templates and contain results and graphics describing one or more chromatograms. Report files can be reviewed and edited within the data system at any time, and can be exported to other applications or computers.

Microsoft Excel is used via DDE links to create reports. For Excel, *.XLS files are employed as report templates. These are loaded into a special subdirectory of the main \Data Ally directory called \REPORTS. Excel can be used directly to edit these *.XLS report template files (See Section 13).

All Data Ally file types can be listed, copied, edited, and deleted using standard DOS and Windows commands. Any type of file can be transferred to another application or another computer. You can translate raw data files and report files into different file formats using the commands in the Utilities screen, if this is desired. Chromatogram files produced by Axxiom 700 Series Data Systems can be easily translated directly into Data Ally formats.

6.3.0 Entering the Data Ally Software

Data Ally can be started from Microsoft Windows. Once the software is initialized, it is loaded into memory and all configured systems are available for running Methods or processing data.

☞ *To start Data Ally from Windows:*

Double-click on the Program Manager icon to see the Program Manager window. Select the Data Ally icon from the Applications menu by double-clicking on that icon. The Data Ally software will load and you will see the Graph screen.

6.4.0 Using the Default System

When you have logged on and obtained entry into Data Ally, you can immediately start collecting chromatograms or for writing programs.

6.4.1 Starting a Run

☞ *To start collecting signal data from any analog detector(s) connected to the current time base*, click the left mouse button with the pointer on the RUN SoftButton, while viewing any screen. If you have not loaded a Method or Sequence file, the default file(s) will be employed as set in configuration. You may wish to click on the METHOD SoftButton to see the DURATION (run time) and RESOLUTION (sampling rate) settings for the current Method before beginning the run. Click inside either of these fields in the Method to change the parameter(s) for the run to be started. You may also wish to change the file name(s) assigned for automatic file saving, or for automatic Integration, by changing the current Method as desired before RUN is executed.

Since Data Ally always logs chromatogram file data directly to the hard disk as it runs, a valid file name must be provided in the default Method file or assigned in any modified Methods before a run start can occur. If the filename(s) currently found in the active detector input entry field(s) in the Method screen is already assigned to any file in the current directory path at the time the RUN command is given, an error message will appear advising you to change the filename(s) or path before a run can be initiated.

If you click on RUN while viewing the Sequence screen, the time base will attempt to run the default Sequence, if any exists.

☞ *To see the incoming signal data traces*, click on the GRAPH SoftButton, if you are not already viewing the Graph screen. The trace(s) for each live detector will be visible on the display, scaled as predefined in the default instructions. Alternatively, you can click on the Graph button at the bottom right corner of the Sequence screen to replace the top-level portion of the Sequence screen with a half-screen Graph window while the Sequence Table remains visible at bottom.

☞ To reset the Display scaling in the Graph, click on Display in the Graph Command Menu bar, and then on Scaling in the submenu. Turn AutoScaling on or off by checking the appropriate check boxes in the scaling dialog box, or click inside the MIN or MAX entry fields for any live detector(s) with AutoScaling off and type in the range value you wish to display. Click on OK to change the display.

☞ To view the Method which is being run (the default Method, unless you have modified it or loaded a new Method), click on the METHOD SoftButton. The top level Method screen will appear. If you wish to determine whether any instrument control functions are in process, click on the EVENTS button to view the Method Events Table display.

☞ To stop the run at any time prior to the programmed DURATION of the current Method, click on the STOP SoftButton. The run will end and the chromatogram(s) will be integrated according to the instructions in the Method. The trace(s) will remain in the Graph screen, already having been saved to disk.

If you allow the current/default Method to run until its DURATION time is completed, the data file(s) acquired will automatically be saved using the filename(s) in the Method screen, and automatic integration will occur for all specified input files using the Integration program(s) referenced in the Method.

☞ To save the Method you have run, click on the METHOD SoftButton to view the Method screen containing the Method in this time base, and select Files/Save As from the Command Menu bar. Enter a new name for your Method and click on OK in the Save dialog box to save the Method for later use.

6.4.2 Analyzing Run Data Manually

If the new data file(s) are not integrated when the Method ends, or if the automatic integration is not completely satisfactory, you can integrate each file in turn manually.

☞ To integrate manually, click the GRAPH SoftButton and using the Graph/Edit/Integrate command. The integration algorithm parameters in the current Integration program are employed, and can be varied as desired.

You can modify the algorithm parameters within the Edit/Integration dialog box as you proceed in order to obtain the desired peak identification and resolution. Retention times will be computed either as absolutes or relative to the first defined peak, depending upon the setting in the current top-level Integration screen.

☞ To manually mark peaks or edit chromatogram baselines, you can use the Graph/Edit/Define Baseline commands, reintegrating after each manual modification.

When you have achieved suitable integration for the chromatogram(s), you can manually produce result reports.

☞ To make a manual report, click on the Graph/Edit/Analyze command. In the Analyze dialog box, select the type of report you want (e.g. Area PerCent) from the Report Library list box, and click on OK to produce the report. Quantitative result reports can be obtained if a Calibration file for the current chromatogram already exists.

6.4.3 Completing the Integration Peak Table

When your chromatogram is properly integrated and edited, and all desired retention times for all peaks are shown in the Graph window, you can easily complete an Integration program suitable for that chromatogram and save it for later use. An Integration Peak Table will be needed to create a Calibration file appropriate for your chromatograms.

➤ To create a suitable integration program, click on the INTEGRATE SoftButton to view the top-level Integration screen - it should now contain the final set of algorithm parameters you have just utilized. Click on the PEAK TABLE button to view the split-screen Peak Table display. The chromatogram will appear in the Graph window at the top of the screen.

➤ To fill the peak retention times into the peak table, click on the RT column header button in the Table, and then on the FILL button at bottom center beneath the table. The AutoFill dialog box will appear. Click on the "Fill From Graph" radio button in the top portion of the dialog, and click on the "Start of Table" and "End of Table" radio buttons in its center section to select filling of the entire peak list. Click on OK - the marked retention times for all the peaks in your chromatogram will be immediately filled into the table, one line or row per peak.

Now, click the pointer inside the top row of the WINDOW column (corresponding to the first peak). Enter a retention window value for the first peak position (a value of 0.1 or 0.2 min should be adequate if your peaks are well resolved). To use the same window value for all other peaks, click on the WINDOW column header button, and then on the "Duplicate" radio button in the AutoFill dialog. Click on OK to copy the first row window value into all other rows in the column. This is the fastest way to set up a new Peak Table; you can always return and set individual peak windows, or other information, later on. You may want to proceed and enter peak names for the identified peaks you will want to include in reports.

If you want to calibrate and quantitate certain peaks in this chromatogram, you will need to enter the type of calibration you wish to use in the CAL column of the new Peak Table. If all peaks are to be calibrated with External Standards, simply click the pointer inside the first row in the CAL column and type "ES". Then click on the FILL button, click on Duplicate, and click OK in the AutoFill dialog to fill the ES entry into all the peak rows.

If you want to run one or more internal standards for quantitation of your chromatogram, you must first identify the internal standard peak(s). To do this, find the first internal standard peak in the table and click the pointer inside its row in the IS column. Type the number "1" in that cell to designate that peak as Internal Standard No. 1. Assign any other internal standard peaks in the same manner by giving those peaks unique numbers (up to 5) in the IS column. Then, make entries of "ISx" for each non-internal standard peak to be computed using any internal standard peak x, either by directly typing those entries into the CAL column (instead of the ES entry for Ext. Standard). Be sure to include the correct Internal Standard Peak reference number, 1 through 5, to link each unknown peak to its corresponding internal standard peak.

You have now defined a basic Peak Table. If you like, you can now edit the table to add peak referencing, comments, and other instructions, depending upon what table options have been configured. You can also click on the INTEG button below the table to return to the top-level Integration screen, to define reports and other options.

➤ To save the Integration program, when both the top-level Integration display and its associated Peak Table are completed, click on the Files/Save As command to save the new Integration program for later use.

6.4.4 Setting Up a Calibration File

Now that you have defined a valid Integration file with Peak Table, you are ready to construct a Calibration file for quantitation.

➤ To create a new Calibration file, click on the Calibrate SoftButton to see the Calibration screen display. Click on Files on its command bar menu to see the Files submenu, and click on Setup. The Setup dialog box will appear.

Click on the arrow selector button to the right of the Reference Integration file entry field to see a listing of current Integration files in memory. The newly-saved Integration program with your new peak table should be shown in the listing. Click on the name of that program to insert it into the field. Now select the number of standard levels you will use to calibrate by clicking on the selector arrow for the No. of Levels entry field, and choosing the desired number - you must have at least one level defined in any Calibration file.

When completed, click on OK. The Calibration Reference table in the window at upper right will be immediately reconstructed to include the peak table from the selected Integration file, and the number of Level columns you have specified. Depending upon the configuration, other columns (FIT, BASIS, etc.) may also be found in this table.

For each peak to be calibrated in the Reference Table, you must now specify the concentration value for the standard(s) to be injected and used to build a Calibration Table. The current Calibration Table at the bottom of the screen will be empty, since no standards have been run yet.

Click anywhere inside the upper right Reference Table to activate it, and use the scroll bars to move to the LEVx columns to begin entering the standard concentration values for all peaks. You may wish to maximize this table using the Maximize command for a better view while editing.

Move the Edit highlight into each cell in the LEVx column(s) corresponding to each peak and type in the assigned standard concentration values. If the values are the same for a number of different peaks in a single LEV column, you can use the FILL function to duplicate the value of a single row into any number of other rows in that column. You can also use the CUT, COPY, and PASTE command buttons beneath the table to move and copy any set of information in the Reference Table for one peak to the row for another peak - with these functions, only the editable columns will be moved or copied, not the fixed columns at left (RT, PeakName, Cal, and IS).

➤ *To save the new Calibration file*, when you have specified one or more LEV values for each row in the Reference Table, use the Files/SaveAs command to assign it a new name. This Calibration is now available for use in performing manual or automated calibrations and for generating quantitative results.

➤ *To perform a manual calibration using a standard chromatogram*, load the standard chromatogram file into the foreground (top left) register of the Graph window, and integrate it as desired, if this has not already been done. Once all the peaks are marked correctly and desired baselines are drawn, click on the Edit command in the Graph command bar menu, and click on Calibrate in the submenu to see the Calibrate dialog. Depending upon the type of standardization desired, you may need to enter internal standard values or other correction factors in the appropriate fields in this dialog before proceeding. When such entries are finished, select the New mode to begin building a new calibration table, and enter the LEV number for the current standard chromatogram. You can begin with any level, and enter calibrations in any level order.

Click on OK to enter the new calibration data into the calibration table of your new Calibration file.

You can verify that the calibration has been accepted by clicking on the Calibrate SoftButton to see the Calibration screen again. The Calibration Table at bottom, which previously had no lines, will now have a single line corresponding to the chromatogram you have just processed manually. You can move the highlight in the Reference Table at upper right up and down to select any desired peak and see that individual peak's calibration table and corresponding Calibration Plot.

➤ *To perform an automated calibration using a Method*, you must first set up the Method screen to run and process a standard. Click the Method SoftButton to see the Method screen, and click inside the filename entry field corresponding to the first active detector input you will use. Enter a valid filename for the standard chromatogram file to be saved from this input. Now, click inside the Integration filename field to the right of the chromatogram filename field, and enter the name of the Integration file which contains the peak table corresponding to your chromatogram. Click on the CAL

selector arrow and select the N(ew) option to specify a new calibration table. Click inside the Calibration filename field to the right of the CAL field and enter the name of the new Calibration file to which the new data will be added. Finally, click on the LEV selector arrow at far right and choose which standard level will be represented by this chromatogram.

If more than one input is to be used, you could repeat the same entry process for the fields corresponding to the second input, to specify both Integration and Calibration files and the mode and level of the calibration to be done. These could be the same file(s) as for the first input or an entirely different set of files.

To validate all your entries, click in any other field on the Method screen. When you are ready to inject or begin the run, click on RUN to start. The run will proceed until the Duration time is reached. At that time, the file logging to disk would be completed, integration and peak identification would occur according to the instructions in the Integration program utilized, and the Calibration table in the referenced Calibration file would be updated using the new response values for each peak found.

You can further automate this process for a large number of individual standard samples or injections by placing one or more Methods programmed in this manner into a Sequence.

Once a Calibration file has been updated with at least one standard chromatogram run, and a Calibration Table exists in that file for each individual peak of interest, unknown samples can be run or processed to determine concentration results. This can be done manually on a stepwise basis using the Analyze command as in Section 6.4.1.2 above, or automatically for single injections using a Method.

➤ *To analyze an unknown sample using a Method*, you can use the same Method used to run standards to acquire and analyze unknown samples automatically. Recall the Method you used for calibrations into the Method screen, and, if necessary, assign a new filename(s) for the active detector inputs. You should leave the Integration and Calibration filename references as they are, but you must change the CAL field entry from a calibration mode code value (New, Replace, or Average) to a "blank" which specifies that the chromatogram(s) to be acquired will not be used as standards.

Before running the Method, click on the Integrate SoftButton to see the Integration screen, and recall the Integration file you are using in your Method if it is not already loaded. Click on the Reports button at the bottom of the top-level screen to see the Reports Selector dialog, and then choose which report format(s) you want to use for reporting the unknown results automatically. Note that you do not need to make any reports at this time, if you do not wish to, since the chromatogram can be reprocessed at any later time for reporting after it has been run and saved. When you have specified the reports you want, click OK, and then resave the modified Integration file by clicking Files on the command menu bar and then Save.

You can now start the new Method by clicking the RUN SoftButton while viewing any screen except the Sequence screen. When the Method finishes, the unknown chromatogram will be logged to disk, integrated, and peak results computed from the current calibration file. Reports will be made and either printed, shown, or saved as you have directed in the Integration program.

6.4.5 Programming a Simple Sequence

➤ *To use your Method to automate processing of a batch of sample runs using an autosampler*, you can create a new Sequence using the Method along with your new Integration program.

First, click on the METHOD SoftButton to view the top-level Method screen. Check to see if the name of the new Integration file you have saved appears in the INTEGRATION name reference field(s) for the active detector(s). If not, click on the arrow adjacent to this field and select the

Integration filename from the listing by clicking on it. You can use the same Integration for more than one detector input, or use different Integration files for each input.

➤ To program your Method for starting via a contact closure from the autosampler, you must edit the Events Table for the Method to include an initial "wait" command using an input flag on the Communications Processor for this time base. Click on the EVENTS button to see the Events Table.

Click inside the TIME column in the first line of the Events Table, and type in 0. Then move the edit highlight by clicking inside the FIN column on the first line. Click on the arrow button which appears at the right of the cell, which will list all the configured input flags for this time base by name. Select the REMOTE START flag by clicking on it, which will automatically enter the code number for this flag into the first line.

Click on the METH button below the table to return to the top-level Method screen. Now, resave the Method with the correct Integration file reference(s) by clicking on the Files/Save command.

➤ To create an appropriate Sequence program for batch processing, click on the SEQUENCE SoftButton to see the Sequence screen. To use this Sequence to acquire new chromatograms, be sure the Acquire Mode button is selected at top left. To automatically fill in a series of twenty lines for 20 samples, click on the METH column cell on the top line of the table, and type in the Method name you wish to use. Then click the FILL button to see the AutoFill dialog, click on the "To End of Table" and "Increment" radio buttons. The "Step Value" field should have a "1" entry - click inside it and change the value to "1" if it does not. Click inside the "No. of Steps" field and enter "20". Click on OK, and a series of blank rows will automatically be added; the ENABLE/line number row at left will show line numbers between 1 and 20.

Once the Method is specified for each sample line in the table, sufficient information is available to run the Sequence. You can make any further edits of the table as desired, using the INSERT, CUT, COPY, and PASTE commands to work within the spreadsheet.

You have now completed a basic Sequence table for your Method.

➤ To start the Sequence, click on the RUN SoftButton while viewing the Sequence screen, with your autosampler loaded with the samples to be run. Before starting, you should verify that the Method DURATION time you are using is sufficiently short so that the next Method will be ready to start before each new autosampler injection.

After clicking RUN, start the autosampler and all the runs will be processed. All the chromatogram raw data files will be saved, chromatograms will be integrated, and any reports you have specified in your Integration program will be printed, displayed, or saved. You can abort the Sequence at any time by clicking on the STOP SoftButton.

6.5.0 Exiting the Data Ally System

Because Data Ally runs in Windows, you do not need to "exit" the system to run any other application(s). To start or select another application, click on the RESTORE button (double headed vertical arrow) at the top right of the current time base window, and then use the mouse pointer to "resize" the time base window so that you can see the Program Manager icon or any other running application icons beneath it. Double click on the Program Manager icon or any other application icon to view a new window and select a new function or application. Data Ally will remain in operation in the background while this is occurring. You can return to Data Ally at any time by using the SWITCH TO command in the Control Application Menu of any other application, or by resizing or moving the other application window(s) to see any Data Ally window, clicking inside the Data Ally window to make it active, and then resizing it as desired.

When you exit Data Ally, you can choose to either stop operation of the Data Ally system while remaining in Windows, without affecting any other applications which may be running in Windows, or

you can exit Windows to DOS. If you exit to DOS, you will end the current Windows session and all other applications which may be running in Windows will be stopped as well.

➡ *To exit Data Ally*, press the UTIL SoftButton on any Data Ally screen to switch to the Utilities screen, and click on Quit in the Command Menu Bar. When the Quit submenu appears, click on Exit to Application. When you exit Data Ally, any runs currently in progress in Data Ally will be terminated automatically.

7.0 Chromatogram Display, Review, and Editing: Graph Screen

The heart of the Data Ally software is the GRAPH screen, which shows a visual image of what is happening in the currently viewed channel. If live data acquisition is in progress, the Graph will show trace(s) of the running chromatograms which are updated in real time. If no live run is being captured, but the Graph is being used to reprocess, review, or edit chromatogram files, traces of those files and any associated baseline codes or annotations will be displayed. As many different windows as desired can be shown in a single Graph screen display, in order to observe different portions of chromatograms at the same time, or to separate individual traces for easier viewing and analysis.

➤ *To see the GRAPH screen for any channel*, switch to that channel (click on its SoftButton) and click the GRAPH SoftButton with the left mouse button, or press the F2 key. The full screen Graph window (Figure 7.1) will appear.




Figure 7.1 Graph Screen

Fully active and functional Graph screens are displayed as part of Data Ally's Method Events Table, Integration Peak Table, and Sequence Table screens, for convenience during graphical editing steps or while live runs are executing. All Graph screens available in a given channel utilize the same configurable settings, menus, and options, and are essentially identical. However, Graph screens from different channels may have unique features and options, depending upon how each channel has been configured.

The main GRAPH screen window, like all Graph windows, has a Command Menu Bar along its top line. Like all the other program screens, its Command Menu is unique to the Graph screen. It is used to manipulate and edit chromatogram data files, to change the display and view parameters defining the current image, to perform mathematical operations on chromatograms, and to print the contents of the Graph portion of the screen.

The visible contents of the GRAPH screen depend upon the current Configuration file. Elements of the display and of many of its associated command menus can easily be modified or "hidden" completely to customize the nature of the GRAPH display for the applications you are performing. See Section 5 for a complete discussion of configuring the GRAPH screen - if you wish to make any changes to the current appearance or features of this screen on a permanent basis, you may want to re-configure this time base to effect such "permanent" changes so that you will not need to reset any GRAPH screen parameters repeatedly.

The examples shown in the following sections illustrate a configured GRAPH screen in which certain elements have been modified.

7.1.0 GRAPH Screen Commands

There are three types of commands which can be executed directly from the Graph Screen: Windows application control commands, the Graph Command Menu Bar commands, and the "screen based" commands which change the graph display using the mouse pointer.

7.1.1 Application and Program Control Menu Commands

These commands are accessed by clicking the left mouse button on the "hyphen bar" button at the top left of either the channel (full screen) window or the Graph program window nested inside the channel window (Figure 7.2). The Maximize, Minimize, and Restore buttons at the top right of either window are also available for resizing and/or restoring that window to its previous state.




Figure 7.2 Graph Application Control Menu

7.1.2 Screen-Based Commands

The Graph screen is always first displayed as a "full screen" window. If Display/Attributes are set in the default Method to show graph axes, a set of axes will be drawn and labeled with the graph units specified in Configuration, and with the display scaling set by Display/Scaling in the default Method.

After chromatogram files are loaded into the Graph screen using the Files commands, the trace for each loaded file will appear in the graph window, and the filename of each chromatogram will be shown inside a "button" in an area below the horizontal graph axis. The display scaling may change after loading a chromatogram if the "auto scaling" option is active in Display/Scaling.

There are four "screen-based" command types that can be directly accessed from the Graph using the mouse pointer.

- **Panning:** The currently-displayed "view" of a loaded chromatogram or chromatograms can be "panned" horizontally or vertically.

➤ *To pan in any direction*, move the pointer cursor with the mouse to the inner left, right, top, or bottom edge of the Graph window until a heavy horizontal or vertical "arrow" appears (Figure 7.3). Then press and hold the left mouse button, and the chromatogram(s) will begin to "move" on the display in the direction of the arrow. Holding the button down momentarily will result in a small shift in position of the chromatogram(s). If the button is held down continuously, higher-speed panning will begin and the contents of the display window will be moved and redrawn continuously.

Panning is very useful for moving about chromatograms with a constant display magnification, in the manner of reviewing a chart recorder plot.




Figure 7.3 Panning Graph Display At Right Edge

- **Zooming:** The current display window can be "zoomed" by using the mouse pointer to select an area of the screen which is to be expanded to fill the full display window. Zooming differs from panning in that it effectively changes the magnification of the display.

➡ To zoom for expanding the current view, move the mouse pointer cursor to any position which defines any corner of the desired new full screen view, press the left mouse button, and move the mouse pointer while holding the left button down (Figure 7.4). A dotted rectangle will be drawn as the pointer moves, which defines the area of the display to be zoomed to fill the full screen. When the desired area is "rubber-banded" inside the dotted line box, release the mouse button, and the display will change to show the "zoomed" view.




Figure 7.4 "Rubber-Band" Display Zooming

The zoom feature can be used repeatedly to expand any area of the screen as much as desired. Data Ally tracks each zoom command and allows you to instantly "unzoom" back to the previous view. To unzoom, press and release the right mouse button. The previous zoomed view will immediately be restored in the active window.

An alternative mode of zooming, "Zoom Center", can be initiated using View/Zoom Center from the command menu bar (Section 7.5.2 below). If Zoom Center is enabled, a new mode is available in which double-clicking the left mouse button anywhere in the Graph window will cause the display to "zoom" by fixed multiplier factors in both the time and signal axes, centered around the current cursor position at the time of the click. The magnification factors

for center zoom can be reset at any time. Zoom Center is most convenient for repeated magnification of the display around a fixed point. Like the normal rubber-band zoom (which can also be used while Zoom Center is enabled), you can reverse repeated zoom center actions by clicking the right mouse button.

- **Stack Position:** Data Ally allows up to sixteen chromatogram files to be loaded into any Graph window. As each new file is loaded, it is placed in a "stack" of current chromatogram files (Figure 7.5). The first file in the stack is termed the "foreground" file. Only the foreground file position in the stack can be used for baseline editing, or for display of baseline or other integration codes and annotations, to prevent confusion. Any loaded file can be selected for editing or displaying other attributes such as baseline codes or peak names, but the desired file must be moved to the "foreground" position in the stack for this to be possible.




Figure 7.5 Chromatograms Loaded into Register Stack

The first chromatogram register(s) in the stack will usually have names associated with detectors whose signals are physically connected to those registers via the Communications Processor. Detector names are assigned in configuration for each channel, and assist you in identifying each trace during live data acquisition. Up to four of the registers can be assigned input names, as well as a fifth register which can be designated as a "Reference" register for automatic loading of a reference chromatogram during live runs. Note that if any detector input register is assigned a display units scheme in Configuration different from the normal voltage units, whenever that input register's contents are shown as the "foreground" file, the signal (y) axis will become labeled with the units appropriate for that input. It is possible to place additional chromatograms into stack registers for concomitant display during live detector signal input, if desired.

Whenever a file is loaded manually, or during live file acquisition in a Method, the filename is listed as an "icon" with a color underline along the bottom of the Graph window. The underline color is used to show the trace corresponding to that file. The first "filename icon" on the left side beneath the horizontal Graph axis always denotes the "foreground" file. Other filename icons are listed horizontally across the screen bottom after the foreground file button in up to two rows. Each filename "icon" can be used to directly manipulate that file in the stack.

➡ To move a non-foreground loaded chromatogram file to the foreground position, simply point to the desired filename icon with the mouse pointer, hold down the left mouse button, and "drag" that filename icon into the foreground button position. When the new filename is in that position, release the button. The new file will now become the foreground file and is

displayed accordingly. All the other loaded files will be shifted in stack position respectively. The same technique can be used to directly move any loaded chromatogram file to any stack position, for example, before performing mathematical operations. The Display/Swap command can also be used to shift file positions in the stack.

- Offset: The filename icons representing loaded chromatograms can also be used to offset any desired files from their normal positions in the Graph window directly.

➡ *To change the offset of any file directly from the Graph screen*, click and hold down the right mouse button with the pointer on the filename icon button corresponding to that file. A special "offset" cursor will appear in the graph window (vertical up/down arrows). The pointer cursor will change to a two-pointed arrow symbol. When the symbol appears, hold down the left mouse button and "drag" the chromatogram to its desired new position on the display, up or down. At the intended new position, release the mouse button and the chromatogram will be anchored in its new position. You can use the Status Box (Section 7.6.4) as a guide to determine how much offset is being applied graphically to any trace.

The Display/Offset command can also be used to offset chromatogram files. Any file offsets created by "dragging" as above will automatically appear in the Display/Offset dialog box.

- Signal Axis: It is possible to configure Data Ally to show plots of Method program information such as HPLC pump gradients. In normal Autoscaling graph mode, the relative graphical scaling of each active signal trace or program profile is determined automatically so that the entire range of each plot is fit on the same graph as all other plots. As a default, the y-axis or signal axis for each graph window is labeled with voltage units regardless of whether or not signals or traces measured in other types of units are present. If pumps are configured in a time base, you can display pump gradient traces using the Display/Pump Traces command (see Section 7.6.3 below).
- Monitor Multiple Time Bases: If the Data Ally system is configured with two time bases, you can use the MONITOR button in the lower right corner of the GRAPH window to split the main GRAPH window into two panes, each of which shows one of the two active time bases. This feature allows simultaneous monitoring of active detector inputs or other information from more than one time base - see Section 7.6.5 below for more information).

If Method parameters such as pump gradient profiles are present in the current Method's Events Table, a set of small "identifier" icons for each such parameter will appear along the left edge of the graph window. Each of these icons will consist of a label describing its meaning (e.g. "%A") and a color code showing the color and style of the plot trace used to show that parameter. The Display/Traces commands (Section 7.7.3 below) determine whether or not each parameter's trace is actually plotted and shown in the current Graph Window.

7.1.3 Graph Command Menu Bar

The Graph screen Menu Bar commands (Files, Edit, Math, Display, View, and Print) control most of Data Ally's chromatogram file display, manipulation, and editing functions (Figure 7.6). Each has its own corresponding pulldown menu, activated by clicking the left mouse button with the pointer on the desired command.




Figure 7.6 Graph Screen Command Menu Bar

The DISPLAY and VIEW Graph commands are accessible both while the channel is running and while it is idle; other commands which are not effective in either mode are "grayed out". The commands in the FILES, EDIT, and MATH submenus are generally not available for use while a Method is being executed.

7.2.0 Working With Chromatogram Files

Graph (chromatogram) files are DOS-compatible files containing the raw data signal versus time data defining one or more chromatograms. Such files normally also include other information describing the current channel, programs in use, date, time, and operator values, all current integration data, and relevant calibration data.

The Command Menu bar Files submenu includes all the commands needed to retrieve and save chromatogram files manually. Automatic file saving can be programmed as part of any Method or Sequence.

7.2.1 Files Commands

Clicking the left mouse button with the pointer on Files shows the Files submenu (Figure 7.7). This submenu controls manual loading, saving, deleting, and clearing of chromatogram files to and from the Graph window, and supports the handling of "group" files.




Figure 7.7 Graph Files Submenu

7.2.1.1 Clearing Chromatograms

The Clear command clears (removes) a single chromatogram or all current chromatograms from the displayed stack.

➡ *To clear a chromatogram file*, click on Files/Clear with the left mouse button to display the Clear Files dialog box (Figure 7.8). A listing of files currently loaded into defined stack register positions will be shown. Use the left mouse button to check or uncheck the file(s) to be cleared, or check the "All" box to clear all displayed files. Once the desired files have been checked, click on "OK" to execute the clear and remove the files.




Figure 7.8 Files/Clear Dialog

If files loaded into registers which are not listed in the dialog box are to be cleared, you can use the filename icons or the Display/Swap command to shift the positions of those files in the stack so that they can be accessed using the Clear dialog box.

When clearing any chromatogram file(s) which have been altered in the Graph screen, a message will appear prompting you to save all such files before they are cleared. This is intended to prevent any unintentional loss of file editing information.

7.2.1.2 Deleting Chromatograms

The Data Ally software does not support direct raw data file deletions from the GRAPH screen menus. To delete a file, use ALT + ESC to access Windows Program Manager, select File Manager, and use the Delete command to delete any desired files.

7.2.1.3 Loading Chromatograms

The Load command recalls a chromatogram file from disk memory into the active Graph window. If several files are to be loaded into the same window, the Load command must be used for each file, or the files can be grouped using the Save Group command (below) for later recall of all the files together in one step.

➡ To load a chromatogram file, click the left mouse button while pointing to the Files/Load command. A menu of files available in the current directory will be shown, with appropriate descriptions if previously entered (Figure 7.9). Select the desired drive and path where the file to be loaded is saved, and then select a file from the directory listing by clicking on its filename. Click on the "radio button" corresponding to the stack register into which the file will be loaded (the next "open" register will automatically be selected). Once the desired file is selected, click on "OK" to recall the file and load it into the current active Graph window.




Figure 7.9 Files/Load Dialog

When you load a chromatogram which has been previously saved after an integration, the last set of baseline marker codes, which is always saved as part of the raw data file, will also be recalled and displayed along with baseline corresponding to those codes. You can use the Graph Display/Attributes settings (Section 7.6.2) to "turn off" baseline and/or code drawing if you prefer.

7.2.1.4 Saving Chromatograms

The Save As command causes a chromatogram currently in the active Graph window to be saved using the existing filename or a new filename. Each file has a unique name and an optional description. If no filename is available, Save As will prompt you to enter a new one. The SaveAll command resaves all loaded chromatograms using their existing names.

➡ To save a chromatogram file, click the left mouse button on the Files/Save As command. A dialog box will appear showing a listing of files in the current directory path and any previous name for the currently selected file (Figure 7.10). Click in the FILENAME field and type in the desired filename. If the filename currently shown is suitable, click on OK to save the current file in the foreground register using that name. A description can also be entered when any file is saved or resaved by clicking in the DESCRIPTION field and typing in the description as needed. Descriptions will appear on the Files/Load dialog box as well and are useful for identifying the source and nature of your saved files.




Figure 7.10 Files/Save As Dialog

The Save All command differs from Save As in that all files in the current Graph window are saved using their existing filenames. If a Save As or Save All command will overwrite an existing chromatogram file, a trap dialog will appear (Figure 7.11).




Figure 7.11 File Overwrite Trap Dialog

Whenever a chromatogram file is saved, all current baseline marker codes and zone codes will also be saved as part of the raw data file.

7.3.0 Editing Graph Files

The Graph screen Command Menu provides many tools which are useful for changing the contents of any chromatogram file, or for generating results from any file.

7.3.1 Editing Commands

The Edit command submenu (Figure 7.12) is used to manually edit and analyze the chromatogram file in the current foreground position in the Graph register stack. It can also be used to "cut" the current Graph image directly into the Windows Clipboard for "pasting" into other documents or applications, or to create a new graphic file using the current image.




Figure 7.12 Graph Edit Command Submenu

➡ To use the Edit commands, point to "Edit" on the Command Menu Bar and click the left mouse button. The Edit submenu will appear.

- **Analyze:** Analyze creates a report for the foreground chromatogram using the report type specified in the current Method file. This report can be displayed on screen, printed, or saved to disk memory as a retrievable report file. The chromatogram must already be integrated in order to be analyzed via this command.

Click the left mouse button with the cursor on Analyze in the Edit submenu to show the Analyze dialog box. This box gives the filename of the foreground chromatogram, and summarizes all options available for processing this file.

- **Calibrate:** Calibrate manually updates the current Calibration file (which can be observed by clicking on the left column CALIBRATE button) using the peak information from the current foreground chromatogram. The chromatogram must already be integrated to perform calibration with this command.
- **Cut:** Cut captures the current active Graph window image, either placing it in the Windows Clipboard or saving it as a new graphic file with a specified format. The image can then be transferred or imported into other documents or applications, which may be running concurrently in Windows.
- **Define Baseline:** The Define Baseline command accesses all Data Ally manual baseline editing functions. You can draw or redraw baseline as desired, place or move individual baseline peak marker codes, or clear codes using the Define Baseline submenus.

To begin baseline editing, click the left mouse button while pointing at Define Baseline in the Edit submenu. The Define Baseline Submenu will appear.

- **Define Zone Codes:** Zone Codes is used to place integration Zone Code markers at any desired position in the current foreground chromatogram, and, if desired, to transfer the codes into the current Method file's Events Table.
- **Integrate:** Integrate causes immediate integration of the chromatogram file in the foreground register.
- **Reintegrate:** Reintegrate causes a recalculation of the results using the baseline as marked. This function is mainly used to update calculations from prior edit sessions,

7.3.2 Integrating Chromatograms Manually

Chromatograms loaded in the active Graph screen window can be integrated manually at any time, or re-integrated if an automatic integration was performed at the time a chromatogram was first

acquired. Only one chromatogram in the Graph window, the one loaded into the foreground register whose filename icon appears at the top left position under the graph time axis, can be manually integrated at one time. This restriction is intended to prevent confusion in visually understanding the integration results on screen. Any desired chromatogram in any register can be swapped into the foreground position for subsequent integration.

Integrate		
Slope Sensitivity	10	uV/s
Slope Interval	5	5
Baseline Drift	5	uV/s
Noise Reject	30	ms
Minimum Peak Height	1000	uV
Minimum Peak Width	5	s
Minimum Peak Area	100	uVs
Dropline-Skim Ratio	10.0	

Buttons: OK, Cancel, Automatic, Help, ☐ View Only

Figure 7.13 Integration Dialog

➤ To integrate the foreground chromatogram, click on Edit in the Graph Command Menu bar, and click the left mouse button on Integrate in the Edit submenu. The Integrate dialog will appear (Figure 7.13), containing several fields with the integration algorithm parameters from the currently-loaded Integration screen (Section 9.1.2). To perform an integration of the entire chromatogram using these parameter values, click OK - the hourglass will appear as the integration begins. When completed, baseline, peak codes, and other information selected by the Display/Attributes setting (Section 7.7.2) will be shown. NOTE: You may click on “**Automatic**” to have the software try to determine the appropriate integration parameters. If you are not satisfied with the results, you may then adjust them as needed.

When you execute an integration, the system will pause for a moment while completing the integration process, and will then redraw the screen with baselines, peak labels, retention times, and other information (e.g. peak names) derived from the current Integration screen, if appropriate. Chromatograms collected with high sampling rates or with many peaks will require a longer integration time than those with low sampling rates or few peaks.

If zone codes have been placed in the chromatogram at the time it was acquired or by previous editing, they will be used by the Data Ally integrator during manual integration unless they are cleared or edited. You can, of course, place additional zone codes or modify the current codes at any time by using the Define Zone Codes command in the Edit submenu.

With every integration, all previous baseline codes associated with the foreground chromatogram are replaced with the new codes from the new integration. However, those codes are not resaved with the chromatogram file until the chromatogram is cleared from the current Graph window or overwritten with a new file.

To change the integration algorithm parameters, click inside any desired parameter entry field in the Integration dialog and type in the new desired value and units (if you wish to change current units). Double-clicking inside any field will allow you to directly type a new entry to replace the old entry. Then click OK again with the left mouse button to integrate using the modified parameters. Each time you modify the algorithm parameters, any changes are automatically duplicated in the corresponding fields in the currently-loaded Integration program screen (Section 9.1.2).

View Only: Checking the View Only box in the Integration dialog will cause integration of only the current view of the foreground chromatogram in the active graph window, instead of the entire chromatogram. This is useful when you wish to check the effects of the current algorithm on only a small region or set of peaks rather than all peaks, since it reduces the time required to complete integration.

7.3.3 Manual Chromatogram Baseline Editing

The chromatogram loaded into the foreground register of the active Graph screen window can be edited at any time to change baseline markers. New peak labels can be placed graphically, existing peak labels from the original integration can be moved, changed, or cleared, and the chromatogram reintegrated to compute new peak areas at any time. Only one chromatogram in the Graph window, the file in the foreground register whose filename icon appears in the top left position beneath the graph axes, can be edited at one time.

7.3.3.1 Entering Baseline Edit Mode

➤ *To begin editing baseline*, with one or more chromatograms visible in the active Graph window, click on Edit in the Command Menu bar to display the edit submenu. Click on Define Baseline to show the Define Baseline submenu (Figure 7.14) from which all the manual editing commands are selected.




Figure 7.14 Baseline Editing Submenu Commands

If the foreground chromatogram has already been integrated, it will already be marked with integration peak code labels. The Display/Attributes command settings determine whether or not these labels and corresponding baseline tics will be displayed. You can use Display/Attributes to select the information you wish to see concerning previous integration. Generally, you will want to have all existing baseline labels visible when manually reviewing and editing the integration.

7.3.3.2 Drawing New Baselines

➤ *To draw baselines under unintegrated peaks*, click on Draw Baseline in the Define Baseline submenu, and then click on Draw Baseline in the poppa menu (Figure 7.15) which appears. A small "B" will now appear on the Graph window above the cursor position with a vertical line beneath it. You can move the cursor so that the vertical line coincides with the beginning of a new peak, and then click the left mouse button to anchor the "B" or Begin Peak code at that position. Now, as you move the cursor to the right, a new baseline segment will "rubber band" under the chromatogram which ends with an "E" and another vertical line. To mark the end point of the peak, move the "E" to coincide with the exact end position and click the left button again to anchor it. The new codes and joining baseline segment will be anchored to the chromatogram, and the newly defined peak will be automatically reintegrated, its new retention time, number, and name (if identified in the current Peak Table) appearing on the display. You can continue to move anywhere in the chromatogram and mark additional peaks in the same manner.




Figure 7.15 Baseline Edit Pop-up Selector

If you use Draw Baseline to mark a rider peak on the leading or trailing edge of an already-marked larger peak, the system will place a "BS" Begin Skim and "ES" End Skim set of codes instead of the usual "B" and "E". Both forward and backward skims can be accommodated, drawing in either direction.

Note that while drawing peaks or skims, the system automatically prevents you from placing baseline marks which would yield "impossible" peak areas, by "forcing" a peak or skim begin or end location forward or backward from the actual position you placed it to the first "acceptable" position. This feature avoids errors resulting from poor code placement techniques or inadequate display resolution.

You can move outside the chromatogram shown in the current Graph window in two ways. First, moving the cursor to any border of the Graph window will cause the display to "pan" in that direction by approximately one third of the current full scale. Panning in this manner can be done while you are drawing a new peak, if necessary, or between defining peaks, without interrupting the Draw Baseline mode.

The second way to change the view while drawing baseline is to click the right mouse button, which redisplay the Draw Baseline pop-up submenu. Clicking on Zoom in this pop-up allows you to temporarily interrupt the peak drawing mode to rubber band a new zoomed view and zoom the display, or unzoom, in the usual manner. You can use a combination of panning and zooming to achieve any level of visual magnification you may need. Of course, you can also terminate the Draw Baseline mode and use the Graph Display and View commands to change the current scaling and/or view of the chromatogram before restarting Draw Baseline and proceeding.

While drawing baselines and skims, each step will be reintegrated as soon as the necessary code placement is completed, giving you an "on-line" display of the results.

When you have finished marking peaks with Draw Baseline, click on the right mouse button to bring up the pop-up submenu, and click on any other option or on End Drawing to leave Draw Baseline mode. It is possible to use any of the Graph Edit, Display, View, or Print commands without exiting Draw Baseline mode - simply click on the desired command(s). Clicking the right mouse button after executing those commands will retrieve the baseline editing pop-up menu.

➡ **To place drop lines on the chromatogram**, select the Drop Lines option on the pop-up menu (Figure 7.15). A "D" character will appear at the current cursor position with a vertical line underneath it. As you move the pointer, the "D" will move with it (Figure 7.16). Place the vertical line at the exact position where you wish to place a drop line and click the left mouse button to anchor the new drop; the effects of the new drop will automatically be reintegrated and the new peak labeled. You can pan

and zoom while manually placing drop lines as described above. To terminate Drop line placement mode, click on the right mouse button and select any other option or End Drawing.




Figure 7.16 Placing a Drop Line With Pointer

7.3.3.3 Reintegration

Reintegration is always automatic, on a stepwise basis, after each new set of baseline codes is placed to define a new peak, or after existing codes are cleared or moved to re-define the existing set of peaks.

Note that reintegration will utilize Peak Table information to name newly-integrated peaks. Only global integration using the Integrate commands, and not Reintegrate, will account for zone codes placed in the graph window. Both global integration and reintegration will affect every other window in which the current foreground chromatogram in the active window is also present. You can adjust the display of baseline codes and other peak information in any window during reintegration by changing each window's Display and View settings.

7.3.3.4 Moving Baseline Codes

➤ To move baseline codes already on the chromatogram, click on the right mouse button to display the pop-up editing menu (Figure 7.15), and click on Move Codes. A "hand" icon will replace the cursor pointer on screen. Place the "hand" on top of the first code to be moved, and click on the left mouse button to "grab" the code. While holding the left button down, moving the hand icon will cause the code to be dragged in any direction, and baseline associated with that code to be redrawn. When the code is in the desired position, release the mouse button to anchor it; the effect of the change will be instantly reintegrated and displayed.

➤ To terminate moving of codes, click on the right mouse button to display the editing pop-up menu and click on any other option or on End Drawing.

7.3.3.5 Clearing/Erasing Baseline Codes

➤ To clear individual baseline codes already placed, click on the right mouse button to display the pop-up editing menu (Figure 7.15), and click on Clear Codes to activate the clear function. A "clear" icon will replace the pointer cursor. Move the icon on top of any code desired, and click on the left mouse button to clear that code. After each Clear step its effect is immediately reintegrated and displayed. Note that clearing one of a pair of B/E codes defining a peak or a skim will automatically clear both codes; clearing a B or E code defining a range of fused peaks will clear some or all of the individual fused peaks, depending upon how the peaks have been separated. You can clear as

many individual codes as desired in this manner. To terminate clearing codes, click the right mouse button to display the editing pop-up, and click on any other option or on End Drawing.

➡ *To erase all codes, or to clear a range of codes in a single step,* use the convenient Erase command in the Define Baseline submenu. Click on Erase, and the submenu (Figure 7.17) will appear containing two commands: All and Range. To remove all codes from the current foreground chromatogram, click on All. To remove only a range of codes, click on Range, then move the pointer into position to draw a "rubber banded" box to define the range of codes to erase. Hold down the left mouse button, move the pointer to draw the rectangular range, and then release the button - all codes in the indicated range will be removed. Automatic reintegration will always occur immediately after any erase command is executed.




Figure 7.17 Erase Baseline Code Commands

7.3.3.6 Undoing Baseline Code Edits

You can "undo" any changes which have been made to the original set of baseline codes for any chromatogram, one change at a time, by using the Undo command in the editing pop-up menu. the manual baseline editor keeps track of every recent series of changes, and undo can be used to "back out" of those changes.

➡ *To "undo" a baseline code change made during manual editing,* click on the right mouse button to show the editing pop-up menu (Figure 7.15), and click the left button on Undo. A number next to the Undo command shows how many baseline modifications are currently registered in memory which can be "undone". Each time you click Undo, the most recent remaining change is undone, and the reintegration resulting from the Undo step is shown. You can continue to click on Undo until the numerical counter reads zero, indicating that all the tracked baseline modifications have been undone. To terminate Undo, click on any other option in the editing pop-up menu or click on End Drawing.

When you have completed all the desired editing functions, you can move another chromatogram into the foreground editing register by dragging its filename icon into the desired top left position, and then proceed to edit the new chromatogram without leaving editing mode. Alternatively, you can exit manual editing mode by clicking the right mouse button to show the editing submenu, and click the left button on End Drawing. When you leave a manual editing cycle, the modified baseline codes you have placed will remain associated with that chromatogram, and will be resaved with the chromatogram raw data automatically when that chromatogram is cleared or another is loaded into its register. This guarantees that the new set of baseline marks will always be present when the chromatogram is next recalled.

7.3.4 Manual Zone Code Editing

Data Ally allows you to enter timed integration event codes, also known as Zone Codes, directly into any chromatogram to adjust and influence the integration algorithm in time-dependent regions as

desired. You can also use Zone Codes in the Method Events Table to automatically insert those codes into newly-acquired chromatogram data files. Once a Zone Code is used to influence the way a chromatogram file is integrated, processed, or displayed, it becomes a permanent part of that file until it is changed or removed.

The menu of available Zone Codes is defined when setting the Integration screen defaults in Configuration. Automatic insertion of zone codes in a Method is possible only if the Zone Codes column has been configured for the Events Table.

Zone Codes differ from baseline marker codes in that they affect the way global integrations are conducted, using the integration algorithm parameters contained in the Integration program screen. After changing or applying zone codes, you must use the Integration command to see their effects. Reintegration after changing baseline marker codes in Define Baseline mode will not be affected in any way by the presence or absence of zone codes.

7.3.4.1 Placing Zone Codes

➡ *To manually place Zone Codes in the foreground chromatogram in the Graph window,* click the left mouse button on Edit in the Command Menu bar, and then click on Define Zone Codes in the submenu. This will display the Zone Codes submenu (Figure 7.18), which works in almost exactly the same manner as the Draw Baseline edit menu.




Figure 7.18 Define Zone Codes Submenu

Now click on DRAW ZONE CODES to pop-up the Zone Codes dialog box (Figure 7.19). The commands for placing, moving, and clearing zone codes work in exactly the same manner as the Draw Baseline commands. See Chapter 9.0 for a complete discussion of all available zone code types.




Figure 7.19 Zone Codes Editor Pop-up

The Zone Codes edit menu shows a listing of the available Zone Codes: Check BCD Input, Integration Off, Group Integrate, Dropline Riders, Skimline Riders, Valley to Valley, and Reset-End Zone.

To defeat integration in portions of the foreground chromatogram, select Integration Off by clicking it with the left mouse button. A "box" pointer containing an "IntOff" code will appear, which can be moved along the time axis. At the first time point where integration is to be disabled, click the left mouse button and an "IntOff" code will be "dropped" on the time axis. The cursor will then change into a Reset code pointer. At the following point where integration is to be resumed, click again and a "Reset" character will appear (Figure 7.20). You can continue to place additional pairs of IntOff/Reset codes as desired - when complete, click the right mouse button again to redisplay the zone codes edit pop-up. **Note: As you place zone codes on the graph, the time location and zone code command will be automatically inserted into the Method Events Table's Time and Zone Columns.**




Figure 7.20 Zone Codes Placed Along Graph Time Axis

To designate timed regions in which the baseline is to be forced into the valley points, click the left mouse button on Valley-to-Valley in the pop-up. A special pointer containing "VtoV" will appear. Move to the first time point where you wish to begin forcing the baseline up into the peak valleys, and click the left mouse button to drop an "VtoV" zone code on the time axis. Then move to the end of the negative peaks region and click again to drop a "Reset" code. Continue to graphically indicate regions for peak quantitation in this manner; when complete, click the right mouse button to redisplay the edit pop-up.

7.3.4.2 Moving Zone Codes

You can change the current chromatogram view without leaving Zone Code editing mode in two ways. Moving the pointer to any edge of the Graph window will cause the chromatogram to "pan" in that direction by approximately one third of the screen scaling. Panning can be done while you are drawing, moving, or clearing codes. Clicking on Zoom in the pop-up edit menu will enable rubber band zooming without leaving editing mode, enabling you to zoom or unzoom to a new view as required. When done zooming, you must click the right button and use the left button to "uncheck" Zoom in the pop-up menu before continuing your editing.

➤ To move Zone Codes which are already present, click the right mouse button to show the edit pop-up menu (Figure 7.19), and click the left button on Move Codes. A "hand" icon will appear instead of the pointer cursor. Move the "hand" on top of the first code you wish to move, and hold down the left mouse button to "grab" that code. Then drag the code to the desired position while holding the button down, and release the button when the code is placed correctly. Repeat this procedure for moving any other codes. To terminate moving codes, click the right mouse button to pop up the edit menu, and click the left button on any other option or on End Drawing.

7.3.4.3 Clearing Zone Codes

➤ To clear Zone Codes which are already present, click the right mouse button to show the edit pop-up menu (Figure 7.19), and click the left button on Clear Codes. The Clear icon will replace the pointer cursor, which can be moved with the mouse to any code. Place the icon on top of the first code to be cleared, and press the left mouse button to clear that code. Repeat the clear process for as many codes as desired. To terminate clearing of individual codes, click the right mouse button to show the pop-up edit menu and click on any other option or on End Drawing.

➤ To clear all Zone Codes or a range of codes rather than individual codes, click on the Erase Codes function in the Define Zone Codes submenu (Figure 7.18, which operates in exactly the same manner as the Erase Baseline Codes functions..

7.3.4.4 Undoing Zone Codes

➤ To undo changes in Zone Codes, click the right mouse button to show the pop-up edit menu (Figure 7.19) and click the left button on Undo. Next to the Undo command is a number indicating the number of recent changes tracked by the Editor which can be undone one at a time. Click on Undo to undo the last change. The counter will decrement by one each time you click Undo. When finished undoing zone code changes, click on any other option in the pop-up menu or on End Drawing.

7.3.5 Manual Chromatogram Analysis

After you have integrated a chromatogram, you can manually create and see any type of result report via the Analyze command.

When Analyze is used, the currently loaded Integration program is used to identify peak components, the Report Template you select from the Report Library (if any) is loaded by the system and "filled in" with the result data from the chromatogram. If the report chosen includes quantitative results, the current Calibration file will be employed to compute and show those results. Reports can be shown on screen, printed out in hard copy, or saved as report files for later recall.

➡ To Analyze an integrated chromatogram to see results or get a report, click on Edit in the Graph Command Menu bar, and then on Analyze. The Analyze dialog box will appear (Figure 7.21). The dialog box includes four types of information: the "Correction Factor" fields as are configured for the currently loaded Method, for entry of sample-specific data relevant to the report; the "Internal Standard" concentration values for up to five peaks as configured; sample information as configured such as name and number; and a Template/Save File Selector allowing a choice of any report format from the current report template library, and setting a filename for any saved report. A set of *.XLS report template files are available for printing via Microsoft Excel, which can be edited (or new templates prepared) for custom reports (see Section 16).




Figure 7.21 Manual Analyze Dialog

If you want or need to enter sample-specific information, such as Sample Name or Dilution factor, click inside the appropriate field and type in or change the entry as desired. Reconfiguring the current Data Ally time base to remove certain types of entry fields will remove those fields from the Analyze dialog as well. Note that only the correction factors, sample information fields, and other items which have been set in configuration for this time base will be active and visible - in order to enable additional items, you must re-configure the current time base before proceeding.

For a complete explanation of the use of correction factors and internal standards in quantitative calculations, see section 10.3 below.

Click on the check boxes to toggle screen display or printing of the report on or off. Any combination of these options can be used. If you wish to save the report, you must enter a filename in the filename entry box.

Click on the arrow at the right side of the "roll down" list box labeled TEMPLATE to see a complete listing of all current report types available. Select the one template to be used for the new report by highlighting it. These templates are all *.XLS files present in the \Data Ally\REPORTS subdirectory. Click on OK when you are finished.

The selected report will now be produced by Microsoft Excel either on screen or on the current printer, or both (Figure 7.22). If you have chosen screen display, you can "page down" through the report sections for review using the Excel "Page" button, or "zoom in" on any section of the report with the "Zoom" button.. **Note: You must have a Printer Driver configured in Windows else a report will not be generated.**




Figure 7.22 Example of Typical Area PerCent Report via Excel

Only one report can be generated via Analyze at a time for a single chromatogram. To manually produce additional reports, repeat the steps above and select a different report template.

To save the report as an Excel .XLS file, type a target directory path and filename for the report file into the FILENAME entry field before exiting the Analyze dialog box with OK. A new .XLS file will be created in the target path including all the information and format for the current report template used.

7.3.6 Manual Calibration

You can manually update a Calibration, if the current foreground chromatogram is for a standard, using the Calibrate command. Based on your entries, Calibrate will recompute the current Calibration file, which can then be reviewed and edited in the Calibration screen (Section 10). Of course, in order to perform valid calibration updates via manual calibration, you must be sure that appropriate Integration and Calibration files are currently loaded.

➡ To perform a manual calibration, with a "standard" chromatogram integrated in the foreground register, click on Edit in the Graph Command Menu bar, and then on Calibrate. The Calibrate dialog box will appear (Figure 7.23).

The screenshot shows the "Calibrate" dialog box. It has a menu bar with "File", "Edit", "Math", "Display", "View", and "Print". The dialog is divided into several sections:

- Correction Factors:** Includes fields for "Dilution" and "Weight".
- Internal Standard Amounts:** Includes a field for "IS1".
- Sample:** Includes fields for "Name" (set to "newsample") and "No.", with a "Descr..." button.
- Calibrate:** Includes a "Level" field (set to 1) and radio buttons for "New", "Average", and "Rep. Lev.".
- Report:** Includes a "Template" dropdown (set to "extstd.xls"), a "Filename" field, and "Print" and "Show" checkboxes.

Buttons for "OK", "Cancel", and "Help" are located on the right side of the dialog.

Figure 7.23 Manual Calibration Dialog

The dialog box contains four types of information: sample-specific factors describing the standard, such as Dilution, as defined in the Sample Information dialog of the current Method screen; Internal Standard concentration entry fields for the number of internal standards currently configured; a Calibration Level selector and Mode radio buttons, for choosing the type of calibration (New Table, Replace Level, Average); and a template selector box, which specifies the calibration report file template, if any, and filename if a report file is to be saved.

Click inside any of the standard data entry fields and type in the appropriate values, if any, or change the default values which appear. If internal standard(s) is called for in the current Peak Table, no entries need to be made in the IS fields provided the Calibration Reference Table (Section 10.1) includes appropriate values for the internal standard peak(s). If entries are made here in the IS fields, they will override any values found in the Reference Calibration Peak Table.

Select the standard level assigned to the current standard for manual calibration. The concentration of this standard for all component peaks should agree with that listed in the current Calibration Peak Reference Table for the selected level. Click on the radio button corresponding to the type of calibration (New, Level Replace, or Average) to be performed - only one option can be marked.

See Section 10.3 below for a complete explanation of calibration functions, including the use of correction factors and internal standards.

Click on the right arrow button in the TEMPLATE field to see a listing of all the available report templates. You may select an appropriate report containing calibration information for printout or display. Of course, calibration data can be included in any Data Ally report, if desired, simply by editing that report's template format (see Section 13).

Click on the check boxes to select calibration report display (Show) or printing. If you wish to save the calibration report as a separate file to disk, you must enter a path/filename for that file in the CAL REPORT entry box. This will result in an .XLS file being saved using the template report format to the target path.

When complete, click on OK to do the calibration update. The current foreground chromatogram will be used to increment the current Calibration file. If the SHOW and/or PRINT boxes are checked, a calibration report will be shown and/or printed (Figure 7.24).

You can click on the CALIBRATE SoftButton to view the current Calibration Table and Plot to verify that the new standard has been accumulated into the table (see Section 10 below).




Figure 7.24 Example of Calibration Report

7.3.7 Manual System Suitability Calculations

System Suitability calculations may be performed on any peaks in a chromatogram. Normally peaks will be baseline resolved to obtain meaningful results.

With Data ally, only selected peaks may have system suitability calculations performed on them. Normally at least two peaks will be selected to enable all parameters to be calculated. In addition the unretained peak should be identified in the Peak Table, if available, or in the Column Parameters section if a specific retention time is required.

➤ To perform a manual System Suitability Calculation

The peaks in the chromatogram should have a corresponding Peak Table relating to them. Select the first peak in the table that is to have System Suitability calculations performed on it and enter an upper case 'S' in the Suit Code column for the peak. Repeat for all other peaks for which you wish to have System Suitability calculations.

If the Unretained or Ko has been integrated, place a 'K' in the Suit Code column corresponding to the peak. If the Unretained peak has not been integrated the retention time and name may be inserted in the Column Parameters dialog on the top level of the Integration Screen.

Re-integrate the chromatogram with the Edit, Integrate Chromatogram option.

The selected peaks will have additional lines drawn on them corresponding to the system suitability calculations.

To produce a report, go to Edit, analyze and select a Suitability report template (A4SUIT.XLS in Europe) to either print or show in the normal way.

The System Suitability report will be produced. Each peak will be displayed separately with the calculations shown next to it. NOTE: For all the calculations to be performed, certain column information will need to be entered in the Column Parameters dialog in the top level of the integration screen.

7.3.8 Cutting the Current Display Image

The Cut command allows you to easily transfer the currently-displayed chromatogram image in the active Graph window into the Windows clipboard or directly into a DOS file which can be passed to another software application or saved for later use. Cut operates on the complete image in the active window, which should be rescaled, zoomed, and otherwise adjusted exactly as desired before using this command.




Figure 7.25 Cut Image Command

➤ To cut a copy of the current Graph image to the Clipboard, begin by making sure the image on display is totally suitable. If not, make whatever changes are needed using the Graph/Edit,

Graph/View, or Graph/Display commands. Then, click on Edit, and then on Cut in the submenu (Figure 7.25). The current image will immediately be copied into the Clipboard, where it will remain available for pasting into any other compatible application. The image will remain in the Clipboard until you either clear it, cut or copy a new image to the Clipboard, or exit Windows.

7.4.0 Math Functions

The Math functions on the Graph Command Menu Bar (Figure 7.26) perform mathematical operations and manipulations upon one or more chromatograms loaded in the current active Graph window. These functions are executed manually using the Math command; they can also be performed automatically as part of a Method if the Auto Math options are utilized. Math functions submenu commands will be "grayed out" and unavailable unless more than one chromatogram file is loaded into the current Graph screen.

Most Math functions are based upon and/or affect only the first two chromatograms in the stack. These chromatograms can be identified as those whose filename icon buttons appear at the top left in the area beneath the horizontal graph axis. You can rearrange the order of chromatograms in the stack by dragging the filename icons into the desired positions or by using the Display/Swap command prior to executing any Math commands.




Figure 7.26 Math Functions Menu

7.4.1 Adding/Subtracting/Ratioing Chromatograms

The three math operation functions use the first two chromatograms in the register stack. You must position the correct chromatograms in those registers via the filename icon buttons or using the Display/Swap function before executing these commands.

- **Add:** Adds the first two chromatograms in the stack (foreground plus next file) and places the result into the first open register in the stack. The result chromatogram is drawn using a different color trace and is labeled "Sum". If the result chromatogram is to be saved, it can be renamed as desired and saved as a new file using the Files/Save As command, or it can be automatically saved via the Result command below.

➡ *To add two chromatograms*, click the left mouse button with the cursor on the Math/Add command. The first two chromatograms in the stack will automatically be added and the result displayed as a new chromatogram named "Sum". If there are no stack register positions available in which to place the result chromatogram, an information box will appear when Add is initiated advising that at least one file must be cleared from the stack before addition is possible.

- **Subtract:** Subtracts the chromatogram in the second register from the chromatogram in the foreground register, and places the difference chromatogram into the next available open

register in the stack. The difference file can then be saved with any desired file name, or used for additional math operations.

➤ *To subtract two chromatograms*, click the left mouse button on the Math/Subtract command. The first two chromatograms in the stack will automatically be subtracted and the resulting difference chromatogram will be displayed as a result named "Diff", in the first available stack position. If there are no open stack positions in which to place the result chromatogram, an information box will appear when Subtract is initiated advising that at least one file must be cleared from a register before subtraction is possible.

- **Ratio:** Ratio computes the ratio of the signals of the first two chromatograms, the foreground chromatogram being divided by the second register file. The result or ratio file is placed into the first open stack register. It can then be saved using any desired filename.

➤ *To ratio two chromatograms*, load the desired files into the first two registers and click on Math/Ratio. The ratio will be computed and will be shown in the first open register with the name "Ratio". If there are no open stack positions to display the ratio file, an information box will appear when you select the Ratio command advising that at least one currently-loaded file must be cleared to open a register before ratioing is possible.

7.5.0 View Commands

The View submenu accesses two special Zoom functions. A "Zoom Box" can be placed in the Graph window as an indicator of the segment of the current chromatogram included in the current view. This feature is convenient when reviewing details within a chromatogram at high magnification. A special "Zoom Center" function is also provided which changes the standard "rubber-band" zoom mode to a "center-zoom" mode which expands the display centered on the cursor position by fixed magnification factors set in Configuration.




Figure 7.27 View Submenu Commands

7.5.1 View Commands

The View submenu commands are as follows:

- **Zoom Box:** Places a special "zoom locator box" on the current Graph window which indicates the current view in relation to the entire chromatogram. Checking the Zoom Box command in the submenu activates the Zoom Box (Figure 7.28), and unchecking it defeats the Zoom Box display. Once displayed, the Zoom Box can be resized and repositioned like any other window. Inside the bottom portion of the Zoom Box are numerical indicators of the current time and signal coordinates range covered by the current view in the Graph window. One Zoom Box can be enabled for each Graph window displayed.

➤ To enable the Zoom Box in the active Graph window, click on Zoom Box. The Zoom Box will appear at the upper right of the Graph window. Inside you will be able to see a miniature image of the current foreground chromatogram file, with the portion of the file "highlighted" in a rectangular area corresponding to the image actually being observed in the Graph. You can now resize and/or reposition the Zoom Box as desired.




Figure 7.28 Zoom Box Enabled in Graph Window

- Zoom Center: Enables "zoom center" zooming mode instead of the normal "rubber-band" zooming. Clicking on the Zoom Center command to check it causes a double-click of the left mouse button within the current Graph window to "zoom" the display in that window by preset factors for both the signal and time axes. Zooming is always centered on the current mouse cursor position. Zoom Center is convenient for repeated zooming at specific points since it requires only one click rather than a rubber-banded definition of an area to zoom, and allows the cursor to remain fixed during zooming..

➤ To enable Zoom Center in addition to rubber-banding zoom, click on Zoom Center to see the subcommands "Enable" and "Set Zoom Scale Factors" (Figure 7.29). Click on "Enable" to activate zoom center. When Enable is checked, a left double click anywhere in the Graph window will cause the display to be zoomed around that point. Clicking the right mouse button will cause "unzooming" by the same factor.




Figure 7.29 Zoom Center Commands

To disable Zoom Center, click on the Zoom Center command again to uncheck it.

➤ To set the zoom center expansion factors, click on Set Zoom Scale Factors. A dialog box will appear (Figure 7.30). Click inside the "Horizontal" or "Vertical" axis fields and enter

the desired expansion factor value. An entry of "10" for either field will cause the display to expand by a factor of ten times in that field's dimension. When the values are set as desired, click OK.




Figure 7.30 Zoom Center Scale Factors Dialog

Zoom Center can be enabled for any specific Graph window(s). It is possible to have some windows set for standard rubber band zooms only and some for zoom center. The rubber band zoom option is always available whether zoom center is active or not.

7.6.0 Changing Graph Display Characteristics

The Display submenu of the Graph Command Menu Bar (Figure 7.31) sets and modifies most aspects of the current active Graph window, including what information is displayed, the colors and relative positions of chromatogram and other traces, the relative placement of multiple chromatograms, and the order of chromatograms in the register stack. Some of the Display functions can also be performed using the filename icon buttons at the bottom of the Graph window.




Figure 7.31 Graph/Display Submenu Commands

7.6.1 Display Attributes

Attributes settings modify the type of graphical information shown in the current Graph window, and provide a wide range of flexibility. As an example, in certain instances, only some types of information may need to be displayed describing the integration of the loaded chromatogram, so it may be desirable to "turn off" baseline markers, peak names, and

retention values. Disabling such selected data from the Graph can be done via the Attributes table.

➡ *To change current Graph Attributes*, click on Display/Attributes to see the Attributes dialog box (Figure 7.32). Any combination or all of the following items can be checked (enabled) to place each checked item on the display whenever it is appropriate (unchecked items will never be displayed):




Figure 7.32 Graph Display/Attributes Dialog

- **Fill Above CG:** This "fills" the display with a color above the foreground chromatogram trace. The color used can be modified with the Display/Colors command below. Fill is often useful for graphically distinguishing peaks by increasing contrast against the screen background. The normal default for Fill Above CG is disabled.
- **Fill Below CG:** This "fills" the display with a color below the foreground chromatogram trace, which can be changed via the Display/Colors command. The normal default is fill off.
- **Fill CG Above BL:** This "fills" the actual areas of all peaks in the foreground chromatogram between the trace and drawn baselines. This Fill Attribute can be very useful if it is enabled during manual baseline drawing, since peaks are filled actively with color as baselines are adjusted. The color can be reset using Display/Colors. The normal default is fill disabled.
- **Show Zone Codes:** This turns on or off display of zone codes. The normal default is for all codes to be displayed.
- **Show Baseline:** This defeats or enables display of baseline drawn under all found peaks. The normal default setting is baseline enabled.
- **Show Baseline Marks:** This Attribute shows the locations of all baseline codes by "tic" marks.
- **Show Baseline Codes:** This causes baseline codes to be displayed or omitted. The normal default for baseline codes is enabled. Note that any combination of codes and their corresponding tic marks can be selected.

- **Show Peak Names:** This causes peak names to be shown or omitted. If enabled, names of all identified and named peaks will always appear in the Graph window for every appropriate peak. The normal default is peak names shown.
- **Show Peak Numbers:** This shows peak numbers, provided they are assigned in the Peak Table. The normal default is for peak numbers shown.
- **Show Retention Times:** This causes absolute retention times to be shown or omitted on the Graph display, if an integration has occurred to assign them.

7.6.2 Screen Display Characteristics Settings

Several of the Display submenu commands change the overall appearance of chromatograms and other types of plots, their order and hierarchy, and their relative positions.

- **Colors:** The Colors command changes all colors assigned to the basic Data Ally Graph screen displays and functions. This includes colors for chromatogram traces, hardware function traces (pump or temperature profiles), graph axes, annotations, chromatogram names, buttons, baseline and zone codes, etc.




Figure 7.33 Display/Colors Dialog

➡ To modify Graph screen colors, click on Colors to see the Colors dialog box (Figure 7.33). The box includes sample color fields for all the components of the screen. To change any color setting, clicking the left mouse button on the sample field for that item - a color "palette" will pop up, showing the selection of possible colors. Select the desired color by clicking on it. The selected color "box" in the palette will have a dark outline drawn around it. Repeat the same process to reset the colors for additional items in the listing. When you are finished, click on OK.

You can change any of the colors assigned to loaded chromatogram traces by clicking on the right arrow in the Traces selection field and highlighting the desired trace name to find and edit its color.

System colors can be changed for each Graph window independently, and can be modified while any window is running or idle. Colors for SoftButtons, program screens, and other Data Ally graphics must be changed in Configuration for each channel, or via the Windows Program Manager/Control Panel.

- **Offset:** The Offset command is used to "move" the apparent position of any chromatogram in the active Graph window in the signal dimension only. You must use the Math/Move command to displace chromatograms in the time dimension. Offset performs the same function which can be done by double-clicking on any chromatogram filename icon and then dragging the chosen trace to a new position directly on the screen.




Figure 7.34 Display/Offset Dialog

➡ To offset one or more chromatograms, click on Offset in the Display submenu. The Offset dialog box will appear (Figure 7.34). It contains offset fields for each defined stack register and for other displayable parameters, such as pump or temperature profiles. Offset entries can be defined differently for the same chromatogram in separate Graph windows.

To specify the amount of offset applied to any chromatogram only, click on its offset entry field, and type in a positive or negative value in appropriate signal units for the desired offset. Then click on OK.

Note that when you offset a chromatogram graphically by using its filename icon button to drag the entire chromatogram to the desired position, the offset created for that chromatogram will appear in its offset field in this dialog. You can refer to either this value or the status box signal value to determine the exact offset applied to a chromatogram.

Offsetting a chromatogram or chromatograms does not change the raw data signal values for that chromatogram but only the apparent display position.

The same technique is utilized for offset control parameter traces also listed in the Offset dialog, depending upon time base configuration. Appropriate units must be entered for each item selected; offsets will, of course, only be relevant if the Traces command (below) is used to enable display of that item in the Graph window.

- **Scaling:** Scaling sets the overall scaling for all types of screen display registers, and enables "Autoscaling" if this is appropriate. The default Scaling mode is usually Autoscaling, which will redraw the active Graph window display axes to contain the dimensions of the largest chromatogram file (signal level or time) in any stack register. Autoscale allows approximately 10% of available screen space for any window to accommodate drawing of baseline codes, retention times, and other labels and annotations, and for easier legibility.

➡ To reset Graph scaling, click on Scaling to see the Scaling dialog box (Figure 7.35). This box contains scale factor boxes for all the named (up to five) input register signal ranges, the time range, and the range for any HPLC or GC control parameter values which are enabled in Configuration. The hardware, input, and time registers also each have separate check

boxes for "AutoScaling", which can be enabled independently for any of these parameter types.




Figure 7.35 Display/Scaling Dialog

To set a fixed range (maximum and minimum value) for the current active graph window for any item, click on the Autoscale check box for that item to defeat Autoscaling, and then click inside that item's corresponding MIN or MAX entry field and type in the desired value. If no MIN is entered, zero value will be assumed. If no MAX is entered, Autoscaling will be enabled automatically. If Autoscaling is checked, Autoscale display will override any simultaneous values in the MIN and MAX fields for each parameter.

Since Autoscaling can be independently set for each Graph axis, it is possible, for example, to set the display for auto-expansion only in the time axis during live runs, or only in the signal axis with fixed-window scrolling in the time axis. The settings for scaling will determine exactly how the display appears during live acquisition runs. Like all the other Display functions, scaling can be changed on a real-time basis while runs are in progress to observe the resulting effects.

All Method Events table parameters which are plottable (pump %, GC oven temp) will appear as separate items in the Scaling dialog for each channel, and can be switched between Autoscaling (the normal mode) and manual scaling in the same manner as the input signal traces.

To set scaling for HPLC pump profiles if these are to be plotted, click Pumps in the Scaling dialog (Figure 7.35) to view the Pumps/Scaling Dialog (Figure 7.36).




Figure 7.36 Pumps/Scaling Dialog

If the Method Events Table in the current time base has a pump control program loaded, this dialog will control scaling of the pump gradient traces superimposed on the signal graph window plot. Enter the desired values for scaling as in Figure 7.35 and click OK when complete.

- Pump Traces: If HPLC pumps have been configured for the current time base, a separate command will appear at the top of the Display submenu, Pump Traces (Figure 7.37).




Figure 7.37 Display Submenu with HPLC Pumps

Click Pump Traces to pop up a listing of all configured pumps (Figure 7.38). Click on each pump related parameter (total flow or individual gradient percents for each pump) to check-mark that parameter and display a plot of that parameter on the GRAPH window along with the signal information.




Figure 7.38 Setting Pump Trace Plotting

Traces will now appear as long as their corresponding parameters are "checked" in this menu - to remove the traces, click on Pump Traces again and click to "uncheck" each item.

- Swap: Swap shifts all loaded chromatograms in the current Graph window into the next forward register, with the previous foreground chromatogram being moved into the "last" register. Repeated use of Swap will place any chromatogram in the stack into the foreground position for editing operations. The alternative to Swap is to drag the chromatogram filename icons into the desired positions.

➡ To swap the second-register chromatogram into the foreground register, click on Swap.

7.6.3 Status Box

The Status box is a powerful and flexible tool which can provide much additional information about the course of a running Method in real time, or assist you with editing functions while the current channel is idle. Using this command will place a "status box" window on the active Graph windows (and in all program windows in the current channel) which contains only the items you select from a configurable listing of many system and computer parameters. The status box can be moved to any position on the display for maximum convenience. When the Status Box is placed in the Graph screen in any channel, it will automatically also appear in all other program window screens in that channel. It is possible to disable the appearance of the Status Box in any or all program windows using the Options/Status command in each such window.

➡ *To set up and enable the status box*, if it is configured, click on Status Box. The Status Box dialog box will appear (Figure 7.39). In the dialog box are a number of items with corresponding check boxes. You can check each item you wish to appear in the status box. When selecting items, you should consider that more items will require a larger box which may obscure some of the active Graph window. It is only possible to have a single status box open in any Graph program window regardless of how many panes you have opened within that window.




Figure 7.39 Display/Status Dialog

Be sure to check the top box for "Enable Status Box", and click on OK. The Status Box will appear as you have defined it in the active window (Figure 7.40). You can now move the box at will.

If the channel is running (you click on RUN), the status box will show a "live" readout of any parameters included, indicating detector input conditions, time, Method program conditions such as pump flows/gradients, flags executing, etc. If the channel is idle (no Method is running), the status box time values will correspond to the current pointer position, signal values will correspond to the signal values of any chromatograms in the selected register(s) at the current pointer position, and Method parameters will correspond to the value of each item in the current Method at the current pointer location. The Status Box can be very useful while editing chromatograms as a time/signal counter, or to track the exact position in a Method where a peak or event appears. If an offset is applied to any trace graphically or by using the Offset command, the signal value shown for that trace in Status will reflect the offset.




Figure 7.40 Status Box Enabled in Graph Window

➡ To disable the status box at any time, click on Status again, and uncheck the Enable Status Box item, then click on OK. The Status Box will disappear until you re-enable or re-define it.

If the Status Box signal display is active during a HOLD condition, or between runs in a Sequence, the signal register(s) will show the current status of each live detector - this can provide valuable information for assessing whether run conditions are suitable for starting the next injection, or if the system has re-equilibrated from a previous run.

Like the other Graph Display parameters, the state and contents of the Status Box can be preset and associated with each individual Method so that it automatically assumes a given "default" state whenever that Method is loaded.

7.6.4 Monitor Time Bases Function

There are two means by which Data Ally allows an operator to monitor the progress of active detector methods or reprocesses in more than one time base at once on the same screen - the MONITOR button at lower right of the GRAPH window (Section 7.1.2) and the Monitor Time Bases command in the Display submenu (Figure 7.31). Both have exactly the same function, and both are only effective if the Data Ally system is configured for two or more time bases.

To use the Monitor Time Bases command to place two windows showing two different active time bases on one display screen, click on the command. The GRAPH window will split horizontally into two panes, each showing one of the two time bases (Figure 7.41).




Figure 7.41 Monitor Time Bases Screen Display

All actions in either time base will now be observed on this screen. You can, of course, resize and reshape each pane as desired using the standard commands and cursors. To restore a single-window view, either simply click the MONITOR button in the lower right hand corner of the screen, or select Display/Monitor Time Bases again from the command menu bar.

7.7.0 Manual Printing of Graph Documents

The Print submenu of the Graph Command Menu Bar (Figure 7.42) allows manual printing of chromatogram files and other information in the current active Graph window. These Print functions can be executed at any time the channel is not running a Method.

Data Ally normally prints via the Windows Print Manager. This means that Windows itself manages the setup of printers and all printing operations, including spooling and print queuing. Printing images from the active Graph window is based on the Windows "what you see is what you get", or WYSIWYG, convention; the printed image will represent exactly those fonts, characters, and graphics shown on the display.

You can use Print Manager spooling by setting its "Low Priority" option to delay print operations to cause minimal slowdown of on-screen displays while Data Ally is running. Alternatively, you can set the "High Priority" option to achieve faster printouts, usually at the expense of screen update times. It is possible and sometimes desirable to use Windows-compatible print spooling or buffering accessories with Data Ally if heavy printing workloads are anticipated, especially if multiple Data Ally Interfaces will be running on a single PC.





Figure 7.42 Print Commands Submenu

Printers are defined for each channel in Configuration/Utilities. You can use the Windows Print Manager at any time while Data Ally is running to change the specifications for a current printer, but you cannot change the printer assigned to a given Data Ally channel without reconfiguring the system.

7.7.1 Printing Chromatograms

➤ *To print various images of chromatograms in the current Graph Window*, click on Print and then on "Chromatograms" to see the Chromatograms Print dialog box (Figure 7.43).



Whole Document: Prints the current active Graph window as a single page, autoscaled for the largest chromatogram loaded in the window.

In Sections: Prints the current active Graph window sized to show the longest loaded chromatogram on several pages, with each page time-scaled to the current scaling. The "From" and "To" entry fields allow selection of only a partial set of pages - the total number of pages to be printed is automatically computed and shown as soon as the option is selected.

If the "Include Printout Info" check box is checked, the informational items chosen using the Printout Info command (see 7.7.2 below) will be included in an area at the top of the printout. Leaving this check box unchecked will cause only the chromatogram file name(s), axis labels, date and time of printout, and page number(s) of printout to be included.

The screenshot displays a Windows 95 desktop environment. In the background, a window titled "Graph - c:\pyramid\pdf\manh-p36.pdf" is visible, showing a technical drawing with dimensions like #9, #10, #12, and #14. Overlaid on this is the "Print Setup" dialog box for the "HP LaserJet III on LPT1:" printer. The dialog box has two tabs: "Printer" and "Paper". The "Printer" tab is selected, showing options for "Default Printer" (currently "Windows Printing System on LPT1:") and "Specific Printer" (selected, "HP LaserJet III on LPT1:"). The "Orientation" section shows "Portrait" selected. The "Paper" section shows "Size: Legal 8 1/2 x 14 in" and "Source: Upper Tray". Buttons for "OK", "Cancel", and "Options..." are on the right. A status bar at the bottom shows coordinates from 10.00 to 22.50.

Figure 7.44 Printer Setup Dialog

Click on the appropriate radio buttons to select either the current printer or select a new printer from the list of Print Manager compatible printers. If you wish to change the orientation for the printout, the page size, or the paper tray used, click these options as desired. Click OK when you are finished.

You can also click on the Options button at center right to make additional modifications in the printout with the Print Options dialog (Figure 7.45).




Figure 7.45 Print Options Dialog

Use the Options shown specific to the current printer driver to change the appearance and handling of printouts as desired. When finished making changes, click on OK or Cancel.

NOTE: THE "PRINT TRUETYPE AS GRAPHICS" CHECK BOX, IF SHOWN, MUST BE CHECKED IN ORDER FOR Data Ally PRINTOUTS TO BE PROCESSED CORRECTLY.

The settings made in the Chromatograms dialogs will remain in effect to control manual chromatogram printing until they are changed again. The same settings will apply to separate-page chromatogram printouts made automatically as part of Method or Sequence execution.

7.7.2 Selecting Information Included On Chromatogram Printouts

The Printout Info command allows a set of informational parameters to be included in any chromatogram printout, whether performed manually using the Print submenu or automatically at the end of a Method or Sequence. Once a group of parameters is chosen using the Print Info dialog, that information is always printed along with chromatogram graphics.

➡ To set up information on chromatogram printouts, click on the Printout Info command to see the Printout Info dialog box (Figure 7.46).




Figure 7.46 Printout Info Dialog

Inside the Printout Info dialog are a number of check boxes corresponding to various system information items. All the items in the Run Info, Chromatogram Info, and Sample Info blocks are relevant to the "foreground" loaded chromatogram file only.

To select the set of desired items, check their respective check boxes, and then click OK. Any or all of the items can be chosen for inclusion.

Once a series of items are checked, they will be printed in a space above the chromatogram printout on each printed page. If multiple-page prints are made using the "In Sections" printing option, all the selected information items will be printed on each page.

The Printout Info parameters will always be included in every chromatogram printout produced manually or automatically via a Method or Sequence. To change the current set of Info parameters, you must use the Printout Info dialog again.

Note that every page of every chromatogram printout will always be numbered at bottom right, whether portrait or landscape orientation, with a page number and the total number of pages (e.g. "Page 1 of 4").

7.7.3 Specifying Report Types

The Report command in the Print submenu allows any type of report for which a template *.XLS file exists to be produced using the foreground chromatogram data currently loaded. To use this function, simply click Report to see the Print Report dialog (Figure 7.47).




Figure 7.47 Print Report Dialog

This dialog allows you to select a report template by name using the entry box with right selector arrow - click on the down arrow button to see a listing of all available types in the reports subdirectory, highlight the one to be used, and click the left mouse button. If you wish to save the new report as an Excel report file, enter a name in the "Save As Filename" entry box. Check either the Print or Show check box, to print out or display the report on screen. When done, click OK to produce the report desired in Excel, or Cancel to abort.

7.8.0 Running From the Graph Screen

Runs can be initiated from the Graph screen, or from any other screen in Data Ally, by double-clicking on the RUN button on the bottom of the left switchbar at any time. Function key F12 on the PC keyboard can also be used to issue the RUN command. When RUN is clicked, the button changes to STOP, which can then be used to terminate the current run at any time. If a manual run is in progress, the STOP button will include the word "RUN" indicating that a single, manual run and not a Sequence is being performed. You must double click on both RUN and STOP to execute those commands.

When a run is executing, the top bar which contains the name of the current channel will include a notation that the channel is running, along with a "three-dot" live indicator which visually prompts that a run is being processed. You can use any of the left switchbar buttons to view the contents of any screen in the current channel, or to switch to any other available channel, during any run.

The following steps occur each time the RUN button is clicked in the Graph screen:

1. Data Ally reviews the Method now present in the current channel. If this Method specifies the use of Integration or Calibration files for processing any data, the referenced files will be loaded for the first designated active detector input.
2. All chromatogram or other files currently loaded into the Graph screen registers corresponding to the detector inputs which are active in the Method are cleared. If those files have not been saved, the software will direct the operator to save those files before proceeding to execute the Method. Chromatogram files not occupying registers assigned to the signal inputs for the Method will remain in the Graph window.
3. The channel will begin to execute the Events Table for the current Method. The "Active" light on the panel of the Communications Processor module for the channel will illuminate while the Method is running. If no Events are included in the table, the Method will simply acquire the data from the specified detectors for the DURATION time and at the RESOLUTIONs specified for each input.

As the data are acquired, the Graph window will either rescale itself automatically, if AutoScaling is specified in the Method, or will assume the fixed scaling directed by the Method. If fixed scaling is used, the active window will "scroll" forward each time the trace approaches 80% of the full screen time value, with the total time span for the window remaining constant. With AutoScaling, the display will resize itself at continual intervals to accommodate the increasing length of the chromatogram traces. Signal Autoscaling can be used with either fixed or auto time scaling. If signal is autoscaled, the signal axis of the active Graph will expand to contain the largest excursion of any chromatogram being acquired.

4. While new chromatograms are being acquired in real time, the operator can pan the display to any previous position, reset the scaling, and use any zoom functions to examine any portion of the already-plotted traces. The Display can be reset between AutoScaling and fixed scaling/scrolling modes. Offset can also be adjusted for any individual trace using either the Display/Offset command or by manipulating the filename icon buttons. Colors of any trace(s) can be changed using the Display/Colors command, and additional graphics such as pump or temperature profile plots can be placed in the active window or removed.

Of course, you can also use the SoftButtons to switch to any of the other program screens in the same channel, or to move to any other channel, without interfering in any way with the run in progress. If more than one channel is configured, you can use the Monitor Timebases display function to show windows for different Timebases on a single screen.

5. When the Method DURATION is finished, the run will stop and the green RUN button will reappear on the left switchbar. The active run indicator will disappear from the top bar. The chromatogram file(s) will be integrated and analyzed automatically if specified in the Method. If multiple chromatograms are to be integrated and analyzed, the first file will be handled with the current Integration and Calibration screens, and then the next file will be integrated after loading its specified Integration file from the Method. Any reports specified in the Integration files being used for the Method will be shown on screen or printed (via Microsoft Excel) following the completion of analysis of each file.

Note that at the completion of each run all chromatogram files are already saved to disk. They will be named according to the instructions in the Method, if any, before the Method ends.

After RUN is pressed to begin a run, the RUN SoftButton will be redrawn as a STOP button. You can click the STOP button at any time to abort the run. If you abort a run before the DURATION time of the currently-loaded Method is completed, you will be prompted as to whether or not you wish to save the "incomplete" data file now written on the disk, and whether you wish to make reports using the partially-acquired chromatogram. If you indicate yes, you will be asked to assign a file name for each captured input signal you wish to save.

8.0 Programming and Executing Runs: The Method Screen

The Data Ally Method program screen manages the definition, setup, and execution of chromatography runs performed with the analytical instruments connected to each time base. You can view the current Method by clicking the left mouse button on the METHOD SoftButton on the left side of the screen at any time.

8.1.0 Method Functions

A Method determines the following parameters:

- The active detector signal inputs and their respective sampling rates.
- The duration time for the run.
- The path and filename for the chromatogram corresponding to each active input, if such files are to be saved.
- The Integration/Peak Identification program which is to be used to process each chromatogram for peak finding and generation of result reports.
- The Calibration program which is to be used to perform or update a calibration, or to compute results.
- Any information pertinent to the sample being injected, including correction factors, internal standard amounts, name, number, and other descriptive items.
- Timed control of relays or contact closures required for system operation.
- Timed control of pumps, GCs, detectors, fraction collectors, or any other programmable instruments.

In other words, every procedural aspect of a run is controlled by the information in a Method.

The Method program consists of and is displayed in two portions: a "top-level" dialog box containing global instructions relevant to the entire Method, and an "Events Table" which contains specific directions for controlling all instrumentation, relay flags, and conditional response mechanisms. The Method screen display should be considered as a "window" showing only one part of the overall Data Ally software on a single screen.

Both the "top-level" screen and the Events Table are saved together comprising a Method file. By loading a saved Method, the user can recall all the parameters associated with automatic operation for a particular application at once. Methods can be loaded manually or automatically as part of a Sequence.

Only one Method program can be present in any time base at any one time. In Configuration, a "default" Method filename is defined, which will always be loaded automatically into the Method screen whenever Data Ally is initialized or when the Method screen is cleared. Configuration also permits many of the settings, display fields, and options in the Method screen to be enabled or disabled, allowing total customization of Method contents for each individual time base. The contents of the default Method file can be changed at any time by reconfiguring. It is recommended that each user set the default file for the most common run conditions, which will minimize the amount of additional programming needed later on.

➡ *To run the currently-loaded Method program at any time*, click on the RUN SoftButton at bottom right. The system will immediately load all the files specified in the current Method, change the Graph display according, and begin to display the signals from all active detector inputs in real time.

After the Method starts, the RUN SoftButton will immediately change to a "STOP" button with an octagonal red "stop" image.

➡ *To stop the currently-running Method*, click on the STOP button at any time. If you do not terminate the Method early, it will continue to run until its DURATION time is completed, and will then process any resulting chromatogram data as specified.

The currently-viewed display can be any screen (Graph, Method, Integrate, Calibrate, or Utilities) except the Sequence screen when starting and stopping Methods with the RUN/STOP SoftButton. If the Sequence screen is being viewed, the system will attempt to start the currently loaded Sequence instead of the current single Method.

8.2.0 The Top Level Method Dialog

When you first click on the METHOD SoftButton, you will see the top-level Method dialog box (Figure 8.1). If the system has just been initialized and a saved Method file has not been loaded, the screen will show all the parameter values associated with the default Method program.




Figure 8.1 Top-Level Method Screen

8.2.1 Method File Identification

The top bar of the Method dialog box shows the current Method filename. A new name can be assigned by using the Files/Save As command to resave the Method with a different name.

To make identification of Method files easier when selecting a Method using the Files/Load command, a description can be entered in the Save As dialog (see section 8.5.3 below) which will appear in the Load dialog box for each Method file.

The remainder of the top-level Method is divided into four "boxed" sections which group items related to the same functions for "Run", "Detectors", "Options", and "Results". At the bottom of the screen is a large center button which switches the display to the Events Table portion of the Method.

8.2.2 Setting Inputs Criteria

The "Detectors" boxed area (Figure 8.2) contains the names of all the active detector inputs which were defined for the current time base in Configuration. Each detector name has an associated check box to its left, which must be checked to enable acquisition from that input, and a Resolution

entry field to its right, which contains the effective sampling rate parameter for that input. Up to four different signal inputs can be set in this area, based on Configuration.




Figure 8.2 Detectors Setup Fields

➡ *To enable detector inputs for a Method*, move the mouse pointer to the check boxes to the left of the available detector names, and click, if necessary, to check the detectors to be active for this Method. If a box is unchecked, that detector signal will not be acquired nor displayed.

Note that unchecking to disable any detector in this box will cause that detector's "Results" fields below to be "grayed out" and inoperable - these fields can only be set for active inputs.

➡ *To set RESOLUTION (sampling rate) for each input*, click the left mouse button with the pointer inside the RESOLUTION field for any selected input. Type in the desired setting in points per second. If a number only is entered with no units code, the system will assume you are using the same units shown for the previous entry. If entering a time, the time should represent the expected width at baseline for the narrowest peak to be found with your application and samples. All RESOLUTION fields for all active detectors must have an entry before a Method can be started. In general, you can set the defaults for RESOLUTION for each channel to be appropriate for the type of procedure (e.g. HPLC) during Configuration, so that you will not need to make frequent adjustments when writing new Methods.

IMPORTANT: NOTE THAT AFTER TYPING ENTRIES INTO FIELDS IN ANY PROGRAM SCREEN, YOU MUST EITHER CLICK INSIDE A DIFFERENT FIELD, CLICK ON ANOTHER SCREEN FUNCTION, OR PRESS ENTER IN ORDER FOR THE LAST FIELD ENTRY TO BECOME "PERMANENT".

RESOLUTION values cannot exceed the equivalent of 400 points/second

8.2.3 Setting Run Criteria

Parameters in the Run boxed region include the total run time for the Method, or DURATION, and any sample-specific information, accessed by clicking the SAMPLE button.

➡ *To set the DURATION or run time for the Method*, click with the pointer inside the DURATION field and type in the desired total run time. The range for run time is determined in Configuration. A time value can be entered in minutes or seconds.

➡ *To set sample-specific information*, click on the SAMPLE button to open Sample Info dialog (Figure 8.3).

Figure 8.3 Sample Information Dialog

This dialog contains only those parameter entry fields which were specified for the Method screen configuration in this time base. You can click inside any of the fields and enter values for each parameter as desired. All this information will be available for inclusion in any reports made with this method - if correction factors such as dilution and weight, or internal standard amounts, are entered, these values can be used to compute reported results.

If Internal Standardization is used, the values for each Internal Standard peak may be defined in the appropriate Calibration file's Reference Peak Table. Those values will be used when a Method is employed to run a standard, if no IS AMT values are entered into the Sample Info dialog. Entries made for any standard Method in the IS AMT fields will override existing values for internal standard concentrations found in the Calibration file. ENTRIES OF IS AMT MUST ALWAYS BE MADE FOR UNKNOWN SAMPLES EITHER HERE IN THE SAMPLE INFO DIALOG OR IN THE SEQUENCE TABLE, or when using the Analyze command.

When all the desired information is entered correctly, click OK to return to the top level Method screen, or click Cancel.

8.2.4 Setting Results Criteria

The boxed area in the center of the screen labeled "Results" (Figure 8.4) instructs the Method to save and integrate chromatogram files which are captured from each detector input, and allows designation of specific calibration schemes if so desired. Default settings and filenames may be present in the fields in this section, entered in Configuration, or they may be blank.

	Filename	Integration	Calibration	Lev	Type
UV-Vis	uvdat	default.int	default.cal		
ECD	ecdat	default.int	default.cal		

Figure 8.4 Result Files Definition Fields

➡ To set names for the chromatogram files captured from any active detector(s), click the left button inside the FILENAME box adjacent to the name of a detector, and type in a path (optional) and filename up to five characters (the last three characters of the eight-character filename are used for automatic series numbering). If the current path is to be used, no path entry needs to be made. These filenames will be assigned to raw chromatogram data files produced when this Method is run. Every active detector input name must be provided with a valid filename here in order to start the

Method. Different files can be saved to different paths or disks, if desired. Note that the final three characters of each filename will ALWAYS be inserted automatically by Data Ally and will overwrite any characters entered in those positions in the file naming fields..

➤ *To set reports to be automatically produced for each detector input at the end of this Method*, click the Reports button at the right of each detector's filename entry field. The Report Set dialog (Figure 8.5) will appear.




Figure 8.5 Report Set Dialog

Up to four different sets/rows of report specifications are shown in this dialog, based on your configuration settings for this time base. To use any combination of the reports shown, click the Enable check boxes at left of each report Template field to checkmark those which are to be used in running and processing this Method. Alternatively, click to uncheck and disable those not needed. Unchecking each row "grays out" the other items in that row which are not required once that report has been disabled.

For each enabled/checked report row, click inside the Template field and enter a template *.XLS filename and path, or click on the adjacent right down arrow selector to see a listing of all report template files present in the current reports subdirectory. Click to highlight and select the one desired. If the report is to be saved as an Excel *.XLS file, click inside the Filename field and type a unique name and path for that report.

To include a separate page printout of the chromatogram information with any report, check the Print CG checkbox in its row. Also click to check (enable) or uncheck (disable) the Print and Show boxes for each report.

Repeat the same steps for each different report desired for this Method. When all reports are entered, click OK to exit, or Cancel to escape with no changes.

➤ *To select the Integration program used to identify peaks for each detector input*, click inside the INTEGRATION field for the first active detector name and type in the desired Integration file, which must exist on disk in the current directory path for this channel. To choose from a selection of all integration program files found in the current path for this file type, click the down selector arrow at right of the field - a list of integration filenames in the currently selected integration program path will appear. Click the desired name to enter it. Repeat this selection for all other inputs which will be active in this Method - chromatograms from any inputs not assigned an Integration program in this manner will not be integrated automatically. The Integration file(s) selected will define what types of reports are created at the end of the Method.

➤ *To specify automatic calibration with this Method, if the sample to be run is a standard*, click with the pointer inside the CAL box for any active input name and type in the desired calibration Type code, which can be either blank (No calibration), N (New calibration), R (Replace calibration for current level), or A (Average into current calibration).

NOTE: The CAL field should only be used (have a non-blank entry) if the sample to be run is a standard which will be used to update an existing Calibration file.

N indicates a new file will be started with this run (all previous data in the calibration file will be replaced), R means the current calibration level will be replaced beginning with this run, and A indicates the new standard's response will be added to the existing calibration table. If you specify a non-blank value in the CAL field, you must also give a level value in the LEV field, which selects the concentration level to be assigned to the standard material being run. Click in the LEV field and type in the desired level number.

If you have configured the Method program in this time base without automatic calibration capability, the Calibration, Lev, and Type fields will be "hidden" and will not be available for programming. You can still perform manual or Sequence-based calibrations unless these functions have also been suppressed in Configuration.

➤ *To select a Calibration file to be loaded automatically for each input*, click inside the CALIBRATION field for that input and type in the path and filename of any current Calibration file. If an invalid Calibration filename is entered, an error message will appear when an attempt is made to start the Method advising that a Calibration file must be initialized and saved. If you intend to use this Method for quantitative result calculations, you should define a Calibration file for each input, whether the sample(s) to be run are standards or unknowns. If the Method will process a standard, the standard at the assigned LEV will be used to increment the assigned CALIBRATION filename automatically as defined in CAL. If the Method will run an unknown, the Calibration filename in CALIBRATION will be automatically loaded and that file will be used to compute the results for the unknown.

You have now provided the Method with all the essential instructions for acquiring and processing new chromatograms. Your Method could be run at this time by simply clicking on the RUN SoftButton at the lower left. There are, of course, additional parameters, settings, and commands you may wish to use, which are accessed by the buttons at the bottom of the screen.

8.3.0 Events Table

➤ *To show the Method Events Table*, click on the EVENTS TABLE button at bottom right. The screen will change to a horizontally split orientation (Figure 8.6) with a working Graph window at the top and a table on the bottom.




Figure 8.6 Method Events Table

8.3.1 Events Table Functions

The Events Table is generally displayed along with a Graph window so that many graphical features for modifying the table entries can be conveniently used. The top Graph window is exactly the same as the full screen Graph accessed by clicking the GRAPH SoftButton, and can be used for both real-

time chromatogram display while a Method is running, or for loading, reviewing, and editing chromatograms when idle. Either the Graph or the Table window can be made active by clicking the pointer anywhere inside either window. The Command Menu bar will change depending upon which window is currently active. Either the Graph or Events window can be resized using Windows conventions or moved around on the display at any time as desired.

The Events Table itself, like all other Data Ally tables, is a "spreadsheet" consisting of an expandable set of columns and rows. At the top of the Events Table is a special header line consisting of a set of "buttons" holding the name of each column entry. The selection of buttons present is defined in Configuration. Each button represents one parameter which can be controlled by the Method from the table, usually on a time basis according to the entry in the TIME column. When each channel is configured, only the types of columns which will actually be relevant to that channel's instruments and applications should be included in the Events Table. The order in the table and column sizes for each column button can be set in the Events Table screen as described below.

Programming a Method which includes control of various instruments or relay flags, time dependent responses to various conditions, or other time dependent events such as Zone Codes is performed by completing rows and columns in the Events Table. Each event which occurs at a specific time in the Method run is normally entered on a separate line labeled with that time value. If more than one type of event or setting is associated with a single time value, several columns may have entries in the single row corresponding to that time.

As a "spreadsheet", the Events Table acts like and can be edited like other types of spreadsheets found in other software applications. A series of buttons along the bottom edge of the table assists in performing spreadsheet operations such as inserting, copying, cutting, and pasting lines or ranges. An extensive set of "Fill" functions can be used to automatically enter lines and/or values into the table, and the interactive graphical editing tools permit fast and simple changes to any table parameter(s).

Whenever the Events Table is displayed, one "cell" is always highlighted, corresponding to the current "editing cell" into which entries can be made. The edit cell can be moved anywhere in the table by moving the cursor pointer to the desired cell and clicking the right mouse button. Alternatively, the <ENTER> or <TAB> key can be used to move the edit cell highlight across each line in the table column by column, from left to right. Pressing <ENTER> at the end of the last existing row in the table will cause a new blank row to be added.

Each row or horizontal line in the Table automatically receives a line number which appears in the column at far left. This column can be used not only as a guide when writing or editing Method programs, but also to disable any individual line(s) during Method execution.

Within each cell, you can place the vertical bar cursor to change the entry at any desired character, and the left and right arrow keys can be used to "step" through characters in any entry one character at a time.

When a Method is running and its Events Table is being executed, the current line being processed, or the last line which was processed, is "highlighted" in the table. The view of the table will automatically be adjusted so that the "highlighted" current line is always approximately two lines from the top line shown in the table window, as the highlight moves down through all the events lines. You can use the vertical scroll bar, or the PgUp/PgDn or CTRL-PgUp/CTRL-PgDn functions to move around in the Events Table while it is running. In order to change a running Method, you must either STOP the Method or wait until its DURATION has elapsed and it stops automatically.

8.3.2 Events Table Column Parameters

The Events Table parameters which may be selected as headers for table columns during Configuration are listed below. In most cases, all available parameters will not be used in any one

Table. Some parameters are appropriate only for certain types of instrumentation. Certain items such as TIME must always be present in the table, and these items are listed as "non-configurable". Note that the actual names for each functional column in the Events Table, along with the widths, colors, and order of those buttons/columns in the table, can be modified in Configuration.

You can easily re-order the columns in the Events Table from left to right or change their individual widths at any time. Both the column order and widths are set as defaults in Configuration. To move a column to another position in the table, simply place the pointer cursor on that column's header button, hold down the left mouse button, and "drag" the header button left or right to the new desired position relative to the other columns. When the pointer is in the desired position, release the mouse button and the entire column will be moved, and the table redrawn accordingly. Modifying the order of columns may make certain types of table editing or review easier by allowing placement of all columns currently being manipulated in view on a single screen - you can easily move "unimportant" column information to the right and off of the displayed table window area.

To change the width of any individual column, move the pointer cursor to the border of its header button and the next column header button in the direction in which you wish to expand or reduce the column width. The cursor will become a "double-sided horizontal arrow" (↔); when this change occurs, immediately hold down the left mouse button and move the arrow cursor in the direction you wish the button edge to expand or contract. When the desired column width is indicated by the position of the cursor, release the mouse button and the column will be resized and the entire table redrawn accordingly. This feature permits you to match the exact sizing of each column to the expected character width of the entries expected in that column - you can also effectively "remove" certain columns without actually eliminating them from the table by reducing their widths to a very small size.

➤ To directly enter a value into any cell in the table, simply click inside that cell so that it becomes highlighted and the vertical entry bar cursor appears. Then type in the desired information. If you wish to delete the entire previous entry, double-click the left mouse button with the entry bar inside the desired cell, so that the entire entry becomes highlighted - typing the first character of the new entry will cause the old entry to be completely deleted.

Available parameters in the Table are:

- **TIME: (Non-Configurable)** This specifies the time for each event or action in each row to occur. Every active row in the table must have a TIME value specified. TIME is entered as a value of minutes or seconds (a numerical value followed by a units designator code). If no units code is given, the units assigned to the previous entry are retained.

➤ To edit TIME, click the left mouse button with the pointer cursor inside any cell in the TIME column and type in the desired entry and units code. This column can also be edited using the AutoFill functions or graphical entry of controllable column parameters (see below).

- **%PUMPA, %PUMPB, %PUMPC, %PUMPD: (Configurable, HPLC Channel only):** The %PumpX columns set relative mixing percentages for from two to four pumps running simultaneously for gradient HPLC. These columns are only included in the Table based on the number of controllable HPLC pumps actually connected in this channel. Entries in each column are in percent, from 0 to 100. The system automatically adjusts entries into each of the columns so that the total per cent for all pumps is always 100%. The entries made in all the %X columns in the table are plotted on the channel's current Graph window if the Graph/Display commands are set accordingly.

When the Events Table is executed during a Method, the software sets all pumps as required at each time point when it finds a command to change any pump. Between the previous command and the next command, the system assumes changes in the pump status/gradient will be linear in nature, and ramps between time dependent states in a linear fashion.

➡ *To edit any of the %PUMPX fields*, click the left mouse button with the pointer cursor inside any cell in any of these columns and type in the desired value in per cent. If certain pumps are not used in certain Methods, their columns can be left blank. These columns cannot be edited using the AutoFill function. As soon as an entry is made and "entered" by pressing ENTER, CTRL-arrow, or TAB, or you click in any other cell in the Table, the entire row immediately updates as the system forces all pump % column entries accordingly to maintain a total of 100%.

- **FLOW:** (Configurable, HPLC or GC) This sets the total flow of solvent in either an HPLC or GC system. FLOW values in any column cells must be within the limits specified in Configuration for the Pump or GC in this channel, and are generally entered in the units specified in Configuration.

When the system encounters FLOW settings, it always assumes that changes from one flow value to the next value encountered are done linearly.

➡ *To edit the FLOW field*, click the left mouse button with the pointer cursor inside any cell in the FLOW column and type in the desired value and units code. Units can be ml/min, ml/sec, etc. If no units code is entered, the units designation for the previous entry is used. This column can also be edited using the AutoFill function (see below).

- **FOUT:** (Configurable) This column sets utilization of any external output relays or contact closures available in the Communications Processor dedicated to this channel, and defined in Configuration. Normally, each configured output relay flag is given a name identifier during Configuration to make it more convenient to program each flag in the Events Table. A flag is entered either by typing its number into the FOUT column. One flag can be used as many times, in as many different rows, as desired in a single Events Table program; any single FOUT column cell can contain multiple flag entries.

➡ *To edit the FOUT field*, click the left mouse button with the pointer cursor inside any cell in the FOUT column and type in the desired flag number.

Each time a flag number code is entered, the flag assumes the "polarity" defined for it during Configuration (High or Low = Closed or Open contact). If polarity is set high for a flag, that flag contact will close when it is executed and remain closed until it is defeated via a second later entry of the same flag in the Events Table. Polarity for flags can only be changed via Configuration. The time duration for each flag execution is adjusted automatically by Data Ally.

You can also use AutoFill editing (below) to change FOUT column entries.

- **FIN:** (Configurable) This column sets utilization of any external input relays or contact closures available in the Communications Processor dedicated to this channel, and defined in Configuration. Normally, each configured input relay flag is given a name identifier during Configuration to make it more convenient to program each flag in the Events Table. A flag is entered by typing its code (S or 1 through 3) into the FIN column. One flag can be used as many times, in as many different rows, as desired in a single Events Table program. The S input flag is usually reserved for remote start of the data system from a contact closure from another device.

➡ *To edit the FIN field*, click the left mouse button with the pointer cursor inside any cell in the FIN column and type in the desired name and/or flag code. As for FOUT, the FIN column can be modified using AutoFill editing (see below).

Each flag entry assumes the "polarity" defined for it during Configuration (High or Low = Closed or Open contact). If polarity is high, a flag will cause the system to "hold" execution if the contact being read is open, and the flag will execute and "unhold" when the contact closes. Note that the time each flag is executed is adjusted automatically.

In the most typical case, you can program the Events Table to use input flag S, Remote Start, to start the Method from an external contact closure by entering "S" into the FIN column at TIME = 1 sec or greater in a Method. Sensing the contact closure will produce the same result as clicking on the RUN SoftButton. A hold-to-inject flag should not be entered at TIME = 0.

- **ZONE:** (Configurable) This permits entry of integration zone codes such as ON/OFF to be executed directly from the timed events table.

➡ To edit the ZONE field, click the left mouse button with the pointer cursor inside any cell in the ZONE column and type in the desired zone code.

- **COMMENT:** (Configurable) This permits entry of row-specific comments as part of the Events Table. It is often useful to include remarks or comments describing each step in the table as an aid to later understanding of how and why the table was created.

➡ To enter text into the COMMENTS column, click inside the column and type the desired entry, or edit the previous entry.

Other columns may be added to the Events Table as required for certain special instrument control interfaces, such as electrochemical detectors (as in Figure 8.6).

8.3.3 Entering Events Table Parameters

There are three ways to enter values for any parameter into the Events Table:

1. Move the pointer into any "cell" in the table and click so that the highlight and the vertical bar entry cursor are moved into that cell. Then, type any desired entry into the cell directly, or, if appropriate, click on the right arrow button to view a listing of all possible entries from which the desired item can be chosen. To move to any other cell, either point to a different cell and click the left mouse button, or press the ENTER key to move the highlight from cell to cell, left to right, and then downward row by row until the desired cell is reached.
2. Copy information from an already entered row into a new row using the COPY and PASTE commands. To copy a row, move the pointer to the row to be copied, and click the left button to place the highlight in that row. Click on the COPY button to make a copy of the highlighted row.

Now, move the pointer to the row just below where you want to insert the copied row, and click anywhere in that row to place the highlight. Click the left mouse button on PASTE, and the copied row will be pasted in as a new line above the highlighted row. You can now make any other changes to the new row by directly entering those changes into each cell.

3. Fill information into any column in the table using the FILL button below the table. The FILL function is a convenient way to fill any table column with a series of entries based on one or more entries already made, or to compute and enter a new series, freeing you from the need to manually type in those entries.

8.3.3.1 Using the FILL Functions

AutoFill operations are performed on one column in the table at a time.

To fill entries into a column, click inside the top cell in that column or the cell containing the entry on which other row values in the column are to be based. Normally a value must be inserted in the column cell as a basis for fill operations. The AutoFill dialog box (Figure 8.7) will appear.




Figure 8.7 AutoFill Dialog for Events Table Time Column

The AutoFill dialog contains three sets of selections: the "Operation" or mechanism of filling on the top portion, the extent or range of filling desired in the middle portion, and special Options for filling in the lower portion.

There are two filling Operations, all of which proceed from the specified Starting Row to the Ending Row:

- "Duplicate" - copy the entry value in the highlighted row of the current column into a range of column cells beginning either after the Current Cell, from the Start of Table, or at any Specified row number, into all the cells in that column until either the Next Entry (next row with a non-blank cell), the End of Table, or any Specified row number. You can also specify the total Number of Steps/cells, set a range for copying based on row time designations via First and Last Times, or overwrite current information in the range of cells specified by setting the appropriate Options.
- "Increment" - increment (or decrement) the value in the highlighted row of the current column by the "Step Value" and copy it into the next row, beginning with the Current Cell or at any Specified row number, and proceeding with consecutive increments until either the "No. of Steps" value is met ("End of Table") or until the next existing entry is found in the column ("Next Entry"). If the "No. of Steps" value requires that new lines be added before the end of the table or the next entry, respectively, new lines will be inserted if the "Allow Insertions" option is checked.

The Increment function can also be used in a "force" mode to automatically compute and fill row values with evenly-spaced increments between two consecutive row entries, or between all sets of consecutive row entries in the entire table. For "forcing", no Step Value is given; this causes the system to find the next consecutive column row entry from the starting point, compute the difference between the highlighted (start) and next entries, divide this difference by the number of empty lines between the two entries, and fill numbers into all intervening rows in even increments.

Note that in most cases the highlighted row must be on the "basis" value to define filling before the AutoFill dialog box is used. If a blank cell is highlighted at the time of fill selection, a blank value (empty) will be filled into the defined range in "Duplicate" mode, and "Increment" mode cannot be used.

Table columns which have numerical entries (e.g. TIME, FLOW) will have the "Increment" function available in their AutoFill dialogs. Columns with non-numerical entries (e.g. FOUT, FIN, COMMENT) will have only the "Duplicate" mode available. Because many instrument control columns include

automatic "ramping" of values linearly between consecutive entries (such as %X, FLOW, or OVEN), filling of intermediate values is not necessary to set gradient functions.

There are four general "ranges" for filling functions:

- (Start Row) "Current Cell"/(End Row) "To Next Entry" - this selection will cause filling in the selected mode based upon the highlighted cell value into existing or newly-inserted blank lines up to the next currently-existing entry in the column.
- (Start Row) "Current Cell"/(End Row) "End of Table" - this selection will cause filling of entries into all remaining rows in the current column, based on the highlighted cell, according to the mode selected and the value set for "No. of Steps".
- (Start Row) "Start of Table"/(End Row) "To Next Entry" - this selection will cause filling in the selected mode from the first row in the current column into existing or newly-inserted blank lines up to the next currently-existing entry in the column. If the highlighted cell is in the top row of the column when this command is used, its effect is identical to that of the "Current Cell" to "To Next Entry" range. This setting is most appropriate when using the "Fill From Graph" mode for Zone Code entries.
- (Start Row) "Start of Table"/(End Row) "End of Table" - this selection will cause filling of entries ranging from the first column row to the end of the table, according to the mode selected and the values set for "No. of Steps". If the highlighted cell is in the top row of the table when this command is used, its effect is identical to that of the "Current Cell" to "End of Table" range. This setting is most suitable when using the "Fill From Graph" mode for Zone Codes.
- (Start Row) "Start of Table"/(End Row) "Specify" - this selection will cause filling of entries from the current highlighted cell (at start of table) to a specified row number.

The "qualifying" commands and fields at the bottom of the AutoFill dialog provide great flexibility for the filling functions. These are:

- "Step Value" - when in Increment mode, this field sets the value, in units corresponding to the current column, by which the highlighted basis cell is to be incremented or decremented. If this value is not given during an increment operation between two existing entries in a column, the system will automatically compute and "force" a series of evenly-incremented values in empty or newly inserted rows between those entries.
- "No. of Steps" - when in Duplicate, Increment, or Fill From Graph modes, this field sets a limit to the number of rows in the current column which are to be filled with new values. If the number of existing empty rows in the specified range exceeds the "No. of Steps" value, only part of that range will be filled; if the "No. of Steps" value exceeds the number of existing empty rows, new rows will be added to the table and filled as specified (provided "Allow Row Insertions" is enabled).
- "First Time" - in Fill From Graph mode, this field allows selection of any starting time for filling a new or existing Events table with Zone Codes, other than the normal time = 0 start. This permits only a partial set of timed zone codes to be automatically filled in any one fill cycle.
- "Last Time" - in Fill From Graph mode, this field allows selection of any ending time for filling a new or existing Events table with Zone Codes, other than the normal "end of table" end. This permits only a partial set of timed zone codes to be automatically filled in any one fill cycle.
- "Allow Insertions" - this checkbox command, if checked, allows automatic insertion and filling of new rows in an existing table if this is necessary to meet the "No. of Steps" criteria in a Duplicate, Increment, or Fill From Graph operation.

- "Allow Progressions" - this checkbox command, if checked, allows values in column cells below the base (start) cell to be employed for filling subsequent rows or for "forcing" values of intermediate rows.
- "Allow Overwrites" - this checkbox command, if checked, allows replacement of existing values in rows in the current table with new values over the selected range in any filling mode.

You can choose various combinations of options in the AutoFill dialog as described below to accomplish the following tasks (click on OK after setting up the Fill dialog as directed):

➤ To fill an existing row's value into all currently empty rows in the column up to the next entry, place the edit highlight on the "basis" row, click FILL to see the AutoFill dialog, select "Duplicate", "Current Cell", and "To Next Entry", with all other fields and options blank/off. If the "basis" row is the first row in the column, you can select "Start of Table" instead of "Current Cell" without highlighting the basis row.

➤ To fill an existing row's value into all currently empty rows in the entire column, place the edit highlight on the "basis" row, click FILL to see the AutoFill dialog, select "Duplicate", "Current Cell", and "End of Table", with all other fields and options blank/off.

➤ To fill an existing row's value into more or less empty rows than the number of empty rows found in the current column up to the next entry, place the edit highlight on the "basis" row, click FILL to see the AutoFill dialog, select "Duplicate", "Current Cell", and "To Next Entry". Type the actual number of empty rows to be filled into the "No. of Steps" column, with all other fields and options blank/off. If the "No. of Steps" value exceeds the number of actual intervening empty rows, sufficient new rows will be added and filled to meet the "No. of Steps" criteria.

➤ To fill all currently empty series of consecutive rows in the column with the values found in the rows prior to each empty group of rows, click FILL to see the AutoFill dialog, select "Duplicate", "Start of Table", "End of Table", and "Allow Progressions" with all other fields and options blank/off.

➤ To overwrite existing values (along with any blank rows) with duplicated values in an entire column, place the edit highlight on the "basis" row for copying, click FILL to see the AutoFill dialog, select "Duplicate", "Current Cell", and "End of Table", and check the "Allow Overwrites" command, with all other fields and options blank/off.

➤ To clear (remove) existing values in an entire column, insert a blank row at the beginning of the table, place the edit highlight on that row, click FILL to see the AutoFill dialog, select "Duplicate", "Current Cell", and "End of Table", and check the "Allow Overwrites" command.

➤ To overwrite existing values in a portion of a column with a duplicated value, place the edit highlight on the "basis" row for copying, click FILL to see the AutoFill dialog, select "Duplicate", "Current Cell", "End of Table", and "Allow Overwrites", and enter a "No. of Steps" value corresponding to the number of existing rows to be filled with duplicates. All other fields and options should be blank/off. If the number of steps exceeds the number of blank rows in the column, the "Allow Insertions" box must be checked to add sufficient rows to meet the "No. of Steps" criteria.

➤ To clear (remove) existing values in a portion of a column, insert a blank row at the beginning of the range to be cleared (or remove an existing row entry), place the edit highlight on that row, click FILL to see the AutoFill dialog, select "Duplicate", "Current Cell", "End of Table", and "Allow Overwrites", and enter a "No. of Steps" value corresponding to the number of existing rows to be cleared (filled with a blank entry). If the number of steps exceeds the number of blank rows in the column, the "Allow Insertions" box must be checked to add sufficient rows to meet the "No. of Steps" criteria.

➤ To fill a column in a new table with incremented values, enter the beginning "basis" value in the first row, place the edit highlight on that row, click FILL to see the AutoFill dialog, select "Increment", "Current Cell", and "End of Table", enter the increment value in the "Increment" field, and enter the number of rows for which incremented values are to be inserted in the "No. of Steps" field.

➤ To fill a column with incremented values in blank rows up to the next existing entry, place the edit highlight on the "basis" row, click FILL to see the AutoFill dialog, select "Increment", "Current Cell", and "To Next Entry", and enter the increment value in the "Increment" field; leave other fields and options blank/off.

➤ To fill an existing row's incremented value into more or less empty rows than the number of empty rows found in the current column up to the next entry, place the edit highlight on the "basis" row, click FILL to see the AutoFill dialog, select "Increment", "Current Cell", and "To Next Entry". Type the increment value into the "Increment" field, and the actual number of empty rows to be filled/incremented into the "No. of Steps" column, with all other fields and options blank/off. If the "No. of Steps" value exceeds the number of actual intervening empty rows, sufficient new rows will be added and filled to meet the "No. of Steps" criteria.

➤ To fill all currently empty series of consecutive rows in the column with incremented values based on the rows prior to each empty group of rows, click FILL to see the AutoFill dialog, select "Increment", "Start of Table", "End of Table", and "Allow Progressions". Type the increment step value for each row into the "Increment" field, with all other fields and options blank/off.

➤ To fill a column with incremented values in all blank rows based on an existing entry, place the edit highlight on the "basis" row, click FILL to see the AutoFill dialog, select "Increment", "Current Cell", and "End of Table", and enter the increment value in the "Increment" field; leave other fields and options blank/off.

➤ To fill a range of rows with evenly-spaced, automatically computed values, with each row's value corresponding to the difference between the basis cell's value and the value of the next row having a non-blank entry in the column divided by the number of intervening rows, place the edit highlight on the "basis" cell (first row in the range), click FILL to see the AutoFill dialog, select "Increment", "Current Cell", and "To Next Entry". If you wish to fill only the existing number of blank intervening rows, leave the "No. of Steps" field and all other commands blank/off. If you wish to fill either only some of the intervening rows, or to add new intervening rows, enter a number in the "No. of Steps" field specifying the total number of intervening rows to be filled. If "No. of Steps" exceeds the number of existing blank lines in the range, you must check the "Allow Insertions" box to automatically add and fill sufficient lines to match the "No. of Steps" entry.

➤ To fill all existing ranges of empty rows in the column with evenly-spaced, automatically computed values, with each filled row's value corresponding to the difference between the value of the cell immediately preceding that range and the value of the next row having a non-blank entry following the range divided by the number of intervening rows, click FILL to see the AutoFill dialog, select "Increment", "Start of Table", "End of Table", and "Allow Progressions". Leave all other selections blank/off.

➤ To overwrite existing values (along with any blank rows) with incremented values in an entire column, place the edit highlight on the "basis" row for incrementing, click FILL to see the AutoFill dialog, select "Increment", "Current Cell", and "End of Table", enter the increment value in the "Increment" field, and check the "Allow Overwrites" command, with all other fields and options blank/off.

➤ To overwrite existing values in a portion of a column with incremented values, place the edit highlight on the "basis" row for copying, click FILL to see the AutoFill dialog, select "Increment", "Current Cell", "End of Table", and "Allow Overwrites", enter the increment value in the "Increment" field, and enter a "No. of Steps" value corresponding to the number of existing rows to be filled with duplicates. All other fields and options should be blank/off. If the number of steps exceeds the

number of blank rows in the column, the "Allow Insertions" box must be checked to add sufficient rows to meet the "No. of Steps" criteria.

Use of AutoFill commands will only overwrite existing cell entries in any table if you have checked the "Allow Overwrites" option checkbox.

After using the automatic fill function, you can of course make any desired changes in any of the filled times before saving the Method program and its Events Table. The FILL feature can be used repeatedly to set up any series of table columns.

8.3.3.2 Method Change Restrictions

Every change made in a current Method screen in any channel is preserved until one of the following occurs:

- Another change is made to that parameter or value by any available means.
- The Method is overwritten by loading a new Method with the Files/Load command.
- The Method is cleared and the default Method loaded with the Files/Clear command.

Before modifying or overwriting a Method, you should save any changes to that Method using the same Method name or as a different Method. Each time you attempt to overwrite a previous Method by resaving it using the same filename, a "trap" message will appear prompting you to confirm that you wish to replace the old Method.

8.4.0 Manual Operation

After you click on RUN to begin a Method, the instructions in the Method program will control the entire operation of the channel while that Method is in process. It is possible to interrupt or override the programmed control in several ways, which permits you to make changes in the executing Method in real time, to manually operate any controllable instrumentation such as HPLC pumps, or to "pause" or "hold" the Method for any length of time until you are ready to proceed.

While a Method is running, you have access to changing any time-dependent steps or events in the Events Table. However, you cannot change the commands and parameters set in the "top-level" Method dialog after a Method has begun to execute.

8.4.1 Hold Function

When a Method starts, you can observe its progress by viewing the Graph screen, with any display criteria desired, or you can select either the "top-level" Method screen or the Method Events Table screen with its companion upper Graph window. If you have configured the top level method to display information such as actual signal values for operating detectors, or current wavelength values, these will be shown in real time in the appropriate fields - if the Status Box is active in the current channel's Graph screen, it can also be viewed while observing the top-level Method. The Events Table screen will show both the Graph window with the chromatogram(s) currently being acquired and Status Box, and the course of execution of any lines and events listed in the table. As the Method proceeds, the table lines will "index" downward with a "highlight" on the currently processed row. At any time during a Method, you can use the scroll bars in the Events Table to change the display to see any line(s), whether or not they have already been executed. The "running" highlight will remain on the current (last executed) line as a constant indicator of the position of that line in the Events Table.

The HOLD button at the bottom of the Events Table can be used to "pause" the Method at any time.

➤ To hold or pause the Method at its current status, click the HOLD button. If the Status Box is active, its registers will remain active, but all its status indicators except for time and signal input values should remain constant during the hold.

➤ To clear a hold and proceed with the current Method program, click on the HOLD button again. The Method will automatically resume its programmed status at the now current time - this may require rapid execution or even skipping of some number of programmed lines in the Events Table to "catch up" to the programmed state of the system for the present time. Data Ally will attempt to execute all delayed rows in the Events Table as quickly as possible after a hold is cleared.

Normally the timer will continue to run during a hold, unless the hold is generated by an Input flag such as "wait for injection". This is because even if you wish to temporarily stop further execution of a Method program, it is important to have a valid record of overall elapsed time once the sample is on the column.

8.5.0 Working with Method Files

The Method program filing system permits you to save and recall complete Methods quickly and easily. The Method screen display, comprising both the top-level dialog box and the Events Table, is configured as a "template" to show information about the current chromatogram acquisition and processing instructions loaded into the channel being observed. This template can be edited using many graphical and other tools to change its contents. The contents of the Method screen display can be saved as a "Method" file, with a unique filename and description, so that it can be recalled and reused at any future time either manually, or as part of an automated Sequence of Methods.

The default Method file, whose name is assigned during Configuration, is always loaded automatically into the Method display screen at the time Data Ally is initialized. You can change any aspects of the default file, including its Events Table, within the active software, except for the features and options which are either "fixed" or "removed" in the Configuration. These items can, of course, also be reset or restored, but the channel must be disabled and a re-Configuration performed to accomplish this.

The Files commands in the Method Command Menu bar (Figure 8.8) are used to load, save, clear, and delete Method program files.




Figure 8.8 Method Files Commands Submenu

8.5.1 Loading Methods

➤ To load a saved Method file from disk memory into the current channel, click the left mouse button on Files in the Command menu bar, and then click on Load. A directory dialog box will appear (Figure 8.9) which lists the Method files (*.MET) available in the current directory path for this channel. If any Description has been entered for individual Methods, this information will appear at the bottom of the directory dialog. Select the desired disk, path, and file by clicking on the desired list items, or by typing the path/filename directly into the entry box. Click on OK to load the new file.




Figure 8.9 Method Files/Load Dialog

If the previous file has been edited and has not been saved or resaved, an information message will appear before the new file is loaded asking if the previous file should be saved. If so, click on YES and assign either the same filename or a new name as desired. If you do not wish to save the edited file, click on NO to proceed. If you are operating in GLP Mode and the previous Method file was used to generate or reprocess chromatograms, you will be forced to save this file if it was not saved before in its current state, before a new Method file can be loaded. When a new file is loaded, it will overwrite the previous Method file.

Load Default File: The Load Default File command causes loading of the Method file defined as the default method in Configuration (see Section 5.2.5.2). Since the default Method file, like all Data Ally defaults, is a "read-only" file, it cannot be changed and resaved, unless the system is reconfigured.

Use the Load Default File command to restore the default method whenever you wish, regardless of what other Methods are saved and/or loaded in the current time base.

8.5.2 Clearing Methods

The Clear command allows all entries from the Method Events Table to be cleared in one step. Global parameters in the top-level Method dialog must be individually deleted.

➤ **To clear the current Method Events Table**, click on File in the Command Menu bar and then click on Clear. A trap message will appear asking you to confirm the clear (Figure 8.10). Click OK, and the Events Table will be cleared. If the previous file has been edited and has not been saved or resaved, an information message will appear before the table is cleared asking if the previous file should be saved. If so, click on YES and assign either the same filename or a new name as desired. If you do not wish to save the edited file, click on NO to proceed.




Figure 8.10 Clearing Method Events Table

8.5.3 Saving Methods

➤ *To save a Method file*, click on Files in the Command Menu bar, and then click on Save. The Method will be resaved under its current name. A new version number will be assigned and the date/time of save and current operator name will be recorded.

➤ *To save a Method with a new name*, click on Files and then Save As. The Save dialog box will appear (Figure 8.11). Click inside the filename entry field and type in the desired new filename. If you wish, click inside the DESCRIPTION field and type a description of the new file, which may help identify this file in the directory listing for future recall. When the correct name and description appear, click on OK to save the file to the current disk path. If another file already exists in that path with the same name, an error message will instruct you to change the filename before saving.




Figure 8.11 Method File/SaveAs Dialog

It is recommended that Method program files saved to disk be "backed up" by simply copying the *.MET files to a floppy disk, removable hard disk, or tape and preserving the copies in a secure place. If modifications are routinely made to the library of *.MET files, backups should occur at routine intervals or even on a daily basis.

8.5.4 Deleting Methods

Use the Windows File Manager Delete function to delete any Data Ally *.met files which are no longer useful.

8.5.5 Printing Methods

➡ *To print a summary of the current Method file*, click on Print Method in the File submenu - the Reporting status window (Figure 8.12) will appear. This window indicates which Method file (the default method is part of the configuration *.cfg file) is currently loaded, the report template file configured for printing the method summary, and any configured saved report name.




Figure 8.12 Reporting Status Window

The method file printing template saved as \reports\met\met.xls can be modified as desired by using Excel (see Section 13).

9.0 Integration Screen: Finding and Identifying Peaks

The Integration program screen, which can be viewed by clicking on the INTEGRATE button on the left SoftButton bar, contains the instructions required for finding peaks during integration of chromatograms, and for identifying individual peaks based on retention times, peak referencing, and other information. It also specifies the type(s) of reports to be prepared after each chromatogram file is processed. You may wish to rename the Integration screen SoftButton "Peak Ident." to clarify its true function for other operators.

Like the other Data Ally programming screens, the Integration screen display should be considered as a "window" showing one portion of the overall software settings utilized in and associated with the current channel. The contents of the Integration screen are in two sections: the so-called "top-level" dialog box, which includes global parameters, and the Peak Table, which comprises a listing of assigned peaks and their identification criteria. Both sections are contained in a single Integration program file.

The values associated with each of the Integration screen fields have defaults which are set during Configuration. Also, various options can be enabled or disabled in both the top-level dialog box and the Peak Table based on the Configuration. Each user can modify all parts of the Integration program screen and save as many modified versions as desired as unique "Integration program files", which have the file extension code ".INT". To distinguish Integration files appropriate for one configured channel from those for a different channel, each channel usually employs a separate directory path for archiving of its program files, that path being assigned during Configuration.

Since the parameters associated with the Integration screen and its Peak Table are essential for identifying peaks, all Data Ally Calibration files are constructed based upon already-defined Integration files. Therefore, a valid Integration file and Peak Table must always be generated for a given chromatogram type before any calibration and quantitations can be done. One Integration file can have an unlimited number of independent and separate Calibration files which use it as a reference, but each Calibration file can have only one Integration file upon which it is based.

9.1.0 Top-Level Integration Screen

The top level Integration screen (Figure 9.1) is used to name the Integration file, to set the global integration algorithm for that file, to specify options for peak finding and calibration, and to select the names and types of reports to be created after processing chromatograms automatically. A large button at the bottom of the screen is used to switch to a view of the Peak Table portion of the Integration file. A second button is used to enter Comments describing the nature of the Integration file.




Figure 9.1 Top Level Integration Screen

9.1.1 File Identification

Like Methods, the filename for the current Integration program file is defined when the file is saved, and is shown in the top window bar. When each channel is first initialized, its Integration screen will display the contents of the default Integration file which has been defined in Configuration using the default file name. If you wish to change any of the information in the Integration screen, you can save the modified Integration file using the same default file name, or as a new file with any desired file name. The extension for filenames of all Integration files is always .INT; only .INT files will be shown in the Integration file directory system and only .INT files can be loaded into the Integration screen. Descriptions of Integration files can be entered at the time the Save As command is used (see Section 9.3.3).

9.1.2 Integration Algorithm

The center portion of the top-level Integration dialog box is portioned into three "boxed" sections. The section on the right labeled "INTEGRATION ALGORITHM" lists a series of entry boxes for the values of each of the integration algorithm settings (Figure 9.2). When the default Integration file is loaded, the values shown in these fields are those assigned as defaults in Configuration. They can, of course, be changed in any way by directly entering new values, or by replacing existing values during manual integration in the Graph screen (which causes the current entries in the ALGORITHM PARAMETERS box to change accordingly).

NOTE: You may click on "**Automatic**" to have the software try to determine the appropriate integration parameters. If you are not satisfied with the results, you may then adjust them as needed.

Integrate		
Slope Sensitivity	10	uV/s
Slope Interval	5	5
Baseline Drift	5	uV/s
Noise Reject	30	ms
Minimum Peak Height	1000	uV
Minimum Peak Width	5	s
Minimum Peak Area	100	uVs
Dropline-Skim Ratio	10.0	

Buttons: OK, Cancel, Automatic, Help, ☐ View Only

Figure 9.2 Integration Algorithm Parameter Fields

9.1.2.1 Algorithm Parameters

The Integration Algorithm parameters are defined as:

Slope Sensitivity: The value of baseline slope (first derivative) at which the integrator determines the start of (positive value on leading edge) or the end of (negative value on trailing edge) a peak. Beginning and end points of peaks are detected as the first point or bunch of points found at which the slope drops below this slope threshold value (Figure 9.4). The default range for slope sensitivity is zero to 1,000,000 microvolts/sec, or zero to 1 Volt/sec. As the value for slope sensitivity decreases, the sensitivity of finding peaks increases, and smaller or broader peaks are detected.




Figure 9.3 Use of Slope Sensitivity Parameter

Slope Interval: The number of consecutive bunched points which are analyzed for each slope/derivative computation when analyzing baseline slope. If the slope interval increases, the sensitivity of finding peaks decreases, since slopes are averaged over larger ranges of signal values. The default range for slope interval is 1 to 100 points. Slope interval can also be expressed in time units (sec or min) if desired.

Baseline Drift: The maximum positive or negative slope possible for baseline drawn under single peaks or groups of co-eluting peaks (Figure 9.4). Baseline drift is often used to limit the integrator in situations where valley-to-valley baselines might be drawn for a series of peaks, in order to force a flatter overall baseline. The default range for baseline drift is zero to 1,000,000 microvolts/sec, or zero to 1 Volt/sec. As baseline drift value decreases, flatter, more horizontal baselines are forced underneath regions of co-eluting peaks.




Figure 9.4 Use of Baseline Drift Parameter

Noise Reject: The noise reject parameter is not functional in v. 1.00.

Minimum Height: The height threshold which must be exceeded for any peak to be detected. The default range for minimum height is zero (no threshold) to 1,000,000 microvolts. This threshold value is expressed as the height of each peak above the apparent baseline. If its value is set at zero, all peaks of any height can be detected and reported.

Minimum Width: The minimum number of points, or minimum time value, for the width of any detected peak at baseline. Peaks which do not exceed this setting will not be reported. The default range for minimum width is 3 points to 1,000,000 points. This value can also be expressed as time in min or sec.

Minimum Area: The area threshold which must be exceeded for a peak to be detected and reported. The default range for minimum area is zero (no threshold) to 10,000,000 microvolt-sec. If this value is set at zero, all peaks of any area can be detected and reported.

Drop Ratio: The drop ratio parameter is not functional in the current version.

9.1.2.2 *How Integration is Performed*

When a chromatogram is integrated manually, or automatically during Method execution, a global integration algorithm is executed which makes use of the algorithm parameters found in the Integration screen dialog box. The global integration process consists of three fundamental steps, in which a chromatogram is scanned, segments containing peaks are identified, those segments are broken down into individual peaks and groups of peaks, and, finally, coeluting peak groups are analyzed and their component peaks separated.

A brief description of the integration process is as follows:

1. The chromatogram to be integrated is scanned and the highest signal value present is found and marked. This point is assigned as the top of the highest peak. The integrator then scans in both directions from the maximum signal point to find the leading and trailing edge thresholds, if they exist, where baseline slopes fall below the Slope Sensitivity value. These points are tentatively marked as begin and end points of the primary peak if they are located.

The integrator then treats the two segments of the chromatogram before and after the borders of the primary peak in the same manner as above, searching for the highest signal value in each segment and then attempting to find borders for each peak based on the Slope Sensitivity threshold. This process is reiterated in as many cycles as necessary until the entire chromatogram has been tentatively classified into one of three categories:

- a. A peak with defined beginning and end thresholds which meet the Slope Sensitivity criteria.
- b. A non-peak "baseline" region which falls between two peaks in category (a).
- c. A "valley" or "saddle" segment which falls between the maxima of two peaks but in which region no leading or trailing edges of peaks can be defined because either the slope sensitivity criteria are not exceeded or because they are not exceeded for a region of more data points than are specified by the Minimum Peak Width value.

When this "stage one" integration is complete, every potential peak has been found.

2. The integrator finds all segments in the chromatogram which correspond to "peaks" bounded by "baseline" on both sides. From left to right, the integrator chooses the first such segment, scans down either side of the "peak", compares the slope/derivative to the Slope Sensitivity threshold, and, when the threshold is passed, uses the Noise Reject value to select a number of data points to test beyond the threshold point to determine if any other peak is encountered.

If this test on either side of the peak finds no unexpected "noise" or sub-peaks, trial baseline is drawn under the peak between the threshold locations. A final test is done to be certain that no "negative area" has been denoted under the trial baseline. If "negative area" is found, the baseline is "backed off" until the negative area is eliminated. The trial baseline is then fixed in that segment, and the process is repeated for all other similar segments.

When this "stage two" integration is complete, all baseline-resolved peaks have been fully identified and marked.

3. For those segments classified as "peak groups", based on the observation that the Slope Sensitivity threshold is not met within the group area but only at its edges, trial baselines are drawn "valley to valley" under each group. The slopes of those valley-to-valley baselines are tested against the Baseline Drift parameter. If Baseline Drift is exceeded by any valley-to-valley baseline, the trial baseline is forced to a lower slope so that its slope does not exceed the drift threshold. All group baselines are tested and redrawn in this manner.

When this "stage three" integration is complete, all co-eluting peak groups have been drawn with baselines, but individual peak components in each group are not yet resolved.

4. For each peak group, the largest (highest maximum signal) peak is found. The integrator scans down from that peak's maximum value along both sides until it finds the first "valley" or location on either side where the Slope Sensitivity threshold is not passed for a number of points exceeding the minimum peak width threshold. At those positions, lines are dropped to the current baseline and the area of the largest peak within those drops is computed. On the trailing side, all the remaining area left in that group of peaks before the next "baseline" segment is computed, and compared to the area of the primary peak. If the ratio of the primary to remaining areas exceeds the Drop Line Ratio value, the drop line is considered valid, and the process is repeated on each side of the large peak using the next peaks found. In any cases where the Drop Ratio is not exceeded, the trailing peaks are skimmed. Leading edge peaks are skimmed only if a Forward Skim Zone Code is present allowing this treatment.

When this phase is completed, all co-eluting peak groups are resolved via either drop lines or skims according to the integration algorithm criteria.

5. For each skimmed peak in each co-eluting group, a final test is applied to determine if a skim or tangential separation should be employed. This test is only performed if the Tangent option is enabled in Configuration - if it is not, the integration finishes at the end of step 4 above.

A test tangent line from the point at which each skim begins to baseline is drawn for each skim line, and its slope is computed. The slope of the current skim line is also computed, and the ratio of the Skim/Tangent slopes is found. If this ratio is less than the Skim Ratio value, the skim line is used. If the ratio exceeds the Skim Ratio value, a tangent is substituted for the existing skim and the integrator repeats the process for the next skimmed peak.

When this phase is completed, all "shoulder" peaks have been classified into either skimmed or tangentially-separated peaks and all final baselines drawn.

After completing step 4 or step 5 (if Tangent drawing is enabled), the integrator displays all retention time values computed for all peaks, draws all tic marks and peak labels, and draws peak names and numbers if available from the peak table. If relative retention times are specified in the Integration screen (see below), the integrator finds the appropriate reference peaks or points before computing and displaying retention times.

It is important to remember that when running Methods, different chromatograms associated with various detector inputs can be assigned different Integration files, which can include algorithm parameter settings most suitable for each unique input. If only one Integration program file is specified for a Method when more than one detector input is active, the same Integration file will be used for automatic integration of all chromatograms from all inputs acquired with that Method.

9.1.2.3 Zone Codes

Zone Codes can be used to create time-dependent modifications in the way the integration algorithm is applied for any specified Method or chromatogram. Any combination of Zone Codes can be placed into a Method Events Table for automatic execution in any chromatogram produced using that Method. Alternatively, Zone Codes can be manually inserted into chromatograms at any time. The can be incorporated into a Method manually or by using the Fill From Graph function (Section 8.3.2) for automatic application to any chromatogram run using the Method. Each time zone codes are applied and utilized during integration, they are included with the chromatogram raw data and other baseline information if the chromatogram is saved or resaved after processing.

A list of all available Zone Codes is as follows:

- **Negative Peaks:** Enables integration of "negative" peaks with areas below baseline in time window until next instance of code is found. All peaks reported with "negative" areas (below baseline) will be marked to indicate this.
- **Integrate Off/Reset:** First code (earliest) defeats integration until second code is encountered. A Reset code enables Integration when it follows an Integrate Off code. Can be used to modify start/end positions for peaks and to ignore solvent fronts, column switch points where baseline shifts, etc.

The Data Ally file storage system saves all baseline codes resulting from each integration operation along with the raw data set for every chromatogram, enabling the results of the last integration performed to be immediately accessible each time a given chromatogram is recalled or reprocessed. If a previously integrated chromatogram is integrated again, the new set of baseline codes from the last integration overwrites the old codes when the data file is resaved.

9.1.3 System Suitability

9.1.3.1 System Suitability Calculations

System Suitability calculations are displayed for each individual peak selected in the Peak Table. The peak itself is displayed with the baseline used and the relevant calculation lines. Calculations where appropriate are calculated according to the USP or EP specifications.

The basic System Suitability calculations performed are:

Name	The Peak Name from the Peak Table
Retention Time	The actual Retention Time at peak top.
Capacity Factor	$\frac{RT_{\text{retainedpeak}} - RT_{\text{unretainedpeak}}}{RT_{\text{unretainedpeak}}}$
Area	The actual area of the peak
Height	The actual height of the peak
PW at 50%	The actual width of the peak at half height (min)
PW at 10%	The actual width of the peak at 10% height (min)
PW at 5%	The actual width of the peak at 5% height (min)

PW at Base The calculated width of the peak at base using tangent method (min)

CFR Asymmetry Factor $\frac{RT_{peakend} - RT}{RT - RT_{peakbegin}}$ at 10% peakheight

Tailing/Symmetry $\frac{RT_{peakend} - RT_{peakbegin}}{2 * RT - RT_{peakbegin}}$ at 5% peakheight

Plates USP = $\frac{16 * (RT^2)}{((RT_{peakend} - RT_{peakbegin})^2)}$ at base

EP = $\frac{5.54 * (RT^2)}{((RT_{peakend} - RT_{peakbegin})^2)}$ at 50% height

Plates/meter USP = $\frac{USP_{Plates}}{ColumnLength}$ EP = $\frac{EP_{Plates}}{ColumnLength}$

Resolution

USP = $\frac{2 * (RT_{Peak1} - RT_{Peak2})}{((RT_{Peak1end} - RT_{Peak1begin}) + (RT_{Peak2end} - RT_{Peak2begin}))}$ at Base

EP = $\frac{1.18 * (RT_{Peak1} - RT_{Peak2})}{((RT_{Peak1end} - RT_{Peak1begin}) + (RT_{Peak2end} - RT_{Peak2begin}))}$ at 50%

HETP mm USP = $\frac{ColumnLength}{USP_{Plates}}$ EP = $\frac{EP_{Plates}}{ColumnLength}$

Note: USP calculations use Tangent Method. EP Calculations use Half Width Method.

9.1.3.2 System Suitability Column Parameters

The Column Parameters options permit input of additional information used in some of the System Suitability calculations. Specific Unretained Peak information may also be input in this dialog.

The Column Parameters are:

Column Length	The actual length of the column in millimeters.
Particle Diameter	The size of the packing material in microns. This figure may not be used in the main report but may be accessed for special calculations.
Flow Rate	The actual flow rate of the solvent through the column. This figure may not be used in the main report but may be accessed for special calculations.

The Unretained Peak Parameters are:

RT	The actual retention time for the Unretained Peak. A value should be entered here if the unretained peak was not integrated in the chromatogram or a specific retention time is required for the calculations.
----	--

Name	The Name of the Unretained Peak if desired will be placed in the system suitability report.
------	---

Note: Do not enter a RT or a Name here if a **K** is used to mark the Unretained Peak in the Peak Table or vice versa.

System Suitability calculations may be performed on any peaks in a chromatogram. Normally peaks will be baseline resolved to obtain meaningful results.

With Data Ally, only selected peaks may have system suitability calculations performed on them. Normally at least two peaks will be selected to enable all parameters to be calculated. In addition the unretained peak should be identified in the Peak Table, if available, or in the Column Parameters section if a specific retention time is required.

To Program System Suitability Codes

Enter a '**S**' in the Suit Code column for each peak to have System Suitability Calculations performed on it. If the Unretained or Ko has been integrated, place a '**K**' in the Suit Code column corresponding to the peak.

If the Unretained peak has not been integrated the retention time and name may be inserted in the Column Parameters dialog on the top level of the Integration Screen.


Note: For all the calculations to be performed certain column information will need to be entered in the Column Parameters dialog in the top level of the integration screen.

Hint: The same Window value may be quickly programmed for all peaks in the Peak table by using The Fill Function.

9.1.4 Peak Identification

Only absolute retention times and windows and "largest area within window" peak logic are available for identifying peaks in the current version.

9.2.0 Integration Peak Table

 To view the Peak Table for the current Integration program, click on the large PEAK TABLE button at bottom center of the top-level Integration dialog. The Peak Table screen (Figure 9.5) will appear.




Figure 9.5 Integration Peak Table Window

The Peak Table is displayed in split-screen fashion beneath a working Graph screen. The Graph screen can be used to display chromatograms during real-time processing or during review, and is "connected" to the Peak Table below to allow graphical updates and modifications of the Peak Table. Click the left mouse button with the pointer anywhere inside either the Graph window or the Peak Table window in this screen, to make that window active.

➤ To return to the top-level Integration dialog box display, click on the INTEG button at the bottom right edge below the Peak Table.

9.2.1 Peak Table Functions

The Peak Table consists of a series of columns defined in Configuration, each of which contains information concerning peak component identification. By definition, every Peak Table must include columns for Retention Time and Window; all other parameters are optional.

Each Peak Table is a "spreadsheet" which operates in much the same fashion as common "spreadsheet" software programs. The Table "spreadsheet" consists of columns and rows, which can be typed into, copied, cleared, deleted, pasted, and appended to. The four buttons at left beneath the table marked Insert, Delete, Copy, and Paste allow rows in the table, and ranges of rows, to be manipulated easily.

Each Peak Table is defaulted to a single horizontal row. Rows are added by either clicking the left mouse button on the INSERT button at bottom, or by repeatedly pressing the ENTER key, which will move to the last column in the last current row and then create a new row below. You must create a new row before entering information into that row.

The INSERT button can be used to insert new blank rows within a series of already-completed rows at any desired position. Move the pointer to any cell in the row below where you wish to add a new row, and click the left button to place the highlight/cursor in that row. Then click the INSERT button to insert the new row above the highlighted row. Repeated clicks on INSERT will insert a series of rows.

The DELETE button is used to permanently delete rows in the table. Move the pointer to the row to be deleted and click the left mouse button to anchor the highlight/pointer in that row. Then click the left button on DELETE to remove that entire row and its contents.

The CUT button can be used to "cut" rows which are already in the table and subsequently "paste" them into a new position. Move the pointer to the row to be cut and click the left mouse button to anchor the highlight/pointer in that row. Then click the left button on CUT to remove that entire row and its contents - the row is placed into an internal memory and can now be recalled using the PASTE command.

Rows which have been cut can be replaced into the table in a different position by using the PASTE command after a CUT. Immediately after cutting a row, move the pointer and click to highlight anywhere in the row just below where you want to re-insert the cut row. Then, click on PASTE and the cut row or rows will be inserted. PASTE cannot be used after rows have been DELETED since deleted rows are permanently discarded. If a row is CUT and not PASTEd back into the table, the row is effectively deleted in the same manner as via the DELETE command.

You can easily re-order the columns in the Peak Table from left to right or change their individual widths at any time. Both the column order and widths are set as defaults in Configuration. To move a column to another position in the table, simply place the pointer cursor on that column's header button, hold down the left mouse button, and "drag" the header button left or right to the new desired position relative to the other columns. When the pointer is in the desired position, release the mouse button and the entire column will be moved, and the table redrawn accordingly. Modifying the order of columns may make certain types of table editing or review easier by allowing placement of all columns currently being manipulated in view on a single screen - you can easily move "unimportant" column information to the right and off of the displayed table window area.

To change the width of any individual column, move the pointer cursor to the border of its header button and the next column header button in the direction in which you wish to expand or reduce the column width. The cursor will become a "double-sided horizontal arrow" (↔); when this change occurs, immediately hold down the left mouse button and move the arrow cursor in the direction you wish the button edge to expand or contract. When the desired column width is indicated by the position of the cursor, release the mouse button and the column will be resized and the entire table redrawn accordingly. This feature permits you to match the exact sizing of each column to the expected character width of the entries expected in that column - you can also effectively "remove" certain columns without actually eliminating them from the table by reducing their widths to a very small size.

9.2.2 Peak Table Components

The contents of the Peak Table columns may be configured for any channel so that some column options are not used, or the order of the columns is different. The column header titles marked "Configurable" are non-essential for running the system and can be omitted from the Peak Table in Configuration if desired. The complete set of options are as follows:

- **ENABLE:** (Not Configurable) This column is both a graphic indicator/selector for whether a row entered in the Peak Table will be used for processing chromatograms and generating reports, and a "peak number" tracking column. If the line number cell in this column is highlighted for a given row, that row is active and will be used - active peaks will be identified as found in reports and calibrated as instructed. If the cell is not highlighted, that row is inactive and will not be used - inactive peaks will still be included in reports, but will not be named and/or calibrated, and no quantitative reports can be generated for such peaks. The ON column is useful for temporarily omitting rows assigned to specific peaks if those rows are not required for a given analysis, without actually deleting the information in that row.

➡ *To manually edit the **ENABLE** column*, move the pointer cursor to the desired cell in that column and double-click the left mouse button to toggle the highlight on or off.

- **RT:** (Not Configurable) This column shows the absolute retention time for the peak component identified in each row. A value for RT is essential for every row which defines a peak. RT is always entered in time units (min or sec).

➡ To manually edit RT, move the pointer cursor to the desired RT cell, click to highlight, and type in the desired value along with a units code, if desired. If no units code is entered, it will be assumed you are entering the value in the units used for the previous entry. The range for RT is -100 to 10,000 minutes. Retention times can be filled automatically from the current foreground chromatogram by using the FILL function (see below).

- **PEAK NAME:** (Configurable) This column shows the name assigned to the peak defined in the current row. Names of up to 30 characters can be given if the column width is set sufficiently large in Configuration.

➡ To manually edit PEAK NAME, move the pointer cursor to the desired cell, click to highlight it, and type in the desired name.

- **WINDOW:** (Non-Configurable) This column defines the retention window for the peak in the current row, as a range of time on either side of the RT value. Windows must be defined for every peak in every row, and must have a value greater than zero. If the top-level Integration dialog specifies Absolute windows, a time value in minutes or seconds is entered. The WINDOW setting must be entered for each row, and can have a unique value for each row. You can use the FILL function below to "fill in" Window settings.

The Data Ally integrator employs the WINDOW along with the RT to distinguish individual peaks when evaluating a chromatogram. The time window defined by each RT value in the current Peak Table, plus or minus the WINDOW value, is examined in each just-integrated chromatogram to determine if any peaks appear inside that region. If so, the largest peak appearing in the retention window is assigned the peak name given in the table. Retention windows are permitted to overlap adjacent windows, the logic rules being applied separately in each window.

➡ To manually edit peak WINDOW, move the pointer cursor to the desired cell, click to highlight it, and type in the desired window value along with a units code, if needed. If no units code is entered, the system will assume that the units of the previous entry are used. WINDOW values can be AutoFilled or graphically edited (see below).

- **RRT Code:** This column is used to designate peaks where the peak windows will be adjusted according to the difference between the expected time and the actual time. The windows for peaks that follow will be shifted by the error amount to compensate for any time skew found. Peaks are designated by placing an upper case R in this column for any peaks to be used for time window adjustments. This can be very helpful when the peaks are taking more or less time that usual to come off the column.
- **Suit Code:** See System Suitability (**section 9.1.3**) for details on producing System Suitability Reports. Codes are "S" or "K".
- **IS:** (Configurable) This column identifies up to five peaks in the table as Internal Standards, which can be applied to any specified peaks to correct for method and sample variations. Entries in the IS column are in the form of numbers, with a maximum range of 1 through 5. Only one peak component in each table can be assigned one particular IS number. The IS row is left blank for peaks which are not Internal Standards.

➡ To manually edit IS, click the pointer cursor inside the desired cell to highlight it, and type in the new value of IS from 1 to 5. An error message will appear if you attempt to designate an internal standard peak number which has already been assigned in the table.

The contents of the IS column will be automatically transferred into a Calibration screen Reference Peak Table when a new Calibration file is set up (Section 10.1). The identities of internal standard peaks can only be edited in the Integration Peak Table, which affects all existing Calibration Reference Tables which reference a given Peak Table.

- **CAL:** (Configurable) This column assigns which peaks/rows in the table can be calibrated, and on what basis. Choices for CAL values are blank (not calibrated), ES (External Standard basis), or IS1 - IS5 (Internal Standard basis with IS peak number). If Internal Standards have been disabled in Configuration, the IS options are excluded. Each peak for which quantitative results are to be determined MUST have an entry in its CAL column.

➡ To manually edit CAL, click the pointer cursor inside the desired cell to highlight it, and type in the new value of CAL. You can also use the FILL function below to fill in CAL entries which are repeated for a series of rows, or graphically select which peaks are to be calibrated (see below)

The contents of the CAL column will be automatically transferred into a Calibration screen Reference Peak Table when a new Calibration file is set up (Section 10.1). The calibration type for each peak can only be edited in the Integration Peak Table, which affects all existing Calibration Reference Tables which reference that Peak Table. Note that only those peaks which have an associated ES or IS code in the CAL column will be accessible to any Calibration file referencing this Peak Table. Therefore, you must determine which peaks are to be calibrated, and the method of calibration desired, before initializing any calibrations based on any peak table.

- **COMMENT:** (Configurable) The COMMENT column is used for entry of text comments describing any line in the Peak Table.

➡ To make a comment entry, click inside the COMMENT cell in the desired row and type in the entry.

9.2.3 Entering Peak Table Values

There are three ways to enter values for any parameter into the Peak Table:

1. Move the pointer into any "cell" in the table and click so that the highlight and the flashing text entry cursor are moved into that cell. Then, type any desired entry into the cell directly. To move to any other cell, either point to a different cell and click the left mouse button, or press the ENTER key to move the highlight from cell to cell, left to right, and then downward row by row until the desired cell is reached.

2. Copy information from an already entered row into a new row using the COPY and PASTE commands. To copy a row, move the pointer to the row to be copied, and click the left button to place the highlight in that row. Click on the COPY button to make a copy of the highlighted row.

Now, move the pointer to the row just below where you want to insert the copied row, and click anywhere in that row to place the highlight. Click the left mouse button on PASTE, and the copied row will be pasted in above the highlighted row. You can now make any other changes to the new rows by directly entering those changes into each cell.

3. Fill information into any column in the table using the FILL button below the table. The FILL function is a convenient way to transfer information from a chromatogram loaded in the Graph window to the table directly, without any typed or transposed entries, or to automatically compute and fill in table values. It works in a manner identical to that of the FILL function for the Method Events Table.

9.2.3.1 FILL Function

The FILL button operates in exactly the same manner as in the Method Events Table and all other Data Ally tables.




Figure 9.6 Automatically Filling Retention Times

The AutoFill dialog box (Figure 9.6) which defines all fill functions is the same dialog as used to fill columns in the Method, Calibration, and Sequence screens - see Section 8.3.3 above for a complete explanation of all terminology. After using the automatic fill function, you can make any desired changes in any of the filled retention times before saving the Integration program and its Peak Table.

In the Peak Table, the "Fill From Graph" function applies only to automatic filling of the RT column with current chromatogram retention times. In order to use this function, a valid integration must already have been completed on the foreground chromatogram in the active Graph window, with all desired peaks and/or peak groups identified by individual RTs. Automatic filling of retention times is an extremely fast and accurate means of completing even long peak tables. When filling from the graph, retentions will be expressed in the current time units shown in the graph window (min, sec).

The AutoFill "Duplicate" mode is available for the PEAKNAME, WINDOW, CAL, and COMMENT columns. Increment mode can be used with the RT, and WINDOW columns.

You can choose various combinations of options in the AutoFill dialog as described below to accomplish the following tasks (click on OK after setting up the Fill dialog as directed):

➤ To fill an existing row's value into all currently empty rows in the column up to the next entry, place the edit highlight on the "basis" row, click FILL to see the AutoFill dialog, select "Duplicate", "Current Cell", and "To Next Entry", with all other fields and options blank/off. If the "basis" row is the first row in the column, you can select "Start of Table" instead of "Current Cell" without highlighting the basis row.

➤ To fill an existing row's value into all currently empty rows in the entire column, place the edit highlight on the "basis" row, click FILL to see the AutoFill dialog, select "Duplicate", "Current Cell", and "End of Table", with all other fields and options blank/off.

➤ To fill an existing row's value into more or less empty rows than the number of empty rows found in the current column up to the next entry, place the edit highlight on the "basis" row, click FILL to see the AutoFill dialog, select "Duplicate", "Current Cell", and "To Next Entry". Type the actual number of empty rows to be filled into the "No. of Steps" column, with all other fields and options blank/off. If the

"No. of Steps" value exceeds the number of actual intervening empty rows, sufficient new rows will be added and filled to meet the "No. of Steps" criteria.

➤ To fill all currently empty series of consecutive rows in the column with the values found in the rows prior to each empty group of rows, click FILL to see the AutoFill dialog, select "Duplicate", "Start of Table", "End of Table", and "Allow Progressions", with all other fields and options blank/off.

➤ To overwrite existing values (along with any blank rows) with duplicated values in an entire column, place the edit highlight on the "basis" row for copying, click FILL to see the AutoFill dialog, select "Duplicate", "Current Cell", and "End of Table", and check the "Allow Overwrites" command, with all other fields and options blank/off.

➤ To clear (remove) existing values in an entire column, insert a blank row at the beginning of the table, place the edit highlight on that row, click FILL to see the AutoFill dialog, select "Duplicate", "Current Cell", and "End of Table", and check the "Allow Overwrites" command.

➤ To overwrite existing values in a portion of a column with a duplicated value, place the edit highlight on the "basis" row for copying, click FILL to see the AutoFill dialog, select "Duplicate", "Current Cell", "End of Table", and "Allow Overwrites", and enter a "No. of Steps" value corresponding to the number of existing rows to be filled with duplicates. All other fields and options should be blank/off. If the number of steps exceeds the number of blank rows in the column, the "Allow Insertions" box must be checked to add sufficient rows to meet the "No. of Steps" criteria.

➤ To clear (remove) existing values in a portion of a column, insert a blank row at the beginning of the range to be cleared (or remove an existing row entry), place the edit highlight on that row, click FILL to see the AutoFill dialog, select "Duplicate", "Current Cell", "End of Table", and "Allow Overwrites Data", and enter a "No. of Steps" value corresponding to the number of existing rows to be cleared (filled with a blank entry). If the number of steps exceeds the number of blank rows in the column, the "Allow Insertions" box must be checked to add sufficient rows to meet the "No. of Steps" criteria.

➤ To fill a column in a new table with incremented values, enter the beginning "basis" value in the first row, place the edit highlight on that row, click FILL to see the AutoFill dialog, select "Increment", "Current Cell", and "End of Table", enter the increment value in the "Increment" field, and enter the number of rows for which incremented values are to be inserted in the "No. of Steps" field.

➤ To fill a column with incremented values in blank rows up to the next existing entry, place the edit highlight on the "basis" row, click FILL to see the AutoFill dialog, select "Increment", "Current Cell", and "To Next Entry", and enter the increment value in the "Increment" field; leave other fields and options blank/off.

➤ To fill an existing row's incremented value into more or less empty rows than the number of empty rows found in the current column up to the next entry, place the edit highlight on the "basis" row, click FILL to see the AutoFill dialog, select "Increment", "Current Cell", and "To Next Entry". Type the increment value into the "Increment" field, and the actual number of empty rows to be filled/incremented into the "No. of Steps" column, with all other fields and options blank/off. If the "No. of Steps" value exceeds the number of actual intervening empty rows, sufficient new rows will be added and filled to meet the "No. of Steps" criteria.

➤ To fill all currently empty series of consecutive rows in the column with incremented values based on the rows prior to each empty group of rows, click FILL to see the AutoFill dialog, select "Increment", "Start of Table", "End of Table", and "Allow Progressions". Type the increment step value for each row into the "Increment" field, with all other fields and options blank/off.

➤ To fill a column with incremented values in all blank rows based on an existing entry, place the edit highlight on the "basis" row, click FILL to see the AutoFill dialog, select "Increment", "Current Cell",

and "End of Table", and enter the increment value in the "Increment" field; leave other fields and options blank/off.

☞ To fill a range of rows with evenly-spaced, automatically computed values, with each row's value corresponding to the difference between the basis cell's value and the value of the next row having a non-blank entry in the column divided by the number of intervening rows, place the edit highlight on the "basis" cell (first row in the range), click FILL to see the AutoFill dialog, select "Increment", "Current Cell", and "To Next Entry". If you wish to fill only the existing number of blank intervening rows, leave the "No. of Steps" field and all other commands blank/off. If you wish to fill either only some of the intervening rows, or to add new intervening rows, enter a number in the "No. of Steps" field specifying the total number of intervening rows to be filled. If "No. of Steps" exceeds the number of existing blank lines in the range, you must check the "Allow Insertions" box to automatically add and fill sufficient lines to match the "No. of Steps" entry.

☞ To fill all existing ranges of empty rows in the column with evenly-spaced, automatically computed values, with each filled row's value corresponding to the difference between the value of the cell immediately preceding that range and the value of the next row having a non-blank entry following the range divided by the number of intervening rows, click FILL to see the AutoFill dialog, select "Increment", "Start of Table", "End of Table", and "Allow Progressions". Leave all other selections blank/off.

☞ To overwrite existing values (along with any blank rows) with incremented values in an entire column, place the edit highlight on the "basis" row for incrementing, click FILL to see the AutoFill dialog, select "Increment", "Current Cell", and "End of Table", enter the increment value in the "Increment" field, and check the "Allow Overwrites Data" command, with all other fields and options blank/off.

☞ To overwrite existing values in a portion of a column with incremented values, place the edit highlight on the "basis" row for copying, click FILL to see the AutoFill dialog, select "Increment", "Current Cell", "End of Table", and "Allow Overwrites Data", enter the increment value in the "Increment" field, and enter a "No. of Steps" value corresponding to the number of existing rows to be filled with duplicates. All other fields and options should be blank/off. If the number of steps exceeds the number of blank rows in the column, the "Allow Row Insertions" box must be checked to add sufficient rows to meet the "No. of Steps" criteria.

☞ To fill Retention Times into the RT column from the current Graph window, click on the RT column header button (or move the edit highlight into that column), click FILL to see the AutoFill dialog (Figure 9.8), select "Fill From Graph", "Start of Table", and "End of Table", leaving all other fields and options blank/off. Lines will be added to the table with RT entries corresponding to the positions of all currently found peaks. If the "Allow Overwrites Data" box is checked, any existing RT column entries will be replaced by newly filled entries at the times assigned to such rows.

☞ To fill Retention Times into a timed region of the RT column from the current Graph window, click on the RT column header button (or move the edit highlight into that column), click FILL to see the AutoFill dialog (Figure 9.8), select "Fill From Graph", "Start of Table", and "End of Table", and enter the time range for filling RTs into the "First Time" and "Last Time" fields, leaving all other fields and options blank/off. Lines will be added to the table with RTs corresponding to the positions of all current peaks, within the specified time limits. If insufficient blank lines exist in the current table to fill all the RTs in the region given, the "Allow Insertions" box must be checked.

☞ To fill Retention Times into a range of rows in the RT column from the current Graph window, click on the RT column header button (or move the edit highlight into that column), click FILL to see the AutoFill dialog, select "Fill From Graph", "Start of Table", and "End of Table", leaving all other fields and options blank/off. Lines will be added to the table with RT entries corresponding to the positions of all current peaks. The number of filled rows can be limited by specifying a "No. of Steps" value less than the total number of peaks in the current graph window.

9.3.0 Working with Integration Files

The Integration program filing system operates in a manner identical to that of every other Data Ally programming screen. The Integration screen display, comprising both the top-level dialog box and the Peak Table, is configured as a "template" to show information about the current integration, peak finding, and reporting instructions loaded into the channel being observed. This template can be edited using many graphical and other tools to change its contents. The contents of the Integration screen display can be saved as a unique "Integration" file, with a unique filename and description, so that it can be recalled and reused at any future time either manually, or as part of an automated Method or a Sequence of Methods.

The default Integration file, whose name is assigned during Configuration, is always loaded automatically into the Integration display screen at the time Data Ally is initialized. You can change any aspects of the default file, including its Peak Table, within the active software, except for the features and options which are either "fixed" or "removed" in the Configuration. These items can, of course, also be reset or restored, but the channel must be disabled and a re-Configuration performed to accomplish this.

The Files commands in the Integration Command Menu bar (Figure 9.7) are used to load, save, clear, and delete Integration program files.




Figure 9.7 Integration Files Command Submenu

9.3.1 Loading Integration Files

➡ To load a saved Integration file from disk memory into the current channel, click the left mouse button on Files in the Command menu bar, and then click on Load. A directory dialog box will appear which lists the Integration files (*.INT) available in the current directory path for this channel (Figure 9.8). Select the desired file by clicking on its name in the directory list, or by typing its name into the entry box. Click on OK to load the new file.




Figure 9.8 Integration Files/Load Dialog

If the previous file has been edited and has not been saved or resaved, an information message will appear before the new file is loaded asking if the previous file should be saved. If so, click on YES and assign either the same filename or a new name as desired. If you do not wish to save the edited file, click on NO to proceed.

When a new file is loaded, it will overwrite the previous Integration file.

To reload the default Integration file programmed during Configuration, use the Load Default command in the Files submenu. This default file can have a complete Peak Table for a "default" chromatogram, if desired - loading it will overwrite any other Integration file currently present in the time base. The default file can be modified at any time for manual operation, but can only be permanently changed in Configuration Mode (Section 5.2.5.3).

9.3.2 Clearing Integration Peak Table

Rather than deleting all lines in an existing peak table when creating a new table, a Clear function is provided to remove all existing values from the current table.

➡ To clear the current Integration Peak Table, click on File in the Command Menu bar and then click on Clear. Then click on Clear Peak Table (Figure 9.9). The Peak Table will be emptied. If the previous file has been edited and has not been saved or resaved, an information message will appear before this file is cleared asking if the previous file should be saved. If so, click on YES and assign either the same filename or a new name as desired. If you do not wish to save the edited file, click on NO to proceed.




Figure 9.9 Clearing Peak Table

9.3.3 Saving Integration Files

➡ *To save an Integration file using the same filename*, click on Files in the Command Menu bar, and then click on Save. A trap message box will appear prompting you to confirm that you wish to overwrite the old file with the same filename. Click on Yes, and the previous version of the file will be overwritten and a new version recorded in the same path.

➡ *To save a modified Integration file with a different filename*, click on Files and on the Save As command in the pulldown. The Save As dialog box will appear (Figure 9.10). Click inside the filename entry field and type in the desired new filename. If you wish, click inside the DESCRIPTION field and modify any existing description of the file, which may help identify this file in the directory listing for future recall. When the correct name and description appear, click on OK to save the file to the current disk path. If another file already exists in that path with the same name, an error message will instruct you to change the filename before saving.




Figure 9.10 Integration Files/Save As Dialog

It is recommended that Integration program files saved to disk be "backed up" by simply copying the *.INT files to a floppy disk, removable hard disk, or tape and preserving the copies in a secure place. If modifications are routinely made to the library of *.INT files, backups should occur at routine intervals or even on a daily basis.

You can use the Files/Integration Description command to directly enter a single or multiline description of the current Integration file. Alternatively, you can enter such a description using the Save As dialog (Figure 9.10) when saving a new or modified Integration file. Any descriptive text entered using the Integration Description command will be shown in the "Description" field when reviewing saved files via the Load dialog (Figure 9.8).

9.3.4 Deleting Integration Files

➡ *To delete an Integration file from the current disk directory*, use the Windows File Manager Delete command.

9.3.5 Printing The Integration File

Use the Print Integration command in the Files submenu to print a summary of the current Integration program file. The template for this printout is saved as \reports\int\int.xls and can be modified as desired using Microsoft Excel (Section 13).

10.0 Reviewing and Editing Calibrations: Calibration Screen

Data Ally Chromatography Manager provides many features for quantitating results to generate calibrations and reports for unknown samples. Any chromatogram acquired or reprocessed in any available mode can be utilized for creating a new calibration or for updating an existing one.

10.1.0 Calibration File Definition

A Data Ally Calibration file is a set of data corresponding to peak components from a specified Integration file's Peak Table, which represents the cumulative responses of the standards already analyzed for determining the concentration of each component in unknown materials. Each time a new standard is injected and its chromatogram is captured, the peak responses in that chromatogram can be applied to any Calibration file(s) associated with the Integration program(s) being executed. Each Calibration file not only holds the actual peak response information, but also contains other items which describe each calibration standard run as predefined in the Configuration for each channel. When the Calibration function is configured, the system manager selects which types of data are to be displayed in the Calibration table, and which options are to be made available for analyzing and interpreting the calibration data.

Calibration files can consist of results from only a single standard chromatogram, and can be as large or contain information from as many different standard runs as desired. Calibration files are updated via instructions from a running Method or Sequence, or from use of the manual Edit/Calibrate command in any Graph window. Only one standard chromatogram can be applied to a Calibration file at one time for updating.

When a calibration is performed manually or as part of a Method or Sequence using a defined Calibration file, one of three options is always selected:

- The standard chromatogram being analyzed is used to establish a new calibration in the file, clearing any previous responses at all levels in that file (the New mode);
- The standard chromatogram is used to clear and update any existing responses in an already-existing calibration in the file for the single calibration level represented by that standard (the Replace Level mode) - all other responses for all other levels in the Calibration file are retained.
- The standard chromatogram is used to update the existing responses in an already-existing calibration in the file by averaging its response at the appropriate level with all other existing responses (the Average mode). This mode is essentially a "moving average" computation since each new standard response is averaged against the previous average response for its unique level.

Regardless of which of the three modes above is used to update a calibration table, every single standard run used for such updates has equal "weighting" of its data compared to all other standard runs, with the overall calibration curve being recalculated at each update from all the available standard results.

Each new calibration file must have an assigned, valid file name, which will end with the automatically-assigned name extension ".CAL". If the currently-loaded Method contains a Calibration filename reference, that calibration file will be automatically loaded when that Method is executed. If the Method does not contain a Calibration filename reference, the Calibration file present in the current channel will be used. If no specific Calibration file has been loaded, the default Calibration filename, if any, assigned in Configuration will be loaded automatically.

Since Calibrations are performed by instructions from a Method, Sequence, or Graph Edit/Calibrate command, actual calibration is always done "behind the scenes". The Calibration file(s) used for real-time or other calibrations need not be loaded into the current Calibration screen display. The

purpose of the Calibration screen, which is viewed by clicking the CALIBRATE SoftButton on the left screen, is to provide a convenient and accessible means of reviewing, editing, and manipulating any desired Calibration file. It is possible to be editing one Calibration file in the Calibration screen of the current channel while that channel is processing standards with one or more different Calibration files in the background.

In order to perform calibrations, either the "default" calibration filename or another user-defined calibration filename must be referenced by the Method or in the manual Calibration command. At any time, only those calibration files actually established and named (regardless of whether they actually contain any calibration data) are accessible for accumulating real calibration data. Data Ally is capable of automatically "resaving" calibration files as updates are made so that the cumulative data is not lost; alternatively, a user can opt for manual-only calibration file saving.

10.2.0 The Calibration Display

➤ To see the Calibration screen (Figure 10.1), click on the CALIBRATE SoftButton. This screen consists of three sections, each of which is a resizable and movable window:

- The Reference Peak Table, at upper right, which captures the table of all calibrated peaks from the Peak Table of the referenced Integration screen, and sets the maximum number of standards to be used and the concentrations of those standards.
- The Calibration Table, at bottom, which lists data describing every standard chromatogram which has been used to construct this calibration file.
- The Calibration Plot window, at upper left, which shows a graphical representation of the calibration data for each individual peak.

All three window sections are interactive, so that changing the presentation in one window will often change the contents or display in the other windows.




Figure 10.1 Calibration Screen

10.2.1 The Reference Peak Table

Every Calibration file, when it is defined, must establish a "connection" with the Peak Table in one Integration program, which is used to identify peaks and set the basic calibration instructions for the

type of chromatogram to be run with a particular Method. Once a reference is entered, each peak in that Peak Table marked for calibration in the Table's CAL column is transferred into the Reference Peak Table window at the upper right of the Calibration screen (Figure 10.2). The referenced Integration file must either be saved to disk in some accessible directory path, or it must be loaded into the current channel, in order to be used properly by the Calibration file.

	RT	Peak Name	Basis	Fit	Unit	Lev1	Lev2	
1	0.546	Peak One	A	L	mg	40.0	80.0	10
2	1.477	Peak Two	A	L	mg	40.0	80.0	10
3	2.164	Peak Three	A	L	mg	40.0	80.0	10
4	2.65	Peak Four	A	L	mg	40.0	80.0	10
5	3.277	Peak Five	A	L	mg	40.0	80.0	10
6	3.692	Peak Six	A	L	mg	40.0	80.0	10
7	4.07	Peak Seven	A	L	mg	40.0	80.0	10
8	4.72	Peak Eight	A	L	mg	40.0	80.0	10
9	5.59	PeakNine	A	L	mg	40.0	80.0	10
10	6.28	Peak Ten	A	L	mg	40.0	80.0	10
11	7.495	PeakEleven	A	L	mg	40.0	80.0	10
12	8.45	Peak 12	A	L	mg	40.0	80.0	10
13	8.83	Peak 13	A	L	mg	40.0	80.0	10
14	9.96	FirstIS	A	L	mg	200.0	80.0	10
15	16.78	Peak21	A	L	mg	40.0	80.0	10
16	12.96	SecondIS	A	L	mg	100.0	80.0	10

Figure 10.2 Reference Peak Table Window

The Reference Peak Table actually specifies the calibration basis, curve fitting for multiple level calibrations, units, number of possible standard concentration levels, and concentration target values for standards at each of those levels. The type of standardization (external or internal standard, internal standard references for each peak) is actually defined in the Peak Table of the referenced Integration screen. This means that a single Integration screen/Peak Table can have multiple associated Calibration files, each with its own specific calibration instructions and responses. In order for a successful calibration to be completed, both the Integration file and Peak Table appropriate for the current standard chromatogram (which identifies all desired peaks correctly and specifies which peaks are standardized and in what modes), and one Calibration file which references that Integration Peak Table must be loaded in memory. When you load a new Calibration file, you will be automatically prompted to load the actual referenced Integration file, if this has not already been done. Alternatively, when you perform a manual or Method-based calibration using a loaded Integration file, Data Ally will automatically make sure you have in fact also loaded an adequate Calibration file which references your Integration file.

A new Calibration file can be created by modifying an existing Calibration file, or the default Calibration file, and then resaving it with a different name using the Save As command (Section 10. below). In order to utilize a Calibration file, it must be saved to disk with a unique filename, and must include a valid Reference Peak table with at least one calibration level defined for at least one peak. If a Calibration file specified when running a Method or Sequence is not valid, or does not contain the proper information, the Method or Sequence will not be permitted to start and an error message will appear advising that an invalid Calibration file was found.

The Reference Peak Table, like other Data Ally tables, is a "spreadsheet" consisting of a series of rows and columns, which can be edited by placing the Edit cell highlight in any row corresponding to a peak, or by using the Fill dialog in the same manner as for other tables. Columns indicating Peak Name, RT, and Calibration Type are automatically filled in with information from the Integration Peak

Table being referenced, to provide the user with a complete description of each peak. Additional columns are provided for FIT, BASIS, UNITS, and LEV for levels from 1 to 16.

You can easily re-order the columns in the Reference Table from left to right or change their individual widths at any time. Both the column order and widths are set as defaults in Configuration. To move a column to another position in the table, simply place the pointer cursor on that column's header button, hold down the left mouse button, and "drag" the header button left or right to the new desired position relative to the other columns. When the pointer is in the desired position, release the mouse button and the entire column will be moved, and the table redrawn accordingly. Modifying the order of columns may make certain types of table editing or review easier by allowing placement of all columns currently being manipulated in view on a single screen - you can easily move "unimportant" column information to the right and off of the displayed table window area.

To change the width of any individual column, move the pointer cursor to the border of its header button and the next column header button in the direction in which you wish to expand or reduce the column width. The cursor will become a "double-sided horizontal arrow" (↔); when this change occurs, immediately hold down the left mouse button and move the arrow cursor in the direction you wish the button edge to expand or contract. When the desired column width is indicated by the position of the cursor, release the mouse button and the column will be resized and the entire table redrawn accordingly. This feature permits you to match the exact sizing of each column to the expected character width of the entries expected in that column - you can also effectively "remove" certain columns without actually eliminating them from the table by reducing their widths to a very small size.

An ENABLE column which numbers lines in the Reference Table is provided as for the other Data Ally tables. This column can be used to enable or disable individual peaks from the Peak Table in order to select whether they are calibrated or not. Normally, each enabled or active row corresponding to one peak will have a colored highlight in this first left hand column. To disable a peak, double-click inside the enable column in its row to remove the highlight. To re-enable it, double click again to restore the highlight.

The CUT/COPY/PASTE command buttons found at the bottom of the Calibration screen can be used with both the Reference Peak Table and, to a limited extent, with the Calibration Data table below. The active table window in this screen will access these buttons. In the Reference Table, because the row contents of this table are fixed by its relationship to a predefined Integration Peak Table, the RT and Peakname columns cannot be edited directly - to change RT or Peakname, simply change the source Integration screen Peak Table, which will automatically modify these parameters in all Calibration Reference Tables based on that Peak Table. However, all other columns in the Reference Table can be cut or copied, pasted, and filled; the FILL button is provided for automatic filling of columns which define Basis, Fit Type, Units, and Levels via the standard AutoFill dialog box.

Every calibrated peak in the Reference Peak Table has its own unique Calibration Table and its own Calibration plot. The Reference Peak table is utilized to select any particular peak, in order to view the Calibration Table and Plot specific to that peak. Both the Table and the Plot will immediately be updated for the newly selected peak. This interactive display mode helps you easily review the current status of calibration of any peak, and modify the calibration quickly if necessary. It may be most convenient to maximize the Peak Reference Table when editing it, by clicking on the Program Control/Maximize button at the window's upper right, especially if a large number of levels is being input. You can use the table's vertical and horizontal scroll bars to move around in the table, as well as the PageDown and PageUp keys (to move up or down a fixed number of rows of peaks), the CTRL-Home or CTRL-End commands (to move to the first or last column on the current row in the table), and the CTRL-PgDn or CTRL-PgUp commands (to move to the beginning or end row/peak in the table, respectively).

The following series of steps must be performed in order to modify an existing or default Calibration Reference Peak Table in order to initialize a new one.

10.2.1.1 Setting Up a New Reference Table

➡ To select the Integration screen which includes the Peak Table for reference for a new calibration file, click on Files in the Calibration Screen Command Menu Bar, to see the Files submenu (Figure 10.3).




Figure 10.3 Calibration Files Submenu

Now, click on Setup to see the Calibration Setup dialog box (Figure 10.4).




Figure 10.4 Calibration Setup Dialog

The Setup dialog defines the Integration File/Peak Table reference for the current Calibration file, which must be an existing Integration file.

➡ To select the referenced Integration file name, click inside the Integration entry field and type in the desired Integration path/filename as reference for constructing the required Peak Table, or click the down-arrow button to the right of the field to see a listing of all Integration files in the current path. To select one of these files as the reference source file, click its name to enter it.

Note that changing the Integration source file name using Setup will cause all data in the current Calibration file to be deleted as the Calibration Reference Table is reconstructed.

➡ To set the total number (maximum number) of different standard concentrations for this calibration, click on the right arrow next to the No. of Standard Levels field, and choosing the desired number of levels (1 to 16). This will cause the selected number of columns to be placed into the current Calibration Reference Table to accommodate all levels. If the Integration source filename is not changed but the number of levels is modified using Setup, existing data in the current Calibration table will not be lost, and the Reference Table will show added levels.

When both the Reference filename and number of levels are shown as desired, click on OK to restructure the Calibration Reference Table. Immediately, the listing of peaks by name and RT in the referenced Integration file will appear in the Calibration Reference window, replacing any previous set of peaks. Also, the Reference Table will immediately update so that it contains the number of level columns specified via Setup. **NOTE THAT EACH TIME SETUP IS USED TO CHANGE THE INTEGRATION FILE REFERENCE, THE EXISTING CALIBRATION TABLE IN THE CURRENT CALIBRATION SCREEN IS CLEARED, AS WELL AS ALL PREVIOUS ENTRIES IN THE REFERENCE TABLE** (since presumably a different set of peaks is being defined). Therefore, you should be certain the previous Calibration file has been saved before using Setup, if it is to be preserved for later use. If you use Setup to change only the number of levels, the existing calibration Reference Table and Calibration Table data will be retained, with the loss or addition of any levels specified. If levels are removed, their data will be lost from the associated Calibration Table.

It is possible to modify the Reference Table of an existing Calibration file without using Setup, if the same set of peaks is to be used for the new table. To change the peaks themselves, the types of calibration (ES, ISx, or none), and the identity or number of Internal Standard peaks (as shown in the IS column), simply load the referenced Integration file, make the desired changes in its Peak Table, and resave the modified Integration with the same filename - the changes made will automatically be echoed into all Calibration files with Reference Tables using that Integration file. To change information other than the peak definitions, you can directly edit the Reference Peak Table at any time.

➡ To use Setup again to change the number of levels in any existing Reference Table, without clearing the entire contents of the Calibration file, bring up the Setup dialog and change the No. of Standard Levels entry. To do this, you must not change the Integration file reference before clicking OK. If the new number of standard levels is greater than the old number, the existing values in all the LEV columns will be preserved, and one or more new blank LEV columns will be added to the Reference Table. If the new number of standard levels is less than the old number, pre-existing LEV columns will be deleted from the old Reference Table along with their contents, beginning with the last LEV column present.

10.2.1.2 Setting Curve Fit Types for Each Peak

The Reference Table defines the manner in which multiple standard levels will be mathematically fit to compute the overall "standard response curve" for each peak. To do this, each peak or row in the Reference Table must have a code entry which determines the type of curve fitting to be done for multilevel calibrations. The FIT column is provided for this function if the number of levels is set to any value greater than one - the default fit type is Linear Regression. Each peak may have its own unique curve fitting type/code, regardless of that used for any other peaks. The Data Ally calibration system treats each peak totally independently of all others.

The following codes are available for programming various types of curve fitting:

L	=	Linear Regression
L0	=	Linear Regression, force through zero
I	=	Interpolated Point to Point
Q	=	Quadratic/Polynomial
Q0	=	Quadratic/Polynomial force through zero

Note that if only one calibration level has been selected, none of these options will apply and only a simple single-level response factor calculation will be made - this is essentially a Linear fit through zero. If only two calibration levels are allowed, only Linear Regression fitting is permitted. If three or more levels are specified, all the fit options are allowed.

FIT codes must be provided for each row in the FIT column. If all peaks will use the same fit, that code needs to be entered only once in the FIT column, in its top row, and the AutoFill function used to fill all subsequent peak rows with the same command.

➤ *To define the curve fit for any peak*, click inside the FIT column cell in that peak's row and either type the code for the desired fit type. You can also use the FILL command to AutoFill Duplicate values of any entered fit code into other rows. To AutoFill the FIT column, place the Edit cell highlight on any desired row, and click on the FILL button to display the AutoFill dialog box. Select Duplicate, select the range desired, and click OK to fill the column.




Figure 10.5 AutoFilling Curve Fit Type

You can change the FIT instructions for any peak(s) at any time after a Reference Table is set up and calibration data have been logged for those peaks. Simply click on the FIT cell for the peak to be modified, and type in the new FIT command code. Press ENTER, CTRL-arrow, or click inside any other cell in the table to make the entry see the results in the Plot window or on Cancel to abort the change. If Cancel is selected, the previous FIT command entry will be automatically restored.

Note that in Configuration each time base can be defaulted to only a limited number of curve fit options or to only a single option. If one fit option only is permitted, usually the FIT column will not be selected for display when configuring.

10.2.1.3 Defining Calibration Basis for Each Peak

Each peak in the Reference Table can be calibrated on either an Area basis (the default) or Height basis. The basis can differ for any individual peak according to the code (A or H) placed in each peak's BASIS column.

➤ *To set the basis for any peak*, click inside the BASIS column cell in that peak's row, and type in the desired basis code (A or H).

If you wish to AutoFill the BASIS column values by Duplicating an existing entry (or clearing a range of rows), click on the BASIS column header button to bring up the FILL dialog box (Figure 10.5). Set the mode and range and click on OK to fill.

Note that the BASIS column can be omitted from the Reference Table when configuring the Calibration screen if the calibration basis is to permanently be locked for either area or height for all peaks. The system defaults to an A (Area) BASIS at the time a new Reference Table is set up.

You can change the BASIS for any peak(s) at any time after a Reference Table is set up and calibration data have been logged for those peaks. Simply click on the BASIS cell for the peak to be modified, and type in the new BASIS command code. Press ENTER, CTRL-arrow, or click inside any other cell in the table to make the entry - a dialog box (Figure 10.6) will appear prompting that you are changing the current calibration and asking for approval. Click OK to make the change for the peak and see the results in the Plot window or on Cancel to abort the change. If Cancel is selected, the previous BASIS entry will be automatically restored.




Figure 10.6 Changing Basis for Existing Peak in Reference Table

10.2.1.4 Setting Concentration/Reporting Units for Each Peak

The UNITS column for each peak defines the label to be applied to all reported results. In Data Ally, units are completely arbitrary for any individual peak. For example, one peak can be reported as mg/μl, with another peak reported as g/ml. It is possible to enter a units label along with any numeric concentration value in any LEVx column in the Reference Table. If such units terms are provided, Data Ally will automatically compensate for differences among units input in different LEV columns for any single peak. However, if an entry is made in the UNITS column, any individual units labels inserted in any LEV column must be interpretable into the type of units specified in the UNITS column. For example, if the UNITS column entry for a peak is "mg/dl", you can make entries of "μg/ml" or "g/l" into any LEVs columns for that peak - however, you cannot make an entry of "L" (liters) or "G" (grams) in any of the LEV columns, since Data Ally has no way to determine the relationship between these entries and the UNITS of mg/dl.

Generally, you will want to use the units designator for each peak component in which the numerical standard concentration values will be entered in the LEVx columns. UNITS entries will be used in every on-screen or printed report, and as part of any saved report files. Note that UNITS values need not be entered in each row in the UNITS column - any entered value will also apply to blank rows in the UNITS column below that value, until the next value in a subsequent row is encountered. For example, if all peaks will be reported with the same units designator, that designator needs to be entered only once in the UNITS column, in its top row corresponding to the first peak.

➤ To set quantitative reporting units for any peak, click inside the peak's UNITS column cell, and type in the desired units. Note that there is no right-button selector available for units entry, since any alphanumeric units designator can be used. You can also use AutoFill to Duplicate or Clear entries in the UNITS column by clicking on the UNITS column header button, and setting the fill mode and range.

There is no UNITS default in the Peak Reference Table - if no entries are made for any peak, no units designator will be printed along with any numerical results. The presence or absence of a units designator will not affect whether or not a peak is actually calibrated and quantitated.

10.2.1.5 Entering Concentration Levels for Each Peak

At the right side of the Reference Peak Table, one or more LEV columns are provided, depending upon the entry made for No. of Standard Levels in Setup. At least one LEV column must always be present in order to define a valid Calibration file. The LEV columns represent up to sixteen different standard samples, with unique or replicate concentration values, to be used to calibrate each peak. It is assumed that when a level is defined in the Reference Table, that level refers to an actual preparation or vial of standard material which will actually include the stated concentrations for each peak component to be calibrated.

LEV values can be entered in any manner desired. At least one LEV concentration value must be entered for at least one peak in a valid Calibration file. Not every LEV column has to contain values for every peak. For any given peak, LEV values consecutive columns can be ascending, descending, mixed, random, or the same. LEV columns can be skipped as desired for any peak. As an example, you can enter LEV1=50, LEV2=100, LEV3 and LEV4 blank, and LEV5=20 for peak A, and LEV3=400, LEV4=400, LEV5 through LEV8 blank, and LEV9=400 for peak B. It is important to make absolutely sure that the concentration values entered for any peak/level actually represent the composition of the standard material to be injected which corresponds to that level.

➤ To enter a LEV concentration value for any peak, click inside the desired peak's LEVX column cell, and type in the numeric concentration value for the standard assigned to that level. You can use AutoFilling to Duplicate or Increment LEV column entries as desired, by clicking on the LEVX column header button and completing the FILL dialog box instructions.

LEV column entries are numerical only and should not have units designators such as g/dl entered as well in addition to a numerical value. The UNITS column should be employed for this purpose.

If you reset No. of Standard Levels to reduce the number of level columns present in an existing Reference Peak Table, by repeating the Setup command, existing LEV column entries in the to-be-deleted columns will be removed and lost.

It is usually advisable to program the Reference Peak Table with sufficient LEV columns to represent the maximum number of different standards you may theoretically wish to use for calibration purposes. Any subset of levels for any combination of peaks can be programmed in the Reference Table at any time.

Once you have assigned all the desired level information for your standards for each peak, you must save your calibration file to prepare it for use in executing calibrations. At least one level concentration value must be defined for at least one peak in order for calibration responses and results to be generated. Note that as many different Calibration files as desired can exist at any time referencing the same Integration Peak Table. However, only one Calibration file can be used in conjunction with any single Integration file at any time. A description can be entered for a Calibration file at the time the Save As command is used (see Section 10.5.3 below).

➤ To save the Calibration file with the newly-edited Reference Peak Table, use the Files/Save As command (Section 10.5.3 below). Once a Calibration file is saved, it can be called automatically in any Method or Sequence for updating or use in results quantitation.

10.2.1.6 Using FILL Functions in the Reference Table

The standard AutoFill functions are accessible for completing entries in the BASIS, FIT, UNITS, and LEVx columns in the Calibration Reference Table. Use the FILL button at bottom right of the Reference Table window to bring up the AutoFill dialog after selecting the desired column to be filled by clicking on its header button. The AutoFill dialog (Figure 10.5) operates in exactly the same manner as that provided in the Events, Peak, and Sequence tables. A complete explanation of the features of the AutoFill dialog is found in Section 8.3.

In the Calibration Reference Table, the "Duplicate" mode is active for the UNITS, BASIS, FIT, and LEVx columns. "Increment" mode operates on the LEVx columns only, and there is no "Fill From Graph" mode.

The following actions can be performed using the AutoFill dialog:

➤ To fill an existing row's value into all currently empty rows in the column up to the next entry, place the edit highlight on the "basis" row, click FILL to see the AutoFill dialog, select "Duplicate", "Current Cell", and "To Next Entry", with all other fields and options blank/off. If the "basis" row is the first row in the column, you can select "Start of Table" instead of "Current Cell" without highlighting the basis row.

➤ To fill an existing row's value into all currently empty rows in the entire column, place the edit highlight on the "basis" row, click FILL to see the AutoFill dialog, select "Duplicate", "Current Cell", and "End of Table", with all other fields and options blank/off.

➤ To fill an existing row's value into more or less empty rows than the number of empty rows found in the current column up to the next entry, place the edit highlight on the "basis" row, click FILL to see the AutoFill dialog, select "Duplicate", "Current Cell", and "To Next Entry". Type the actual number of empty rows to be filled into the "No. of Steps" column, with all other fields and options blank/off. If the "No. of Steps" value exceeds the number of actual intervening empty rows, sufficient new rows will be added and filled to meet the "No. of Steps" criteria.

➤ To fill all currently empty series of consecutive rows in the column with the values found in the rows prior to each empty group of rows, click FILL to see the AutoFill dialog, select "Duplicate", "Start of Table", "End of Table", and "Allow Progressions", with all other fields and options blank/off.

➤ To overwrite existing values (along with any blank rows) with duplicated values in an entire column, place the edit highlight on the "basis" row for copying, click FILL to see the AutoFill dialog, select "Duplicate", "Current Cell", and "End of Table", and check the "Allow Overwrites" command, with all other fields and options blank/off.

➤ To clear (remove) existing values in an entire column, insert a blank row at the beginning of the table, place the edit highlight on that row, click FILL to see the AutoFill dialog, select "Duplicate", "Current Cell", and "End of Table", and check the "Allow Overwrites" command.

➤ To overwrite existing values in a portion of a column with a duplicated value, place the edit highlight on the "basis" row for copying, click FILL to see the AutoFill dialog, select "Duplicate", "Current Cell", "End of Table", and "Allow Overwrites", and enter a "No. of Steps" value corresponding to the number of existing rows to be filled with duplicates. All other fields and options should be blank/off. If the number of steps exceeds the number of blank rows in the column, the "Allow Insertions" box must be checked to add sufficient rows to meet the "No. of Steps" criteria.

➤ To clear (remove) existing values in a portion of a column, insert a blank row at the beginning of the range to be cleared (or remove an existing row entry), place the edit highlight on that row, click FILL to see the AutoFill dialog, select "Duplicate", "Current Cell", "End of Table", and "Allow Overwrites", and enter a "No. of Steps" value corresponding to the number of existing rows to be cleared (filled with a blank entry). If the number of steps exceeds the number of blank rows in the column, the "Allow Insertions" box must be checked to add sufficient rows to meet the "No. of Steps" criteria.

➤ To fill a column in a new table with incremented values, enter the beginning "basis" value in the first row, place the edit highlight on that row, click FILL to see the AutoFill dialog, select "Increment", "Current Cell", and "End of Table", enter the increment value in the "Increment" field, and enter the number of rows for which incremented values are to be inserted in the "No. of Steps" field.

➤ To fill a column with incremented values in blank rows up to the next existing entry, place the edit highlight on the "basis" row, click FILL to see the AutoFill dialog, select "Increment", "Current Cell", and "To Next Entry", and enter the increment value in the "Increment" field; leave other fields and options blank/off.

➤ To fill an existing row's incremented value into more or less empty rows than the number of empty rows found in the current column up to the next entry, place the edit highlight on the "basis" row, click FILL to see the AutoFill dialog, select "Increment", "Current Cell", and "To Next Entry". Type the increment value into the "Increment" field, and the actual number of empty rows to be filled/incremented into the "No. of Steps" column, with all other fields and options blank/off. If the "No. of Steps" value exceeds the number of actual intervening empty rows, sufficient new rows will be added and filled to meet the "No. of Steps" criteria.

➤ To fill all currently empty series of consecutive rows in the column with incremented values based on the rows prior to each empty group of rows, click FILL to see the AutoFill dialog, select "Increment", "Start of Table", "End of Table", and "Allow Progressions". Type the increment step value for each row into the "Increment" field, with all other fields and options blank/off.

➤ To fill a column with incremented values in all blank rows based on an existing entry, place the edit highlight on the "basis" row, click FILL to see the AutoFill dialog, select "Increment", "Current Cell", and "End of Table", and enter the increment value in the "Increment" field; leave other fields and options blank/off.

➤ To fill a range of rows with evenly-spaced, automatically computed values, with each row's value corresponding to the difference between the basis cell's value and the value of the next row having a non-blank entry in the column divided by the number of intervening rows, place the edit highlight on the "basis" cell (first row in the range), click FILL to see the AutoFill dialog, select "Increment", "Current Cell", and "To Next Entry". If you wish to fill only the existing number of blank intervening rows, leave the "No. of Steps" field and all other commands blank/off. If you wish to fill either only some of the intervening rows, or to add new intervening rows, enter a number in the "No. of Steps" field specifying the total number of intervening rows to be filled. If "No. of Steps" exceeds the number of existing blank lines in the range, you must check the "Allow Insertions" box to automatically add and fill sufficient lines to match the "No. of Steps" entry.

➤ To fill all existing ranges of empty rows in the column with evenly-spaced, automatically computed values, with each filled row's value corresponding to the difference between the value of the cell immediately preceding that range and the value of the next row having a non-blank entry following the range divided by the number of intervening rows, click FILL to see the AutoFill dialog, select "Increment", "Start of Table", and "End of Table". Leave all other selections blank/off.

➤ To overwrite existing values (along with any blank rows) with incremented values in an entire column, place the edit highlight on the "basis" row for incrementing, click FILL to see the AutoFill dialog, select "Increment", "Current Cell", and "End of Table", enter the increment value in the

"Increment" field, and check the "Allow Overwrites" command, with all other fields and options blank/off.

➤ *To overwrite existing values in a portion of a column with incremented values*, place the edit highlight on the "basis" row for copying, click FILL to see the AutoFill dialog, select "Increment", "Current Cell", "End of Table", and "Allow Overwrites", enter the increment value in the "Increment" field, and enter a "No. of Steps" value corresponding to the number of existing rows to be filled with duplicates. All other fields and options should be blank/off. If the number of steps exceeds the number of blank rows in the column, the "Allow Insertions" box must be checked to add sufficient rows to meet the "No. of Steps" criteria.

10.2.2 The Calibration Point Table

The Calibration Point Table, like the Method Events and Peak Tables, is a "spreadsheet" with columns and rows (Figure 10.7). It is intended as a simple means of displaying all the desired information about each standard calibration run which is currently accumulated in the Calibration file. The columns listed in the Table are configurable, so that only the types of information needed can be included in the Table, in the order desired. When you view the calibration screen, you are seeing the Point Table (and Plot) corresponding to the currently-selected peak in the Peak Reference table; highlighting a different peak in the reference table will cause that peak's Calibration Point Table and Plot to be displayed instantly. The name of the currently-selected peak will always appear on the top bars of both the Calibration Point Table and Plot windows.

Because you can select which columns/types of information will appear in this table, and change the column order and individual widths, the Calibration Point Table is an extremely versatile way to display data regarding standard chromatograms. The table is also fully interactive with the Plot window (below) to assist in reviewing the nature of this data. Of course, the Calibration Point Table window can be moved or resized as desired at any time.

➤ *To view the Calibration Point Table and Plot for any peak in the Reference Table*, move the Reference Table highlight to the desired peak's row by clicking anywhere in that row, or by using the scroll arrows. The Point Table and Plot will immediately change to show the current data for the new peak - the top bars of both these windows will show the name assigned to the peak.

	Level	Conc.	Dil	Int File	Cg Filename	Area
HPLC1	1	4	10.0	c:\pyramid\caltest.int	c:\pyramid\pdfcal10is.pdf	1524.79937
Timebase	2	2	80.0	c:\pyramid\caltest.int	c:\calfiles\cal80is.pdf	12198.2331
HPLC1	3	3	100.0	c:\pyramid\caltest.int	c:\calfiles\cal100is.pdf	15247.7914
Edit	4	2	80.0	4.0	c:\pyramid\caltest.int	3049.55829
TB2						
Timebase						
TB2						
Edit						

Figure 10.6 Calibration Table

Each time a new standard is run and its responses applied to a component peak at any defined level, a new line is added to the Calibration Point Table describing that standard. Each line includes information on the Level, Concentration, Area Response, Height Response, Response Factor, and Chromatogram filename for that standard, as well as dilution or weighting factors (see Section 5.2.5.4 for a complete list of configurable column items for the Calibration Point Table). Lines describing standard runs are listed in the order which corresponds to the position of columns in the table. Entries are sorted based first on the contents of the leftmost column in each row (in ascending order), then on the contents of the second column, third column, etc.

Like all other Data Ally tables, the Calibration Point Table features a dual-function line numbering ENABLE column at far left, which is always visible regardless of the current view of the table

selected. Each time a new line is added to the table it automatically receives the next available line number. Normally, all new lines have the ENABLE column highlighted, the data in each line being enabled and used to compute the calibration response or curve fit. Double-clicking on the ENABLE column cell for any row in the table causes the highlight to disappear and that row to become disabled, thus removing it from the computation and effectively "failing" that standard from the calibration of that peak. This is an extremely convenient tool for testing the effect of removing outlier runs from certain calibrations.

You can use the scroll bars at the bottom and right of the table to move through all of the rows and columns, as well as the <PgUp>/<PgDn> keys (to move up and down a fixed number of rows), the CTRL-Home and CTRL-End commands (to move to the leftmost and rightmost columns), the <TAB> or <ENTER> keys (to move one cell to the right), and the CTRL-PgUp and CTRL-PgDn commands (to move to the first row or last row). At any time, one line in the table will be highlighted. The highlight can be moved to a different line by clicking the left mouse button on any desired column in the row to be highlighted. As in the other Data Ally tables, the highlight represents the "editing point" inside the table. Unlike the other Data Ally tables, the Calibration Point Table is mainly a reviewing tool rather than a programming tool and thus permits very limited editing.

Like all the other Data Ally tables, you can resize any table columns at any time by using the special header column button double-arrow cursor to "drag" column header borders to desired positions. You can also "drag" entire columns to any relative position in the table in order to achieve the most useful view of the table contents at any time.

Note that the COPY/CUT/PASTE/FILL command buttons at the bottom of the Calibration screen are not active for the Calibration Point Table window. This is because the values in the Calibration Point Table are derived from actual standard runs, and are not intended for manual editing. The only ways to remove standard run data which has been logged into the Calibration Table via manual or automatic calibration is to either use Setup again to restructure and reinitialize the table or to use the Clear Calibration/Clear Point Table command (Section 10 below).

Every line in the Table, representing a single standard injection with all its corresponding data, has a corresponding "data point" on the graphical plot in the Calibration plot window above. If standards at more than one concentration level have been run, a "multilevel calibration" has been produced, and a calibration curve is fit and plotted to the points in the Plot Window representing all the data.

The current "highlighted" position in the Calibration Point Table is dynamically linked to the Plot window, so that the plotted data point corresponding to that line is always itself highlighted (circled) in the plot. Moving the highlight to a different point in the plot automatically moves the table highlight to indicate the data for that new point, and vice versa.

10.2.3 The Calibration Plot

The Calibration Plot (Figure 10.7) is a visual tool for analyzing the quality and characteristics of a Calibration file for any peak which is being used to generate quantitative results for unknowns. All standard response data in the Point Table will always be plotted automatically, in a manner selected via Calibration Options (see Section 10.3.6 below). Usually, the default plot setup is for a plot of standard area versus concentration; the nature and axes units for this plot can be modified. Each time a new calibration run is completed and its result added to the calibration table, a corresponding point representing that result is added to the plot. Points corresponding to individual standard runs which are currently "active" (included in the Calibration Curve computation) are shown in one color in the Plot, while points which are "disabled" (via double-clicking the ENABLE/line number column) and not included in the current Calibration computation are shown in red. The Plot window normally uses AutoScaling to size its axes to contain the highest-value points from the current calibration. It is also normally defaulted to plot standard concentration (corrected concentration) versus basis (area or height).




Figure 10.8 Calibration Plot Window

Superimposed over the scatter plot of all the standard data points is the current "calibration curve" which is computed using the curve fit method currently selected in the Reference Peak Table for the current peak. The fit method can be changed easily for replotting. For close visual inspection of a calibration plot, it may be convenient to maximize the Plot window by clicking on its Program Control Maximize button at upper right. The Plot window can, of course, also be resized and moved to any desired location or shape on screen.

The current curve fitting coefficients and FIT mode are displayed in the upper left corner of the Plot window at all times for the peak whose data is being viewed.

The Plot window is normally autoscaled so that the scaling in both axes represents the range between zero and the maximum actual response and concentration values logged into the Calibration table. If no data are present in the Calibration table (i.e. no standards have actually been run yet), the scaling of the plot window represents the full range of the calibration levels defined in the Reference Peak Table. Each time you change the Fit type, Basis, standardization type (external or internal), or other relevant criteria defined in the Reference Table for any peak, the Plot will immediately be modified and replotted accordingly. You can use this feature to graphically and quickly test the effects of, for example, changing a curve fit for any peak, or deleting certain standard points from a calibration curve as a trial.

It is possible to use "rubber-band" zooming in the Plot Window to see an expanded view of the plotted information. This is done in exactly the same manner as normal zooming in the Graph windows. Zooming is useful to get a clearer view of regions in which several replicate points may be nearly overlapping.

➤ To execute a zoom inside the calibration plot, move the pointer cursor to the position in the plot where you want the zoomed view to begin, hold down the left mouse button, and rubber-band the desired view - when the dotted line on the plot shows this view, release the mouse button and the new view will be drawn. You can zoom repeatedly as many times as you like until the limit of resolution is reached. To unzoom, click the right mouse button to see the previous view. Repeated clicking of the right button will always return to the original autoscaled plot.

Note that when a peak is first selected via the Reference Table, one point in the Plot window is always highlighted (surrounded with a small circle), which corresponds to the top calibration run row highlighted in the Calibration Point Table. Double-clicking the left mouse button on a different point in the plot will cause the highlight in the Table to move instantly to the run represented by that point. By

the same token, clicking the edit cell anywhere in any row of the Calibration Point Table will cause the circle highlight in the Plot to move to the point represented by that row's standard injection. Thus, the Plot and Point Table windows are completely interactive.

If any rows are present in the Calibration Point Table, and/or any points are present in the Plot, for any peak, that peak has been calibrated and a result can be computed, since a response factor is present for the peak in memory. An unlimited number of additional calibrations (lines and points) can be added for any given peak(s). In typical operating mode, Data Ally will continue to capture cumulative calibration data into an open Calibration file until the user decides to either clear and restart a new calibration in that file (via the New calibration command option) or the user selects the Files/Clear command to remove all existing calibration information from the current file.

As soon as a single point appears for any peak in the Calibration Point Table, a plotted line will be drawn in the plot window between that point's position and the plot's origin. This line represents the relationship between concentration and response, and will be updated and re-averaged as additional points are added. Use of the Linear (L) and Quadratic (Q) curve fit options will result in a plot of multiple standard level data ranging from the highest standard run to the x or y axis of the plot - the y-intercept value will always be included in the plot. If the fit through zero option codes are used (L0 or Q0), the plotted line will always pass through the origin. For Interpolated fit options, the plot will only span the actual range of standards run and shown in the Calibration table. An Interpolated fit plot will only pass through the origin if one or more standards run has been designated as a zero concentration level (blank) or if the I0 code has been used for FIT to extrapolate through the origin.

If the current Calibration includes internal standardization, or correction factors for dilution, weight, injection volume, or user-defined factors, the Plot Window will generally show a plot of "corrected" concentration and/or basis (area or height) values, with the plot labeling indicating the correction. For example, if weight factors are used for all standards and unknowns, the standard curve in the plot window will actually be a graph of area/height versus a "corrected" concentration which equals the given concentration for each standard level in the Reference Table multiplied by its given Weight factor.

10.3.0 Theory of Calibration

Data Ally calibrations are performed using standard models for computing response factors and multiple level curve fits, with or without internal standards and the use of other correction factors to account for sample dilutions, injected volumes, and measured weights.

There are two basic modes of calibration available. In External Standard mode, standards consisting of varying known amounts of mixtures of the peak components to be measured in all unknowns are injected, and response factors or curves generated from the known responses of the standards. In Internal Standard mode, each peak to be quantitated is assigned one Internal Standard peak as a reference, and the responses computed via external standards are corrected on a peak-by-peak basis to account for differences in recovery of the Internal Standard peak(s) from run to run. Any individual peak in a Peak Table and any Calibration Reference Table using that Peak Table can be calibrated in External (ES) or Internal (IS) Standard mode, provided suitable Internal Standard peaks are incorporated in each run and identified as references in the Peak Table's IS column. All curve fitting, basis, and correction factor options are available, regardless of the mode selected for any peak.

10.3.1 External Standard Calibration

In External Standard calibrations, an ES code is placed in the CAL column of the Peak Table for all peaks to be calibrated in this mode. Up to sixteen different concentration levels for external standards are specified for each individual peak in the Calibration Reference Table, which must correspond to the exact concentration values expected for each external standard mixture injected.

For single point External Standard mode, the Response Factor is defined as

$$RF = \frac{\text{Corrected Concentration}}{\text{Basis (Area or Height)}}$$

where

$$\text{Corrected Concentration} = \frac{\text{Nominal Concentration (from Ref. Table LEV value)}(\text{Weight})(\text{Inj. Vol.})}{\text{Dilution}}$$

In most cases, the nominal concentration given for each level in the Reference Table should represent the actual concentration of the standard in the injected aliquot. If the injected volume will vary among standards and unknowns, an Injection Volume can be input for each sample. If a predilution of stock standards is to be injected, Dilution factors can be entered to compensate for the dilution of the injected material. If standards and unknowns are to be weighed out for injection, a Weight factor can be input for each sample. Note that if Injection Volume or Weight factors are input for any individual standard, they must be input for all standards and all unknowns for results to be obtained. If such a factor is not input, an error condition will result when attempting to perform or update a calibration or generate quantitative results for unknowns. Dilution factors can be entered as needed on an individual basis for standards or unknowns at any time - if no specific dilution factor is provided, Data Ally assumes the dilution factor is one.

The standard definition of RF above, the change in concentration versus change in basis, is represented by the default plotting of Corrected Concentration as the ordinate and Basis (area or height response) as the abscissa in the Calibration Plot window.

10.3.1.1 Single Level Calibration

For single level calibration, the results for each unknown are computed by the formula

$$\text{Actual Concentration} = \frac{(RF)(\text{Basis of unknown})(\text{Dilution})}{(\text{Weight})(\text{Injection Volume})}$$

Single-level calibration can include as many replicates of the same corrected concentration as desired, all replicates being averaged together with equal weighting. As long as at least one standard run appears in the Calibration Point Table, a single level calibration will be plotted in the Plot window with the theoretical response curve drawn from the plot origin through the logged response value.

For multiple level external standard calibrations, corrected concentrations for each standard are compared to measured peak areas or heights via regression fitting, utilizing either linear, quadratic, exponential, or interpolated methods, to determine the relationship between concentration and response over a range of values.

10.3.1.2 Multiple Level Calibration - Linear Regression

If linear regression curve fit is selected, regression is used to determine the coefficients which best fit the linear relationship

Corrected Concentration = KB (Basis) + KC, where

$$KB = \frac{n\sum C_x B_x - (\sum B_x)(\sum C_x)}{n\sum B_x^2 - (\sum B_x)^2} \quad (\text{slope})$$

$$KC = \frac{\sum C_x}{n} - \frac{(KB)(\sum B_x)}{n} \quad (\text{concentration intercept})$$

with C_x = the corrected concentration of component x
 B_x = basis (area or height) of component x
 n = number of standards run

Once the corrected concentration is found for a peak's basis, the actual concentration for that peak is computed by

$$\text{Actual Concentration} = \frac{(\text{Corrected concentration})(\text{Dilution})}{(\text{Weight})(\text{Injection Volume})}$$

For linear fits, a correlation coefficient (r) can also be computed as a measure of goodness of fit:

$$r = \frac{n(\sum C_x B_x) - (\sum B_x)(\sum C_x)}{\sqrt{[n\sum B_x^2 - (\sum B_x)^2][n\sum C_x^2 - (\sum C_x)^2]}}$$

In most instances, when linear regression fitting is used, the concentration intercept term KC will not be zero. If you wish to cause the KC term to be "forced" to go through zero, in order to equate a peak basis (area or height) of zero with a concentration of zero, you can use the "L0" curve fit option instead of the normal "L" fit. The L0 fit code causes the regression to be computed in the usual manner, as above, and then changes the KC term to zero regardless of its computed value to force the curve through the origin. It is advisable to evaluate the true responses of low-concentration standards before using the L0 option, in order to avoid inappropriate alterations of unknown results. With either L or L0 curve fits, the Calibration Plot will always draw the plotted theoretical response curve from the concentration intercept value through the maximum standard concentration/basis value found in the Calibration Point Table.

10.3.1.3 Multiple Level Calibration - Quadratic Regression

If quadratic (polynomial) regression fitting is selected, the existing standards will be used to determine the coefficients which best fit the polynomial curve

Corrected Concentration = KA (Basis²) + KB (Basis) + KC, where

KA =

$$KA = \frac{n[(\sum B_x^2)(\sum B_x^2 C_x) - (\sum B_x C_x)(\sum B_x^3)] + \sum B_x[(\sum C_x^2)(\sum B_x^3) + (\sum B_x C_x)(\sum B_x^2)] - (\sum B_x^2 C_x)(\sum B_x)^2 - (\sum C_x^2)(\sum B_x^2)}{n[(\sum B_x^2)(\sum B_x^4) - (\sum B_x^3)^2] + 2(\sum B_x)(\sum B_x^2)(\sum B_x^3) - (\sum B_x)^2(\sum B_x^4) - (\sum B_x^2)^3}$$

KB =

$$KB = \frac{n[(\sum B_x C_x)(\sum B_x^4) - (\sum B_x^3)(\sum B_x^2 C_x)] + \sum B_x[(\sum B_x^2)(\sum B_x^2 C_x) - (\sum C_x^2)(\sum B_x^4)] + \sum B_x^2[(\sum C_x^2)(\sum B_x^3) - (\sum B_x^2)(\sum B_x C_x)]}{n[(\sum B_x^2)(\sum B_x^4) - (\sum B_x^3)^2] + 2(\sum B_x)(\sum B_x^2)(\sum B_x^3) - (\sum B_x)^2(\sum B_x^4) - (\sum B_x^2)^3}$$

$$KC = \frac{\sum C_x^2[(\sum B_x^2)(\sum B_x^4) - (\sum B_x^3)^2] + \sum B_x C_x\{(\sum B_x^2)(\sum B_x^3) - (\sum B_x)(\sum B_x^4)\} + \sum B_x^2 C_x[(\sum B_x)(\sum B_x^3) - (\sum B_x^2)^2]}{n[(\sum B_x^2)(\sum B_x^4) - (\sum B_x^3)^2] + 2(\sum B_x)(\sum B_x^2)(\sum B_x^3) - (\sum B_x)^2(\sum B_x^4) - (\sum B_x^2)^3}$$

After determining the corrected concentration for any unknown peak from the standard curve, its actual concentration is computed from

$$\text{Actual Concentration} = \frac{(\text{Corrected concentration})(\text{Dilution})}{(\text{Weight})(\text{Injection Volume})}$$

If the "Q" curve fit option is used in the Calibration Reference Table, the concentration intercept KC will usually not equal 1. The "Q0" fit option can be employed to force a quadratic curve fit to go through the origin so that KC = 0. For the Q or Q0 fit options, the Calibration Plot window will always draw the theoretical best-fit response curve plotted from the concentration intercept value through the maximum standard concentration/basis found in the current Calibration Point Table.

10.3.1.4 Multiple Level Calibration - Interpolation

To use the Interpolated curve fitting option, select the "I" code in the Reference Table FIT column for any peak(s). Interpolated fitting computes linear curve fit segments between each corrected concentration level found in the current Calibration Point Table, with all replicates at any single level cumulatively averaged to find the response for that level. Within each segment, the same linear regression curve fitting procedure is applied as described above in Section 10.3.1.1, to determine a slope (KB) and intercept (KC) specific for that segment. The curve fit for each segment is used only to compute results for unknowns whose response basis falls within that particular segment. As many different segments are calculated as needed to represent the entire response curve in the Calibration Point Table. Although a maximum of 16 nominal concentration levels can appear in the Calibration Reference Table, the number of corrected concentration levels can exceed 16 if various correction coefficients such as dilution, weight, injection volume, or user defined factors are utilized.

With Interpolated fits, the Calibration Plot window always plots only the theoretical curve between each of the interpolated segments. The I0 "fit through zero" code can be used to extend an interpolated fit from the lowest standard response value in the Calibration Point Table to the origin. Unknown peaks which have a response basis below the lowest value in the plot or above the highest value will not be computed, and will be labeled "beyond interpolation limits" in any quantitative reports produced.

10.3.2 Internal Standard Calibration

In order to compensate for inter-run or inter-sample differences, internal standard components can be added to all standard and unknown samples and employed to correct for such variations. Although no single component peak in any chromatogram can have more than one associated internal standard peak, it is possible to have up to five different sets of peaks associated with five different internal standard components in a single chromatogram. It is also possible for some peaks in a chromatogram to be calibrated only with external standards, while others are referenced to internal standards.

Note that when Internal Standards are used, all calibration standards and unknowns must have values supplied by the operator indicating the concentrations of every Internal Standard peak. Normally, when running standards during calibration, the Internal Standard concentrations for each standard level will be listed in the Calibration Reference Table LEV columns for each IS peak, and will be employed automatically for the required calculations each time a standard is run and used to update the Calibration Point Table. Alternatively, the operator can key in exact values for each IS peak at each standard level manually or as part of a Method or Sequence program. If all required IS peak concentration values are not supplied when running a standard, those peaks which are calibrated against IS peaks lacking values will not have their current Calibration Point Tables updated. If all required IS peak concentration values are not supplied when running any unknown sample, those peaks calibrated against IS peaks lacking values will not be reported.

Components added to a sample mixture as internal standards must normally not be expected to occur in unknowns to be analyzed. It is helpful if internal standard peaks are well resolved from other

peaks, so that their areas can be clearly determined in each run. Use of differing Internal Standard peak concentrations is allowed for different standard levels during calibration, or for different unknowns. However, most internal standard calibration techniques involve the "spiking" of the same known amount of each internal standard component into all standards and unknowns, making programming and computational requirements very simple.

For a single-level internal standard calibration for any peak x, the "relative response factor" is defined as

$$RF_{Rel} = \frac{(\text{Corrected Concentration of } x)(\text{Basis of IS})}{(\text{Basis of } x)(\text{Concentration of IS})}$$

where IS is the internal standard peak used for peak x. By this definition, peaks used as internal standards will always have a defined $RF_{Rel} = 1$. Whenever any internal standards are used in any Calibration file's Reference Table, it is assumed that all calculated and reported RF values are RF_{Rel} .

Note that internal standard concentrations are never corrected in any manner for dilution, weight, or injection volume, since it is always assumed that the stated value for the amount/concentration of internal standard injected represents the exact injected value. In most cases, if internal standards are used, factors correcting for dilution, weight, injection volume, or other sample-specific differences are not necessary since the recovery of the internal standard itself will automatically compensate for any of these elements and conditions. The value input as IS AMT for any internal standard in the Graph/Edit/Calibrate or Analyze dialogs, the Method Sample Info dialog, or the Sequence table will always be returned in any result reports as the actual target value, regardless of any other correction factors used.

Correction factors such as Dilution, Weight, Injection Volume, and User-Defined factors may still be used with Internal Standardization, but such factors will be applied only to the final result after any internal standard correction is made. To avoid confusion and potential result errors, it is strongly recommended that such factors should not be employed when internal standard calibrations are performed.

10.3.2.1 Single Level Calibration With Internal Standard

For single level calibration, the results for each unknown are computed by the formula

$$\text{Concentration} = \frac{(RF_{Rel})(\text{Rel. Basis})(\text{Dilution})}{(\text{Weight})(\text{Injection Volume})}$$

where

$$\text{Rel. Basis} = \frac{(\text{Basis (Area or Height) of Unknown})(\text{Conc. of Internal Standard})}{(\text{Basis of Internal Standard})}$$

As for external standard calibrations, single-level internal standard calibration can include as many replicates of the same corrected concentration as desired. The basis used for each Internal Standard peak can be either area or height regardless of the basis of the unknown peak to be measured, as long as the same basis is used for any single Internal Standard peak in all standards and unknowns.

For multiple level internal standard calibrations, the corrected concentration for each standard level is compared to its measured relative basis (area or height) against its internal standard peak, via regression fitting, utilizing either linear, quadratic, exponential, or interpolated methods. Relative basis is computed for each unknown peak in each unknown run, and then employed to convert the relative RF obtained from the curve fit to the final result.

10.3.2.2 Multiple Level Internal Standard Calibration - Linear Regression

If linear regression curve fit is selected, regression is used to determine the coefficients which best fit the linear relationship

Corrected Concentration = KB (Rel. Basis) + KC, where

$$KB = \frac{n\sum C_x R_x - (\sum R_x)(\sum C_x)}{n\sum R_x^2 - (\sum R_x)^2} \quad (\text{slope})$$

$$KC = \frac{\sum C_x}{n} - \frac{(KB)(\sum R_x)}{n} \quad (\text{concentration intercept})$$

with C_x = the corrected concentration of component x

$$R_x = \frac{(\text{basis of component x})(\text{concent. of IS})}{(\text{basis of IS})}$$

n = number of standards run

The actual concentration of an unknown peak is then found as

$$\text{Actual Concentration} = \frac{(\text{Corrected concentration})(\text{Basis of IS})(\text{Dilution})}{(\text{Concent. of IS})(\text{Weight})(\text{Injection Volume})}$$

For linear fits, a correlation coefficient (r) can also be computed as a measure of goodness of fit:

$$r = \frac{n(\sum C_x B_x) - (\sum B_x)(\sum C_x)}{\sqrt{[n\sum B_x^2 - (\sum B_x)^2][n\sum C_x^2 - (\sum C_x)^2]}}$$

As with external standard calibrations, the "L0" curve fit option can be used instead of the normal "L" fit to change the KC term to zero regardless of its computed value, to force the curve through the origin. With either L or L0 curve fits, the Calibration Plot will always draw the plotted theoretical response curve from the concentration intercept value through the maximum standard concentration/basis value found in the Calibration Point Table.

10.3.2.3 Multiple Level Internal Standard Calibration - Quadratic Regression

If quadratic (polynomial) regression fitting is selected, the existing standards will be used to determine the coefficients which best fit the polynomial curve

Corrected Concentration = KA (Basis²) + KB (Basis) + KC, where

KA =

$$\frac{n[(\sum B_x^2)(\sum B_x^2 C_x) - (\sum B_x C_x)(\sum B_x^3)] + \sum B_x[(\sum C_x^2)(\sum B_x^3) + (\sum B_x C_x)(\sum B_x^2)] - (\sum B_x^2 C_x)(\sum B_x)^2 - (\sum C_x^2)(\sum B_x^2)}{n[(\sum B_x^2)(\sum B_x^4) - (\sum B_x^3)^2] + 2(\sum B_x)(\sum B_x^2)(\sum B_x^3) - (\sum B_x)^2(\sum B_x^4) - (\sum B_x^2)^3}$$

KB =

$$n[(\sum B_x C_x)(\sum B_x^4) - (\sum B_x^3)(\sum B_x^2 C_x)] + \sum B_x[(\sum B_x^2)(\sum B_x^2 C_x) - (\sum C_x^2)(\sum B_x^4)] + \sum B_x^2[(\sum C_x^2)(\sum B_x^3) - (\sum B_x^2)(\sum B_x C_x)]$$

$$KC = \frac{n[(\Sigma B_x^2)(\Sigma B_x^4) - (\Sigma B_x^3)^2] + 2(\Sigma B_x)(\Sigma B_x^2)(\Sigma B_x^3) - (\Sigma B_x)^2(\Sigma B_x^4) - (\Sigma B_x^2)^3}{\Sigma C_x^2[(\Sigma B_x^2)(\Sigma B_x^4) - (\Sigma B_x^3)^2] + \Sigma B_x C_x[(\Sigma B_x^2)(\Sigma B_x^3) - (\Sigma B_x)(\Sigma B_x^4)] + \Sigma B_x^2 C_x[(\Sigma B_x)(\Sigma B_x^3) - (\Sigma B_x^2)^2]}$$

$$n[(\Sigma B_x^2)(\Sigma B_x^4) - (\Sigma B_x^3)^2] + 2(\Sigma B_x)(\Sigma B_x^2)(\Sigma B_x^3) - (\Sigma B_x)^2(\Sigma B_x^4) - (\Sigma B_x^2)^3$$

The actual concentration of each unknown peak is then found as

$$\text{Actual Concentration} = \frac{(\text{Corrected concentration})(\text{Basis of IS})(\text{Dilution})}{(\text{Concent. of IS})(\text{Weight})(\text{Injection Volume})}$$

If the "Q" curve fit option is used in the Calibration Reference Peak Table, the concentration intercept KC will usually not equal 1. The "Q0" fit option can be employed to force a quadratic curve fit to go through the origin so that KC = 0. For the Q or Q0 fit options, the Calibration Plot window will always draw the theoretical best-fit response curve plotted from the concentration intercept value through the maximum standard concentration/basis found in the current Calibration Table.

10.3.2.4 Multiple Level Internal Standard Calibration - Interpolation

To use the Interpolated curve fitting option, select the "I" or "I0" code in the Reference Table FIT column for any peak(s). Interpolated fitting operates in exactly the same manner as for external standard calibrations; the only distinction is that Relative Basis is used to fit each linear response curve segment between adjacent concentration levels. Each segment has a unique slope (KB) and intercept (KC), which is used only to find results for relative responses falling within that segment.

As for external standardization, the Calibration Plot window plots only the theoretical curve between each of the interpolated segments, unless the I0 code is used to force the lowest standard response value in the Calibration Point Table to be extrapolated through zero level. With the I code, unknown peaks which have a response basis below the lowest value in the plot or above the highest value will not be computed, and will be labeled "beyond interpolation limits" in any quantitative reports produced.

10.4.0 Building, Reviewing and Modifying Calibrations

10.3.1 Building a New Calibration Table

Posting the results of standard chromatograms into a Calibration Point Table can be done manually or automatically in a Method or Sequence of Methods. Any existing Calibration file can be updated at any time, as long as that file's Reference Peak Table matches the Integration Peak Table used to find and identify peaks. NOTE THAT ANY STANDARD CHROMATOGRAM WHICH WILL BE USED TO PERFORM A CALIBRATION (UPDATE THE CALIBRATION POINT TABLE) MUST BE INTEGRATED USING THE MATCHING INTEGRATION PEAK TABLE SO THAT ALL PEAKS TO BE CALIBRATED ARE PROPERLY IDENTIFIED. To assure that this occurs, in some cases already-integrated chromatograms must be re-integrated with the desired Integration screen/Peak Table loaded in the current time base.

Each time a "calibration" is performed, one line is added to the Calibration Point Table, and its corresponding response point to the Calibration Plot, for each identified peak. Only the found peaks which match the Reference Peak Table specifications, and which have a concentration value assigned for the standard level being processed, can actually have results accumulated into the calibration. For example, if a certain series of peaks are not found in a new standard chromatogram, the Calibration Point Tables for those peaks will not and cannot have a new line appended representing the results for those peaks in that chromatogram.

If Internal Standardization is called for in the Integration Peak Table being used to identify peaks, valid entry values for the concentrations of all the assigned internal standard peaks must be present in the existing Calibration Reference Table, or must be entered by the operator in the Graph Edit/Calibrate dialog, Method Sample Info dialog, or Sequence table, as appropriate. The Calibration Reference Table window includes a column identifying which peaks represent internal standards, for convenience when assigning and reviewing such standards.

If Dilution or Weight correction factors are entered for any standard sample, they will be automatically accounted for in the calibration computations, and listed in the Calibration Table if suitable columns have been configured. It is recommended that if such correction factors are frequently used, appropriate columns should be included in the Calibration Table to assist operators in recognizing the effects of these factors upon calibrations.

To perform a manual calibration and add a standard response for each peak to the current Calibration file's Point Table, load the standard chromatogram desired into the Graph window and be sure the correct, matching Integration file is loaded in the current time base. Check to verify that the standard chromatogram has already been integrated and all its peaks are identified properly - if this is not the case, use the Graph Edit/Integrate command to integrate again. Then select Edit/Calibrate (Figure 10.9) and complete the required entries for correction factors such as Dilution, Weight, or Sample Volume (if configured), Internal Standard Amounts/Concentrations, Level No., and Mode (New, Average, Replace Level). Click OK to proceed with the manual calibration.




Figure 10.9 Using Manual Calibrate Function to Process a Standard Chromatogram

Click the Calibration SoftButton to view the Calibration screen - the previous Calibration Point Table and Plot will be updated with a new row/point corresponding to the newly processed standard chromatogram. All its relevant data will be shown in the Point Table.

10.3.2 Reviewing the Calibration Point by Point

Since the Table and Plot are interactive, it is possible to examine the nature of each individual plotted point in detail, by graphic selection.

First, click inside the Reference Peak Table window to make this the active window, and use the scroll bars to select the peak of interest with the highlight, and to update the Table and Plot windows.

➤ To select a point in the Plot, make the Plot Window active by clicking anywhere inside it. Now, move the pointer cursor to the point about which you would like information, and double-click the left mouse button on that point.

If the Plot window is too small and the points too closely placed for individual point selections, you can expand the Plot window to any desired size by clicking the expansion arrows for its sides or corners and "dragging" its borders to the new position. You can also rubber-band zoom inside the Plot window to see an expanded view. It is suggested that you leave the first two lines of the Table showing on the display at the bottom, so that you can see the information on each peak.

When you double-click on a single point inside the Plot window, the Table will automatically scroll backward or forward so that the descriptive line corresponding to that point will appear highlighted on the top line of the table (Figure 10.10). In this manner, you can check the detailed characteristics of any points in the Plot quickly and easily.




Figure 10.10 Automatic Indexing of Calibration Table After Plot Point Selection

➤ To select a row/standard in the Table, use the pointer and table scroll bars to move to any row in the table corresponding to the standard chromatogram whose corresponding point you would like to locate on the plot. Click anywhere on that row, and the point representing that line will be highlighted in the Plot window. You can zoom the plot display for closer study if you like.

You can use this mechanism to review, check, and verify any standard chromatogram or replicate included in the Calibration file, for any peak.

10.3.3 Editing Individual Standard Points

In addition to locating and studying individual points, you can also test the effect of each point on the current Calibration curve fit. To perform point-by-point editing, you must have configured the ENABLE column to be present in the Calibration table in this channel.

The ENABLE column is identical to that found in the other Data Ally program tables. Each cell in the column contains either a colored highlight, indicating that cell's row is enabled or active in the current calibration, or is blank, indicating that cell's row is disabled or inactive in the calibration. If this column is not present, you will not be able to edit individual replicate points.

➤ To disable a single point and remove it from the current Calibration curve, enter the Plot window and double-click on the desired point. The Table will immediately index to the line corresponding to that point. Now double-click on the ENABLE cell of the highlighted row in the Table - the ENABLE

highlight will be toggled off. The point in the Plot window will change to a red "x", and a new calibration curve fit will appear superimposed over the old fit (Figure 10.11).




Figure 10.11 Failing a Calibration Data Point

Alternatively, you can directly click on any row in the Table to be removed from the calibration, and double-click its ENABLE cell to toggle that standard replicate on or off.

➡ **To re-enable a point, which has been disabled**, repeat the above procedure and double-click on the un-highlighted ENABLE cell of any disabled row to re-highlight it and re-insert it into the curve fit.

Note that it is not possible to completely "delete" calibration values from the table under normal circumstances.

When you have finished designating all points/rows to be disabled for the current calibration for all peaks, you can either Analyze chromatograms manually in any Graph window using the modified calibration, or you can save the modified calibration with the same or with a different filename and then run new Methods or Sequences to process unknowns using the new calibration. If you have edited calibration points, you will be prompted to save the modified Calibration file before you run a new Method or Sequence or load a new Calibration. SAVING A CALIBRATION FILE WITH ROWS DISABLED FOR ANY PEAK(S) WILL PRESERVE THE DEFEAT OF THOSE ROWS UNTIL THEY ARE RE-ENABLED AND THE FILE RESAVED AGAIN.

10.3.4 Changing the Calibration Curve Fit

You can change the type of calibration curve fit performed for any peak at any time using the Peak Reference Table.

➡ **To change the curve fit type for a peak**, click inside the Reference Table Window to make it active, click inside the desired peak's FIT column cell, and change the entry as desired. Press ENTER, TAB, or CTRL-arrow, or click inside any other cell, to make the new entry. A dialog box will appear advising that the calibration is being changed and prompting for approval. Click OK to proceed with the change - the curve fit will immediately be recomputed and the new Plot drawn. Clicking Cancel will circumvent the change and restore the previous entry.

You can change the fit type as often as you like in this manner. When the desired fit is determined, you can manually process chromatograms using Analyze, or save the modified Method file and run new Methods or Sequences to process unknown samples.

The FIT parameter only changes the curve fit mathematics used if more than one level of standards have been run and included in the Table. If only a single level is present, the FIT column button will be non-functional until additional levels are entered.

10.3.5 Changing the Calibration Basis

Each peak in the Calibration file can have its own individual "basis" for calibration, either area or height. You can change the basis for any peak by modifying the BASIS column code for that peak in the Reference Peak Table; this can be done in conjunction with changing the current plotting/weighting for any peak, or its curve fit type.

➡ *To modify the basis for a given peak*, first click on that peak's BASIS column cell in the Reference Peak Table to select it, and type in the desired code (A or H). Press ENTER, TAB, or CTRL-arrow, or click inside any other cell, to make the new entry. A dialog box will appear advising that the calibration is being changed and prompting for approval. Click OK to proceed with the change - the Plot window will instantly change to show the effects of using the new basis. You can use this switching to quickly assess the impact of changing basis from area to height or vice versa on any peak. Clicking Cancel will disable the command and restore the previous entry.

10.5.0 Working With Calibration Files

The Calibration filing system permits you to save and recall calibrations quickly and easily. The contents of the current Calibration table and its Reference Peak Table, as well as all the options and plot settings, can be saved as a "Calibration" file, with a unique filename and description, so that it can be recalled and reused at any future time either manually, in an automated Method, or in a Sequence of Methods.

You can specify a default Calibration file by name during Configuration, which is always loaded automatically into the calibration display screen at the time Data Ally is initialized. You can change any aspects of the default file within the active software, except for the features and options which are either "fixed" or "removed" in the Configuration, but you can only save such changes to the default file in Configuration Mode (Section 5.2.5.4). If you wish, it is possible to run a series of standard chromatograms, use them to update any calibration file you create, and then use that file to "assign" parameters to the default Calibration file in Configuration mode.

The Files commands in the Calibration Command Menu bar are used to load, save, clear, and delete Calibration files, and can also provide a calibration version historical record.

10.5.1 Loading Calibration Files

➡ *To load a saved Calibration file from disk memory into the current channel*, click the left mouse button on Files in the Command menu bar, and then click on Load. A directory dialog box will appear (Figure 10.12) which lists the Calibration files (*.CAL) available in the current directory path for this channel. Select the desired file by clicking on its name in the directory list, or by typing its name into the entry box. Click on OK to load the new file.




Figure 10.12 Calibration Files/Load Dialog

If the previous file has been edited and has not been saved or resaved, an information message will appear before the new file is loaded asking if the previous file should be saved. If so, click on YES and assign either the same filename or a new name as desired. If you do not wish to save the edited file, click on NO to proceed.

When a new file is loaded, it will overwrite the previous Calibration file.

At any time, you can re-load the default Calibration File which is set in Configuration (Section 5.2.5.4). This file is unique in that it can only be modified in Configuration mode - it can contain a complete set of calibration data, if desired. To load the default file, use the Files/Load Default command, which will show the name of the default calibration. When this command is executed the default file will overwrite any existing calibration file present in the current time base.

10.5.2 Clearing Calibrations

The Clear function allows individual clearing of either the Calibration Point Table (the standard data) or the Reference Table (the Level values, plus fit, basis, and units parameters) without using Setup or deleting an existing file, or both.

➡ To clear any part of the current Calibration file, click on File in the Command Menu bar and then click on Clear. The Clear Select submenu will appear (Figure 10.13).




Figure 10.13 Clear Submenu for Calibration Files

Click on the desired command to clear either or both tables. If the previous file has been edited and has not been saved or resaved, an information message will appear before this file is cleared asking if the previous file should be saved. If so, click on YES and assign either the same filename or a new name as desired. If you do not wish to save the edited file, click on NO to proceed.

10.5.3 Saving Calibrations

➤ To save a modified Calibration file using the same filename, click on Files in the Command Menu bar, and then click on Save. A trap message will appear prompting for confirmation that file should be resaved with the current name, overwriting the previous file. Click on OK to proceed and resave - the file's internal version number will be updated.

➤ To save a modified Calibration file using a new name, click on the SaveAs item in the Files menu. The SaveAs dialog box will appear (Figure 10.14). Click inside the filename entry field and type in the desired new filename. If you wish, click inside the DESCRIPTION field and type a description of the new file, which may help identify this file in the directory listing for future recall. When the correct name and description appear, click on OK to save the file to the current disk path. If another file already exists in that path with the same name, an error message will instruct you to change the filename before saving.




Figure 10.14 Calibration Files/SaveAs Dialog

➤ To apply a description to the current Calibration file, click on Calibration Description in the Files submenu. The Description dialog box (figure 10.15) will appear. You can type in a single or multi-line description which will be saved with the file and will appear in the Description box in the Load directory (Figure 10.12). Alternatively, you can enter the description directly into the Description box in the Save As dialog (Figure 10.14) at the time a new or modified Calibration file is saved.




Figure 10.15 Set Calibration Description Dialog

10.5.4 Deleting Calibrations

➡ To delete a Calibration file from the current disk directory, use the Windows File Manager Delete function.

10.5.5. Printing Calibrations

Use the Files/Print Calibration command to print a complete summary of the current calibration file. Normally, each single page in such a printout will summarize the tabular data in the Calibration Point Table for one peak over the complete "history" of the current Calibration file. You can modify the format of this report by using Excel to edit the \reports\cal\cal.xls template file (see Section 13).

11.0 Running and Reprocessing Batches: The Sequence Screen

Data Ally can operate in conjunction with and control virtually any automatic sampler for the capture and reprocessing of batches of chromatogram files. All batch processing operations are programmed in the Sequence screen, which is viewed by clicking on the SEQUENCE SoftButton on the left screen.

The Sequence screen is useful if more than one sample or chromatogram is to be run automatically or an autosampler is being controlled by Data Ally. If only single injections are required, the entire process of acquiring, saving, analyzing, calibrating, and reporting chromatograms can be automated via the Method program screen alone. Effectively, the Sequence function allows you to run a series of Method programs as a single operation. It is possible to configure each Data Ally channel so that the Sequence function is not accessible to users, if it will not be needed in that channel.

11.1.0 Sequence Screen Functions

Click on the SEQUENCE SoftButton to display the Sequence program screen. This screen consists of two sections: a "top-level" dialog box in the upper screen, which includes settings and parameters affecting the entire Sequence, Autosampler Dialog access via the "Autosampler" button and the "Sequence Table" in the lower screen, which is a "spreadsheet" containing a list of samples and Methods to be processed. The Table portion of the screen operates in a very similar fashion to that of the Method Events and Integration Peak Tables.

The default Sequence file selected in configuration will be automatically loaded into the Sequence screen display whenever each channel is initialized, or when a different Sequence is cleared. It is recommended that the default Sequence file name be programmed for the most common running conditions, to insure the least possible amount of subsequent programming required for routine users.

If you wish to "split" the screen with a Graph window on top and the Sequence Table below it, click on the GRAPH button at bottom right. This type of display is useful during the processing of live Sequences, since it allows viewing of chromatograms developing in real time as well as the current injection position in the sample table.

11.2.0 Top-Level Sequence

The "top-level" Sequence dialog identifies global sequence processing parameters (Figure 11.1).




Figure 11.1 Sequence Screen


➤ To program a new Sequence or modify a Sequence program already loaded, click the left mouse button with the pointer cursor inside any of the entry fields, and enter the desired value.

Like other program files, Sequence files are assigned names and descriptions at the time they are saved using the Files/SaveAs command (Section 11.5.3 below).

11.2.1 Operating Modes

There are two Modes available for Sequences, listed in the very top left of the dialog (Figure 11.1)..

- **Acquire Mode:** For acquiring new chromatograms from live detector inputs as specified in the Method program files associated with each line in the Table. This mode is not available in any Edit Time Base.
- **Reprocess Mode:** For reprocessing chromatograms already acquired in previous runs and saved to disk. Reprocess loads each chromatogram file specified from a given disk location, integrates and analyzes the file as instructed by the elected Method, and creates reports according to the selected Integration file(s). Reprocess will perform the above operations as quickly as possible and in a continuous manner. Reprocess mode can be accessed in both "live" and Edit time bases.
- **Summarize Mode:** For producing tabulated Summary Reports from an entire sequence of runs which have existing integration results. **NOTE:** A Summary Report can also be automatically produced at the end of an Acquire Mode sequence by enabling the Summary Report and choosing a Summary Template via the Summary Dialog box. Click on the Summary button located in the Reports sections of the top level sequence screen to display the Summary Dialog box. See example of a Summary Reports dialog box below:



Summary Reports Dialog

➤ To run a Sequence in any mode, click the RUN button **while the Sequence screen is displayed** - IF THE SEQUENCE SCREEN IS NOT VIEWED, A SINGLE METHOD WILL BE RUN. Only one mode can be selected at any one time, and only those modes which are supported by the current Configuration will be available; the mode cannot be changed after a Sequence has been started. If the current channel is to be used for data acquisition only, only the Acquire mode should be allowed in Configuration. For an Edit-Only channel, only the Reprocess modes are available.

Select the desired mode by clicking on its button to enable it before using RUN to start.

➤ To stop a running Sequence, click the STOP SoftButton in any Data Ally display screen. Both the currently executing Method and the Sequence will be halted.

The top-level Sequence also contains fields for START and END rows in the Sequence Table. These operate interactively with the ENABLE column (see below) to define which lines in a preprogrammed Sequence table will actually be executed each time a Sequence is run.

➤ To limit Sequence processing to any consecutive series of rows in a programmed Sequence Table, click inside the START entry box and type the number of the first row to be processed, then

click inside the END box and type the number of the last row to be processed. Only those rows inclusive of the START and END row number will actually be executed in the Sequence.

You can also select whether or not the reports associated with the Methods to be run as part of a Sequence will be printed or shown on screen. In order to print, the Print checkbox in top center of the screen must be checked - in order to show programmed reports on screen, the Show checkbox must be checked. Leaving both checkboxes unchecked during a Sequence run or reprocess will cause all reports, either printed or shown, to be skipped during the entire Sequence. Note that you cannot check or uncheck these boxes after a Sequence has been started in order to change the reporting mode during a Sequence. The most common utilization of the check boxes for reporting is to defeat reporting without changing the report specifications for any of the programmed Methods, for example, during a reprocess Sequence in which a recalibration is performed.

At top right of the Sequence screen is another group of check boxes inside a boxed area marked Options. These check boxes are for Validation, Suppression of Integration, and Suppression of Calibration during the Sequence.

Checking the Validation box causes the Sequence to verify that all programmed files (Method, Integration, Calibration, chromatogram raw data) called for during all rows in the Sequence table are in fact present in the defined directory paths. This is a useful feature to avoid errors during long term, unattended Sequence processing - the drawback to checking the Validation check box is that a certain amount of time is required before actual Sequence processing begins in order for Data Ally to perform all the necessary file checks. It is recommended that the Validation feature be used if multiple-row Sequences will be processed without operator attention such as during overnight sessions.

Checking the Suppress check box for either Integration or Calibration disables any integration or calibration commands which may be included in the Methods incorporated into the Sequence table. These commands are most useful if you wish to recall previously-integrated (and manually edited) chromatograms for reprocessing purposes, or if you wish to use only the current Calibration file data sets in order to produce reports. Suppressing Integration or Calibration by checking the respective boxes will avoid any new integrations, which could change the quantitation of identified peaks, or calibrations, which could inadvertently modify the calibration response factors.

11.3.0 Sequence Table

The Sequence Table (Figure 11.1) is a queuing list for processing a series of samples using any selected group of Methods. Each line or row in the "spreadsheet" corresponds to a single vial on the typical autosampler, from which one or more injections can be made. Processing of an entire programmed Sequence can be done continuously, under Data Ally control, or the series of runs and samples programmed can be interrupted at will by the operator. You are free to disable or modify any lines or commands in the Sequence Table which have not yet been executed, including the number of samples or injections, the Methods used, or other information. Sequence Tables can be of any desired length, but will usually correspond to the maximum tray size assigned for the current autosampler, if this is chosen in Configuration.

11.3.1 Sequence Table Functions

Unlike the Method and Integration tables, the Sequence Table has no "graphical" editing functions, since it is basically a worksheet list of runs to be performed and has no relationship to graphic events. The Sequence Table can be edited in a number of ways, and a set of FILL functions are available to make programming new Sequences less tedious.

The simplest way to make any entry in the Table is to move the highlighted edit cell inside the desired cell (column and row), click to anchor the highlight there, and type in the desired value. You can use the <ENTER>, <TAB> and PgDn/PgUp keys to move around within the Table when reviewing or editing, as well as the usual CTRL-Home, CTRL-End, CTRL-PgUp, CTRL-PgDn, and CTRL-arrow

key commands. The vertical and horizontal scroll bars at the bottom and right of the Table can be used in conjunction with the mouse to quickly change your position in the Table when editing.

Each time a new line is added to the Table, a new number will automatically be assigned in the ENABLE/line number column at far left. This column is always visible regardless of the view of the table to make it easier to locate and edit particular items or entries.

The "spreadsheet" commands, INSERT, DELETE, COPY, PASTE, and CUT, at the bottom of the Table operate in exactly the same way as for the Events, Peak, and Calibration Reference Tables.

➡ To INSERT new lines into the Table, position the edit highlight anywhere in the line/row below the intended insertion point, and click to anchor it. Then click the INSERT button to insert one or more empty lines above that position. You can also insert new lines at the end of the Table by moving the highlight to the last cell in the table and pressing ENTER.

➡ To delete a row, move the edit highlight into any column in the row to be deleted, and click on DELETE. The line will be removed. The DELETE command differs from CUT in that a row's contents cannot be PASTEd back into a different position after being DELETED.

➡ To copy a row, move the edit highlight to the desired row to be copied, and click to anchor. Then click on the COPY button to copy an image of that row.

➡ To cut (move) a row from one position in the table to another, click the edit highlight in the desired row to be moved, and click on the CUT button.

➡ To paste a copied or cut image of a row into the Table, move the edit highlight anywhere in the row below the position for the row to be pasted. Then click on the PASTE button to paste the copied or cut row above the highlighted row.

If you INSERT or CUT rows from a Sequence, the line number listing in the ENABLE column will automatically be altered to renumber rows accordingly.

You can re-order the columns in the Sequence Table from left to right or change their individual widths at any time. Both the column order and widths are set as defaults in Configuration. To move a column to another position in the table, simply place the pointer cursor on that column's header button, hold down the left mouse button, and "drag" the header button left or right to the new desired position relative to the other columns. When the pointer is in the desired position, release the mouse button and the entire column will be moved, and the table redrawn accordingly. Modifying the order of columns may make certain types of table editing or review easier by allowing placement of all columns currently being manipulated in view on a single screen - you can easily move "unimportant" column information to the right and off of the displayed table window area.

To change the width of any individual column, move the pointer cursor to the border of its header button and the next column header button in the direction in which you wish to expand or reduce the column width. The cursor will become a "double-sided horizontal arrow" (↔); when this change occurs, immediately hold down the left mouse button and move the arrow cursor in the direction you wish the button edge to expand or contract. When the desired column width is indicated by the position of the cursor, release the mouse button and the column will be resized and the entire table redrawn accordingly. This feature permits you to match the exact sizing of each column to the expected character width of the entries expected in that column - you can also effectively "remove" certain columns without actually eliminating them from the table by reducing their widths to a very small size.

11.3.2 Sequence Table Columns

The following is a listing of possible columns which can be included in a Sequence Table. Typically, not all items are present in any single Table. It is indicated whether or not each individual item is configurable or must always be present in the Table.

- **ENABLE:** This field determines whether any individual lines/rows in the Sequence Table will be executed once the Sequence is started. Each row in this column can be either highlighted (enabled) or not highlighted (disabled) - disabled rows will not be executed. Double-clicking the left mouse button inside any cell in this column will "toggle" the highlight on or off. Disabling or enabling of rows in this manner can be done while a Sequence is executing if desired.
- **NO: (Configurable)** This field specifies VIAL NUMBER for an autosampler tray if autosampler control is included and configured (this column should only be configured if such control is provided with the Data Ally system for the current time base).
- **# INJS: (Configurable)** This field sets the variable number of injections **(from 1 to 9)** for the current row/sample in the Sequence table.
- **SMP. NAME: (Configurable)** This field defines the sample name to be run on any line/row in the Table. The sample name can have up to 30 characters, depending upon the column width for SAMPLE NAME set in Configuration.
- **SMP. DESC.: (Configurable)** This field can be used to enter a sample description of up to 30 characters in length.
- **CG FILENAME (For each configured detector input): (Configurable)** This field defines the chromatogram raw data path\filename to be created/saved for the incoming data from the first detector configured. Although the name can be specified in the Method file, this row entry permits a different filename to be used on a row by row basis, if desired, from the name specified in the selected Method file - any entry in the CG FILENAME fields will override the filename entry for this detector signal input in the corresponding Method file. sample number for the line/row, with up to 30 characters depending upon the column width value set in Configuration. Note that if more than one detector input is available in this time base, each detector can have its own chromatogram filename set independently of each other and of the Method file specifications. The ability to enter unique chromatogram filenames for each row also allows the creation of "reprocessing" sequences which handle any group of data files from any source and directory.

Note: Exclusive of the Drive and Path portion of the filename, the actual filename is limited to 5 characters maximum because the right most 3 character positions will be used to construct line number/injection number for each file. XXXXX011.PDF for example would be the first injection from the first sequence line. XXXXX021.PDF would be the first injection from sequence line number two. The PDF extension is automatically added by Data Ally.

- **INT. FILENAME (For each configured detector input): (Configurable)** Like the CG FILENAME field, this column permits Integration files to be assigned in the Sequence independently of those selected in any included Methods, on a sample or row-specific basis if desired. This field is used to enter a path/filename for integration and peak identification - a separate column will be available for configuration for each detector input assigned in each time base. Entries in the INT. FILENAME columns will override any Integration filename entries in any Method in the same row in the Sequence table. This feature provides maximum flexibility for either modifying how a Sequence will be processed without actually changing and resaving any Method files, or for reprocessing a series of chromatogram files requiring vastly different integration criteria.
- **CAL. FILENAME (For each configured detector input): (Configurable)** As with the CG and INT. FILENAME fields, this column permits Calibration files to be used as desired in the Sequence regardless of those selected in any included Methods, on a sample or row-specific basis if

desired. This field is used to enter a path/filename for calibration information - a separate column will be available for configuration for each detector input assigned in each time base. Entries in the CAL. FILENAME columns will override any Calibration filename entries in any Method in the same row in the Sequence table.

- **RPT. TEMPL.** (For each configured detector input): (Configurable) This column permits Excel report template files to be used in the Sequence independently of those selected in any included Methods, on a sample or row-specific basis if desired. This field is used to enter a path/filename for the *.XLS file describing the desired report. In the current version, only one such report template can be entered per input in a single Sequence. These entries will override any Integration filename entries in any Method in the same row in the Sequence table.
- **RPT. FILE** (For each configured detector input): (Configurable) Like the previous fields, this column permits independent specification of saved report file names for each defined detector input, regardless of those names selected in any included Methods, on a sample or row-specific basis if desired. This field is used to enter a path/filename for the report(s) defined either as part of called Methods in the Sequence or by the RPT. TEMPL. column entry. Entries in this column will override any report save filename entries in any Method in the same row in the Sequence table. Again, this offers ultimate flexibility in determining specific report generation for each Sequence without having to review or modify any saved Methods you wish to use.
- **CAL. HIST. TEMPL.** (For each configured detector input): (Configurable) This column permits row by row definition of the Calibration History reports to be prepared describing the current calibration for each configured detector input. These entries will override any Method-based entries for calibration reports. The CAL. HIST. column identifies a path/filename for the *.XLS Excel template file used to create each desired calibration report.
- **WEIGHT:** (Configurable) This factor allows corrections for weighed samples, as configured in the Method default set screen.
- **DILUTION:** (Configurable) This factor allows corrections for sample dilutions, as configured in the Method default screen.
- **METH:** (Configurable) This field determines which Method program is loaded from disk when this Sequence Table is executed. A different Method can be associated with each line in the Table, in any order or combination. Note that the Method references Integration and Calibration files in the normal manner. If the METH column is disabled in Configuration, the system will always load the default Method file for all lines in every Sequence Table. Entries in the METH column rows are always valid Method filenames, with or without path names. You can double-click on any row in the METH column to see a pop-up listing of all available Method files in the current path.

Method Name must be set for all rows in the table, but it is also possible to modify most of the parameters in each Method on a sample or row-specific basis by configuring and using the series of sample-specific entry rows for chromatogram filename, integration and calibration filename, report template and saved filenames, and calibration report filename described above. Any entries in any of these columns will override any entries of the same parameters in the referenced Method file on the same row in the Sequence table. Note that it is thus possible to create a Sequence table using one basic source Method file and to modify the actual processing of each row in the Sequence in a nearly unlimited manner without incorporating many different Method files.

- **CAL. CODE:** (Configurable) This field identifies the current line/row as a standard sample which is to be used to update the selected Calibration file(s). Entries in the CAL. CODE field can be either **N**, start new calibration file, **R**, replace existing level, or **A**, average new result into current

calibration. Entries in this field will override any calibration instructions contained in any Method files in the Sequence on a row specific basis.

- **CAL. LEVEL:** (Configurable) This field selects the calibration level for the standard defined in the current line/row. The range of level entries depends upon the number of levels programmed for the currently referenced Integration file(s). Entries in this field will override any calibration instructions contained in any Method files in the Sequence on a row specific basis.
- **IS1 AMT - IS5 AMT:** (Configurable) These five fields set concentration values for the Internal Standard peaks number 1 through 5, if such peaks have been defined in the Integration Peak Table. These columns are only configured for as many Internal Standards as are to be used.
- **PRT. CG:** (Configurable) This column defines whether or not the chromatogram generated during processing of each row in the Sequence Table will be printed according to its respective Method instructions. Entries in this column are either Y for Yes or N for No. Any entry in this column, if configured, will override any printing instructions contained in any Method specified for any sample/row.
- **COMMENT:** (Configurable) The COMMENT field for each line/row can be used to make notes on a line-by-line basis which will be saved with the Sequence Table. Comments can be up to 30 characters depending upon the column width setting in Configuration.

Note that many of the Sequence Table column parameters can be defined either here or in the Method screen Sample Information button dialog box. If a Method is designated in a Sequence Table which has Sample Information associated with it, those values will be automatically used for processing or reporting as long as the corresponding cells in the Sequence Table remain blank for all lines/rows employing that Method. **ANY AND ALL SAMPLE-SPECIFIC ENTRIES MADE IN THE SEQUENCE TABLE WILL ALWAYS OVERRIDE ANY CORRESPONDING ENTRIES MADE IN ANY METHOD FILES.**

11.3.3 Fill Functions

Information can be entered into the Sequence Table by using the automatic FILL functions associated with certain columns in the Table. The FILL button beneath the Table can be used in a similar manner to that in the Method Events or Integration Peak Tables to "fill" series of new values into blank column cells. The AutoFill dialog box (Figure 11.2) used to set all filling functions is of the same format and works in exactly the same manner as that for the other Data Ally tables (see Section 8.3 for a detailed explanation of its features and terminology).

The "Duplicate" filling mode is useable for the # INJS, METH, CAL CODE, CAL LEVEL, ISAMT, WT, DIL, SAMPLE NAME, SAMPLE DESC., SAMPLE NO., and COMMENT columns, as well as all configured filename and report specification fields.

The "Increment" filling mode is operable for the # INJS, CAL LEVEL, ISAMT, WT, DIL, and SAMPLE NO., columns.

In Data Ally Sequence programming, incrementing filling is normally used to begin by setting up a series of METH entries to define the length of the Sequence table, one line/row per sample. Other information can then be entered using fill functions or via manual entries. .




Figure 11.2 Auto Filling the METH Column

The following AutoFill operations can be performed for the Sequence table:

➤ To fill an existing row's value into all currently empty rows in the column up to the next entry, place the edit highlight on the "basis" row, click FILL to see the AutoFill dialog, select "Duplicate", "Current Cell", and "To Next Entry", with all other fields and options blank/off. If the "basis" row is the first row in the column, you can select "Start of Table" instead of "Current Cell" without highlighting the basis row.

➤ To fill an existing row's value into all currently empty rows in the entire column, place the edit highlight on the "basis" row, click FILL to see the AutoFill dialog, select "Duplicate", "Current Cell", and "End of Table", with all other fields and options blank/off.

➤ To fill an existing row's value into more or less empty rows than the number of empty rows found in the current column up to the next entry, place the edit highlight on the "basis" row, click FILL to see the AutoFill dialog, select "Duplicate", "Current Cell", and "To Next Entry". Type the actual number of empty rows to be filled into the "No. of Steps" column, with all other fields and options blank/off. If the "No. of Steps" value exceeds the number of actual intervening empty rows, sufficient new rows will be added and filled to meet the "No. of Steps" criteria.

➤ To fill all currently empty series of consecutive rows in the column with the values found in the rows prior to each empty group of rows, click FILL to see the AutoFill dialog, select "Duplicate", "Start of Table", "End of Table", and "Allow Progressions", with all other fields and options blank/off.

➤ To overwrite existing values (along with any blank rows) with duplicated values in an entire column, place the edit highlight on the "basis" row for copying, click FILL to see the AutoFill dialog, select "Duplicate", "Current Cell", and "End of Table", and check the "Allow Overwrites" command, with all other fields and options blank/off.

➤ To clear (remove) existing values in an entire column, insert a blank row at the beginning of the table, place the edit highlight on that row, click FILL to see the AutoFill dialog, select "Duplicate", "Current Cell", and "End of Table", and check the "Allow Overwrites" command.

➤ To overwrite existing values in a portion of a column with a duplicated value, place the edit highlight on the "basis" row for copying, click FILL to see the AutoFill dialog, select "Duplicate", "Current Cell", "End of Table", and "Allow Overwrites", and enter a "No. of Steps" value corresponding to the number of existing rows to be filled with duplicates. All other fields and options should be blank/off. If the

number of steps exceeds the number of blank rows in the column, the "Allow Insertions" box must be checked to add sufficient rows to meet the "No. of Steps" criteria.

➤ To clear (remove) existing values in a portion of a column, insert a blank row at the beginning of the range to be cleared (or remove an existing row entry), place the edit highlight on that row, click FILL to see the AutoFill dialog, select "Duplicate", "Current Cell", "End of Table", and "Overwrite Existing Data", and enter a "No. of Steps" value corresponding to the number of existing rows to be cleared (filled with a blank entry). If the number of steps exceeds the number of blank rows in the column, the "Allow Insertions" box must be checked to add sufficient rows to meet the "No. of Steps" criteria.

➤ To fill a column in a new table with incremented values, enter the beginning "basis" value in the first row, place the edit highlight on that row, click FILL to see the AutoFill dialog, select "Increment", "Current Cell", and "End of Table", enter the increment value in the "Increment" field, and enter the number of rows for which incremented values are to be inserted in the "No. of Steps" field.

➤ To fill a column with incremented values in blank rows up to the next existing entry, place the edit highlight on the "basis" row, click FILL to see the AutoFill dialog, select "Increment", "Current Cell", and "To Next Entry", and enter the increment value in the "Increment" field; leave other fields and options blank/off.

➤ To fill an existing row's incremented value into more or less empty rows than the number of empty rows found in the current column up to the next entry, place the edit highlight on the "basis" row, click FILL to see the AutoFill dialog, select "Increment", "Current Cell", and "To Next Entry". Type the increment value into the "Increment" field, and the actual number of empty rows to be filled/incremented into the "No. of Steps" column, with all other fields and options blank/off. If the "No. of Steps" value exceeds the number of actual intervening empty rows, sufficient new rows will be added and filled to meet the "No. of Steps" criteria.

➤ To fill all currently empty series of consecutive rows in the column with incremented values based on the rows prior to each empty group of rows, click FILL to see the AutoFill dialog, select "Increment", "Start of Table", "End of Table", and "Allow Progressions". Type the increment step value for each row into the "Increment" field, with all other fields and options blank/off.

➤ To fill a column with incremented values in all blank rows based on an existing entry, place the edit highlight on the "basis" row, click FILL to see the AutoFill dialog, select "Increment", "Current Cell", and "End of Table", and enter the increment value in the "Increment" field; leave other fields and options blank/off.

➤ To fill a range of rows with evenly-spaced, automatically computed values, with each row's value corresponding to the difference between the basis cell's value and the value of the next row having a non-blank entry in the column divided by the number of intervening rows, place the edit highlight on the "basis" cell (first row in the range), click FILL to see the AutoFill dialog, select "Increment", "Current Cell", and "To Next Entry". If you wish to fill only the existing number of blank intervening rows, leave the "No. of Steps" field and all other commands blank/off. If you wish to fill either only some of the intervening rows, or to add new intervening rows, enter a number in the "No. of Steps" field specifying the total number of intervening rows to be filled. If "No. of Steps" exceeds the number of existing blank lines in the range, you must check the "Allow Insertions" box to automatically add and fill sufficient lines to match the "No. of Steps" entry.

➤ To fill all existing ranges of empty rows in the column with evenly-spaced, automatically computed values, with each filled row's value corresponding to the difference between the value of the cell immediately preceding that range and the value of the next row having a non-blank entry following the range divided by the number of intervening rows, click FILL to see the AutoFill dialog, select "Increment", "Start of Table", "End of Table", and "Allow Progressions". Leave all other selections blank/off.

➡ To overwrite existing values (along with any blank rows) with incremented values in an entire column, place the edit highlight on the "basis" row for incrementing, click FILL to see the AutoFill dialog, select "Increment", "Current Cell", and "End of Table", enter the increment value in the "Increment" field, and check the "Allow Overwrites" command, with all other fields and options blank/off.

➡ To overwrite existing values in a portion of a column with incremented values, place the edit highlight on the "basis" row for copying, click FILL to see the AutoFill dialog, select "Increment", "Current Cell", "End of Table", and "Allow Overwrites", enter the increment value in the "Increment" field, and enter a "No. of Steps" value corresponding to the number of existing rows to be filled with duplicates. All other fields and options should be blank/off. If the number of steps exceeds the number of blank rows in the column, the "Allow Insertions" box must be checked to add sufficient rows to meet the "No. of Steps" criteria.

➡ To fill information from the current Method into a column for the entire table, click on the desired column header button (or move the edit highlight into that column), click FILL to see the AutoFill dialog, select "Fill From Method", "Start of Table", and "End of Table", leaving all other fields and options blank/off. If the "Allow Overwrites" box is checked, any existing values in the selected column will be replaced by newly filled values from the current Method.

➡ To fill information from the current Method into a portion of a column, click on the desired column header button with the edit highlight in the first row to be filled in that column), click FILL to see the AutoFill dialog, select "Fill From Method", "Current Cell", and "End of Table", and set the number of rows to be filled with the value from the current Method as "No. of Steps". If existing row values are to be overwritten, you must also check the "Allow Overwrites" box. If the "No. of Steps" entry exceeds the number of existing blank rows to be filled, the "Allow Insertions" box must be checked in order to automatically add sufficient new lines to meet the "No. of Steps" filling criteria.

Alternatively, you can fill a range of blank rows in any column from the edit highlight cell up to the next entry in that column by using the "Fill From Method", "Current Cell", and "To Next Entry" commands, leaving all other fields and options blank/off.

11.4.0 Running a Sequence

A properly-programmed Sequence can be started simply by clicking on the RUN SoftButton at lower left while the Sequence screen is in the active window. This will cause the first Method in the Table to be loaded along with its referenced Integration and Calibration files, if any, and will start that Method's execution. Any specific file references programmed in the Sequence table for that row will be used in place of those specified in the Method, if different. After the first Method is complete, the second row Method will be loaded and executed, and all other rows in the Table processed in turn.

The Sequence can be aborted at any time by clicking on the STOP/RUN SoftButton. If a Sequence is aborted before completion, it can be restarted at the Table row where the termination step occurred by double-clicking the ENABLE column cells in all the rows for which execution was completed to disable them and then restarting the Sequence by clicking RUN.

As the Sequence runs, the table will be "indexed" upward so that the currently-in process row (highlighted) which corresponds to the currently-running Method/injection always appears two lines from the top of the table window. You can use the horizontal and vertical scroll bars while running to examine any row in the table, even those rows which have already been processed. The SoftButtons can be used to jump to any other screen or channel and back to the running Sequence display. You can also press the GRAPH button at bottom right to select a split-screen view of the Sequence table while it is executing and an active top Graph screen (Figure 11.3). Click on the SEQ button at bottom right of this display to switch back to the complete Sequence dialog screen.




Figure 11.3 Split-Screen Sequence Table with Graph

When a Sequence is executing, the Status Box can be placed on screen to monitor the progress of individual Methods and the entire Sequence.

While between Methods, or while holding at the start of each new Method, the input signal Status registers will be "live" and will show the actual signals being read in real time, to help the operator judge baseline level and stability, even though the signal data are not being saved during such periods because the system timer is not running.

11.4.1 Modifying a Sequence While Executing

Like all other Data Ally tables, it is possible to modify a running Sequence Table while the Sequence is executing. You can use the INSERT, DELETE, COPY, PASTE, and CUT commands to add new rows or delete existing rows, change command settings for any row, or even use the FILL functions as described above to complete programming of new rows. The only limitation to on-line changes in the Table is that only those lines which have not yet been executed can be modified. One very convenient tool for on-line Sequence changes is the ENABLE column, which can be employed to quickly disable selected groups of vials without rewriting the Table.

If the Reprocess mode is in use, you will observe each chromatogram loaded from disk using the filename(s) in the selected Method(s) in any Graph window as it is redrawn and integrated. In Browse mode, the reprocessing operations will be on a stepwise basis after starting. The first chromatogram will be displayed and integrated with the currently referenced Integration, and the results will be shown in all Graph window. When all steps are complete for the programmed reprocessing of the first chromatogram file, the next will be loaded and these steps repeated, and so on.

11.5.0 Working with Sequence Files

The Sequence program filing system operates in a manner identical to that of every other Data Ally programming screen. The Sequence screen display, comprising both the top-level dialog box and the Sequence Table, is configured as a "template" to show information about the current sequence instructions loaded into the channel being observed. The contents of the Sequence screen display can be saved as a unique "Sequence" file, with a unique filename and description, so that it can be recalled and reused at any future time.

The default Sequence file, whose name is assigned during Configuration, is always loaded automatically into the Sequence display screen at the time Data Ally is initialized. You can change any aspects of the default file within the active software, except for the features and options which

are either "fixed" or "removed" in the Configuration. These items can, of course, also be reset or restored, but the channel must be disabled and a re-Configuration performed to accomplish this.

The Files commands in the Sequence Command Menu bar (Figure 11.4) are used to load, save, clear, and delete Sequence program files.




Figure 11.4 Sequence Files Menu

11.5.1 Loading Sequence Files

To load a saved Sequence file from disk memory into the current channel, click the left mouse button on Files in the Command menu bar, and then click on Load. A directory dialog box will appear which lists the Sequence files (*.SEQ) available in the current directory path for this channel (Figure 11.5). Select the desired file by clicking on its name in the directory list, or by typing its name into the entry box. Click on OK to load the new file.




Figure 11.5 Sequence Files/Load Dialog

If the previous file has been edited and has not been saved or resaved, an information message will appear before the new file is loaded asking if the previous file should be saved. If so, click on YES and assign either the same filename or a new name as desired. If you do not wish to save the edited file, click on NO to proceed.

To load the default Sequence file programmed during Configuration, click the Load Default command in the Files submenu - the default sequence will be loaded and will overwrite the previous Sequence file in the current time base. The default file can be modified only by performing Configuration.

11.5.2 Clearing Sequence Files

A Clear function allows you to remove all entries from the Sequence Table without deleting the lines, or deleting the entire file.

To clear the current Sequence Table, click on File in the Command Menu bar and then click on Clear. Click Clear Sequence Table to clear the entire table (Figure 11.6). If the previous file has been edited and has not been saved or resaved, an information message will appear before this file is cleared asking if the previous file should be saved. If so, click on YES and assign either the same filename or a new name as desired. If you do not wish to save the edited file, click on NO to proceed.




Figure 11.6 Sequence Clear Function

11.5.3 Saving Sequence Files

To save a Sequence file using the same name, click on Files in the Command Menu bar, and then click on Save. A trap message box will appear prompting confirmation of whether the existing file with the current name should be overwritten - click OK if you wish to resave the file. The previous version will be overwritten with the current contents of the Sequence screen.

To save a modified Sequence under a new name, click on Save As in the Files menu. The Save As dialog box will appear (Figure 11.7). Click inside the filename entry field and type in the desired new filename. If you wish, click inside the DESCRIPTION field and type a description of the new file, which may help identify this file in the directory listing for future recall. When the correct name and description appear, click on OK to save the file to the current disk path. If another file already exists in that path with the same name, an error message will instruct you to change the filename before saving.




Figure 11.7 Sequence Files/Save As Dialog

It is recommended that Sequence program files saved to disk be "backed up" by copying the *.SEQ files to a floppy disk, removable hard disk, or tape and preserving the copies in a secure place. If modifications are routinely made to the library of *.SEQ files, backups should occur at routine intervals or even on a daily basis.

You can enter a text description for the current Sequence file, which will appear in the "Description" field in the Load dialog box when scanning saved files, either by directly typing a description into this field when using the Save As dialog above, or by selecting the Sequence Description command in the

Files submenu, which permits you to type a multi-line text entry into the description field which will be saved with the Sequence file.

11.5.4 Deleting Sequence Files

To delete a Sequence file from the current disk directory, use the Windows File Manager Delete command.

11.5.5 Printing Sequence Files

Use the Print Sequence command in the Files submenu to print a summary report for the current Sequence file. The Print function utilizes an Excel template file which is saved as reports\seq\seq.xls, and which can be modified like any other report template using Excel (see Section 13).

11.5.6 Automatically Starting Sequence Files

A special feature built into the Data Ally System is the ability to automatically run the default Sequence file at boot up time. You may specify which time bases' sequence files are started via the command line switches. The following is an example of a command line switch to start time base TB1 and TB1 EDIT. This command could be placed in the program group/program item command line area of your Windows Data Ally Group. If you wanted to start automatically from the DOS prompt you would need to precede the command with WIN and a space or whatever command and path to start your particular Windows program.

C:\Data Ally\Data Ally.EXE /start/1/2

12.0 Utilities

The Data Ally Utilities screen (Figure 12.1), displayed via the UTIL SoftButton in any channel, includes the Configuration functions accessible only to the System Manager, and the Exit command which terminates Data Ally operation.

None of the Utilities functions are essential to routine acquisition, processing, and reporting of chromatography data. Any function in the Utilities Command Menu can be used at any time in any channel, whether or not that channel is actively running a Method or Sequence. In many cases, the Utilities screen and its contents will not be required at all by certain system operators.



Figure 12.1 UTILITIES Screen

12.1.0 Utilities Files Commands

The Utilities screen contains a set of file-related functions accessible by clicking the Files command on the menu bar to see the Files submenu (Figure 12.2).




Figure 12.2 Utilities Files Submenu

The two commands available in the submenu allow you to load a new configuration file to reset the overall system parameters, without entering Configuration mode, or to translate existing data system raw data files from other formats into Data Ally-readable format.

12.1.1 Resetting Configuration

To change the configuration file which controls the current setup of the entire system, click on Load Configuration in the submenu - the Load Configuration dialog will appear (Figure 12.3). Select the desired Configuration file from any path or directory, and click OK, or CANCEL to abort re-configuring. Clicking OK causes the selected configuration to be initialized - the system may "time out" briefly and show the hourglass cursor while this is happening.




Figure 12.3 Load Configuration Dialog

Use of the Load Configuration command provides instant and easy access to any alternative programmed configuration. For example, if you wish to set two different configurations specifying completely different detector inputs and controlled instrumentation, representing two totally different HPLC setups, creating and saving two specific configuration files which select each of the optional system configurations allows you to switch the Data Ally software to match the analytical components and functions as needed, at any time.

12.1.2 Translate Functions

The Translate command can convert files saved using Axxiom Chromatography 700 Series data systems to Data Ally compatible formats, or can convert Data Ally-generated data files to ASCII and other output formats for use with other software or devices. Click Translate in the Files submenu to see the Translate submenu (Figure 12.4).

NOTE: Your version of Data Ally may contain other translator utilities that have been added since the user manual was published. Contact your distributor for translator information.




Figure 12.4 Translate Submenu Commands

➡ To convert Axxiom 700 Series chromatogram data files to Data Ally data file (*.pdf) format, click "Axxiom 700 to Data Ally" to see the 700 Series conversion dialog (Figure 12.5).




Figure 12.5 Axiom 700 Series File Conversion Dialog

This function will find all files in the typical 700 series data file format, with file names having the extensions *.0*, in any disk path, and will allow you to convert any or all of the files to Data Ally *.pdf format. You can also type in the name of a known 700 Series data file (path and filename) to convert single files on a one-at-a-time basis.

To convert one data file only, click inside the "File to Convert" entry box at upper left and type in any valid path/filename for a 700 Series raw data file. The name of the new converted file will appear in the box marked "Data Ally PDF Files" at right, which will have the same root filename as the existing 700 series file and will appear in the same path by default. If you wish to edit the new filename or path location, simply click on the PDF file name shown in the right hand box and edit its name and path as desired. Then click the Convert Single File button at bottom right to perform the file conversion. Both a copy of the existing 700 Series data file and the new converted file will remain on disk after each conversion event.

To convert a series of files, select the disk path and directory using the selector dialog in the center box for the 700 series files desired. All the names of the relevant data files in the selected path will appear in the box at left. Now, click to highlight each of the desired filenames in the box for conversion - their respective converted *.pdf filenames will appear in the box at right after each file is chosen. When all files are selected and highlighted, click the Convert All Files button at bottom to perform all the conversions as a batch.

When all conversions are completed, click the Quit button to quit the conversion function.

➡ To export Data Ally *.pdf data files, click the Export Data Ally To ... command in the Translate submenu. The Export dialog will appear (Figure 12.6).




Figure 12.6 File Export Dialog

This dialog allows reformatting of single *.pdf files from any disk path into ASCII format, or reconversion to native Data Ally or "binary" format. This dialog can also be used to "scan" files to evaluate their data structures and determine the range of signal values in each file.

To select a data file, click File in the dialog's command menu bar to view the Files Recall dialog (Figure 12.7). Use this dialog to choose a Data Ally chromatogram file in any disk path - when the desired file is shown, click OK, or on CANCEL to abort.




Figure 12.7 Export/Files Dialog

After clicking OK, the Files dialog will disappear and the name of the selected file will appear in the filename field of the File Export Dialog. This file will be scanned, requiring a few seconds, and when the scan is finished the number of data points, minimum and maximum signal values, and maximum signal delta value (difference between highest and lowest point values) will be computed and will appear in their respective boxes. You can scan through the file on a point by point basis by clicking the Next Point button at the bottom - each time this is clicked the next data point in the file, beginning with the first point, will be scanned individually and data describing that point will appear instantly in the fields in the upper right portion of the dialog. Click on the Do All Points button to repeat the scan of the entire file. You can click the STOP button at any time to immediately halt the scanning process.

To convert to ASCII format, click the radio button marked ASCII at lower right and then click the Do All Points button. The file will be reprocessed point by point and resaved as a new file in ASCII format with the same name in the same directory path - this file will no longer be readable as a Data Ally *.pdf binary file. However, you can reconvert to Data Ally ("binary") format by selecting an ASCII format file and choosing the radio button marked Binary, then repeating the Do All Points command to rescan and reformat the file.

12.2.0 Configure

The Configure command in Utilities accesses the Configuration Mode (Section 5). This is the only point in the Data Ally software at which Configuration Mode can be entered. Instructions for using the Configure submenu are found in the section on Configuration.

12.2.1 Exiting the Data Ally System

You can quit the Data Ally software only from the Utilities screen via the Quit command.

➡ To quit running Data Ally and remain in Windows, click on Quit in the Command Menu bar, then click on Exit to Application. A trap message dialog box will appear to verify that you wish to quit. Click on YES to quit. All active Data Ally windows and icons will be closed and any other software running in Windows will continue unaffected.

13.0 Data Ally Help Functions

Like other Windows applications, Data Ally provides a substantial amount of built-in "Help" which is instantly accessible from anywhere in the software. The yellow HELP SoftButton at the top left of each Data Ally screen can be used to display Help dialog boxes at any time.

13.1.0 Using the HELP System

➡ *To get HELP at any time while running Data Ally*, click the HELP SoftButton or click the "Help" command in the current dialog menu bar, if one is present - a HELP window will appear (Figure 13.1).




Figure 13.1 Data Ally Help Window

The overall Data Ally HELP system is organized into a series of topics, any of which can be reviewed each time you select HELP. The Help topics are indexed by both topic categories and classes and by "key words". You can move forward or backward to scan through the entire listing of Help topics until you find what you need, or you can Search through all topics for a key word or phrase describing the feature or action desired.

Key words in the Help dialogs can be marked with either a dotted or a solid underline. If you click the left mouse button with the pointer on a word having a solid underline, the Help system will "jump" to a new Help topic related to that word. Clicking on a word or phrase having a dotted underline will show a definition of that word or phrase.

All Data Ally Help screens are context sensitive in that they will appear with specific topics keyed to the current program screen or function being used. Figure 13.1 shows the initial Help display keyed to the Graph screen, which deals with explaining chromatogram manipulation and other graphical functions. You can of course select any Help topics at any time once any Help screen is being displayed in an open window.

➡ *To see a listing of main Help topics*, click the Contents button on the Help window command menu. The Contents listing window will appear (Figure 13.2) showing all the primary Help topics available.




Figure 13.2 Help Contents Index Window

➡ To search for a specific desired topic in Help, click the Search button on the Help window command bar menu. The Help Search dialog (Figure 13.3) will appear.




Figure 3.3 Help Search Dialog

The top portion of this dialog contains a selection box which pages down through all the listed Help topic headings in the Help system alphabetically. Scan down through the listing to find the desired heading and click on it to highlight - it will appear in the upper entry box. Alternatively, click inside this box and directly type in any topic description in which you are interested. To see the list of associated items with each topic, click the Show Topics button, and such a listing will appear in the lower list box in the window. You can now highlight any item in the sublisting (these are the specific items in the Help system) and then click the Go To button to see that specific Help window with a detailed explanation. Click Close if you wish to close the dialog and escape.

To step back to the last Help window/topic viewed before the currently-displayed topic, click the Back button on the top.

To view a "history" of all recent Help topic windows you have viewed, in order, click the History button - a Help History window (Figure 13.4) will appear. This window lists in descending chronological order all the recent screens and items you have viewed.




Figure 13.4 Help History Window

You can also use the ">>" and "<<" buttons to move, one screen at a time, through all the available Help topics in the Data Ally Help system. Stepping through the available screens in this manner is a good way to get an overview of all the topics which are covered.

All the other Help functions, as well as the overall organization of Data Ally Help, are standard Windows functions and operate the same way as for Windows applications of all types in general. You can get Help which describes how to use the Data Ally Help system by selecting the Help command on the Command Menu bar of the Help dialog box, and then clicking on How to Use Help in the submenu. Clicking the "Always on Top" item in this submenu to checkmark it will maintain the current Help window on top of any other Data Ally display window to make using Help easier.

Any Help dialog box can be moved or resized in any position on the display. You may wish to change the location of the Help box in order to see something it otherwise obscures underneath. Any Help dialog can be "minimized" to an icon by clicking the down arrow/minimize button at the upper right corner of the Help dialog, or by using Minimize in the Application Control Menu at upper left. After iconizing, the Help window will remain on screen as an icon, and can be redisplayed at any time by double-clicking on that icon.

➡ To close Help dialog boxes, click on Files/Exit in the Command Bar menu or use the Application Control Menu Close command. Whenever you re-open Help, the Help window will retain its last size and location.

Note that most of the topics and entries in the Data Ally Help system are indexed with a section number reference for this Operators Manual, to assist you in quickly finding the most appropriate explanation in the manual for any item you are viewing in Help.

13.2.0 Marking Help Topics

You can use the Help "bookmark" function to indicate specific reference points within the Help system to which you want to return often. After a "bookmark" has been inserted in Help, you can return instantly to that place by using the Bookmark menu in the Help window.

➤ *To insert a "bookmark" at a Help topic*, first show the desired topic in the Help window, and then click on Define in the Bookmark menu. The Bookmark Define dialog (Figure 13.5) will appear. Click on OK to use the title of the Help topic as the bookmark name, or click inside the name field and type a new bookmark name, then click on OK. This will now cause the new bookmark name to appear on the Help Bookmark menu.




Figure 13.5 Bookmark Define Dialog

➤ *To move to a "bookmarked" topic*, click on Bookmark in the Help window Command Menu bar, and select the name of the bookmark to be viewed. You can also press the underlined number assigned to the desired bookmark instead of clicking on its name. The bookmarked topic will appear.

➤ *To remove a bookmark*, click on Bookmark to see the submenu, then click on Define. Click on the name of the bookmark to be removed, and click on Delete. Click on OK to execute the deletion.

13.3.0 Annotating Help Topics

It is possible to add your own customized annotations and notes to any of the Data Ally Help topics. Annotations can be up to 2000 characters long and can include any type of characters.

➤ *To add an annotation to a Help topic*, with that topic in the current Help window, click on Edit and then on Annotate in the submenu - the Annotation dialog (Figure 13.6) will appear. Type in the annotation text you wish to add. The text will automatically "wrap" from line to line as you type. When you have finished the annotation, click on OK to complete it.




Figure 13.6 Help Annotation Dialog

These steps can be repeated to add more text after existing annotation(s). However, only the sum of all annotations made can be erased (below).

➡ To see a Help annotation, view the Help topic with the annotation, then click on the "paper clip" symbol. Alternatively, you can press the TAB key to highlight the paper clip symbol, and then press ENTER. The annotation will appear.

Click on CANCEL to clear the annotation.

➡ To erase an annotation, click the paper clip symbol on the current Help topic to see the annotation. Click on the DELETE button to remove the annotation. Use of DELETE will remove all the contents of the annotation window, even if the window includes several "annotations" added at different times.

14.0 Reporting with Microsoft Excel

Data Ally utilizes Microsoft Excel, v. 5.0 or higher, to format, prepare, display, and print any type of customized report desired. **Italicized** words in this section are intended to be consistent with terms referred to in Excel documentation. Information is passed from the Data Ally software to Excel via Direct Data Exchange (DDE) and Visual Basic for Applications (VBA) functions. Desired information fields and records are placed into positions on a report *worksheet* contained in an Excel .XLS *workbook* file, called a template by Data Ally, residing in a subdirectory under ..\REPORTS. Any number of template .XLS files can be reviewed and edited for use in reporting with Data Ally.

Each time the Data Ally software is initialized, a copy of Excel is also automatically loaded, so that Excel is available for creating any run, calibration, or program summary reports whenever these are desired. Excel generates each report from the information in Data Ally, and actually manages the display and/or printing of the report, supplying such functions as *print preview* and view expansion. The Data Ally data system should not be running, however, when modifying or creating new templates.

14.1.0 Using Excel to Review and Edit Report Templates

Microsoft Excel is generally loaded from its own directory on the computer's hard drive. It can be accessed while the Data Ally software is running via Windows Program Manager directly or through File Manager.

14.1.1 Loading Excel

To load Excel from Program Manager, activate the Program Manager window and select the Excel program icon from its own application window or from an "Applications" icon/group, depending upon how it has been installed. Double click with the left mouse button on the Excel icon to load Excel. The Excel spreadsheet will appear.

To load Excel from File Manager, find the Excel directory in the master directory listing and click on its directory icon to open the files. Locate the "excel.exe" file and click on its icon to execute it and load Excel.

14.1.2 Opening An Existing Report Template File

To open a Data Ally report template, .xls file, for any report type, left-mouse click on the Files command in the main Excel worksheet screen and click on Open. The Open dialog will appear (Figure 16-1).

Double-click on the "C:\\" drive icon, and then double-click on the ..\Data Ally directory icon to open it and expose its subdirectories. Then double-click on the ..\REPORTS directory icon to open it and expose its subdirectories. A listing of files contained in the ..\Data Ally\REPORTS directory will appear. This directory will contain files like pyramid.xls and pyrrpt.xll which communicate with the data system and should not be touched.

UNDER NO CIRCUMSTANCES SHOULD ANY OF THE FILES PRESENT IN THIS ..\REPORTS DIRECTORY BE CHANGED OR DELETED! DOING SO CAN CAUSE FAILURE OF Data Ally REPORTING MODES.

Double-click on the ..\REPORTS\ directory icon to open it and expose its subdirectories:

- ..\REPORTS\CAL
- ..\REPORTS\INT
- ..\REPORTS\MET
- ..\REPORTS\SEQ
- ..\REPORTS\PDF
- ..\REPORTS\SUM.

Now, double-click on one of these directories to open it and expose its template files. To load one of these templates, select the specific template file, then PRESS AND HOLD DOWN THE SHIFT KEY and

simultaneously click on OK (or double-click on the template filename.) The template file selected will appear showing the Excel worksheet (Figure 14.2). Note that it is important to hold the SHIFT key down while opening template files so that the Auto_Open macro, which normally generates the report during the running of Data Ally, will not execute now and thereby corrupt the template you expect to modify. If you accidentally open a template without using the shift key, simply close the template answering **no** to the question: “Do you want to save it?”. Then open the template as described above. Any of the .XLS template files in these directories can be recalled, viewed, and edited.




Figure 14.1 Excel Files/Open Dialog Box




Figure 14.2 Examining an Excel Worksheet Report Template

14.1.3 Elements of a Report Template File

An Excel report template file for Data Ally usually contains several elements:

- A common header section with information about the user, the analytical system, the method, etc.
- One or more sections containing information about the files used (names and descriptions, date/time of changes, etc.), the operator (name, position, etc.), or other data concerning a chromatogram run.
- A “table” section in which run data (peak numbers, names, areas, heights, concentrations, other parameters) is summarized in a tabular format, with a top title line and a sum line at its bottom giving totals for certain parameters in the table.

- An *embedded* graphic(s) or chromatogram(s) occupying a specified range of cells in the worksheet. The image of one or more chromatograms can be directly imported from Data Ally into Excel via DDE, also certain other information (such as peak data) could be “graphed” as a bar chart, pie chart, or x-y plot in a saved report with the aid of user programming.
- Individual or groups of “calculated fields” which utilize various data elements from Data Ally for performing some calculation or manipulation in the report.

A single report can include any combination of these elements in any desired order.

When an existing report template file is opened in Excel, any fixed text it contains will be shown in the currently selected font and style, and in the correct position. Fields which will contain data captured from Data Ally via DDE links are shown with a proprietary notation that contains a leading exclamation point indicating that a DDE reference will be created for the field in that position. These DDE fields will be activated when they are listed in a defined name, *Dde_Cells* (See section 14.2.6).

You can use the scroll bars at the right side and bottom of the Excel worksheet, or the PageDown/PageUp keys, or the up/down/left/right arrow keys to move through the worksheet and to advance to subsequent pages, if more than one page is formatted.

14.2.0 Excel Basics

PLEASE REFER TO YOUR MICROSOFT EXCEL USERS MANUAL AND DOCUMENTATION FOR A MORE DETAILED EXPLANATION OF EXCEL’S FEATURES AND CAPABILITIES. THE FOLLOWING SECTION IS A BRIEF INTRODUCTION TO THE BASIC ELEMENTS AND FUNCTIONS IN EXCEL NEEDED IN FORMATTING REPORTS WITH Data Ally.

The Excel worksheet consists of a series of numbered rows and columns making up a large number of individual *cells*, which can be assigned values or formulas and can be edited or moved. A complete set of tools for formatting the appearance of the worksheet (*fonts, borders, shading*) and for performing graphical and mathematical functions is available. These tools are usually accessed via the *tool bars* located above the top of the worksheet. A *formula bar* appears immediately below the toolbars showing the value entered in the currently selected cell (*text or formula*). To the far left of the formula bar is the *name box* which contains the currently-selected worksheet cell reference number in its left corner (e.g. R10C2 for Row 10, Column 2). If the current selection has a *name* associated with it, the name will show up in the name box instead of the cell reference. *Defined names* are employed for a number of Data Ally reporting options. Defined names can be used to reference a cell, a range of cells or even a Visual Basic procedure. Open the *name box* by clicking the arrow to the right of it to see the names that correspond to cell references (See section 14.2.4). To see all the names, it is necessary to choose open the *Define Name Dialog Box* (See section 14.2.5).

The Excel worksheet has a “what you see is what you get” display format, so that the actual appearance of text, graphics, and other elements in the worksheet window reflects the appearance of these items when printed. A *print preview* function is available in the File submenu to allow examination of the final format of reports without having to print them.

When the *mouse pointer* is inside an Excel worksheet, it looks like a “white cross”. It can be used to select any cell or cell range in the worksheet. You can move the mouse pointer over any cell and click the left mouse button to place the *insertion point* in the cell which permits editing its contents. A range or series of cells can be selected by clicking the left mouse and holding it down as the pointer is dragged across a series of rows and/or columns to highlight a range. Excel draws a solid border line around the selected cell range to make it visible. This temporary border will not be visible on the printed report. Menu item **Format-Borders-And-Shading** is used to create borders and patterns that will print.

14.2.1 Changing Column Widths and Row Heights

The width of columns and height of rows can be changed in two ways:

- Move the mouse pointer in the worksheet over the top column “button” (with the column number) or onto the left-edge row “button” (with the row number) which is to be resized. When you move the pointer over either edge of the column or row button, the pointer will change to a split horizontal (column) or split vertical (row) arrow. Hold down the left mouse button and drag the edge of the button in either direction to make the button (and the attached row or column) larger or smaller. The entire column or row border will be redrawn to follow the button dragging. When the desired size is obtained, release the mouse button and the worksheet will be redrawn.
- Select the desired column or row or range with the mouse pointer. Click on **Format** in the Excel Command Menu bar to see the Format submenu. Click on either Row Height or Column Width to bring up the dialog box for changing that parameter. Type in the new desired height or width, or specify the “standard” value, and click on OK. The selected range will be resized.

14.2.2 Changing Contents of a Single Cell

In Data Ally report template .XLS files, each cell may contain either no information, a value of fixed alphabetic or numeric data of any type, an Excel *formula*, a *function* call or a Data Ally DDE field.

Once the *insertion point* is set in a cell, the following changes can be made to it:

- Text or numbers can be directly typed into the report cell.
- An Excel formula can be directly typed into a report cell. Formulas must always begin with the “=” character. Read your Excel User’s Guide and Excel Function Reference to find out all you need to know about available formulas, functions and their use.
- A proprietary Data Ally DDE field can be directly entered into a report cell. Data Ally DDE references all contain an exclamation point as the first character. The name is taken from the list of available DDE fields as listed in Section 14.4. It is also required that each DDE field be individually listed in a Dde_Cells defined name as described below.
- You can clear any existing or previous DDE field entry by selecting the cell where that field appears, clicking on the right button to see the Edit pop-up, and clicking on *Clear Contents*
- You can directly type fixed text characters into any cell(s) adjacent to a DDE field cell as labels or explanations for the DDE field parameter. Cells containing descriptive text of this nature can be moved or copied from other locations in the existing worksheet using the Cut or Copy commands in Excel as described in section 14.2.3.

14.2.3 Moving Data From One Cell To Another

- *Cut* or move the value, formula or DDE field from an already-existing cell in the current worksheet by selecting that cell, clicking the right mouse button to see the Edit pop-up menu, and clicking on Cut. Move the mouse pointer to the desired new cell position, click on that cell to highlight it, and click the right button for the Edit pop-up menu. Then select Paste to place the contents in the new cell.
- *Copy* the value, formula or DDE field from an already-existing cell in the current worksheet by selecting that cell, clicking the right mouse button to see the Edit pop-up menu, and clicking on Copy. Move the pointer to the new cell position, click on that cell to highlight it, and click the right button for the Edit pop-up menu. Then select Paste to place the contents in the new cell.

14.2.4 Using The Name Box

- The *name box* appears in the top left corner of the top screen message line and has a down arrow to the right of it. When a cell is selected, the name box usually displays row/column description for the cell.

When a cell has a name associated with it, however, the name shows up in the name box instead. The down arrow to the right opens the name box. Then selecting a particular name causes the cells named to be highlighted. This is very useful for verifying that defined names have the correct syntax and thus point to all and only the desired cells. The name box can also be used to define names. However, since the Dde_Cells name described below requires cells to be listed with comma separators, it is necessary to define them using the menu bar to open the *Define Name Dialog Box* instead.

14.2.5 Using the Define Name Dialog Box

- Make a pre-selection and open the Define Name Dialog Box by clicking menu items **Insert** then **Names** then **Define**. At this point, one can select a predefined name or enter a new name to be defined in the box at the top. The RefersTo: box at the bottom of the dialog contains the references for the name. In the case of a new name, it contains the reference to the current selection. See sections below regarding the appropriate methods for pre-selection of cells in each case. If the contents of the RefersTo: box, at the bottom, is not correct, it can be manually edited. As with all formulas, the first character must be the character =. When the RefersTo: box is correct, click the Add button and OK to exit.

14.2.5.1 Defining A Named Range

- Preselect the range before opening the dialog. To preselect a range, click on the top left cell and press and hold the left mouse button while dragging the mouse to highlight the desired range. Open the Define Name Dialog Box and the range as necessary (See section 14.2.5). Examples of defined names that require ranges are: Table_Origin, Summary_Origin (See section 14.2.6).

14.2.5.2 Defining A Named Row Reference

- Preselect the row(s) before opening the dialog. To preselect a row or rows, click on the *row heading* and draw the mouse down to select the row(s) to highlight the desired rows. Open the Define Name Dialog Box and edit the row reference if necessary (See section: 14.2.5) Examples of defined names that require row references are Print_Titles and Summary_Bound (See section 14.2.6).

14.2.5.3 Defining A Named Multiple Selection Of Cells

Make a *multiple selection* of cells before opening the dialog. To make a multiple selection of cells, left mouse click on the first cell in the list, press and hold down the control key, clicking on each successive cell in the list. Open the Define Name Dialog Box and edit the list of cells if necessary (See section 14.2.5). An example of a Defined Name that requires a list of cells is Dde_Cells (See section 14.2.6). Though the process is more arduous, the RefersTo: box at the bottom could be edited by hand without preselecting cells. A sample entry might be “=R5C2,R6C2,R7C2,R8C4”. Colon characters are not permitted in multiple selections for use with Dde_Cells.

14.2.6 Inspecting Defined Names Used in Report Templates

- **Dde_Cells** Every cell in the report template containing a DDE name **must be separately listed** in a Dde_Cells separated by a comma with the exception of cells which are contained within the ranges of Table_Origin or within Summary_Origin. **This is essential for proper operation of DDE** when reports are produced. Regardless of where or how many times any DDE field is placed, all its cell references MUST be included in Dde_Cells. It makes no difference in what order these fields are listed but each DDE field must be separated from each other with a comma. Since selection of cells by dragging the cursor results in ranges (includes a colon), the procedure in section 14.2.5.3 must be followed instead.
- **Dde_Cells1 through Dde_Cells9** Since defined names have room for a maximum number of fifteen or so cell references, additional defined names may be used to list all the DDE fields. The defined names reserved for this purpose are Dde_Cells1 through Dde_Cells9. It makes no difference which of the Dde_Cells names are used for which DDE field. It is convenient if all the DDE fields listed in a name are grouped together. If the particular Dde_Cells name is clicked on in the name box, the spreadsheet will

scroll to the first DDE field and all of the DDE fields listed will be highlighted for inspection. These are also comma delimited multiple selections of cells.

- **Print_Titles** This name is used to define headings for tables which will print at the top of each page. This must be a row(s) reference which is created as described in section 14.2.5.2.
- **Summary_Bound** This name, when it is used, describes the row below Summary_Origin. It is only used when Summary_Origin defines multiple rows and is then only required when using the special Visual Basic Routines for manipulating tuples with multiple rows. This must be a row(s) reference which is created as described in section 14.2.5.2.
- **Summary_Origin** This name is used to describe a range on a Summary Report that contains summary data. This must be a Named Range which is created as described in section 14.2.5.1.
- **Table_Origin** This name is used to describe a table containing DDE names on a report template. This must be a Named Range which is created as described in section 14.2.5.1. When more DDE fields are required than will fit on the page, the report can be adjusted for landscape mode to add the extra fields or the DDE fields can be divided among multiple tables. Additional tables can be included on the report by defining their rectangles using of DDE name Table_Origin1, Table_Origin2 under Table_Origin.

14.3.0 Modifying Existing Report Template Styles

It is easy to change Excel-based report template files to create exactly the report formats desired using information that will be sent from Data Ally. To get started, the best approach is usually to change an already-existing working report format and then resave it as a new template file. Modifying and saving templates should be done without Data Ally running and the template opened with the shift key (See Section 14.1.2).

14.3.1 Modifying A Template Paper Size

Data Ally templates are designed for three pauper sizes: A4, Letter and Legal. They are designed so that the template can be converted from one paper-size to another with a modicum of effort.

- **Setting Default PaperSize.** The program files: Calibration, Integration, Method and Sequence each have a default template that is used whenever the **Files-Print** menu item is used. This system default template can be modified in the Data Ally.ini file. Under the [ini.rpt] section, define: PaperSize=Letter or PaperSize=Legal or PaperSize=A4. If this definition is missing, the Data Ally data system assumes PaperSize=A4.
- **Required Program Filenames.** The following table summarizes the program templates that will be opened depending upon the PaperSize option.

Path:	.\reports\cal	.\reports\met	.\reports\int	.\reports\seq
In Europe:				
PaperSize=A4	a4cal.xlsa4met.xls	a4int.xls	a4seq.xls	
In the U.S.A.				
PaperSize=Let	letcal.xlsletmet.xls	letint.xls	letseq.xls	
PaperSize=Legal	lgcal.xls	lgmet.xls	lgint.xlslgseq.xls	

- **Hints for Designing A Template For More Than One Paper Size.** The *Summary Info Dialog Box* can be opened with menu item **File-Summary Info....** For program files, the comment box contains information about the Top Margin and the Bottom Margin. This comment should be consistent with the actual Top and Bottom Margins set in the *Page Setup Dialog Box* opened by **File_Page Setup....** In the Page Tab, select the paper size; in the Margins Tab, select the Top and Bottom margins. These are the only parameters that need to be changed to convert from one Paper size to another once you figure out what the correct Top and Bottom margins are for the different page sizes. If you want to be able to create a

report from scratch that converts easily this way, a landscape report should be created based upon letter paper because it has the shortest length. A portrait report should be created based upon A4 paper because it has the narrowest width.

14.3.2 Placing And Formatting Report Headers

A report header is a set of information and/or graphics that prints once at the beginning of a report. In most cases, the report header will include some identification of the type of report, the institution or laboratory, the date and time of the report, the operator name, the data file(s) used for the report, the analytical instrumentation or system employed for the run, sample-specific information, and a description of the procedural files used to generate the run and report. Generally, report header information is global in nature.

When DDE fields are placed on the report, they must be preceded by an exclamation point and the cell containing the DDE reference must be listed in one of the Dde_Cells defined names as described in Section 14.3.2. DDE names available for reports are summarized in section 14.4.

14.3.3 Placing And Formatting Page Headers Using Print_Titles

A page header is a set of information that prints at the top of each page. If you wish some DDE fields to be included in the header to be repeated on each page of the report they must be included in a Dde_Cells name. Also, you must include any rows containing desired header fields, labels, and graphics in the range of rows that define the name *Print_Titles*. This name assignment will cause those rows to be printed at the top of each page, along with any table column titles which may be present in the named range of rows. Usually, the column headings for the table will be the last row or two of the *Print_Titles* range. Indeed, *Print_Titles* usually only contains column headings.

14.3.4 Placing and Formatting Data/Result Tables

Tables in Data Ally reports are of three fundamental types:

Peak information tables based on the set of “enabled” peaks from the current Integration Peak Table, with selected result and/or calibration parameters listed for each peak in a series.

Summary peak tabulation tables based on “enabled” peaks found in a Sequence of runs, with selected parameters listed for the same peak over a series of runs or injections.

Summary run tabulation tables, with selected parameters listed for each chromatogram filename over a series or group of runs. For these tables,

The last two types of tables can only be placed into “summary reports” executed as part of Sequences. Only one table can be placed into any single report.

14.3.4.1 Defining The Table Position

The position of each table in any report is defined by the range of the first data row being named “Table_Origin”. Any table will be printed only if a cell with this name appears on the worksheet. In an existing report template, select “Table_Origin” in the name box (See section 14.2.4) to see the present location of the table origin cell, if any. In order to define a new table position, use the Define Name Dialog Box to define Table_Origin (See section 14.2.4.1) as a rectangular region. You can click inside the location box in the Define Names dialog and edit the current cell reference, or delete it and then use the cross pointer to select the new desired cell for the table origin. Alternatively, the named cell can be moved to a different location by using the Cut/Paste function.

14.3.4.2 Tables With Folded Rows Or Tuples

For most Data Ally Templates, all the printed data for a data row (called a tuple) fits on one printed row on the report. The first tuple in this case is simply the first printed row of the table and Table_Origin is a range one

row high. When all the data for one tuple will not fit on one row, the data may be folded onto several rows and Table_Origin may be many rows high. That is the case for the Summary and System Suitability reports. In that case, the bottom row of Table_Origin usually does not contain any data and is reduced in height to form a tuple separator. This allows the tuples to be visibly separate on the report instead of blending in with all of the rest of the rows. Tables with tuples of more than one row have a restriction that all rows in a given column must be similar. That is, a column defined in Table_Origin may contain only DDE fields or text labels, but not both.

14.3.4.3 Defining Table Contents/Fields

Once the position of the table is set by Table_Origin, the actual contents of the table must be defined. To do this, individual fields containing the information reference formulas must be entered into their positions in the first tuple of the table.

Type the name of each table field representing the items to be included in the table in a set of consecutive columns within Table_Origin, in the same format as for the DDE fields described above. The table can include as many field entries as desired. To make all the table data fit on one row, if the column widths and number of field columns exceed the total width available on a single page, you can either select a page width appropriate for printing the table on more than one page, or change column widths used to decrease the total table width, or change from portrait to landscape orientation for printing the page. If you would rather have multiple report rows per data row (or tuple), you must follow the restrictions described in section 14.3.3.2 for folded rows.

Section 14.4 contains the list of fields available for a peak information table (from a single run).

Note that column width, which can be set for each table column, will effectively determine how many digits or characters can be shown or printed for each field. You can also format numerical results in the table using the Excel Format/Number command, and align text or numerical column values in each column using the Format/Align commands.

To insert one or more blank columns between any other columns in a table, simply leave them blank.

14.4.0 Data Ally DDE Names For Use On Reports

The easiest way to create a user specific report is to modify an existing report template that starts out working perfectly. Test this original template yourself with the data system before you make any modifications or you might end up starting all over. Make a copy of this template and step-wise modify and test it. Make small modifications, saving and testing each of them by running the data system, before you change more things. For instance, it is reasonable to add many DDE fields to a report header, edit the appropriate Dde_Cells for them and then test the whole lot. It is not reasonable to also start changing column widths and modifying table entries at the same time and expect it all to work the first time. Don't forget that you must not have the data system running when you step-wise modify a template. Also don't save a template if you forgot to hold down the Shift key when you opened it.

The following sections list the proprietary DDE fields available for placement on reports. The column labeled as Templates Using This Field include some of the report Templates that now use the field. It may be helpful to you to examine those templates and notice how the DDE field is entered and how the cell is included in a Dde_Cells name.

14.4.1 Current Context Sensitive DDE Names

DDE names specific to the current context: Method, Integration, Sequence or Calibration

<u>Proprietary DDE Field</u>	<u>Range</u>	<u>Templates Using This Field</u>
!hdr.Description		Cal, Int, Met, Seq
!hdr.Filename		Cal, Int, Met, Seq, Full
!hdr.Type		
!hdr.Version		Cal, Int, Met, Seq
!hdr.stamp.Date		
!hdr.stamp.Time		
!signal.Cols		
!signal.Description		
!signal.Name		Met, Seq
!sum.pdf.Description		
!sum.pdf.Filename		
!sys.Description		
!sys.Name		Full
!sys.Number		Cal, Int, Met, Seq, Full
!sys.SerialNumber		
!sys.box.Id		
!sys.box.slave.Enable		
!sys.box.slave.Id		

14.4.2 Calibration DDE Names

For Report Heading: DDE names referring to fields stored in Calibration files (.cal)

<u>Proprietary DDE Field</u>	<u>Range</u>	<u>Templates Using This Field</u>
!cal.int.Filename		
!cal.is.Peaks:1	(to 5 is)	
!cal.levels.Cols		
!cal.peaks.Rows		Cal
!cal.points.Rows		Cal

For Tables: DDE names referring to fields stored in Calibration file tables (.cal)

!cal.peaks.enable.Row::<"*,*>	Enable	
!cal.peaks.enable.Row::<#,#>	Line No.	
!cal.peaks.table.Basis		Cal
!cal.peaks.table.cal.Mode		Cal
!cal.peaks.table.Conc:1	(to 16 levels)	Cal
!cal.peaks.table.Fit		Cal
!cal.peaks.table.is.Number		Cal
!cal.peaks.table.Name		Cal
!cal.peaks.table.RT		Cal
!cal.peaks.table.Unit		Cal
!cal.points.enable.Row::<"*,*>	Enable	Cal
!cal.points.enable.Row::<#,#>	Line No.	Cal
!cal.points.table.cal.Level		Cal
!cal.points.table.correction.Dilution		Cal
!cal.points.table.correction.InjectionVolume		Cal
!cal.points.table.correction.Weight		Cal
!cal.points.table.int.Filename		Cal
!cal.points.table.is.Area		Cal

!cal.points.table.is.Conc		Cal
!cal.points.table.is.Height		Cal
!cal.points.table.met.Filename		Cal
!cal.points.table.pdf.Filename		Cal
!cal.points.table.process.Date		Cal
!cal.points.table.process.Time		Cal
!cal.points.table.results.AdjConc		
!cal.points.table.results.AdjResp		
!cal.points.table.results.Area		Cal
!cal.points.table.results.AreaPct		Cal
!cal.points.table.results.Coefficients:1	(to 3 coeff.)	Cal
!cal.points.table.results.Conc		Cal
!cal.points.table.results.Height		Cal
!cal.points.table.results.HeightPct		Cal
!cal.points.table.results.RF		Cal
!cal.points.table.results.RT		Cal
!cal.points.table.run.Date		Cal
!cal.points.table.run.Time		Cal
!cal.points.table.sample.Description		Cal
!cal.points.table.sample.Name		Cal
!cal.points.table.sample.Number		Cal

For Tables: Column heading or field label text for calibration file tables (.cal)

!view.cal.peaks.enable.Row.ShortTitle		
!view.cal.peaks.table.Basis.ShortTitle		Cal
!view.cal.peaks.table.cal.Mode.ShortTitle		Cal
!view.cal.peaks.table.Conc.ShortTitle:1	(to 16)	
!view.cal.peaks.table.Fit.ShortTitle		Cal
!view.cal.peaks.table.is.Number.ShortTitle		Cal
!view.cal.peaks.table.Name.ShortTitle		Cal
!view.cal.peaks.table.RT.ShortTitle		Cal
!view.cal.peaks.table.Unit.ShortTitle		Cal
!view.cal.points.enable.Row.ShortTitle		Cal
!view.cal.points.table.cal.Level.ShortTitle		Cal
!view.cal.points.table.correction.Dilution.ShortTitle		Cal
!view.cal.points.table.correction.InjectionVolume.ShortTitle		Cal
!view.cal.points.table.correction.Weight.ShortTitle		Cal
!view.cal.points.table.int.Filename.ShortTitle		
!view.cal.points.table.is.Area.ShortTitle		Cal
!view.cal.points.table.is.Conc.ShortTitle		Cal
!view.cal.points.table.is.Height.ShortTitle		Cal
!view.cal.points.table.met.Filename.ShortTitle		Cal
!view.cal.points.table.pdf.Filename.ShortTitle		
!view.cal.points.table.process.Date.ShortTitle		Cal
!view.cal.points.table.process.Time.ShortTitle		Cal
!view.cal.points.table.results.Area.ShortTitle		Cal
!view.cal.points.table.results.AreaPct.ShortTitle		Cal
!view.cal.points.table.results.Coefficients.ShortTitle:1	(to 3)	Cal
!view.cal.points.table.results.Conc.ShortTitle		Cal
!view.cal.points.table.results.Height.ShortTitle		Cal
!view.cal.points.table.results.HeightPct.ShortTitle		Cal
!view.cal.points.table.results.RF.ShortTitle		Cal
!view.cal.points.table.results.RT.ShortTitle		Cal
!view.cal.points.table.run.Date.ShortTitle		Cal

!view.cal.points.table.run.Time.ShortTitle	Cal
!view.cal.points.table.sample.Description.ShortTitle	Cal
!view.cal.points.table.sample.Name.ShortTitle	Cal
!view.cal.points.table.sample.Number.ShortTitle	Cal

14.4.3 Integration DDE Names

For Report Heading: DDE names referring to fields stored in Integration files (.int)

<u>Proprietary DDE Field</u>	<u>Range</u>	<u>Templates Using This Field</u>
!int.column.Flow		
!int.column.Length		
!int.column.ParticleDiameter		
!int.parm.BaselineDrift		Int
!int.parm.DroplineRatio		Int
!int.parm.MinArea		Int
!int.parm.MinHeight		Int
!int.parm.MinWidth		Int
!int.parm.NoiseReject		Int
!int.parm.SlopeInterval		Int
!int.parm.SlopeSensitivity		Int
!int.peaks.Rows		Int
!int.retention.Time		Int
!int.suit.Enable		
!int.unretained.Name		
!int.unretained.RT		
!int.window.Logic		Int

For Tables: DDE names referring to fields stored in Integration file tables (.int)

!int.peaks.enable.Row:<“ “*>	Enable	Int
!int.peaks.enable.Row:<#,#>	Line No.	Int
!int.peaks.table.cal.Mode		Int
!int.peaks.table.Comment		Int
!int.peaks.table.is.Number		Int
!int.peaks.table.Name		Int
!int.peaks.table.rel.Code		
!int.peaks.table.RT		Int
!int.peaks.table.suit.Code		Int
!int.peaks.table.Window		Int

For Tables: Column heading or field label text for integration file tables (.int)

!view.int.peaks.enable.Row.ShortTitle	Int
!view.int.peaks.table.cal.Mode.ShortTitle	Int
!view.int.peaks.table.Comment.ShortTitle	Int
!view.int.peaks.table.is.Number.ShortTitle	Int
!view.int.peaks.table.Name.ShortTitle	Int
!view.int.peaks.table.rel.Code.ShortTitle	
!view.int.peaks.table.RT.ShortTitle	Int
!view.int.peaks.table.suit.Code.ShortTitle	
!view.int.peaks.table.Window.ShortTitle	Int

14.4.4 Method DDE Names

For Report Heading: DDE names referring to fields stored in Method files (.met)

<u>Proprietary DDE Field</u>	<u>Range</u>	<u>Templates Using This Field</u>
!met.cal.Code:1	(to 2 signals)	Met
!met.cal.Filename:1	(to 2 signals)	Met
!met.cal.Level:1	(to 2 signals)	Met
!met.cal.rpt.Enable:1	(to 2 signals)	
!met.cal.rpt.Template:1	(to 2 signals)	
!met.correction.Dilution		Met
!met.correction.InjectionVolume		Met
!met.correction.Weight		Met
!met.events.Rows		Met
!met.int.Filename:1	(to 2 signals)	Met
!met.is.Conc:1	(to 1 is)	Met
!met.pdf.Filename:<1,1>	(to 2 sig,4 files)	Met
!met.pdf.rpt.Chromatogram:<1,1>:"Y","N">	(to 2 sig,4 files)	Met
!met.pdf.rpt.Enable:<1,1>:<"Y","N">	(to 2 sig,4 files)	Met
!met.pdf.rpt.Filename:<1,1>	(to 2 sig,4 files)	Met
!met.pdf.rpt.Print:<1,1>:<"Y","N">	(to 2 sig,4 files)	Met
!met.pdf.rpt.Show:<1,1>:<"Y","N">	(to 2 sig,4 files)	Met
!met.pdf.rpt.Template:<1,1>	(to 2 sig,4 files)	Met
!met.run.Duration		Met
!met.sample.Description		Met
!met.sample.Name		Met
!met.sample.Number		Met
!met.signal.Enable:1	(to 2 signals)	Met
!met.signal.Resolution:1	(to 2 signals)	Met

For Tables: DDE names referring to fields stored in Method file tables (.met)

!met.events.enable.Row:<" ",*>	Enable	Met
!met.events.enable.Row:<#,#>	Line No.	Met
!met.events.table.Comment		Met
!met.events.table.Time		Met
!met.events.table.det.EC.Command		Met
!met.events.table.det.EC.dc.Command:1	(to 4 ecdc)	Met
!met.events.table.det.EC.dc.guard.Command		Met
!met.events.table.det.EC.pulse.Command		Met
!met.events.table.det.EC.pulse.guard.Command		Met
!met.events.table.det.WL.Command:1	(to 2)	Met
!met.events.table.det.WL.Lambda:1	(to 2)	Met
!met.events.table.flag.in.Code		Met
!met.events.table.flag.out.Code		Met
!met.events.table.oven.Temp		Met
!met.events.table.pump.Flow		Met
!met.events.table.pump.Pct:1	(to 4 pumps)	Met
!met.events.table.zone.Code		Met

For Tables: Column heading or field label text for method file tables (.met)

!view.met.events.enable.Row.ShortTitle	Met
!view.met.events.table.Comment.ShortTitle	Met
!view.met.events.table.det.EC.Command.ShortTitle	Met
!view.met.events.table.det.EC.dc.Command.ShortTitle:1 (to 4)	Met
!view.met.events.table.det.EC.dc.guard.Command.ShortTitle	Met
!view.met.events.table.det.EC.pulse.Command.ShortTitle	Met
!view.met.events.table.det.EC.pulse.guard.Command.ShortTitle	Met
!view.met.events.table.det.WL.Command.SuperTitle	Met

!view.met.events.table.det.WL.Lambda	
!view.met.events.table.flag.in.Code.ShortTitle	Met
!view.met.events.table.flag.out.Code.ShortTitle	Met
!view.met.events.table.oven.Temp.ShortTitle	
!view.met.events.table.pump.Flow.ShortTitle	Met
!view.met.events.table.pump.Pct.ShortTitle	Met
!view.met.events.table.Time.ShortTitle	Met
!view.met.events.table.zone.Code.ShortTitle	Met

14.4.5 PDF DDE Names

DDE names referring to fields stored in PDF files (.pdf)

<u>Proprietary DDE Field</u>	<u>Range</u>	<u>Templates Using This Field</u>
!pdf.cal.Code		Full
!pdf.cal.Description		Full
!pdf.cal.Filename		Full
!pdf.cal.Level		Full
!pdf.correction.Dilution		Full
!pdf.correction.InjectionVolume		Full
!pdf.correciton.Weight		Full
!pdf.int.column.Flow		
!pdf.int.column.Length		
!pdf.int.column.ParticleDiameter		
!pdf.int.Description		Full
!pdf.int.Filename		Full
!pdf.int.unretained.Name		
!pdf.int.unretained.RT		
!pdf.is.Conc:1	(to 5 is)	Full
!pdf.met.Description		Full
!pdf.met.Filename		Full
!pdf.process.Date		
!pdf.process.Time		
!pdf.results.Rows		
!pdf.run.Date		Full
!pdf.run.Time		Full
!pdf.sample.BCD		
!pdf.sample.Description		Full
!pdf.sample.Injection		Full
!pdf.sample.Name		Full
!pdf.sample.Number		Full
!pdf.sample.Vial		
!pdf.seq.Description		Full
!pdf.seq.Filename		Full
!pdf.totals.Area		Full Area%
!pdf.totals.AreaPct		Area%
!pdf.totals.Conc		Full
!pdf.totals.ConcPct		Full
!pdf.totals.Height		Full
!pdf.totals.HeightPct		Full

For Tables: DDE names referring to fields stored in PDF files (.pdf)

<u>Proprietary DDE Field</u>	<u>Range</u>	<u>Templates Using This Field</u>
!pdf.results.table.Area		Full, Area%, RFKC
!pdf.results.table.AreaPct		Area%

!pdf.results.table.Coefficients:1	(to 3 coeffs.)	RFKC
!pdf.results.table.Conc		Full, RFKC
!pdf.results.table.ConcPct		Full, RFKC
!pdf.results.table.Height		Full
!pdf.results.table.HeightPct		Full
!pdf.results.table.int.RT		
!pdf.results.table.Name		Full, Area%, RFKC
!pdf.results.table.Number		Full, Area%, RFKC
!pdf.results.table.RF		Full, RFKC
!pdf.results.table.RT		Full, Area%, RFKC
!pdf.results.table.suit.Code		
!pdf.results.table.suit.begin.Signal:1	(to 5 suitability cols)	
!pdf.results.table.suit.begin.Time:1	(to 5 suitability cols)	
!pdf.results.table.suit.end.Signal:1	(to 5 suitability cols)	
!pdf.results.table.suit.end.Time:1	(to 5 suitability cols)	
!pdf.results.table.Unit		Full, RFKC

14.4.6 Sequence DDE Names

For Report Heading: DDE names referring to fields stored in Sequence files (.cal)

<u>Proprietary DDE Field</u>	<u>Range</u>	<u>Templates Using This Field</u>
!seq.auto.Enable		
!seq.auto.FillTime		
!seq.auto.FlushCount		
!seq.auto.InjectionVolume		
!seq.auto.NeedleWash		
!seq.auto.PreflushVolume		
!seq.cal.Suppress		
!seq.int.Suppress		
!seq.process.First		Seq
!seq.process.Last		Seq
!seq.process.Mode		Seq
!seq.process.Validate::<"Y","N">		Seq
!seq.rpt.Print::<"Y","N">		Seq
!seq.rpt.Show::<"Y","N">		Seq
!seq.sequence.Rows		Seq
!seq.sum.rpt.Enable:1	(to 2 signals)	
!seq.sum.rpt.Filename:1	(to 2 signals)	
!seq.sum.rpt.Print:1	(to 2 signals)	
!seq.sum.rpt.Show:1	(to 2 signals)	
!seq.sum.rpt.Template:1	(to 2 signals)	

For Tables: DDE names referring to fields stored in Sequence file tables (.cal)

!seq.sequence.enable.Row::<"*,*>	Enable	Seq
!seq.sequence.enable.Row::<#,#>	Line No.	Seq
!seq.sequence.table.auto.FillTime		
!seq.sequence.table.auto.FlushCount		
!seq.sequence.table.auto.InjectionVolume		Seq
!seq.sequence.table.cal.Code		Seq
!seq.sequence.table.cal.Filename:1	(to 2 signals)	Seq
!seq.sequence.table.cal.Level		Seq
!seq.sequence.table.cal.rpt.Template:1	(to 2 signals)	Seq
!seq.sequence.table.Comment		Seq

!seq.sequence.table.correction.Dilution		Seq
!seq.sequence.table.correction.InjectionVolume		Seq
!seq.sequence.table.correction.Weight		Seq
!seq.sequence.table.Injections		Seq
!seq.sequence.table.int.FileName:1	(to 2 signals)	Seq
!seq.sequence.table.is.Conc:1	(to 5 is)	Seq
!seq.sequence.table.met.FileName		Seq
!seq.sequence.table.pdf.FileName:1	(to 2 signals)	Seq
!seq.sequence.table.pdf.rpt.FileName:1	(to 2 signals)	Seq
!seq.sequence.table.pdf.rpt.Template:1	(to 2 signals)	Seq
!seq.sequence.table.rpt.Chromatogram:1:<"Y","N">	(to 2)	Seq
!seq.sequence.table.sample.Description		Seq
!seq.sequence.table.sample.Name		Seq
!seq.sequence.table.sample.Number		Seq
!seq.sequence.table.sample.Vial		Seq

For Tables: Column heading or field label text for sequence file tables (.seq)

!view.seq.sequence.enable.Row.ShortTitle		Seq
!view.seq.sequence.table.auto.FillTime.ShortTitle		
!view.seq.sequence.table.auto.FlushCount.ShortTitle		
!view.seq.sequence.table.auto.InjectionVolume.ShortTitle		Seq
!view.seq.sequence.table.cal.Code.ShortTitle		Seq
!view.seq.sequence.table.cal.FileName.ShortTitle		Seq
!view.seq.sequence.table.cal.Level.ShortTitle		Seq
!view.seq.sequence.table.cal.rpt.Template.ShortTitle		Seq
!view.seq.sequence.table.Comment.ShortTitle		Seq
!view.seq.sequence.table.correction.Dilution.ShortTitle		Seq
!view.seq.sequence.table.correction.InjectionVolume.ShortTitle		Seq
!view.seq.sequence.table.correction.Weight.ShortTitle		Seq
!view.seq.sequence.table.Injections.ShortTitle		Seq
!view.seq.sequence.table.int.FileName.ShortTitle		Seq
!view.seq.sequence.table.is.Conc.ShortTitle:1	(to 5 is)	Seq
!view.seq.sequence.table.met.FileName.ShortTitle		Seq
!view.seq.sequence.table.pdf.FileName.ShortTitle		Seq
!view.seq.sequence.table.pdf.rpt.FileName.ShortTitle		Seq
!view.seq.sequence.table.pdf.rpt.Template.ShortTitle		Seq
!view.seq.sequence.table.rpt.Chromatogram.ShortTitle		
!view.seq.sequence.table.sample.Description.ShortTitle		Seq
!view.seq.sequence.table.sample.Name.ShortTitle		Seq
!view.seq.sequence.table.sample.Number.ShortTitle		Seq
!view.seq.sequence.table.sample.Vial.ShortTitle		

15 Hardware Specific Interface Information

15.1 Autosamplers

15.1.1 Marathon Autosampler by Spark-Holland

Three Marathon Autosampler models are currently controlled by Data Ally. They are the Marathon Basic™, the Marathon-XT™, and the Marathon Basic-Plus™. The Marathon Basic does not support variable injection volume. The user is advised to study the Marathon user's manual to become familiar with the operational details of the specific model in use.

15.1.1.1 Working with Spark Marathon Autosamplers

Data Ally provides control over the Marathon Autosampler via the Sequence mode of operation. In order to control the Marathon, the Marathon control must be configured in by entering the configuration mode and then adding the Marathon autosampler model to be controlled. Next while still in configuration, you must go to the Sequence screen and click on the “Autosampler button”. This will open the Autosampler setup dialog box, where default values can be entered if desired. The “Enable” checkbox must be checked for autosampler control to be enabled.




figure 15.1 Marathon Autosampler Setup Dialog

While still in configuration and on the Sequence screen, enable the columns in the Sequence Table that are desired for control of the autosampler. Typically the Vial Number column, Number of injections per vial and volume(for the XT model) are enabled. See section 5 for general configuration details. Figure 15.2 shows a typical Sequence Table. Note Figure 15.2 does not show an injection volume column because the Marathon Basic model was configured in and it does not support variable injection volume.




Figure 15.2 Typical Sequence Table For Autosampler Operation

15.1.2 Parameters which can be controlled on all Models

- Loop Fill Time *Range is 1 to 999 seconds ***
- Pre Flush Volume *Range is 50* to 999 mL ***
- Needle Wash *On or Off*

*Preflush range is minimum of 50 mL when using variable injection volume, else 30 mL.

**** Only one of these parameters should be used!** If you are using the "Dispenser controlled Injection" only fill in the Pre Flush Volume. (see discussion of modes below)

15.1.3 Parameters which apply to XT and Basic-Plus Models Only

- Variable Injection Volume *10 to 490 mL in 10 mL Steps (maximum of 50% of loop volume)*

15.1.4 Connections Between Data Ally and Marathon

1. Locate the Serial Control Cable provided and connect the DB25 connector end to the Marathon RS232 plug. Connect the other end to the autosampler control port of the Data Ally System. For a Data Ally P1/P2 Interface Module, connect it to the port labeled SER 2. For a Data Ally Personal System, connect it to the forth serial port connection. This is the top port when the connector bracket is installed vertically in the computer chassis.
2. Locate the Inject Flag cable and connect the DB15 end to the Marathon XT I/O Control port. Connect the other end to the S input on the Data Ally Interface Module or the S input of the Data Ally Personal flag connector.
3. If the Marathon does not have a DB15 it may have a terminal strip. In the case of the terminal strip type I/O control connector, use a standard Data Ally Flag cable and connect the leads to the N.O. and COMM. positions of the CONTACT section of the Injection Marker.

15.1.5 Placing the Marathon in Remote Control Mode

- You must always manually place the autosampler into the Remote (Serial Control) mode by first turning on the autosampler and then pressing the keys "F,4" The autosampler display should show "Serial Mode"
- NOTE The Autosampler should always be ready prior to booting up the Data Ally software. If this is forgotten, initialization commands will be missed. Should this happen, you may reset the Data Ally Interface Module or press the reset button located on the rear edge of the Data Ally Personal main PCB. This will cause the initialization commands to be re-transmitted to the autosampler.

15.1.6 Autosampler Sequence Control Issues

- It is necessary to configure the Data Ally Control Setup for either Marathon or Marathon-XT autosampler control. Choose the Marathon-XT for either the XT or the Basic-Plus Model. These models support variable injection volume.
- While in Data Ally Configuration, you may enable the Vial No. column so that you can perform random vial number selection during a run. Also enable the number of injections column and enable the autosampler enable check box in the top Sequence screen after clicking the Autosampler button. **NOTE: The maximum number of injections per vial in Data Ally is limited to 9. This is due to the way file names are constructed from the sequence line and injection numbers.**
- When a Sequence is started, Data Ally waits for the autosampler to communicate and start before the Method file is allowed to run. Should the autosampler send an error message to Data Ally, you will be informed that the autosampler is not responding and to check the connections. This error can be a result of invalid parameters. See the ranges listed above.

- Usually the method events table will be programmed to contain an S input flag at 0.10 minutes. This allows any pump control to begin and hold until an injection contact closure is sent from the autosampler.

15.1.7 Injection Modes

There are two mode of filling the injection loop. The first mode is the Dispenser Controlled Injection mode. This mode uses the Digital Dispenser System pump. The second mode is the Pressure Controlled Injection Mode. The latter mode requires Air pressure to operate.

15.1.8 Dispenser Controlled Injection Mode

1. Dispenser aspirates approximately 10 mL air bubble in the sample needle for more efficient removal of the previous sample.
With the injection valve in inject position, the coaxial needle pair is inserted into the sample. Air pressure, applied to the sample through the outer air needle, ensures that no air or vapor bubbles are formed during sample withdrawal.
2. Dispenser withdraws "Pre Flush Volume" from sample vial to fill the sample line with sample and to remove previous sample. (50--999 mL, programmed via the sequence screen)
3. Injection valve switches into load position; dispenser transports 2 or 3 x loop volume through the loop to fill the loop quantitatively, in case of flushed loop injections.
4. Injection valve switches into inject position. Sample loop is now part of the HPLC mobile phase flow path: sample is transported to the column. Analysis time starts.

15.1.9 Pressure Controlled Injection

Requires Air pressure and uses the Loop Fill Time parameter to determine how long the loop is filled. See the Marathon user's manual section 2 for complete details.

Should you need further assistance, contact your distributor.

15.1.9 Marathon Error Codes

The table below list the Error Code and the Marathon Function associated with them as an aid to solving any Errors reported during control the Autosampler via Data Ally.

Code	Function	Notes
107	Loop Volume	maximum of 50% of loop volume
108	First Vial	
110	Flush Time	
111	PreFlush Volume	If you are using the "Dispenser controlled Injection" only fill in the Pre Flush Volume
150	Vial Number	

15.2 Kontron 460 Autosampler

The Kontron 460 family of autosamplers can be controlled by Data Ally. Parameters which can be controlled are injection volume, needle wash, vial number selection and the number of injections per vial (1 to 9)

15.2.1 Setting up the Kontron for Remote Control

In order for the Kontron autosampler to be controlled remotely by Data Ally, the baud rate setting must be set for 9600 baud. This is accomplished by setting DIP SWITCH Position #8 to the ON position inside the autosampler. Please review the instructions in the autosampler instruction manual for details on how to set the DIP SWITCH.

Failure to set the baud rate properly will produce an error #3 on the autosampler display when control is attempted.

15.2.2 Connections Between Data Ally and Kontron 460

1. Locate the Serial Control Cable provided and connect the DIN connector end to the Kontron RS232 plug. Connect the other end to the autosampler control port of the Data Ally System. For a Data Ally P1/P2 Interface Module, connect it to the port labeled SER 2. For a Data Ally Personal System, connect it to the forth serial port connection. This is the top port when the connector bracket is installed vertically in the computer chassis.
2. Locate a Flag cable and attach a proper connector to the other end to connect to the connect the DB9 RELAY CONTROL port located on the autosampler. Use the REMOTE CONTINUE contact positions #1 & #6. Connect the other end to the S input on the Data Ally Interface Module or the S input of the Data Ally Personal flag connector.

15.2.3 Working with Kontron 460 Series Autosamplers

Data Ally provides control over the Kontron Autosampler via the Sequence mode of operation. In order to control the Kontron, the Kontron control must be configured in by entering the configuration mode and then adding the Kontron autosampler model to be controlled. Next while still in configuration, you must go to the Sequence screen and click on the "Autosampler button". This will open the Autosampler setup dialog box, where default values can be entered if desired. The "Enable" checkbox must be checked for autosampler control to be enabled.




figure 15.3 Kontron Autosampler Setup Dialog

While still in configuration and on the Sequence screen, enable the columns in the Sequence Table that are desired for control of the autosampler. Typically the Vial Number column, Number of injections per vial and volume are enabled. See section 5 for general configuration details. Figure 15.4 shows a typical Sequence Table.




figure 15.4 Typical Sequence Table For Autosampler Operation

15.2.4 Kontron Autosampler Sequence Control Issues

- It is necessary to configure the Data Ally Control Setup for Kontron 460 autosampler control. Choose the Kontron 460 Autosampler Model from the ADD dialog while in Utilities/Configuration/Setup mode.
- While in Data Ally Configuration, you may enable the Vial No. column so that you can perform random vial number selection during a run. Also enable the number of injections column and enable the autosampler enable check box in the top Sequence screen after clicking the Autosampler button. **NOTE: The maximum number of injections per vial in Data Ally is limited to 9. This is due to the way file names are constructed from the sequence line and injection numbers.**
- When a Sequence is started, Data Ally waits for the autosampler to communicate and start before the Method file is allowed to run. Should the autosampler send an error message to Data Ally, you will be informed that the autosampler is not responding and to check the connections. This error can be a result of invalid parameters.
- Usually the method events table will be programmed to contain an S input flag at 0.10 minutes. This allows any pump control to begin and hold until an injection contact closure is sent from the autosampler. This mode of synchronization requires a flag cable between the Data Ally S input connection and the inject relay contact connector of the autosampler.

15.3 Detectors

15.3.1 General Instructions for connection of detectors to Data Ally

Detectors with analog outputs are connected to the PCP via the Detector Signals INPUT 1 and INPUT 2 on the rear of the unit. On Data Ally Personal the connectors are located on the rear of the Data Ally Personal card in the computer system.

Two cables with mini-phono plugs on one end and 'spade connectors' on the other end are supplied with each PCP or Data Ally Personal board.

Tip: To Connect a Detector connect the spade connectors to the + and - integrator output of the detector. This output may be marked 'Integrator' or '1 Volt'. Where possible avoid using the 10mv chart recorder output. Plug the RCA Phono plugs into INPUT 1 or INPUT 2 on the rear of the PCP or board. Repeat for an additional detector on the same time base or second signal from the same detector.

Warning: The voltage input range of Data Ally is 0 - 1 Volt. Do not connect a 10 Volt signal to the PCP or board using standard cables. If a 10 Volt output is required then contact Axxiom Chromatography or your dealer to obtain 0 - 10 Volt cables.

15.3.2 UV and Fluorescence Detector Control

Most UV and Fluorescence detectors are controlled in the same way regardless of the brand or model. The number of possible control columns may vary depending on the detector.

Tip: To Program UV Detectors, In the '**Time**' column of the Events Table enter the first time for detector control followed by the Enter key. In the 'Command' column for the specific detector enter the desired command and value.

See details on specific detectors for commands. These are:

Detector Command	Where x =	Meaning
ZERO x	0, 5, 10, 50, 100	Zero to preset mv level
RESP x	0, 1, 2	Fast, Standard, Slow
RANGE x	0, 13	

In the **UV-Vis** column enter the required Wavelength for detection for a UV detector.

Fluorescence detectors have a **Em Wvlgth** and **Ex Wvlgth** column for programming emission and excitation wavelengths. If a wavelength change is required during the run, enter the time for the next step and repeat step Repeat until the required commands are complete.

15.3.3 Controlling the Antec Decade Electrochemical Detector

Both Pulse and DC modes are supported.

- The detector is connected to **SER1** on the PCP via a multi-colored cable to the RS232 port on the rear of the unit.
- Check the Decade EPROM version for RS232 support option. (See the Decade manuals for details) Contact your Decade dealer if you do not already have the RS232 support option installed.
- The Decade Detector must be placed in the RS232 mode by pressing the **F5** key on the Decade front panel.
- The Decade may be run in DC Mode or Pulse Mode. The choice of mode is made in Data Ally System Configuration and may only be changed if Data Ally is re-started.

Since most applications require only one mode, the decision on which mode to use may be made at installation after which time it will remain the same.

The control of the detector is made through the Method Screen Events Table in which a number of dedicated columns are displayed.

In order to control the Decade detector the Decade driver must be configured in by entering the configuration mode in Data Ally. Choose either the DC or Pulse driver and then, while still in configuration mode, go to the Method screen events table and enable the columns desired to control the Decade detector. (See section 5 in the Data Ally User's manual for general configuration details.

Below are the commands for the two modes:

Global Commands for either mode

Commands

INJECTx (Valve option)
 RESETx (Resets time)
 TEMPx (Sets oven temp)
 MARK (Crashes comm!)
 AUXx (Various outputs)
 TCx (Time Constant Clamp)
 KBx (Keyboard lock)

Value

1=load, 2=inj.
 1
 00=off, else 15-60
 1* Dangerous!
 0, ?, 1111-2222
 1=on, 2=off
 1=locked, 2=not locked

Cell Commands DC mode

OFFx (Offset %)
 MVx (DC Cell Pot.)
 FILTERx (Time Const.)
 CELLx (On/Off)
 PAx (Range)
 NAx "
 UAx "

Value

10-15(1=+), 20-25(2=-)
 1000-1200(1=+), 2000-2200(2=-)
 1-6 = 0.1-5
 1=on, 2= off)
 11-19= 10pA-5nA
 21-29 = .1nA-50nA
 31-39 = 10uA-5.0uA

Cell Commands Pulse Mode

OFFx (Offset %)
 P1Tx (Pulse Time1)
 P2Tx (Pulse Time2)
 P3Tx (Pulse Time3)
 STIMEx (Sample Time)
 CELLx (On/Off)
 PAx (Range)
 NAx "
 UAx "
 P1MVx (Pulse Pot.)
 P2MVx (Pulse Pot.)
 P3MVx (Pulse Pot.)

Value

10-15(1=+), 20-25(2=-)
 10-200
 10-200
 10-200
 1-5 = 20, 40, 60, 80, 100 ms.
 1=on, 2= off)
 11-19= 10pA-5nA
 21-29 = .1nA-50nA
 31-39 = 10uA-5.0uA
 1000-1200(1=+), 2000-2200(2=-)
 1000-1200(1=+), 2000-2200(2=-)
 1000-1200(1=+), 2000-2200(2=-)

15.3.4 Controlling the ESA Electrochemical Detector

The Coulochem is connected to COM1 on the PCP via a multi-colored cable to the RS232 port on the rear of the unit. Before it can be used with Data Ally it is necessary to switch the detector to Remote Communication Mode. To Switch the Coulochem to Remote Mode

1. Press the MODE button on the front panel of the detector.
2. Press the UP arrow key 4 times until 'SYSTEM MENUS' is shown in the display.
3. Press the Enter key.
4. Press the UP arrow key 4 times until 'ENTER

REMOTE COMMUNICATIONS MODE' is shown in the display.

5. Press the Enter key.

The detector will now be in Remote Mode and the green **REMOTE** light will show.

15.3.4.1 General Notes about Coulochem II Operation with Data Ally

Coulochem II may be run in DC Mode or Pulse Mode. The choice of mode is made in Data Ally System Configuration and may only be changed if Data Ally is re-started. Since most applications

require only one mode, the decision on which mode to use may be made at installation after which time it will remain the same.

The control of the detector is made through the Method Screen Events Table in which a number of dedicated columns are displayed.

These are:

EC Cmd Sends global commands which effect the whole detector, e.g. Autozero, Cell on etc.

Guard Sends commands which effect only the guard cell, e.g. Potential setting.

Anal #1 - #4 Send commands to each of 4 possible analytical electrodes.

Note: Only one command may be input on a single time line of the Events Table for the Coulochem. If multiple events are programmed on a single line the detector may not be able to accept all of the necessary codes sent to it.

15.3.4.2 Global Coulochem Commands

Global commands are programmed through the EC Cmd column in the Events Table. Global commands may be used in both DC and Pulse detector modes.

The possible commands are:

CELL Switches Analytical Cells in or out. Signals may only be measured when cells are placed in.

Command Where x = Effect

CELLx 0 Cells Out

1 Cells In

LOCAL Returns detector to local control. Once in local mode, the detector cannot be returned to remote control except through the front panel.

Command Effect: **LOCAL** Return to local

ZERO Sends an Autozero command to detector. No other commands will be accepted by the detector during Autozero.

Command Effect: **ZERO** Sends Autozero

RDCx Sends a run command to the detector using the currently defined DC parameters. This command must be sent after all cell parameters have been set. **NOTE:** Failure to send and RDC will result in the detector not running correctly.

Command Where x = Effect

RDC0 Puts detector into run.

15.3.4.3 Analytical Cell Commands - DC Mode

Analytical Cell commands are programmed through the ANAL#1 to ANAL#4 columns in the Events Table.

The DC mode of operation is set through the system configuration. The possible commands are:

MV Sets the analytical cell potential.

Command Where x = Effect

MVx - 2000 to Sets -2000mV

+2000 (any integer between) Sets +2000mV

CU Sets the full scale current.

	Command	Where x =	Effect
CUx	0	100pA	
	1	200pA	
	2	500pA	
	3	1nA	
	4	2nA	
	5	5nA	
	6	10nA	
	7	20nA	
	8	50nA	
	9	100nA	
	10	200nA	
	11	500nA	
	12	1 μ A	
	13	2 μ A	
	14	5 μ A	
	15	10 μ A	
	16	20 μ A	
	17	50 μ A	
	18	100 μ A	

FILTER Sets the filter time constant

	Command	Where x =	Effect
FILTERx	0	0.2sec	
	1	0.5sec	
	2	1.0sec	
	3	2.0sec	
	4	5.0sec	
	5	10.0sec	

OFF Sets the channel offset in %.

	Command	Where x =	Effect
OFFx	-100 to	Sets -100%	
	+100 (and any integer between)	sets +100%	

15.3.4.5 Analytical Cell Commands - Pulse Mode

In Pulse Mode, Analytical Cell commands are programmed through the ANAL#1 column in the Events Table. The Pulse mode of operation is set through the system configuration. The Possible Commands are:

PMV Sets the analytical pulse potential E1.

Command	Where x =	Effect
---------	-----------	--------

PMVx - 2000 to Sets -2000mV
 +2000 (any integer between) Sets +2000mV

PMS Sets the analytical pulse duration T1.

Command	Where x =	Effect
---------	-----------	--------

PMSx 100 to (Sets 100msec)+1000(any integer between) Sets 1000msec

PACQ Sets the analytical pulse acquisition delay. The acquisition delay must be between 60 and the value of the T1 parameter -10. E.g. if T1 above = 200msec then PACQ may be between 60 and 190msec.

Command	Where x =	Effect
---------	-----------	--------

PACQx 60 to (Sets 60msec) T1-10(any integer between) Sets T1-10msec

CMV Sets the cleaning pulse potential E2.

Command	Where x =	Effect
---------	-----------	--------

CMVx - 2000 to Sets -2000mV
 +2000 (any integer between) Sets +2000mV

CMS Sets the cleaning pulse duration T2

Command	Where x =	Effect
---------	-----------	--------

CMSx 10 to Sets 10 msec
 +1000(any integer between) Sets 1000msec

EMV Sets the equilibration pulse potential E3

Command	Where x =	Effect
---------	-----------	--------

EMVx - 2000 to Sets -2000mV
 +2000 (any integer between) Sets +2000mV

EMS Sets the equilibration pulse duration T3.

Command	Where x =	Effect
---------	-----------	--------

EMSx 10 to Sets 10msec
 +1000(any integer between) Sets 1000msec

OFF Sets the channel offset in %.

Command	Where x =	Effect
---------	-----------	--------

OFFx -100 to Sets -100%
 +100 (and any integer between) sets +100%

15.3.4.6 Guard Cell Commands - DC Mode

Guard Cell commands are programmed through the Guard column in the Events Table.

The DC mode of operation is set through the system configuration. The possible commands are:

MV	Sets the Guard Cell potential.
Command	Where x = Effect
MVx	- 2000 to Sets -2000mV
	+2000 (any integer between) Sets +2000mV

15.3.4.7 Guard Cell Commands - Pulse Mode

Guard Cell commands are programmed through the Guard column in the Events Table.

The Pulse mode of operation is set through the system configuration. The possible commands are:

PMV	Sets the Guard Cell potential.
Command	Where x = Effect
PMVx	- 2000 to Sets -2000mV
	+2000 (any integer between) Sets +2000mV

15.4 Photo Diode Array Detectors

Photo Diode Array Detectors and Rapid Scanning Detectors that emulate a Photo Diode Array detector are controlled by Data Ally in the following manner.

1. The entire spectra scan is captured in a master PDF file.
2. One to four individual wavelengths are captured and saved in real time and are also displayed on the upper window of the PDA graph screen as shown in Figure 15.5
3. The spectra plot is displayed in real time in the lower window of the PDA graph screen
4. The **A4Purty.XLS template** is provided for Peak Purity Reporting.
5. Special RS232 Windows drivers are employed, along with a 16550 compatible RS232 port to ensure reliable data capture in Microsoft Windows.




Figure 15.5 Typical Data Ally PDA Graph Screen

15.4.1 General PDA Formulas used for Spectra, Similarity, and Purity Plots

In order to plot or graph the spectra and Peak Purity and to report the Similarity Index (0 to 1) the following formulas are employed.

15.4.1.1 Spectra Display

Spectra are displayed at any given time in the lower graph window using the following formulas:

$$S_n = S_R - S_0$$

Where: S_n = Spectrum at time n

S_R = Raw spectrum

S_0 = Spectrum at time zero

15.4.1.2 Similarity Index

The similarity index (SI) compares two spectra across the entire wavelength range of the captured file. The search is performed using the following formulas:

$$(a(\lambda_1), a(\lambda_2), \dots, a(\lambda_n))$$

Where: $a(\lambda_1)$ is the absorbance at λ_1

The representative vector is one spectrum per vector as follows:

$$\vec{S} = (a(\lambda_1), a(\lambda_2), \dots, a(\lambda_n))$$

Since there are two spectra involved, \vec{S}^1 corresponds to Spectrum 1, and \vec{S}^2 corresponds to Spectrum 2..

$$\vec{S}^1 = (a1(I1), a1(I2), \dots, a1(In))$$

$$\vec{S}^2 = (a2(I1), a2(I2), \dots, a2(In))$$

The Similarity between two spectra is calculated as follows:

SI=

$$\vec{S}^1 \times \vec{S}^2$$

$$|\vec{S}^1| |\vec{S}^2| = \cos \theta$$

Peak Purity is determined using the average of Similarity of the Apex vs. the Upslope spectra and the Apex vs. Downslope spectra. The spectra used to calculate the purity index are as shown below:

Purity =

$$\frac{SI_{[S_A-S_B]} + SI_{[S_A-S_B][S_D-S_B]}}{2}$$

Where: S_A = Apex spectrum

S_B = Background spectrum

S_C = Upslope spectrum

S_D = Downslope spectrum

The Purity index value will range from zero to 1, with values nearer to 1 being more pure.

15.4.1.3 Setting the Number of Peak Spectra

To set the number of Peak Spectra go the Integration Peak Table and enter a value indicating the desired number of Peak Spectra (choice of 3,5,7 or 9). See Figure 15.6 below.




Figure 15.6 Integration Peak Table with Peak Spectra Column

15.4.1.4 Setting the Number of Real-Time Single Wavelength Channels

The number of individual wavelength channels may be set between one and four channels by enabling the desired channels in the top level method screen. See Figure 15.7.




Figure 15.7 Top Level Method Screen for PDA with four individual channels available

15.5 LabAlliance Model 506 or Linear 206 UV/Vis Control Installation Instructions

Data Ally Software Installation

1. Install the 3-disk set as usual.
2. Verify that the wscddrv.386 file is in the Data Ally directory.

3. Add the following lines to the [386Enh] section of the Windows System.ini file
device=c:\data ally\wcscdrv.386
WCSCInBufLen=65535
WCSCOutBufLen=1024
WCSCMaxPorts=10
4. Check the RS232 port number in the following section of the Data Ally.CFG file:
[sys.com]
Base=(1016, 760, 1000, 744)
IRQ=(4,3,4,3)
Port=2

NOTE: port=2 means comm2, port=1 means comm1, etc.

NOTE: The computer should have a 16550A or compatible serial port. Use MSD diagnostic to verify the computer's RS232 ports if you are not sure.

16 Data Ally System Security & Log On Feature

16.1.0 System Log On

When you first install Data Ally and boot up for the first time, the “Administrator” is prompted to enter his name and password. Please make a note of the name and password used as it is not visible after it is entered into the system.




Figure 1 System Security Startup Greeting Dialog

16.2.0 Adding Users

Once the Administrator has entered his or her name and password, proceed to the Utilities screen and Logon as the Administrator. Then choose Edit Users. The Set Users dialog box will appear as shown in Figure 3. Click on the Add button to add a new user and set the new user’s access levels.




Figure 2 System Security menu items Edit Users & Logon.




Figure 3 Set Users Dialog

When the Edit Users Dialog appears (see Figure 4) , you may fill in the User Name, Password, and other information as appropriate. Note that the default access right codes are shown in the lower text window of the Edit Users Dialog box. If you wish to change the access rights for the user you are editing, you would do so by clicking on the Access button with the left mouse button.




Figure 4 Edit Users Dialog




Figure 5 User Access Rights Dialog

16.2.2 Setting Users Role

Choose an Access Domain by clicking on the item in the Domain List on the left. Next choose an Access category for that Domain by clicking on an item in the Roles List on the Right. (Figure 6)




Figure 6 Set Access Role Dialog

Access Roles

The User has an assigned access role for each domain. This list is shown in order of minimum access to maximum access role. Access Roles are as follows:

- **No Access** The user has no access in this domain. The Softbutton action for that domain will be disallowed.
- **Reader** The user can only view documents or files in this domain.
- **Executor** The user can specify documents in this domain.

Specific Domains

- Chromatogram
 - The user can reprocess an existing chromatogram subject to other access roles.
- Method
 - The user can execute (start) an existing method.
- Integration
 - The user can integrate (if they have appropriate chromatogram access).
- Calibrate Template
 - The user can calibrate.
- Calibrate Curve
 - The user can analyze unknowns based on the curve.
- Sequence
 - The user can execute (start) existing sequences.
- Reports
 - The user can print/show reports through excel.
- **Author**
 - The user can create new documents in this domain.
 - Specific Domains:
 - Chromatogram
 - The user can acquire new data files.
 - Calibration Templates
 - The user can create new calibration files.
 - Calibration Curve
 - The user can create a new curve and/or add calibrators to an existing curve.
 - Reports
 - The user can create .XLS files.
- **Editor**
 - The user can edit documents, subject to future auditing.
 - Specific Domains
 - Calibration Curve
 - The user can modify an existing curve, including enabling/disabling points.
- **Manager**
 - The user can delete existing documents in this domain.
 - The user can control/circumvent (future) auditing of this domain.
 - Specific Domains
 - Configuration
 - The user can manage the user list.
- **Example Users**
 - **Administrator**

- **All Domains: Manager Access Role**
- **Lab Manager**
 - **Configuration: Editor Access Role**
 - **All Other Domains: Manager Access Role**
- **Methods Developer**
 - **Configuration: No Access Role**
 - **Chromatogram: Editor Access Role**
 - **All Other Domains: Manager Access Role**
- **Lab Tech**
 - **Configuration: No Access Role**
 - **Chromatogram: Author Access Role**
 - **Calibration Curve: Editor Access Role**
 - **All Other Domains: Executor**

Two additional features have been added to the System Security and Logon Features. The new features are listed below:

1. Option to enable and disable requiring password and logon.
2. Default Button to set a particular user as the default user. What ever access rights the default user has been provided by the Administrator, those are the access rights allowed while the default user is logged on.

17 Installation and Routine Troubleshooting

If difficulties are encountered during the installation process, check the following section for guidelines to resolving the most common types of problems. If none of the suggestions below applies or helps you complete a successful installation, contact your distributor for assistance.

17.1.0 Serial Mismatch Errors

A "serial mismatch" error occurs when the error message box appears during the initial bootup period for the Data Ally software.

Serial mismatch indicates that one of the following situations has occurred:

- The software verifies that the Data Ally Data Ally Card installed in the PC has a different encoded serial number and is thus not compatible.
- The software cannot for some reason locate the Data Ally Card even though it is physically installed in the PC, so that no serial number match verification is possible.
- The PC has not found or recognized the special command **device=c:\Data Ally\vpdtd.386** which must be present in the Windows SYSTEM.INI file before Data Ally can be loaded.

If a Serial Mismatch error occurs, you must first confirm that both the Data Ally Card and Windows have been installed correctly.

The Data Ally Card installed must not be damaged in any way, must be completely and properly inserted into a 16-bit slot in the PC bus, and must have a serial number which matches the serial number printed on your Data Ally software diskettes. To ensure that the card is mounted correctly, turn the PC power OFF, open the computer and remove the Host Card from its position in the bus, recheck the serial number if you have not already done so, and replace the card in a different slot position, if one is available, taking care to firmly and fully insert the connectors on the card's bottom into the slot. If another slot is not available, re-insert the card into the same slot position used previously.

Before attempting to restart Data Ally, check the Windows SYSTEM.INI file as described in Section 17.1.4 above, if you have not already done so. You should also check the Virtual Memory specification for the Windows swap file, as described in Section 17.1.4, since this may also affect the PC's ability to load Data Ally. Once the correct Windows setup is verified, restart the computer as described above and reinitialize Windows and Data Ally.

If the Serial Mismatch repeats, it is possible there may be an incompatibility between the I/O interrupt settings and DMA channel settings used by other, non-Data Ally cards or interfaces in the PC and those required by the Data Ally Data Ally Card. Data Ally uses the following settings:

Interrupts IRQ 10 and 11 (defaults)
DMA Channels 1 and 3

Two lines in the **System.INI** file in the **[386Enh]** section set the **IRQs** for Data Ally. Their default values are:

VPTDTimerIRQ=10
VPTDCommIRQ=11

Should you need to change the IRQs, you may reposition the IRQ jumpers on the Host Adapter card or Personal card to the new IRQ setting and edit the System.Ini file as needed.

In case of a serial mismatch not resolved as above, it is recommended that all cards, interfaces, and subsystems on the PC be checked to verify that the interrupt and DMA channel numbers above are not being utilized elsewhere. Reviewing the documentation for the PC and for other installed cards and interfaces may be necessary to confirm such settings, which can also usually be determined by observing the positions of jumpers or switches physically located on the PC's mother board or other cards. In some instances it will be possible to disable the interrupt or DMA settings used by other cards or devices, if they conflict with those used by Data Ally, or to reset them differently. **YOU SHOULD NOT RESET THE DMA SETTINGS AND JUMPERS ON THE Data Ally DATA ALLY CARD.** If there appears to be a conflict here, contact your distributor for assistance.

Another possible cause of Serial Mismatch is a "speed incompatibility" between the PC and the Data Ally Data Ally Card. In certain cases, it may be possible to change the computer's AT-Buss wait states to compensate for such an effect. Contact your distributor for instructions before attempting to make any such adjustments.

17.1.2 Other Startup Error Messages

It is possible that certain error messages other than Serial Mismatch will appear when attempting to initialize the Data Ally system.

- General Protection Faults at start up are usually caused by faulty Graphic drivers or SCSI hard disk set up problems.
- For Graphic driver problems, always try using the default Microsoft VGA or SVGA driver to see if that clears up the problem.
- For **SCSI** hard disk problems check with the **SCSI** interface manufacturer to be sure that their installation instructions have been followed and that you have their latest SCSI Windows driver file installed. Usually there will be MS-DOS related setups as well for SCSI drives.

Note: In the USA you can check Microsoft Download Service at (206) 936-6735, using a modem, to search for any driver file updates. This service is free except for the cost of the phone call. You can search for categories such as printer drivers, graphic card drivers, etc.

This same service is available on the Internet at www.microsoft.com.

NOTE:: Do not assume that your computer does not suffer from these setup problems just because you have not seen a problem. You must be sure that everything is correct, or your computer will exhibit strange behavior when larger applications such as Data Ally are run.

17.1.3 Errors During Data Ally Operation

There are three types of errors which can be encountered during routine Data Ally operation. Each type of error will elicit the appearance of a message box of a certain type on screen, regardless of whether Data Ally is running in full screen mode, in a window, or "iconized".

Run Errors: A run error occurs when an unexpected condition occurs during operation. After a run error, a message box will appear stating the type of error. **Figure 17.9** shows an example of a "**Communications Error**" dialog box which occurs when an attempt is made to start a Data Ally Method when the Communications Processor connected to that time base is not powered or its cable is disconnected.

NOTE: To clear such messages, press the reset button located on the rear edge of the circuit board. This will cause the software to re-synch with the hardware.

Generally, these errors are "non-fatal" and operation will resume immediately after the message is cleared - instructions are given in the message box concerning the probable cause of the error. However, in some cases, the same message box will be repeated indefinitely, signifying a "semi-fatal" error which will prevent further normal operation. To clear this condition, exit the Data Ally software in the normal manner, or reboot the computer using the CTRL+ALT+DEL command.




Figure 17.9 Communications Error Information Message

Page Faults: A "page fault" error occurs during certain types of operations related to DOS or Windows functions, such as file logging or printing. Page Fault errors are usually "fatal" in that they cause a permanent halt to normal Data Ally operation which cannot be remedied by simply clearing the displayed error message box. If Page Fault errors persist, the DOS and Windows setup commands detailed in the above sections should be carefully reviewed to determine if all recommendations are being followed.

Another potential cause of Page Faults or General Protection Faults occurs on certain computers not compatible with driver commands, associated with Western Digital disk controllers, often embedded in the Windows SYSTEM.INI file. These commands are related to the "Use 32 Bit Disk Access" check box command in the Virtual Memory setup. If Page Fault errors recur while Data Ally is running, check the SYSTEM.INI file using Windows Notepad to determine if the two lines "device=*int13" and "device=*wdctrl" are present (Figure 17.10).




Figure 17.10 Removing Driver Lines in SYSTEM.INI File

If present, disable both lines as in Figure 17.10 by typing a semicolon (;) character at the beginning of each line. When finished, resave the modified SYSTEM.INI file before exiting.

After disabling these driver lines, use the Windows Control Panel subprogram to "uncheck" the "Use 32 Bit Disk Access" check box.

Reboot the computer after making these changes and initialize Data Ally to determine if the Page Fault errors still persist. Other possible causes of Page Faults are incorrect network software driver setup and improper CONFIG.SYS files.

General Protection Faults: These faults appear in the form of a message box when a condition occurs in Windows/DOS which prevents the system from proceeding. In the event of a General Protection error, the Data Ally software will become nonfunctional, and the CTRL+ALT+DEL command must be used to quit Data Ally. This command will usually not affect other Windows or DOS applications which may be running concurrently with Data Ally.

General Protection faults may be caused by the same factors which elicit Page Fault errors, by timing problems in Windows, or by excessive memory use by any application.

Printer Problems: These are almost always caused either by using the improper printer driver in Windows or by faulty printer drivers. Emulating other printers by choosing the wrong driver file in Windows can cause graphics to drop out. The proper driver for the brand and model of printer should always be used. The latest version of the driver should be obtained either from the printer manufacturer or Microsoft Download or the Internet. Using older drivers with bugs can cause many problems, including General Protection faults.